

Super Sani-Cloth®

CLEANING AND DISINFECTANT

WIPE



Medical Equipment Compatibility Guide

BE THE DIFFERENCE®



Super Sani-Cloth®

CLEANING AND DISINFECTANT

WIPE



The compatibility between medical equipment and the disinfectants used to clean them is crucial.

Disinfecting equipment with products not designed for their surface materials can cause surfaces to become unsightly, brittle and cracked (Environmental Stress Cracking). That can result in surfaces becoming increasingly hard to disinfect and possibly harbour harmful bacteria. Eventually this could render the equipment unusable or in need of replacement.

Unfortunately, there is no one disinfectant that is compatible with all surfaces, fabrics and equipment in your facility.

The following document provides a list of medical equipment which is compatible with Super Sani-Cloth®. The information contained in this guide is based upon and representative of manufacturers' IFUs and recommendations. It does not indicate testing conducted by PDI.

At PDI International, we are constantly working with equipment manufacturers to add more equipment to our list. If your equipment does not appear in this guide, please contact us.



Compatibility legend/key:

	Compatible
	Not compatible
	Action required, no information – contact equipment manufacturer for recommendations

Product Details:	Item Code	Wipe Count	Wipe Size	Qty Per Case
100ct Canister	PDIESC1621E	100	200mm x 200mm	6



Contents

4
4
4
4
4
5
5
6
6
6-7
7
7
7
8
9
9
9
9-10
10
10
10
11
11
11
11-12
12
12
13
13
13
13
13
13
13
14
14
14-22
23-24
24
24



Anesthesia		
GE Healthcare	Aestiva Aspire Aisys Aestiva Avance	>
Philips	Anesthetic Gas Monitor	>
Bed / Mattress / Pads / Vinyls		
Posey	Bed 8060 Bed 8070 Floor Cushions	>
Stryker	BS3 Stretcher	▶
ArjoHuntleigh	Breeze	▶
Omnova Solutions	Boltaflex PreFixx	▶
Bladder Volume Instrument		
Verathon	BladderScan BVI 3000 BladderScan BVI 9400	▶
Verathon	BladderScan BVI 6400	>
Verathon	BladderScan BVI 6100	▶

Blood Pressure Monitors / Cuffs		
GE Healthcare	Critikon	▶
GE Healthcare	Tonoport V Blood Pressure System	▶
Nihon Kohden	BP Cuffs and Telemetry	>
Panasonic	Upper Arm BP Monitor Wrist Blood Pressure Monitor	■
Spacelabs	Tri-Cuff	■
Technicuff	Blood Pressure Cuff Model UA-1745L	▶
Welch Allyn	767 Wall and Mobile Aneroids	■
Welch Allyn	Connex ProBP 3400	■
Welch Allyn	BP Cuffs	■
Blood/Fluid Warming Units		
Arizant Healthcare	Bair Ranger	▶
Arizant Healthcare	Ranger Model 245 Ranger Model 247 Bair Paws Model 850 Bair Paws Model 875	▶
Philips	InnerCool RTx Models 23-05 InnerCool RTx Models 23-06	▶



Compression System		
ArjoHuntleigh	Flowtron ACS900	▶
ArjoHuntleigh	Flowtron Excel	▶
ArjoHuntleigh	Flowtron Trio	▶
ArjoHuntleigh	Flowtron Universal	>
Covidien	Kendall SCD 700	▶
Covidien	SCD 700	▶
Covidien	SCD Express	■
Microtek Medical	Venodyne System 610 Venodyne System 710 Venodyne System 810	>

BCT / X-Ray / Radiography		
AccuVein	AV300	▶
AccuVein	AV400	▶
ACIST Medical Systems	Angiographic Contrast Injection System	▶
AGFA	DR 10s and DR 14s	▶
AGFA	DR 400 and DR 600 Systems	▶
AGFA	DX-D 10C, DX-D 10G, DX-D 20C and DX-D 20G	>
AGFA	DX-D 30C and DX-D 35C	>
AGFA	DX-D 40C, DX-D 40G, DX-D 45C, DX-D 45G, DX-D 60C and DX-D 60G	▶
AGFA	DX-D System	>
Bard	Sapiens Tip Confirmation System, Sherlock Tip 3CG, Sherlock Tip Location System	▶
Carestream	DRX-1 System Detector 7050, 7066, 8820, 8821	>
GE Healthcare	BrightSpeed 16 BrightSpeed 4 BrightSpeed 8 BrightSpeed Elite	>
GE Healthcare	Discovery PET/CT 590 RT	>
GE Healthcare	Discovery PET/CT 600 Discovery PET/CT 690 Infinia Precision 500D	>
GE Healthcare	Optima CT580 (RT and W)	>
GE Healthcare	Panda iRes Warmer	■



Defibrillator		
Philips	HeartStart FR3, HeartStart MRx	▶
Philips	HeartStart M5066A	>
Physio-Control	Lifepak 12, Lifepak 15, Lifepak CR Plus, Lifepak Express	▶
Welch Allyn	AED 10	>
Zoll Medical	AED Plus, AED Pro, AuthPulse Model 100, E Series, M Series, R Series	▶
Zoll Medical	Propaq MD, X Series	■
Diagnostics		
Cepheid	GeneXpert	>

ECG / Cardiograph / Telemetry / Cables / Leads		
GE Healthcare	ApexPro Telemetry	■
GE Healthcare	ECG Cables & Guidewires Multilink Cables (E9008n), MAC 1600ECG Analysis System, MAC 3500, MAC 5500, MAC 5500HD, MAC 800 Resting ECD Analysis System	■
GE Healthcare	Mac 500 EKG Machine	▶
Masimo	LNCS Sensors	>
Nihon Kohden	ZM-520 Transmitters, ZM-530 Transmitters	▶
Philips	ECG Cables Leads	>
Philips	ECG Lumb Clamp Electrodes 40494E, PageWriter Xli, Signal-Averaged ECG, StressVue Exercise Stress Testing System	>
Philips	PageWriter 10/10i, Agilent M2662A, PageWriter 100 Cardiograph M1772A, PageWriter 200&300pi M1771A/1770A, PageWriter TC50, PageWriter TC70, PageWriter Touch, PageWriter Trim I, PageWriter Trim II, PageWriter Trim III, PageWriter Trim Rx	▶
Philips	PageWriter TC	▶



ECG / Cardiograph / Telemetry / Cables / Leads		
Spacelabs	Cables and Lead Wires, Cables and Transmitters, Sensors and Analyzers	■
Welch Allyn	CardioPerfect Workstation	>
Welch Allyn	CP 100	▶
Welch Allyn	CP 150, CP 50	■
Enteral Feeding Pump		
Abbott Healthcare	Flexiflo Quantum	■
Abbott Healthcare	Patrol ePump	▶
Covidien	Kangaroo ePump	▶
Moog	EnteraLite ePump, Infinity Orange ePump	>
Fetal Monitor		
Philips	Avalon CL Fetal Monitor Base Station, Avalon Fetal Monitor FM20/30, Avalon Fetal Monitor FM40/50	▶
Philips	Fetal Monitor FM-2 (M2922A), Fetal Telemetry System Series 50 A (M1351A), Fetal Telemetry System Series 50 IP-2 (M1353A), Fetal Telemetry System Series 50 XM (M1350B), Fetal Telemetry System Series 50 XMO (M1350C)	▶
Philips	Fetal Telemetry System Series 50 T (M1310A)	>

Furniture / Upholstery / Apparel		
3M	M-Series Headgear	>
Hill-Rom	Surgical and Bariatric Beds	>
ICS Jail Supplies Inc.	Shower Shoes	▶
Medline Industries, Inc.	Mattresses and Mattress Pads	>
Micromedical	InView Goggles	>
Midmark	Barrier Free Examination Table, Barrier Free Universal Procedures Table, Model 109 Pediatric Table, Model 222 Exam Table, Model 223 Exam Table, Model 230 Universal Procedures Table, Model 244 Bariatric Power Treatment Table, Model 309 Pediatric Table	■
Midmark	Model 409 Pediatric Table, Model 491 Exam Chair, Model 604 Pediatric Table, Model 622 Exam Table, Model 623 Exam Table, Model 630 Universal Procedures Table, Model 640 Pediatric Table, Model 641 Power Procedures Chair, Model 646 Procedures Chair, Model 647 Procedure Chair, Model 712 Surgery Table, Model 95 Treatment Table, Otolaryngology Chair	>
MTI	Furniture	>



Glucometers / Blood Monitoring Systems		
Abbott Healthcare	FreeStyle Flash, FreeStyle Freedom, FreeStyle Freedom Lite, FreeStyle Lite, I-Stat	▶
Abbott Healthcare	Freestyle InsuLinx	>
Abbott Healthcare	OptiumEZ, Xceed Optium	▶
Abbott Healthcare	Precision Xceed Pro	▶
Abbott Healthcare	Precision Xtra	▶
AgaMatrix	Blood Glucose Meter	▶
Alere	PT-INR Monitoring System	▶
Arizant Healthcare	SpotOn Temp	▶
Arkray	Assure 4	▶
Arkray	Assure Brilliance	▶
Arkray	Assure Platinum	▶
Arkray	Assure Pro	▶
Bayer	Breeze 2	■
Bayer	Contour	■
Bayer	Contour Next, Contour Next EZ	>
Cambridge Sensors	Microdot	▶
Cambridge Sensors	Microdot Xtra	▶
Dexcom	Dexcom G4 Platinum	>

Glucometers / Blood Monitoring Systems		
Glucocard Vital	Assure	▶
Lifescan	Onetouch Ping, Onetouch VerioIQ	>
Lifescan	Onetouch Select, Onetouch Ultra, Onetouch Ultra2, Onetouch Ultramini	>
Lifescan	SureStep Meters, SureStepFlexx	▶
McKesson	TRUEtrack	▶
Nova Biomedical	Max, Max Link, Plus	▶
Nova Biomedical	StatSensor Creat, StatStrip Glu, StatStrip Glu/Ket, StatStrip LAC	▶
Oak Tree International	EasyMax LTC Meter, EasyMax NG Meter	▶
Prodigy Diabetes Care	Prodigy AutoCode, Prodigy Pocket, Prodigy Preferred, Prodigy Voice	▶
PSS World Medical	Quintet Meter	>
Roche	Accu-Check Inform I	▶
Roche	Accu-Check Inform II, Accutrend Plus	▶
Roche	Coagu-Chek XS-System	■
Wavesense	Keynote Pro	■



Handheld Scanner / Phone / Device		
Ascom	d41, d61, d62, d81, i62	>
Ascom	Myco	▶
Honeywell	Xenon™ 1900h, Xenon™ 1902h	▶
Motorola	DS6707 Scanner	▶
Hemodialysis		
Baxter	Prismaflex Dialysis Machine	>
Fresenius	2008K, Crit-Line Monitor III	>
Fresenius	GranuFlo Concentrate Dissolution	>
Fresenius	Liberty Cyclor, Newton Cyclor	■
Gambro	Phoenix	>
Incubators / Chambers		
Draeger	Air-Shields Isolette ModelC400QT	>
SECHRIST Industries	Hyperbaric Chambers	>

Infusion Pump / Systems		
Baxter	Auto Syringe AS50	■
Baxter	Colleague	■
Baxter	Flo-Gard 6200, Flo-Gard 6201, Flo-Gard 6301 Dual Channel, Syndeo	■
Baxter	Ipump	▶
Baxter	PCA II	■
Baxter	Sigma Spectrum	▶
BBraun	Horizon Nxt	■
BBraun	Outlook 100, Outlook 200	>
BBraun	Outlook ES	>
BBraun	Vista Pumps	>
Cardinal Health	Alaris 8000, Alaris 8015	▶
Carefusion	Alaris EtCo2 Module, Alaris Syringe Module, Alaris System Signature Edition Gold	▶
Carefusion	Alaris PC Unit, Alaris Pump Module	▶
Carefusion	Alaris PCA Module, Alaris SE Module, Alaris SpO2 Module	▶
Carefusion	Alaris Pump - IUI Connector (only)	■



Infusion Pump / Systems		
Hospira	Acclaim encore, OmniFlow 4000 Plus, Plum 1.6, Plum 360 Infuser - compatible with Hospira MedNet Software 30010-04,30011-04, 30012-04.30010-65, 30011-65,30012, Plum XL, Plum XL3 and Plum XL3 Micro/Macro™	>
Hospira	GemStar, LifeCare PCA, Plum A+, Plum A+ 3 w/Hospira MedNet, Plum A+ w/Hospira MedNet, Plum A+3	>
Hospira	LifeCare PCA 3, LifeCare PCA Plus II	>
Hospira	Sapphire Pump	>
Hospira	Symbiq Infusion System	>
Moog	4000 CMS, 4000 Plus, 6000 CMS, PainSmart IOD	▶
RMS Medical Products	Freedom 60 Syringe Infusion System	■
Smiths Medical	CADD Legacy PCA	▶
Smiths Medical	Medfusion 3500	■
Smiths Medical	Medfusion 4000	■

Laryngoscope		
Ambu-King Vision	King Vision	▶
Karl Storz	Karl Storz C-MAC Video Laryngoscopes	>
Verathon	GlideScope Cobalt, GlideScope GVL, Ranger, Ranger Single Use	▶
Verathon	GlideScope Cobalt AVL	▶
Lifts / Gurney / Harness		
ArjoHuntleigh	BHM V4i	▶
ArjoHuntleigh	Ergolift-600	▶
ArjoHuntleigh	Graduate, IR System, Maxi Sky 600, Maxi Sky 600 ECS, V10, V5 Duo, Voyager, Wipedown Standing Sling	▶
ArjoHuntleigh	Maxi Move, Sara 3000	>
ArjoHuntleigh	RotoRest Delta	>
Lighting		
Midmark	Model 354 Lighting System, Model 355 Lighting System	■
Steris	Amsco SQ240 Surgical Lighting Systems, Amsco SQ240 SurgiVision Surgical Lighting and Video, Harmony LA Surgical Lighting and Visualization System, Polaris Surgical Lighting	▶
Welch Allyn	Green Series Exam Lights	■
Welch Allyn	GS 777 Wall Transformer	>



Mammography		
GE Healthcare	Senographe 2000D, Senographe DMR, Senographe DRM+, Senographe DS, Senographe Essential	▶
Hologic	Mammography	▶
Siemens	Mammography Machines	>
Siemens	Somatom	■
Mobile Computing / Health IT		
Artromick	Aci, TX10, TX11, TX20	■
Artromick	LTC Autholock, NX15	>
Capsa Solutions	LX10, LX15, LX5	>
Capsa Solutions	VX25	▶
Enovate	Enovate	■
GE Healthcare	Case	▶
GE Healthcare	CIC Pro Clinical Info Center	■
Honeywell	Dolphin 9700hc mobile computer, Dolphin 99EXhc mobile computer, Dolphin7800hc mobile computer	▶
KSI (Key Source Intl.)	San-a-Key	■
Midmark	Workstation	■
Motorola	CA50 Device, MC9090 Mobile Computer	▶
Panasonic	Toughbook	▶
Socket Mobile	SOMO	>
Vecna	Vkiosk V7 Freestanding Kiosk	>

Parts		
Welch Allyn	3.5 V Lithium Ion Rechargeable Handle	■
Patient Monitors		
Edwards Lifesciences	OM-2 Optical Module	▶
Fukuda Denshi	DS-7100	▶
GE Healthcare	B30	■
GE Healthcare	Carescape V-100, Solar 7000, Solar 8000, Solar 9500	>
GE Healthcare	Corometrics 170 series	▶
GE Healthcare	Corometrics 250 series, Dash 2000, Dash 2500, Dash 3000, Dash 4000, Dash 5000, Dinamap ProCare	■
GE Healthcare	Dinamap Pro 1000, Dinamap Pro 100V2, Dinamap Pro 200V2, Dinamap Pro 300V2, Dinamap Pro 400V2, Solar 8000M	■
GE Healthcare	Giraffe Omnibed	■
GE Healthcare	Vital View	■
Hewlett Packard	AiO Monitor	■
Mindray Healthcare	Accutorr 3, Accutorr 7, DPM3	>
Mindray Healthcare	DPM-VS Monitor	>
NDS	Radiance 19" Medical Display	>
NEC	NEC Monitors	■
Nihon Kohden	BSM-6000, BSM-9000	■



Patient Monitors		
Philips	865299 15" LCD Color Display, 865299 17" LCD Color Display, 865299 19" LCD Color Display, Intellivue Multi-Measurement Module, Intellivue Patient Monitor MP 20/30, Intellivue Patient Monitor MP 40/50, Intellivue Patient Monitor MP 60/70/80/90, Intellivue Patient Monitor MP2, Intellivue Patient Monitor MP5, Intellivue Patient Monitor MP5SC, Intellivue Patient Monitor MP5T, Intellivue Patient Monitor MX600, Intellivue Patient Monitor MX700, Intellivue Patient Monitor MX800, Intellivue Transcutaneous Gas Measurement, M8031B 15" LCD Color Display, M8031B 17" LCD Color Display, M8031B 19" LCD Color Display, M8033C 15" LCD Color Display, M8033C 17" LCD Color Display, M8033C 19" LCD Color Display	▶
Philips	CMS Patient Monitoring System M1165-66-67-75-76-77A	▶
Philips	Criterion 60 Pressure-Oxygen Monitor	▶▶
Philips	IntelliVue MX40 Wearable Patient Monitor	▶▶
Philips	Intellivue Patient Monitor MP 30, Intellivue Patient Monitor MP 80, Intellivue Patient Monitor MP 90	▶
Philips	Intellivue Patient Monitor MP 70	>
Philips	Patient Monitor M3 & M4, Patient Monitor System M1205A V24 and V26	▶
Philips	SureSigns VS2	■
Philips	SureSigns VS3	>
Philips	SureSigns VS4	>

Patient Monitors		
Spacelabs	Patient Monitors	■
SunTech	Blood Pressure Monitor Model 247B	▶
Welch Allyn	Atlas, Consulting Pneumatic Otoscope, Digital MacroView	>
Welch Allyn	Connex Spot Vital Signs	▶
Welch Allyn	Connex Vital Signs Monitor 6000 Series	■
Welch Allyn	Propaq CS	■
Welch Allyn	Propaq Encore	■
Welch Allyn	Spot Vital Signs Lxi	■
Pulse Oximeter		
Masimo	Radical-7	>
Mindray Healthcare	DPM1	>
Nonin Healthcare	8500 Oximeters, 8500M Oximeters	■
Radiation Therapy		
Elekta	Gamma Knife	■



Respirators / Respiratory Therapy		
3M	PAPR GVP Powered Air Purifying Respirator	▶
Manufacturer	PAPR Powered Air Purifying Respirator, Reusable Respirator FF-400 Series, Reusable Respirator 6000 Series, RR Reusable Respirators, Versaflo TR-300 PAPR	>
Manufacturer	Reusable Respirator 7000 Series	▶
Prodigy Diabetes Care	Prodigy Mini-Mist	>
VapoTherm	Precision Flow	>
RFID / Security Tags / Badges		
Accutech	ES 2100, ES 2200, IS 3200, BR 4200	▶
Cardinal Health	Pyxis Bio ID	>
CenTrak	RFID Tags	▶
RF Technologies	Safe Place Infant Transmitter	>
Versus Technology	Badge	■

Stethoscope		
3M	Littman Cardiology III, Littman Cardiology STC, Littman Classic II Pedia & Infant, Littman Classic II S.E., Littman Master Cardiology, Littman Master Classic II, Littman Model 3000, Littman Model 3100, Littman Model 3200, Littman Model 4000, Littman Model 4100, Littman Select	>
Welch Allyn	Elite Electronic Stethoscope	>
Suction Unit		
Ohio Medical	Push-To-Set™ Intermittent Suction Unit PTS-ISU	>
Surfaces		
Corian	Corian	■
Formica Corporation	Formica	▶
Telemetry		
Cardiac Science	Quinton Exlipse Premier ECG	>
Fukuda Denshi	Dynascope Telemetry Transmitter LX-7120	■
Nihon Kohden	NTX ZM-540, NTX ZM-541PA	■
Philips	Intellivue Patient Monitor M2636B, Transmitters	>
Philips	Telemetry System	▶
Philips	TRx Tranceiver	>



Temperature Management Unit		
Arizant Healthcare	Bair Hugger, Bair Paws	▶
Arizant Healthcare	Bair Hugger Model 505, Bair Hugger Model 750, Bair Hugger Model 775	▶
Thermometers		
Bbraun	Thermoscan Ear Pro	■
Exergen	TAT-2000C	▶
Welch Allyn	Sure Temp Plus 690 and 692	■

Transducers / Probes		
Aloka	ProSound SSD-3500, SSD-4000, SSD-Alpha 10, SSD-Alpha 5, SSD-Alpha 6, SSD-Alpha 7	■
BK Medical	Transducer Model 2050, Transducer Model 2052, Transducer Model 8560, Transducer Model 8660, Transducer Model 8666-RF, Transducer Model 8667, Transducer Model 8670, Transducer Model 8802, Transducer Model 8808, Transducer Model 8808e, Transducer Model 8809, Transducer Model 8811, Transducer Model 8814, Transducer Model 8815, Transducer Model 8816, Transducer Model 8818, Transducer Model 8819, Transducer Model 8820e, Transducer Model 8822, Transducer Model 8823, Transducer Model 8824, Transducer Model 8826, Transducer Model 8827, Transducer Model 8830, Transducer Model 8836, Transducer Model 8837, Transducer Model 8838, Transducer Model 8848, Transducer Model 8862, Transducer Model 8863, Transducer Model 8870, Transducer Model UA1237	➤



Transducers / Probes		
Fujifilm SonoSite	X-Porte Ultrasound System C11xp Transducer	■
Fujifilm SonoSite	X-Porte Ultrasound System C35xp Transducer	■
Fujifilm SonoSite	X-Porte Ultrasound System D2xp Transducer, X-Porte Ultrasound System HFL38xp Transducer, X-Porte Ultrasound System HFL50xp Transducer, X-Porte Ultrasound System ICTxp Transducer, X-Porte Ultrasound System L25xp Transducer, X-Porte Ultrasound System L38xp Transducer, X-Porte Ultrasound System P10xp Transducer, X-Porte Ultrasound System P21xp Transducer	■
Fujifilm SonoSite	X-Porte Ultrasound System HSL25xp Transducer	■

Transducers / Probes		
GE Healthcare	Logiq Probe Model 3S-D, Logiq Probe Model 4C-A, Logiq Probe Model 4D10L, Logiq Probe Model 4D16L, Logiq Probe Model PA2-5P, Logiq Probe Model PCW4.0, Logiq Probe Model RAB6-D, Logiq Probe Model RM14L, Logiq Probe Model RM6C, Logiq Probe Model RRE6-10, Logiq Probe Model RRE6-10-D, Logiq Probe Model RSM5-14, Logiq Probe Model RSP5-12, Logiq Probe Model RSP6-12, Logiq Probe Model RSP6-16, Logiq Probe Model RSP6-16-D, Logiq Probe Model SCW2.0, Logiq Probe Model SP10-16, Logiq Probe Model SP10-16-D, Logiq Probe Model SP4-10	■
GE Healthcare	TEE Probe, Vivid	>
GE Healthcare	Logiq Probe Model 10L, Logiq Probe Model 10LB, Logiq Probe Model 10Lb-RS, Logiq Probe Model 10S, Logiq Probe Model 10S-RS, Logiq Probe Model 11L, Logiq Probe Model 11L-D, Logiq Probe Model 12L, Logiq Probe Model 12L-RS, Logiq Probe Model 12L-SC, Logiq Probe Model 12S-D, Logiq Probe Model 12S-RS, Logiq Probe Model 16L-RS, Logiq Probe Model 3.5C (continued on next page)	>



Transducers / Probes	
GE Healthcare	Logiq Probe Model 3.5CS, Logiq Probe Model 348C, Logiq Probe Model 3C, Logiq Probe Model 3C-RS, Logiq Probe Model 3Cb, Logiq Probe Model 3CRF, Logiq Probe Model 3CRF-D, Logiq Probe Model 3S, Logiq Probe Model 3S-LC, Logiq Probe Model 3S-RC, Logiq Probe Model 3S-RS, Logiq Probe Model 3S-SC, Logiq Probe Model 3Sc-RS, Logiq Probe Model 3SR, Logiq Probe Model 3V, Logiq Probe Model 3V-D, Logiq Probe Model 4C, Logiq Probe Model 4C-D, Logiq Probe Model 4C-RS, Logiq Probe Model 4C-SC, Logiq Probe Model 4S, Logiq Probe Model 4V-D, Logiq Probe Model 546L, Logiq Probe Model 5C, Logiq Probe Model 5CS, Logiq Probe Model 5S, Logiq Probe Model 5S-LC, Logiq Probe Model 5S-RS, Logiq Probe Model 6S, Logiq Probe Model 6S-D, Logiq Probe Model 6S-RS, Logiq Probe Model 739L, Logiq Probe Model 7L, Logiq Probe Model 7S, Logiq Probe Model 7S-RS, Logiq Probe Model 8C, Logiq Probe Model 8C-RS (continued on next page)



Transducers / Probes	
GE Healthcare	Logiq Probe Model 8L, Logiq Probe Model 8L-RS, Logiq Probe Model 9L, Logiq Probe Model 9L-D, Logiq Probe Model 9L-RS, Logiq Probe Model BE9C, Logiq Probe Model BE9CS, Logiq Probe Model C1-5-D, Logiq Probe Model C1-5-RS, Logiq Probe Model C1-6-D, Logiq Probe Model C2-5-RC, Logiq Probe Model C2-9-D, Logiq Probe Model C31, Logiq Probe Model C358, Logiq Probe Model C36, Logiq Probe Model C364, Logiq Probe Model C386, Logiq Probe Model C4-8-D, Logiq Probe Model C55, Logiq Probe Model C551, Logiq Probe Model C721, Logiq Probe Model CAE, Logiq Probe Model CBF, Logiq Probe Model CS, Logiq Probe Model CZB, Logiq Probe Model E72, Logiq Probe Model E721, Logiq Probe Model E8C, Logiq Probe Model E8C-RS, Logiq Probe Model E8CS, Logiq Probe Model E8CS-RC, Logiq Probe Model E8CS-SC, Logiq Probe Model G3S, Logiq Probe Model i12L, Logiq Probe Model i12L-RS, Logiq Probe Model I739, Logiq Probe Model I739-RS (continued on next page)





Transducers / Probes		
GE Healthcare	Logiq Probe Model 4D3C-L, Logiq Probe Model 4D5C-L, Logiq Probe Model 4D8C, Logiq Probe Model 4DE7C, Logiq Probe Model AB2-5, Logiq Probe Model AB2-7, Logiq Probe Model AB2-7-D, Logiq Probe Model AB4-8, Logiq Probe Model AC2-5, Logiq Probe Model IC5-9, Logiq Probe Model PA2-5, Logiq Probe Model PA4-7, Logiq Probe Model PA6-8, Logiq Probe Model PA6-8D, Logiq Probe Model RAB2-5, Logiq Probe Model RAB2-5-D, Logiq Probe Model RAB2-5L, Logiq Probe Model RAB4-8-D, Logiq Probe Model RAB4-8L, Logiq Probe Model RAB4-8P, Logiq Probe Model RAM3-8, Logiq Probe Model RIC5-9, Logiq Probe Model RIC5-9-D, Logiq Probe Model RIC5-9H, Logiq Probe Model RIC5-9W, Logiq Probe Model RIC6-12-D, Logiq Probe Model RNA5-9, Logiq Probe Model RNA5-9D, Logiq Probe Model RRE5-10-D, Logiq Probe Model SP6-12	▶
MediStim	Doppler X-Plore	>
MediStim	Flow probes, PA-Probes, PB-Probes, PQ-Probes, PR-Probes	▶

Transducers / Probes		
Philips	Affiniti Series Transducers 3D9-3v, BP10-5ec, C9-4v, C10-3v, C10-4ec, Affiniti Series Transducers C5-1; C6-2; C8-5; C9-2, Affiniti Series Transducers D2tcd; D2cwc; D5cwc, Affiniti Series Transducers L12-3; L12-4; L12-5 50 Compact; L15-7io; L18-5, Affiniti Series Transducers S4-2; S5-1; S8-3; S12-4, Affiniti Series Transducers V6-2; VL13-5, CX30 Ultrasound System, Philips IE33 Transducers C5-1, C5-2, C8-5, C9-4, Philips IE33 Transducers D2cwc, D5cwc, D2tcd, Philips IE33 Transducers L11-3, Philips IE33 Transducers L15-7io, Philips IE33 Transducers L8-4, L9-3, Philips IE33 Transducers S5-1, S8-3, S12-4, Philips IE33 Transducers X3-1, X5-1, X7-2	>
Philips	ClearVue 350 and 550 Ultrasound Systems Transducer	>
Philips	ClearVue 580 Ultrasound System Transducer, ClearVue 650 Ultrasound System Transducer	>
Philips	iE33	>



Transducers / Probes		
Philips	Affiniti Series Transducers S7-3t, Affiniti Series Transducers X7-2t, Philips IE33 Transducers S7-2t omni,X7-2t	>
Philips	iU22	■
Philips	Philips IE33 Transducers S7-3t,S8-3t	>
Philips	Avalon CTS Cordless Fetal Transducer System M2720A	▶
Philips	EnVisor, HD11, HD11 XE, HD15, HD7, HD9, Sparq	▶
Siemens	Acuson S2000 Transducer Model 6C1 HD, Acuson S3000 Transducer Model 6C1 HD, Acuson SC2000 Transducer Model 6C1HD, Acuson X700 Transducer Model VF 12-4	■
Siemens	Acuson X300 PE Transducer Model BP9-4	>
Siemens	Acuson S2000 Transducer Model 18L6 HD, Acuson S3000 Transducer Model 18L6 HD, Acuson Sequoia Transducer Model 15L8W	>

Transducers / Probes		
Siemens	Acuson S1000 Transducer Model EV-8C4, Acuson S1000 Transducer Model V7M, Acuson S2000 ABVS Transducer Model 14L5 BV, Acuson S2000 Transducer Model EV-8C4, Acuson S3000 Transducer Model 8C3 HD, Acuson S3000 Transducer Model EV-8C4, Acuson S3000 Transducer Model V7M, Acuson Sequoia Transducer Model 17L5 HD, Acuson Sequoia Transducer Model EC-10C5, Acuson Sequoia Transducer Model EV-8C4, Acuson Sequoia Transducer Model V7B, Acuson Sequoia Transducer Model V7M, Acuson X150 Transducer Model L9-5, Acuson X500 Transducer Model BE9-4, Acuson X500 Transducer Model EC9-4, Sonoline G20 Transducer Model EC9-4, Sonoline G40 Transducer Model L9-5, Sonoline G50 Transducer Model BE9-4, Sonoline G50 Transducer Model Endo V 3D (continued on next page)	>



Transducers / Probes		
Siemens	Sonoline G60 Transducer Model BE9-4, Sonoline G60 Transducer Model EC9-4	>
Siemens	Acuson Antares Transducer Model C5-2, Acuson Antares Transducer Model C5F1, Acuson Antares Transducer Model C7F2, Acuson Antares Transducer Model CH4-1, Acuson Antares Transducer Model CH6-2, Acuson Antares Transducer Model CW2, Acuson Antares Transducer Model CW5, Acuson Antares Transducer Model CX5-2, Acuson Antares Transducer Model EC9-4, Acuson Antares Transducer Model EV9F4, Acuson Antares Transducer Model P10-4, Acuson Antares Transducer Model PH4-1, Acuson Antares Transducer Model PX4-1, Acuson Antares Transducer Model VF-X9-4, Acuson Antares Transducer Model VF10-5, Acuson Antares Transducer Model VF13-5 (continued on next page)	▶

Transducers / Probes		
Siemens	Acuson Antares Transducer Model VF13-5 SP, Acuson Antares Transducer Model VF7-3, Acuson Antares Transducer Model VFX13-5, Acuson S2000 Transducer Model 7CF2	▶
Siemens	Acuson X300 Transducer Model VF13-5, Acuson X500 Transducer Model EV9-4	▶
Siemens	Acuson SC2000 Transducer Model V7M	▶
Siemens	Acuson CV70 Transducer Model 5.0L45, Acuson CV70 Transducer Model 5.0P10, Acuson CV70 Transducer Model C6-2, Acuson CV70 Transducer Model CW2, Acuson CV70 Transducer Model CW5, Acuson CV70 Transducer Model L10-5, Acuson CV70 Transducer Model P4-2, Acuson CV70 Transducer Model V5Ms, Acuson Cypress Transducer Model 7V3C, Acuson Cypress Transducer Model Aux CW, Acuson Cypress Transducer Model V5Ms, Acuson S1000 Transducer Model 4V1, Acuson S1000 Transducer Model 7CF1 (continued on next page)	▶



Transducers / Probes	
Siemens	<p>Acuson S1000 Transducer Model 7CF2, Acuson S1000 Transducer Model 9EVF4, Acuson S1000 Transducer Model CW2, Acuson S1000 Transducer Model CW5, Acuson S1000 Transducer Model EC9-4, Acuson S1000 Transducer Model MC9-4, Acuson S1000 Transducer Model V5Ms, Acuson S2000 Transducer Model 7CF1, Acuson S2000 Transducer Model 9EVF4, Acuson S2000 Transducer Model CW2, Acuson S2000 Transducer Model CW5, Acuson S2000 Transducer Model EC9-4, Acuson S2000 Transducer Model V5Ms, Acuson S3000 Transducer Model 7CF1, Acuson S3000 Transducer Model 7CF2, Acuson S3000 Transducer Model 9EVF4, Acuson S3000 Transducer Model CW2, Acuson S3000 Transducer Model CW5, Acuson S3000 Transducer Model EC9-4, Acuson S3000 Transducer Model MC9-4, Acuson S3000 Transducer Model V5Ms, Acuson SC2000 Transducer Model 4Z1c, Acuson SC2000 Transducer Model AUX CW, Acuson SC2000 Transducer Model CW2, Acuson SC2000 Transducer Model CW5, Acuson SC2000 Transducer Model V5Ms, Acuson Sequoia Transducer Model 15L8, Acuson Sequoia Transducer Model 3V2c, Acuson Sequoia Transducer Model 5V2c, Acuson Sequoia Transducer Model 6L3, Acuson Sequoia Transducer Model 7V3c, Acuson Sequoia Transducer Model 8C4, Acuson Sequoia Transducer Model 8L5, Acuson Sequoia Transducer Model 8V5, Acuson Sequoia Transducer Model Aux CW, Acuson Sequoia Transducer Model V5Ms, Acuson X150 Transducer Model EC9-4 (continued on next page)</p>



Transducers / Probes	
Siemens	<p>Acuson X150 Transducer Model EV9-4, Acuson X300 PE Transducer Model C6-2, Acuson X300 PE Transducer Model C7F2, Acuson X300 PE Transducer Model C8-5, Acuson X300 PE Transducer Model CH5-2, Acuson X300 PE Transducer Model CW2, Acuson X300 PE Transducer Model CW5, Acuson X300 PE Transducer Model EC9-4, Acuson X300 PE Transducer Model EV9-4, Acuson X300 PE Transducer Model EV9F4, Acuson X300 PE Transducer Model P5-1, Acuson X300 PE Transducer Model P9-4, Acuson X300 PE Transducer Model V5Ms, Acuson X300 PE Transducer Model VF8-3, Acuson X300 Transducer Model CW2, Acuson X300 Transducer Model CW5, Acuson X300 Transducer Model EC9-4, Acuson X300 Transducer Model EV9-4, Acuson X300 Transducer Model VF8-3, Acuson X500 Transducer Model C6-2, Acuson X500 Transducer Model C7F2, Acuson X500 Transducer Model C8-5, Acuson X500 Transducer Model CW2, Acuson X500 Transducer Model CW5, Acuson X500 Transducer Model L10-5, Acuson X500 Transducer Model MPT7-4, Acuson X500 Transducer Model P4-2, Acuson X500 Transducer Model P9-4, Acuson X700 Transducer Model C7F2, Acuson X700 Transducer Model CW2, Acuson X700 Transducer Model CW5, Acuson X700 Transducer Model EC9-4w, Acuson X700 Transducer Model EV9F4, Acuson X700 Transducer Model V5Ms, Sonoline G20 Transducer Model 7.5L75S, Sonoline G20 Transducer Model C5-2, Sonoline G20 Transducer Model EV9-4 (continued on next page)</p>





Transducers / Probes		
Siemens	<p>Sonoline G20 Transducer Model L10-5, Sonoline G40 Transducer Model CH5-2, Sonoline G40 Transducer Model EC9-4, Sonoline G40 Transducer Model EV9-4, Sonoline G40 Transducer Model P4-2, Sonoline G40 Transducer Model VF10-5, Sonoline G50 Transducer Model, Sonoline G50 Transducer Model 5.0 C50+, Sonoline G50 Transducer Model 5.0P10, Sonoline G50 Transducer Model 7.5L70, Sonoline G50 Transducer Model C5-2, Sonoline G50 Transducer Model C6-2, Sonoline G50 Transducer Model C6-3 3D, Sonoline G50 Transducer Model C8-5, Sonoline G50 Transducer Model CW2, Sonoline G50 Transducer Model CW5, Sonoline G50 Transducer Model Endo VII, Sonoline G50 Transducer Model EV9-4, Sonoline G50 Transducer Model L10-5, Sonoline G50 Transducer Model MPT7-4, Sonoline G50 Transducer Model VF13-5, Sonoline G50 Transducer Model VF13-5 SP, Sonoline G60 Transducer Model C5-2, Sonoline G60 Transducer Model VF13-5 SP, Sonoline G60 Transducer Model 5.0 C50+, Sonoline G60 Transducer Model 5.0L45, Sonoline G60 Transducer Model 5.0P10, Sonoline G60 Transducer Model 7.5L70, Sonoline G60 Transducer Model C6-2, Sonoline G60 Transducer Model C6-3 3D, Sonoline G60 Transducer Model C8-5, Sonoline G60 Transducer Model CW2, Sonoline G60 Transducer Model CW5, Sonoline G60 Transducer Model Endo V 3D, Sonoline G60 Transducer Model Endo VII, Sonoline G60 Transducer Model EV9-4, Sonoline G60 Transducer Model L10-5 (continued on next page)</p>	▶

Transducers / Probes		
Siemens	<p>Sonoline G60 Transducer Model MPT7-4, Sonoline G60 Transducer Model VF13-5</p>	▶
Siemens	<p>Acuson CV70 Transducer Model VF13-5 SP, Acuson Cypress Transducer Model 3V2C, Acuson Cypress Transducer Model 4C1, Acuson S1000 Transducer Model 10V4, Acuson S1000 Transducer Model 14L5, Acuson S1000 Transducer Model 14L5 SP, Acuson S1000 Transducer Model 4C1, Acuson S1000 Transducer Model 4P1, Acuson S1000 Transducer Model 4V1c, Acuson S1000 Transducer Model 6C2, Acuson S1000 Transducer Model 8V3, Acuson S1000 Transducer Model 9L4, Acuson S2000 Transducer Model 10V4, Acuson S2000 Transducer Model 14L5, Acuson S2000 Transducer Model 14L5 SP, Acuson S2000 Transducer Model 4C1, Acuson S2000 Transducer Model 4P1, Acuson S2000 Transducer Model 4V1, Acuson S2000 Transducer Model 4V1c, Acuson S2000 Transducer Model 6C2, Acuson S2000 Transducer Model 8V3, Acuson S2000 Transducer Model 9L4, Acuson S3000 Transducer Model 10V4, Acuson S3000 Transducer Model 14L5, Acuson S3000 Transducer Model 14L5 SP, Acuson S3000 Transducer Model 4C1, Acuson S3000 Transducer Model 4P1, Acuson S3000 Transducer Model 4V1, Acuson S3000 Transducer Model 4V1c, Acuson S3000 Transducer Model 6C2, Acuson S3000 Transducer Model 8V3, Acuson S3000 Transducer Model 9L4, Acuson SC2000 Transducer Model 10V4, Acuson SC2000 Transducer Model 4V1c, Acuson SC2000 Transducer Model 8V3 (continued on next page)</p>	▶



Transducers / Probes		
Siemens	Acuson SC2000 Transducer Model 9L4, Acuson Sequoia Transducer Model 10V4, Acuson Sequoia Transducer Model 4C1, Acuson Sequoia Transducer Model 4V1/4V1c, Acuson Sequoia Transducer Model 6C2, Acuson Sequoia Transducer Model 8L5T, Acuson Sequoia Transducer Model 8V3c, Acuson Sequoia Transducer Model 9L4, Acuson X150 Transducer Model CH5-2, Acuson X150 Transducer Model P4-2, Acuson X150 Transducer Model P8-4, Acuson X150 Transducer Model VF10-5, Acuson X150 Transducer Model VF13-5, Acuson X300 PE Transducer Model P4-2, Acuson X300 PE Transducer Model P8-4, Acuson X300 PE Transducer Model VF10-5, Acuson X300 PE Transducer Model VF13-5, Acuson X300 PE Transducer Model VF13-5 SP, Acuson X300 Transducer Model CH5-2, Acuson X300 Transducer Model P4-2, Acuson X300 Transducer Model P8-4, Acuson X300 Transducer Model VF10-5, Acuson X300 Transducer Model VF13-5, Acuson X300 Transducer Model VF13-5 SP, Acuson X500 Transducer Model VF13-5, Acuson X500 Transducer Model VF13-5 SP, Acuson X700 Transducer Model 4C1, Acuson X700 Transducer Model 4V1c, Acuson X700 Transducer Model 6C2, Acuson X700 Transducer Model VF10-5	▶
Sonosite	L25n	■
Sonosite	L25n, L38n, L38xi, L52n Vet Only, P21n	■
Sonosite	L25, L25e, L25x, L38, L38e, L38x, L52 Vet Only, L52e Vet Only, L52x Vet Only, P10, P10x, P17, P21x	■

Transducers / Probes		
Sonosite	D2, D2x, HFL38, HFL38x	➤
Sonosite	C11n, C60n, ICTn	➤
Sonosite	C11, C11 Case Urethane, C11e, C11x, C15, C15 Cast Urethane, C15e, C60, C60e, C60x, C8, C8e, C8x, HFL50x, ICT, ICTx, SLA, SLAx, SLT	➤
Sonosite	P11x	▶
Spacelabs	Transducers	■
Toshiba	Toshiba Transducers	■
Toshiba	Toshiba Transducers	■
Toshiba	Toshiba Transducers	▶
Zonare	Model	■
Manufacturer	Transducer	➤



Ultrasound Systems		
Bard	Site-Rite 8	▶
Bard	Site Rite Vision	▶
Bard	Site Rite 6, Site Rite Prevue	▶
Echosens	FibroScan 502	■
Fujifilm SonoSite	X-Porte Ultrasound System Footswitch	■
Fujifilm SonoSite	M-Turbo System, X-Porte Ultrasound System ECG Cable	■
Fujifilm SonoSite	X-Porte Ultrasound System	■
Fujifilm SonoSite	X-Porte Ultrasound System ACDC Power Supply	■
Fujifilm SonoSite	180PLUS	■
Fujifilm SonoSite	iViz Ultrasound Protective Case	■
Fujifilm SonoSite	iViz Ultrasound P21V	■
Fujifilm SonoSite	X-Porte Ultrasound System PowerPark Docking Station	■
Fujifilm SonoSite	iViz Ultrasound Carry Case, X-Porte Ultrasound System Docking Stand, X-Porte Ultrasound System X-Porte TTC Accessory	■
Fujifilm SonoSite	Edge System	▶
Fujifilm SonoSite	iViz Ultrasound	▶

Ultrasound Systems		
GE Healthcare	Logiq E9	■
GE Healthcare	LOGIQ A5, Logiq Book XP Enhanced, Logiq E, Logiq E9 with XDclear, Logiq I, Logiq P3, Logiq P5, Logiq P6, Logiq S7, Somo V Abus, Somo V Abus Platinum, Vivid 7, Vivid E, Vivid E9 with XDclear, Vivid i, Vivid q, Vivid S5, Vivid S6, Voluson, Voluson 730, Voluson e, Voluson E6, Voluson E8 / E8 Expert, Voluson i, Voluson P8, Voluson S6, Voluson S8, Vscan Pocket Ultrasound	>
GE Healthcare	Logiq S8	▶
GE Healthcare	Venue 40	▶
Philips	CX50	■
Philips	EPIQ 5 Ultrasound System, EPIQ 7 Ultrasound System, iE33 Ultrasound System	■
Philips	Affiniti Series Ultrasound System, Philips iE33 Ultrasound System	■
Philips	ImagePoint, OptiGo, Sonos 5500, Sonos 7500	▶
Philips	Apogee, HDI	▶
Philips	Ultramark	▶



Ultrasound Systems		
Rich Mar Corp	Therasound 3 Ultrasound	>
Siemens	Acuson X300	■
Siemens	Acuson Sequioa	>
Siemens	Acuson Antares, Acuson CV70, Acuson Cypress, Acuson S1000, Acuson S3000, Acuson SC2000, Acuson X150, Acuson X300 PE, Acuson X500, Acuson X700, Sonoline G20, Sonoline G40, Sonoline G50, Sonoline G60	▶
Sonosite	TITAN	■
Sonosite	iLook15, iLook25	■
Sonosite	MicroMaxx	■
Sonosite	NanoMaxx	■
Sonosite	S-Cath, S-FAST, S-ICU, S-MSK, S-Nerve, S-Women's Health	■

Ventilator		
Carefusion	Avea Ventilator Systems	>
Covidien	Puritan-Bennet 840 Ventilator, Puritan-Bennet 980 Ventilator	▶
Draeger	Evita Infinity V500, Evita XL	>
Philips	BiPAPFocus Ventilator System, Respironics Trilogy02, Respironics Trilogy100, Respironics Trilogy200, Respironics Trilogy202	▶
Philips	Respironics V200	>
Philips	Respironics V60 Ventilator	■
Women's Health		
Ameda	Platinum Breast Pump	■
Welch Allyn	KleenSpec Cordless Illumination System	■
Welch Allyn	KleenSpec	■

Authorised distributor

In Australia:

For customer service, call 1300-735-292
To email an order, ordersanz@thermofisher.com
To order online: thermofisher.com

In New Zealand:

For customer service, call 0800-933-966
To email an order, ordersanz@thermofisher.com
To order online: thermofisher.com

ThermoFisher
SCIENTIFIC

Find out more about how our Layered Approach to Infection Prevention can help you.

PDI Ltd, Pywell Road, Willowbrook Est. Industrial Estate, Corby - Northants, NN17 5XJ, UK

11.09.2023 V1 THERMO FISHER SCIENTIFIC



▶	Compatible	■	Not compatible	>	Action required, no information – contact equipment manufacturer for recommendations
---	------------	---	----------------	---	--