



promed®

Latex Powder Free

Promed® Examination Gloves

- » Natural Rubber Latex
- » Non-Sterile
- » Single-Use
- » Powder Free
- » Micro Textured Surface

Comfort

Natural Latex Rubber:

- Natural elasticity provides a soft, comfortable fit and reduces hand fatigue; provides optimal sensitivity

Low Protein Level <50µg/g:

- Reduces risk of allergic reactions

Micro Textured Surface:

- Provides non-slip grip in wet and dry conditions

Ambidextrous:

- Greater comfort, fit and feel

Beaded Cuff:

- Ease of donning and durability during use

Cuff Length:

- Greater comfort, fit and feel

Protection

Suitable Use:

- Medical examination - non-invasive diagnostic and therapeutic procedures

Barrier Protection:

- Effective against potentially infective materials and other contaminants

Strength and Durability:

- Resistant to puncture and tears; resistant to punctures with re-seal qualities

Quality Standards:

- CE Mark - Medical 93/42/EEC, ASTM3578-00, Exceeds AS4011 for pinholes, AQL1.5, ARTG No 98364



Latex Powder Free

Promed® examination gloves

Product	Size	Catalogue No.	Packaging
Latex Powder Free examination gloves; 100 gloves	X-small	PMD148060	10/Carton
Latex Powder Free examination gloves; 100 gloves	Small	PMD148070	10/Carton
Latex Powder Free examination gloves; 100 gloves	Medium	PMD148080	10/Carton
Latex Powder Free examination gloves; 100 gloves	Large	PMD148090	10/Carton
Latex Powder Free examination gloves; 100 gloves	X-large	PMD148095	10/Carton

Glove Holder

One, two and three tier

Product	Size	Catalogue No.	Packaging
Glove holder, wall mounted	One tier	APD3210	Each
Glove holder, wall mounted	Two tier	APD3450	Each
Glove holder, wall mounted	Three tier	APD3640	Each



Glove Box Colours

Size	Colour
X-Small	Yellow
Small	Light Blue
Medium	Pink
Large	Blue
X-Large	Purple

ThermoFisher
SCIENTIFIC

For customer service, call 1300-735-292
To fax an order, use 1300-731-771
Visit us online: www.thermofisher.com.au