



promed.

Nitrile Soft

Promed® Examination Gloves

- » Latex Free
- » Non-Sterile
- » Single-Use
- » Powder Free
- » 200 gloves/box for efficient storage, less waste and fewer dispenser changes
- » Blue coloured glove for easy identification

Comfort

Latex Free and Powder Free:

- Reduces risk of allergic reactions and prevents type 1 allergic reactions

Improved Conformability:

- Similar flexibility to natural rubber, conforms to the wearer's hand during use and allows optimal sensitivity.

Micro Textured Surface:

- Provides non-slip grip in wet and dry conditions

Ambidextrous:

- Greater comfort, fit and feel

Beaded Cuff:

- Ease of donning and durability during use

Cuff Length:

- Available in standard 240mm

Protection

Suitable Use:

- Medical examination - non-invasive diagnostic and therapeutic procedures

Barrier Protection:

- Effective against potentially infective materials and other contaminants

Strength and Durability:

- Resistant to puncture and tears; resistant to development of holes

Quality Standards:

- CE Mark - Medical 93/42/EEC, ASTM D 6319, Exceeds AS4011 for pinholes, AQL1.5, ARTG No 100975

Nitrile Soft

Promed® examination gloves

Product

Size

Catalogue No.

Packaging

Nitrile Soft examination gloves; 200 gloves	X-small	PMD173060	10/Carton
Nitrile Soft examination gloves; 200 gloves	Small	PMD173070	10/Carton
Nitrile Soft examination gloves; 200 gloves	Medium	PMD173080	10/Carton
Nitrile Soft examination gloves; 200 gloves	Large	PMD173090	10/Carton
Nitrile Soft examination gloves; 200 gloves	X-large	PMD173095	10/Carton

Glove Holder

One, two and three tier

Product

Size

Catalogue No.

Packaging

Glove holder, wall mounted	One tier	APD3210	Each
Glove holder, wall mounted	Two tier	APD3450	Each
Glove holder, wall mounted	Three tier	APD3640	Each

Glove Box Colours

Size

Colour

X-Small	Yellow
Small	Light Blue
Medium	Pink
Large	Blue
X-Large	Purple



ThermoFisher
SCIENTIFIC

For customer service, call 1300-735-292
To fax an order, use 1300-731-771
Visit us online: www.thermofisher.com.au