

Custom Standard DNA Oligos

Custom oligos are also available in multi-well plates. Orders of at least 24 oligos can be delivered in 96-well plates and at least 192 oligos in 384-...

[Plate Oligos ordering tool](#)

1

Our Oligo Inquiries team is here should you have questions or require assistance placing your order. Please contact our technical team at

OligoInquiries@thermofisher.com

2

Use the Custom Standard DNA Oligo Tool below to add your oligo sequences or benefit from our Bulk Upload option to import up to 200 oligo sequences at one time.

[Bulk Upload](#)

3

Project 1

[Set project preferences](#)

4

1. Oligo Name

Researcher Name

My Oligo

Dr. Doe

[See sample label](#)

5

Oligo Sequence (5' to 3')

ACG TAC GTA CGT ACG TAC GT

[Add mixed bases](#)

Base Count: 20

Tm: 53.4 (50mM), 75.0 (1M) GC%: 50 Minimum Yield: 55.9 (µg), 9.1 (nmol)

6

Synthesis Scale

25 nmole

Purification

Desalted

5' Mod

None

3' Mod

None

7

Formulation

Dry

[Show more options](#)

8

Standard Shipping

Same Day Shipping

USD 4

Ship by **May 10**. Order within **1 hr 44 min**.

1 Oligo

USD 4

[Add to Cart](#)

Custom DNA oligo online ordering guide

Ordering custom DNA and modified oligos is easy with these features and tips

1. Use the Plate Oligos ordering tool if you would like to receive your oligos in 96- or 384-well plates instead of tubes.
2. Get help from our dedicated support team if you have any questions.
3. Bulk upload is a quick alternative to entering oligos individually. You can use the tool to access a Microsoft[™] Excel[™] template and upload your desired sequences and oligo configuration.
4. Set default synthesis scale, purification, and formulation options for the entire project to save time entering your oligos.
5. Cut and paste or type in your sequence and then press “Enter,” or start selecting configuration options to refresh the oligo price. Use the “Add mixed bases” button to insert phosphorothioates or degenerate bases at specific positions in the oligo.
6. The minimum yield will update in real time as changes are made to the selected scale and purification options.
7. Selection options are organized to be completed from left to right based on dependencies. If you are unable to select your desired purification or modification options because they appear grayed out, try increasing your synthesis scale.
8. Same-day shipping is available for a subset of oligo types and scales based on your local time when placing an order.

Order now at thermofisher.com/ordercustomoligos

invitrogen