

ThermoFisher
SCIENTIFIC

IL-17A

IL-6

IFN

TNF- α

VEGF

Choose high-quality Invitrogen ELISA kits for your protein quantitation

PROMO CODE
RELI25

Valid till 30 Sep 2022

Purchase any **2 ELISA kits** to get **25% off**

AND redeem a **FREE**

**FUEL YOUR
RESEARCH** **PLAY**




Find out more at thermofisher.com/elisapromo2022

ThermoFisher
SCIENTIFIC

Over 2000 highly-referenced ELISA kits you can trust

| | Uncoated ELISA | Coated ELISA | Instant ELISA | ProQuantum |
|-----------------------------------|-------------------------------|-------------------------------|-------------------------------|---------------------|
| Overnight coating required | Yes | No, pre-coated plate | No, pre-coated plate | No coating required |
| Analytical sensitivity | <5 pg/ml | <5 pg/ml | <5 pg/ml | <1 pg/mL |
| Dynamic range | <7.8-500 pg/ml | <7.8-500 pg/ml | <7.8-500 pg/ml | <1.6-5,000 pg/mL |
| Sample volume | 10-100 µL | 10-100 µL | 25-50 µL | 2-5 µL |
| Hands-on time | 1 hr 45 mins | 1 hr 20 mins | 20 mins | 30 mins |
| Incubation time | Overnight | 2.5-4 hr | 4 hr | 1 hr |
| Instrumentation needed | Microplate reader, absorbance | Microplate reader, absorbance | Microplate reader, absorbance | qPCR instrument |

Experience our one-wash Instant ELISA kit! Switch now!



- **Reduce steps and hands—on time by 60%**
- **Maximum accuracy**—reduced handling means less error and more consistent results
- **Better value**—generate standard curve data in parallel with additional well strips provided to enable use of all 96 wells for your samples



Biomarker Guide to Protein and RNA Quantitation

Navigate the different tools available for protein and RNA quantitation in oncology, inflammation research, immunology, metabolism, apoptosis, cell signaling, neurology.



Download Now

This promotion is open only to customers in Singapore. Customers who purchase any 2 ELISA kits are entitled to 25% off the list price and a free Fuel Your Research Play. To enjoy the discount and redeem your free play, please present the promo code RELI25 at time of order no later than Sep 30, 2022. For the full list of eligible items, please visit thermofisher.com/elisapromo2022. The QR code for the free play will be sent to the email address from which the order is sent from. Due to the redemption process, please allow 2 weeks for the QR code to be sent. Offer valid till 30th Sept 2022, or until promotional supplies are depleted, whichever comes first. There is no limit to the number of redemptions a lab may make. To find out where our Fuel Your Research Kiosks are, please visit thermofisher.com/FuelYourResearch. Cannot be combined with other discounts or promotions. Limited to Research Use Only products. By participating in the promotion, you warrant that you are not prohibited by employment, contract, or law from accepting a gift from Thermo Fisher Scientific. Healthcare professionals may not participate in this promotion. Offer void where prohibited, licensed or restricted by local laws or regulation or agency/institutional policy. Other restrictions may apply. Thermo Fisher Scientific reserves the right to change these terms and conditions at any time without prior notice.