

製品に関する重要なお知らせ

2023年7月

お客様各位

東京都港区芝浦4-2-8
サーモフィッシャーサイエンティフィック
ライフテクノロジーズジャパン株式会社

抗体製品の防腐剤に関するお知らせ

拝啓 時下ますますご清栄のこととお慶び申し上げます。平素は弊社製品につきまして格別のお引き立てを賜り、厚く御礼申し上げます。

この度、弊社抗体製品^{※注}におきまして、製造ロット毎に異なる防腐剤が含有するという事案が発生しております。含有の可能性のある防腐剤は、チメロサル（含有量0.01%）、アジ化ナトリウム（含有量は0.1%未満）、プロクリン300（含有量は0.1%未満）の3種類となります。製品番号が同一の場合でも、製造ロット番号が異なる場合は防腐剤の種類も変更となる可能性があります。

0.01%チメロサルを含有する当該ロット製品につきましては、毒物及び劇物取締法における医薬用外劇物に該当します。加えてチメロサル（エチル水銀チオサリチル酸ナトリウム）は有機水銀化合物であり、水銀使用製品となります。保管・輸送の際には適切なお取り扱いをお願い申し上げます。また、水銀使用製品を廃棄する際には特別管理廃棄物として管理を行い、適切な廃棄業者へ依頼していただくように合わせてお願い申し上げます。

なお、対象製品は、受注発注品につきご注文後のキャンセル、特定の防腐剤あるいは製造ロットの選択、納品後の返品はできかねますので予めご了承下さい。お客様には大変ご迷惑お掛けしておりますことを深くお詫び申し上げます。

※注 製品番号の頭文字が「MA」または「PA」で始まる Invitrogen 社製の抗体製品。対象製品の一覧につきましては別表1をご参照下さい。

製品番号	製品名
MA534962	HDAC1 Recombinant Rabbit Monoclonal Antibody
MA535044	MMP14 Recombinant Rabbit Monoclonal Antibody (ARC0211)
MA535046	CDK6 Recombinant Rabbit Monoclonal Antibody (ARC0224)
MA535047	CCR7 Recombinant Rabbit Monoclonal Antibody (ARC0231)
MA535048	MKP3 Recombinant Rabbit Monoclonal Antibody (ARC0237)
MA535049	HSC70 Recombinant Rabbit Monoclonal Antibody (ARC0258)
MA535050	HSP60 Recombinant Rabbit Monoclonal Antibody (ARC0260)
MA535051	IFNAR1 Recombinant Rabbit Monoclonal Antibody (ARC0262)
MA535052	LAMP2 Recombinant Rabbit Monoclonal Antibody (ARC0274)
MA535053	LRP1 Recombinant Rabbit Monoclonal Antibody (ARC0275)
MA535054	Lysozyme Recombinant Rabbit Monoclonal Antibody (ARC0276)
MA535055	GARS Recombinant Rabbit Monoclonal Antibody (ARC0514)
MA535056	CPT2 Recombinant Rabbit Monoclonal Antibody (ARC0516)
MA535057	CD8a Recombinant Rabbit Monoclonal Antibody (ARC0329)
MA535058	hnRNP K Recombinant Rabbit Monoclonal Antibody (ARC0512)
MA535059	SR-BI Recombinant Rabbit Monoclonal Antibody (ARC0334)
MA535060	IGF1 Recombinant Rabbit Monoclonal Antibody (ARC0505)
MA535061	Fibrillarin Recombinant Rabbit Monoclonal Antibody (ARC0506)
MA535062	ISL1 Recombinant Rabbit Monoclonal Antibody (ARC0511)
MA535064	HLA-DRA Recombinant Rabbit Monoclonal Antibody (ARC0518)
MA535066	AP2M1 Recombinant Rabbit Monoclonal Antibody (ARC0522)
MA535067	MMP11 Recombinant Rabbit Monoclonal Antibody (ARC0523)
MA535068	IL1RA Recombinant Rabbit Monoclonal Antibody (ARC0524)
MA535069	IRF5 Recombinant Rabbit Monoclonal Antibody (ARC0525)
MA535072	HSP27 Recombinant Rabbit Monoclonal Antibody (ARC0531)
MA535073	CD3z Recombinant Rabbit Monoclonal Antibody (ARC0533)
MA535074	MBP Recombinant Rabbit Monoclonal Antibody (ARC0535)
MA535075	IDE Recombinant Rabbit Monoclonal Antibody (ARC0537)
MA535076	CD206 Recombinant Rabbit Monoclonal Antibody (ARC0538)
MA535077	Complement C3 Recombinant Rabbit Monoclonal Antibody (ARC0541)
MA535078	Galectin 3 Recombinant Rabbit Monoclonal Antibody (ARC0542)
MA535080	FABP1 Recombinant Rabbit Monoclonal Antibody (ARC0545)
MA535082	ATP5A1 Recombinant Rabbit Monoclonal Antibody (ARC0549)
MA535083	Ku70 Recombinant Rabbit Monoclonal Antibody (ARC0551)
MA535085	G6PD Recombinant Rabbit Monoclonal Antibody (ARC0553)
MA535086	Annexin A2 Recombinant Rabbit Monoclonal Antibody (ARC0554)
MA535087	APE1 Recombinant Rabbit Monoclonal Antibody (ARC0556)
MA535088	Ferritin Light Chain Recombinant Rabbit Monoclonal Antibody (ARC0557)
MA535089	GPX4 Recombinant Rabbit Monoclonal Antibody (ARC0558)
MA535090	SERPINC1 Recombinant Rabbit Monoclonal Antibody (ARC0559)
MA535092	HSPA9 Recombinant Rabbit Monoclonal Antibody (ARC0566)
MA535093	CD55 Recombinant Rabbit Monoclonal Antibody (ARC0568)
MA535094	C1QBP Recombinant Rabbit Monoclonal Antibody (ARC0569)
MA535095	Angiopoietin 2 Recombinant Rabbit Monoclonal Antibody (ARC0571)
MA535096	GluR2 Recombinant Rabbit Monoclonal Antibody (ARC0572)
MA535097	MCM7 Recombinant Rabbit Monoclonal Antibody (ARC0573)
MA535098	Ephrin B2 Recombinant Rabbit Monoclonal Antibody (ARC0576)
MA535099	ACE Recombinant Rabbit Monoclonal Antibody (ARC0577)
MA535101	Myoglobin Recombinant Rabbit Monoclonal Antibody (ARC0582)
MA535102	CD47 Recombinant Rabbit Monoclonal Antibody (ARC0584)
MA535103	Cytokeratin 5 Recombinant Rabbit Monoclonal Antibody (ARC0585)
MA535104	delta Catenin Recombinant Rabbit Monoclonal Antibody (ARC0586)
MA535105	BMP4 Recombinant Rabbit Monoclonal Antibody (ARC0587)
MA535106	HLA-A Recombinant Rabbit Monoclonal Antibody (ARC0588)
MA535107	GM130 Recombinant Rabbit Monoclonal Antibody (ARC0589)
MA535109	DNA Ligase IV Recombinant Rabbit Monoclonal Antibody (ARC0595)
MA535110	CD23 Recombinant Rabbit Monoclonal Antibody (ARC0596)

MA535111	ENO1 Recombinant Rabbit Monoclonal Antibody (ARC0599)
MA535112	ADAR1 Recombinant Rabbit Monoclonal Antibody (ARC0602)
MA535113	MAD2 Recombinant Rabbit Monoclonal Antibody (ARC0603)
MA535114	GYPC Recombinant Rabbit Monoclonal Antibody (ARC0605)
MA535115	MCM3 Recombinant Rabbit Monoclonal Antibody (ARC0607)
MA535116	MATH1 Recombinant Rabbit Monoclonal Antibody (ARC0609)
MA535118	Mast Cell Chymase Recombinant Rabbit Monoclonal Antibody (ARC0614)
MA535119	ARF6 Recombinant Rabbit Monoclonal Antibody (ARC0617)
MA535120	FGF2 Recombinant Rabbit Monoclonal Antibody (8W6M5)
MA535121	CD41 Recombinant Rabbit Monoclonal Antibody (ARC0620)
MA535122	ALDH2 Recombinant Rabbit Monoclonal Antibody (ARC0623)
MA535123	LOX Recombinant Rabbit Monoclonal Antibody (ARC0624)
MA535126	BACE1 Recombinant Rabbit Monoclonal Antibody (ARC0696)
MA535128	RhoGDI Recombinant Rabbit Monoclonal Antibody (ARC0629)
MA535130	GABRB2 Recombinant Rabbit Monoclonal Antibody (ARC0631)
MA535131	Calreticulin Recombinant Rabbit Monoclonal Antibody (ARC0632)
MA535132	hnRNP A1 Recombinant Rabbit Monoclonal Antibody (ARC0633)
MA535133	FAP Recombinant Rabbit Monoclonal Antibody (ARC0634)
MA535134	eIF4A1 Recombinant Rabbit Monoclonal Antibody (ARC0636)
MA535136	MC1R Recombinant Rabbit Monoclonal Antibody (ARC0638)
MA535137	DIAPH1 Recombinant Rabbit Monoclonal Antibody (ARC0639)
MA535138	Monoamine Oxidase B Recombinant Rabbit Monoclonal Antibody (ARC0641)
MA535139	CtBP1 Recombinant Rabbit Monoclonal Antibody (ARC0642)
MA535140	CBS Recombinant Rabbit Monoclonal Antibody (ARC0643)
MA535141	DDB2 Recombinant Rabbit Monoclonal Antibody (ARC0644)
MA535142	Albumin Recombinant Rabbit Monoclonal Antibody (ARC0647)
MA535144	Kininogen 1 Recombinant Rabbit Monoclonal Antibody (ARC0653)
MA535145	Collagen X Recombinant Rabbit Monoclonal Antibody (ARC0659)
MA535146	MiTF Recombinant Rabbit Monoclonal Antibody (ARC0660)
MA535148	HBEGF Recombinant Rabbit Monoclonal Antibody (ARC0663)
MA535149	FBP1 Recombinant Rabbit Monoclonal Antibody (ARC0664)
MA535150	RHAMM Recombinant Rabbit Monoclonal Antibody (ARC0667)
MA535151	IGFBP-1 Recombinant Rabbit Monoclonal Antibody (ARC0671)
MA535152	CD127 Recombinant Rabbit Monoclonal Antibody (ARC0672)
MA535153	AZGP1 Recombinant Rabbit Monoclonal Antibody (ARC0673)
MA535154	Glypican 3 Recombinant Rabbit Monoclonal Antibody (ARC0675)
MA535194	Cathepsin B Recombinant Rabbit Monoclonal Antibody (ARC0395)
MA535155	SERCA2 ATPase Recombinant Rabbit Monoclonal Antibody (ARC0679)
MA535157	Phospho-NR4A1 (Ser351) Recombinant Rabbit Monoclonal Antibody (ARC0686)
MA535159	MERTK Recombinant Rabbit Monoclonal Antibody (ARC0229)
MA535158	beta Tubulin Recombinant Rabbit Monoclonal Antibody (ARC0203)
MA535160	VPS34 Recombinant Rabbit Monoclonal Antibody (ARC0286)
MA535161	RAB5 Recombinant Rabbit Monoclonal Antibody (ARC0297)
MA535162	GLUT2 Recombinant Rabbit Monoclonal Antibody (ARC0305)
MA535165	LGR5 Recombinant Rabbit Monoclonal Antibody (ARC0321)
MA535166	Dot1L Recombinant Rabbit Monoclonal Antibody (ARC0688)
MA535167	ACSS2 Recombinant Rabbit Monoclonal Antibody (ARC0690)
MA535168	NSD3 Recombinant Rabbit Monoclonal Antibody (ARC0707)
MA535169	TOP1 Recombinant Rabbit Monoclonal Antibody (ARC0708)
MA535170	AGR2 Recombinant Rabbit Monoclonal Antibody (ARC0709)
MA535172	MEK5 Recombinant Rabbit Monoclonal Antibody (ARC0711)
MA535173	PGK1 Recombinant Rabbit Monoclonal Antibody (ARC0700)
MA535174	PEBP1 Recombinant Rabbit Monoclonal Antibody (ARC0704)
MA535175	MKK7 Recombinant Rabbit Monoclonal Antibody (ARC0712)
MA535176	HDAC7 Recombinant Rabbit Monoclonal Antibody (ARC0713)
MA535177	HSPB8 Recombinant Rabbit Monoclonal Antibody (ARC0715)
MA535178	VWF Recombinant Rabbit Monoclonal Antibody (ARC0716)
MA535179	TK1 Recombinant Rabbit Monoclonal Antibody (ARC0719)
MA535180	S100A10 Recombinant Rabbit Monoclonal Antibody (ARC0720)

MA535181	PDHA1 Recombinant Rabbit Monoclonal Antibody (ARC0722)
MA535182	SDHA Recombinant Rabbit Monoclonal Antibody (ARC0726)
MA535185	ABCG1 Recombinant Rabbit Monoclonal Antibody (ARC0336)
MA535188	Protein APC Recombinant Rabbit Monoclonal Antibody (ARC0346)
MA535189	beta-3 Tubulin Recombinant Rabbit Monoclonal Antibody (ARC0456)
MA535190	BMI-1 Recombinant Rabbit Monoclonal Antibody (ARC0376)
MA535191	SDHB Recombinant Rabbit Monoclonal Antibody (ARC0733)
MA535192	BTK Recombinant Rabbit Monoclonal Antibody (ARC0485)
MA535193	Catenin alpha-1 Recombinant Rabbit Monoclonal Antibody (ARC0348)
MA535197	CD105 Recombinant Rabbit Monoclonal Antibody (ARC0446)
MA535198	CD11a Recombinant Rabbit Monoclonal Antibody (ARC0344)
MA535199	CD11b Recombinant Rabbit Monoclonal Antibody (ARC0414)
MA535200	CD14 Recombinant Rabbit Monoclonal Antibody (ARC0353)
MA535201	CD18 Recombinant Rabbit Monoclonal Antibody (ARC0400)
MA535202	CD34 Recombinant Rabbit Monoclonal Antibody (ARC0219)
MA535203	CD36 Recombinant Rabbit Monoclonal Antibody (ARC0461)
MA535205	CD154 (CD40 Ligand) Recombinant Rabbit Monoclonal Antibody (ARC0399)
MA535207	CD46 Recombinant Rabbit Monoclonal Antibody (ARC0428)
MA535209	CD79a Recombinant Rabbit Monoclonal Antibody (ARC0482)
MA535210	CD80 Recombinant Rabbit Monoclonal Antibody (ARC0384)
MA535211	CD86 Recombinant Rabbit Monoclonal Antibody (ARC0333)
MA535213	CDC42 Recombinant Rabbit Monoclonal Antibody (ARC0484)
MA535214	CDK5 Recombinant Rabbit Monoclonal Antibody (ARC0223)
MA535215	CDX2 Recombinant Rabbit Monoclonal Antibody (ARC0450)
MA535216	Cullin 1 Recombinant Rabbit Monoclonal Antibody (ARC0436)
MA535219	Cyclin B1 Recombinant Rabbit Monoclonal Antibody (ARC0474)
MA535217	CXCR4 Recombinant Rabbit Monoclonal Antibody (ARC0381)
MA535218	Cyclin A2 Recombinant Rabbit Monoclonal Antibody (ARC0359)
MA535220	Cytokeratin 14 Recombinant Rabbit Monoclonal Antibody (ARC0351)
MA535222	Cytokeratin 20 Recombinant Rabbit Monoclonal Antibody (ARC0288)
MA535223	Doublecortin Recombinant Rabbit Monoclonal Antibody (ARC0380)
MA535224	TRAIL-R2 (DR5) Recombinant Rabbit Monoclonal Antibody (ARC0406)
MA535225	eIF4E Recombinant Rabbit Monoclonal Antibody (ARC0447)
MA535226	4EBP1 Recombinant Rabbit Monoclonal Antibody (ARC0374)
MA535229	Ezrin Recombinant Rabbit Monoclonal Antibody (ARC0392)
MA535230	FADD Recombinant Rabbit Monoclonal Antibody (ARC0415)
MA535231	FGFR2 Recombinant Rabbit Monoclonal Antibody (ARC0480)
MA535232	FGFR3 Recombinant Rabbit Monoclonal Antibody (ARC0398)
MA535234	GAP43 Recombinant Rabbit Monoclonal Antibody (ARC0246)
MA535236	GBA Recombinant Rabbit Monoclonal Antibody (ARC0500)
MA535238	GRB2 Recombinant Rabbit Monoclonal Antibody (ARC0430)
MA535239	GSK3 alpha Recombinant Rabbit Monoclonal Antibody (ARC0462)
MA535240	GSTP1 Recombinant Rabbit Monoclonal Antibody (ARC0439)
MA535243	ITGA2 Recombinant Rabbit Monoclonal Antibody (ARC0457)
MA535245	CD61 (Integrin beta 3) Recombinant Rabbit Monoclonal Antibody (ARC0460)
MA535246	JNK3 Recombinant Rabbit Monoclonal Antibody (ARC0366)
MA535247	Leptin Receptor Recombinant Rabbit Monoclonal Antibody (ARC0454)
MA535248	METTL3 Recombinant Rabbit Monoclonal Antibody (ARC0487)
MA535249	MMP2 Recombinant Rabbit Monoclonal Antibody (ARC0432)
MA535252	NEFM Recombinant Rabbit Monoclonal Antibody (ARC0396)
MA535254	NGFR Recombinant Rabbit Monoclonal Antibody (ARC0284)
MA535255	NMDAR2A Recombinant Rabbit Monoclonal Antibody (ARC0410)
MA535256	NOTCH1 Recombinant Rabbit Monoclonal Antibody (ARC0285)
MA535257	P glycoprotein Recombinant Rabbit Monoclonal Antibody (ARC0470)
MA535258	DJ-1 Recombinant Rabbit Monoclonal Antibody (ARC0354)
MA535259	Parvalbumin Recombinant Rabbit Monoclonal Antibody (ARC0385)
MA535260	PAX6 Recombinant Rabbit Monoclonal Antibody (ARC0367)
MA535261	Paxillin Recombinant Rabbit Monoclonal Antibody (ARC0368)
MA535265	RXRA Recombinant Rabbit Monoclonal Antibody (ARC0468)

MA535266	RhoA Recombinant Rabbit Monoclonal Antibody (ARC0372)
MA535267	S100A4 Recombinant Rabbit Monoclonal Antibody (ARC0475)
MA535270	SMAD5 Recombinant Rabbit Monoclonal Antibody (ARC0448)
MA535271	SOX2 Recombinant Rabbit Monoclonal Antibody (ARC0449)
MA535272	SUMO1 Recombinant Rabbit Monoclonal Antibody (ARC0215)
MA535273	TDP-43 Recombinant Rabbit Monoclonal Antibody (ARC0492)
MA535274	TLR2 Recombinant Rabbit Monoclonal Antibody (ARC0433)
MA535275	TLR7 Recombinant Rabbit Monoclonal Antibody (ARC0401)
MA535276	TNFR2 Recombinant Rabbit Monoclonal Antibody (ARC0397)
MA535277	TNFAIP3 Recombinant Rabbit Monoclonal Antibody (ARC0355)
MA535278	Transferrin Recombinant Rabbit Monoclonal Antibody (ARC0338)
MA535279	VCAM-1 (CD106) Recombinant Rabbit Monoclonal Antibody (ARC0312)
MA535280	VEGF Receptor 1 Recombinant Rabbit Monoclonal Antibody (ARC0363)
MA535281	WNT5A Recombinant Rabbit Monoclonal Antibody (ARC0405)
MA535283	EpCAM Recombinant Rabbit Monoclonal Antibody (ARC0277)
MA535285	H4K5ac Recombinant Rabbit Monoclonal Antibody (ARC0002)
MA535286	HMGB1 Recombinant Rabbit Monoclonal Antibody (ARC0001)
MA535287	Prohibitin Recombinant Rabbit Monoclonal Antibody (ARC0007)
MA535288	PDGFRB Recombinant Rabbit Monoclonal Antibody (ARC0009)
MA535290	AIF Recombinant Rabbit Monoclonal Antibody (ARC0015)
MA535289	HIF-1 beta Recombinant Rabbit Monoclonal Antibody (ARC0010)
MA535291	C/EBP beta Recombinant Rabbit Monoclonal Antibody (ARC0017)
MA535293	TSC2 Recombinant Rabbit Monoclonal Antibody (ARC0019)
MA535294	CHK2 Recombinant Rabbit Monoclonal Antibody (ARC0022)
MA535295	TCF7L2 Recombinant Rabbit Monoclonal Antibody (ARC0027)
MA535296	Caspase 8 Recombinant Rabbit Monoclonal Antibody (ARC0028)
MA535297	RARA Recombinant Rabbit Monoclonal Antibody (ARC0030)
MA535298	BRG1 Recombinant Rabbit Monoclonal Antibody (ARC0035)
MA535299	Caspase 6 Recombinant Rabbit Monoclonal Antibody (ARC0038)
MA535300	STAT1 Recombinant Rabbit Monoclonal Antibody (ARC0042)
MA535301	MEK1 Recombinant Rabbit Monoclonal Antibody (ARC0044)
MA535303	TRIM28 Recombinant Rabbit Monoclonal Antibody (ARC0047)
MA535304	HSP90 beta Recombinant Rabbit Monoclonal Antibody (ARC0053)
MA535305	SOD2 Recombinant Rabbit Monoclonal Antibody (ARC0055)
MA535306	E2F1 Recombinant Rabbit Monoclonal Antibody (ARC0058)
MA535308	FAS Recombinant Rabbit Monoclonal Antibody (ARC0061)
MA535309	Glucocorticoid Receptor Recombinant Rabbit Monoclonal Antibody (ARC0062)
MA535310	NQO1 Recombinant Rabbit Monoclonal Antibody (ARC0065)
MA535311	CTCF Recombinant Rabbit Monoclonal Antibody (ARC0067)
MA535312	ATG3 Recombinant Rabbit Monoclonal Antibody (ARC0073)
MA535313	BAD Recombinant Rabbit Monoclonal Antibody (ARC0074)
MA535314	NFATC1 Recombinant Rabbit Monoclonal Antibody (ARC0076)
MA535315	Drosha Recombinant Rabbit Monoclonal Antibody (ARC0077)
MA535316	ETS1 Recombinant Rabbit Monoclonal Antibody (ARC0082)
MA535317	ATG7 Recombinant Rabbit Monoclonal Antibody (ARC0083)
MA535318	NFkB p52/p100 Recombinant Rabbit Monoclonal Antibody (ARC0084)
MA535321	PAK1 Recombinant Rabbit Monoclonal Antibody (ARC0087)
MA535322	PARP1 (cleaved Asp214) Recombinant Rabbit Monoclonal Antibody (ARC0091)
MA535324	HuR Recombinant Rabbit Monoclonal Antibody (ARC0101)
MA535328	GATA3 Recombinant Rabbit Monoclonal Antibody (ARC0115)
MA535329	H3R2me1 Recombinant Rabbit Monoclonal Antibody (ARC0124)
MA535330	PML Recombinant Rabbit Monoclonal Antibody (ARC0125)
MA535331	SP1 Recombinant Rabbit Monoclonal Antibody (ARC0128)
MA535334	DNMT3A Recombinant Rabbit Monoclonal Antibody (ARC0138)
MA535335	Caspase 3 p12 Recombinant Rabbit Monoclonal Antibody (ARC0143)
MA535338	PPAR gamma Recombinant Rabbit Monoclonal Antibody (ARC0155)
MA535339	ATG5 Recombinant Rabbit Monoclonal Antibody (ARC0156)
MA535341	Cathepsin D Recombinant Rabbit Monoclonal Antibody (ARC0160)
MA535343	cIAP1 Recombinant Rabbit Monoclonal Antibody (ARC0168)

MA535344	FANCD2 Recombinant Rabbit Monoclonal Antibody (ARC0172)
MA535345	BCL-2 Recombinant Rabbit Monoclonal Antibody (ARC0173)
MA535346	IKK alpha Recombinant Rabbit Monoclonal Antibody (ARC0174)
MA535347	PIM1 Recombinant Rabbit Monoclonal Antibody (ARC0175)
MA535349	VDAC Recombinant Rabbit Monoclonal Antibody (ARC0187)
MA535351	HDAC9 Recombinant Rabbit Monoclonal Antibody (ARC0735)
MA535352	N-WASP Recombinant Rabbit Monoclonal Antibody (ARC0737)
MA535354	HEC1 Recombinant Rabbit Monoclonal Antibody (ARC0741)
MA535356	ROCK2 Recombinant Rabbit Monoclonal Antibody (ARC0744)
MA535357	PRDX3 Recombinant Rabbit Monoclonal Antibody (ARC0747)
MA535358	TGFBI Recombinant Rabbit Monoclonal Antibody (ARC0757)
MA535359	POLB Recombinant Rabbit Monoclonal Antibody (ARC0752)
MA535360	SLC7A11 Recombinant Rabbit Monoclonal Antibody (ARC0753)
MA535361	REST Recombinant Rabbit Monoclonal Antibody (ARC0755)
MA535362	gamma Synuclein Recombinant Rabbit Monoclonal Antibody (ARC0762)
MA535363	RAD51 Recombinant Rabbit Monoclonal Antibody (ARC0764)
MA535366	NRF1 Recombinant Rabbit Monoclonal Antibody (ARC0768)
MA535367	RPA70 Recombinant Rabbit Monoclonal Antibody (ARC0773)
MA535368	RRM2 Recombinant Rabbit Monoclonal Antibody (ARC0775)
MA535369	USP7 Recombinant Rabbit Monoclonal Antibody (ARC0777)
MA535370	PSMA1 Recombinant Rabbit Monoclonal Antibody (ARC0780)
MA535371	IKK epsilon Recombinant Rabbit Monoclonal Antibody (ARC0781)
MA535372	PEDF Recombinant Rabbit Monoclonal Antibody (ARC0784)
MA535373	Coilin Recombinant Rabbit Monoclonal Antibody (ARC0785)
MA535374	PTBP1 Recombinant Rabbit Monoclonal Antibody (ARC0783)
MA535375	SFPQ Recombinant Rabbit Monoclonal Antibody (ARC0788)
MA535376	OGT Recombinant Rabbit Monoclonal Antibody (ARC0790)
MA535377	TLE1 Recombinant Rabbit Monoclonal Antibody (ARC0793)
MA535378	SMARCA5 Recombinant Rabbit Monoclonal Antibody (ARC0795)
MA535379	RNF2 Recombinant Rabbit Monoclonal Antibody (ARC0802)
MA535380	IKAROS Recombinant Rabbit Monoclonal Antibody (ARC0803)
MA535381	SIRT6 Recombinant Rabbit Monoclonal Antibody (ARC0808)
MA535382	Ub-K48 Recombinant Rabbit Monoclonal Antibody (ARC0811)
MA535383	Rb Recombinant Rabbit Monoclonal Antibody (ARC0810)
MA535384	RbAp48 Recombinant Rabbit Monoclonal Antibody (ARC0813)
MA535385	NOX4 Recombinant Rabbit Monoclonal Antibody (ARC0815)
MA535386	PSMD4 Recombinant Rabbit Monoclonal Antibody (ARC0818)
MA535387	PFKM Recombinant Rabbit Monoclonal Antibody (ARC0819)
MA535388	RPS19 Recombinant Rabbit Monoclonal Antibody (ARC0820)
MA535389	SNRPA Recombinant Rabbit Monoclonal Antibody (ARC0823)
MA535390	FBXO32 Recombinant Rabbit Monoclonal Antibody (ARC0830)
MA535391	HIF1AN Recombinant Rabbit Monoclonal Antibody (ARC0826)
MA535393	NCF2 Recombinant Rabbit Monoclonal Antibody (ARC0828)
MA535394	P2Y6 Recombinant Rabbit Monoclonal Antibody (ARC0831)
MA535395	PHD1 Recombinant Rabbit Monoclonal Antibody (ARC0230)
MA535396	Filamin A Recombinant Rabbit Monoclonal Antibody (ARC0242)
MA535397	HP1 alpha Recombinant Rabbit Monoclonal Antibody (ARC0244)
MA535398	PUMA alpha Recombinant Rabbit Monoclonal Antibody (ARC0247)
MA535399	Angiopoietin 1 Recombinant Rabbit Monoclonal Antibody (ARC0248)
MA535400	M6PR Recombinant Rabbit Monoclonal Antibody (ARC0263)
MA535401	Cytokeratin 17 Recombinant Rabbit Monoclonal Antibody (ARC0271)
MA535402	NDUFA13 Recombinant Rabbit Monoclonal Antibody (ARC0833)
MA535403	Nudel Recombinant Rabbit Monoclonal Antibody (ARC0835)
MA535404	NONO Recombinant Rabbit Monoclonal Antibody (ARC0836)
MA535405	RACK1 Recombinant Rabbit Monoclonal Antibody (ARC0837)
MA535406	XPO5 Recombinant Rabbit Monoclonal Antibody (ARC0838)
MA535407	BAF57 Recombinant Rabbit Monoclonal Antibody (ARC0839)
MA535408	TREX1 Recombinant Rabbit Monoclonal Antibody (ARC0841)
MA535409	MYLK Recombinant Rabbit Monoclonal Antibody (ARC0842)

MA535410	MTCO2 Recombinant Rabbit Monoclonal Antibody (ARC0844)
MA535412	RAD50 Recombinant Rabbit Monoclonal Antibody (ARC0854)
MA535413	AKR1C3 Recombinant Rabbit Monoclonal Antibody (ARC0857)
MA535414	SAM68 Recombinant Rabbit Monoclonal Antibody (ARC0858)
MA535416	DKK3 Recombinant Rabbit Monoclonal Antibody (ARC0861)
MA535417	HTRA2 Recombinant Rabbit Monoclonal Antibody (ARC0864)
MA535420	SGK1 Recombinant Rabbit Monoclonal Antibody (ARC0868)
MA535421	MTH1 Recombinant Rabbit Monoclonal Antibody (ARC0869)
MA535422	ERK5 Recombinant Rabbit Monoclonal Antibody (ARC0871)
MA535423	Nogo Receptor Recombinant Rabbit Monoclonal Antibody (ARC0872)
MA535424	PRDX5 Recombinant Rabbit Monoclonal Antibody (ARC0873)
MA535425	G3BP1 Recombinant Rabbit Monoclonal Antibody (ARC0875)
MA535426	TRAP1 Recombinant Rabbit Monoclonal Antibody (ARC0876)
MA535428	PI3K p55 gamma Recombinant Rabbit Monoclonal Antibody (ARC0320)
MA535429	Cyclin D3 Recombinant Rabbit Monoclonal Antibody (ARC0325)
MA535432	TRPC1 Recombinant Rabbit Monoclonal Antibody (ARC0884)
MA535433	SKP2 Recombinant Rabbit Monoclonal Antibody (ARC0886)
MA535434	NDRG1 Recombinant Rabbit Monoclonal Antibody (ARC0887)
MA535435	CD49d (Integrin alpha 4) Recombinant Rabbit Monoclonal Antibody (ARC0888)
MA535436	Calpain 2 Recombinant Rabbit Monoclonal Antibody (ARC0891)
MA535437	Cyclin H Recombinant Rabbit Monoclonal Antibody (ARC0893)
MA535438	RBPJ Recombinant Rabbit Monoclonal Antibody (ARC0894)
MA535439	JAB1 Recombinant Rabbit Monoclonal Antibody (ARC0895)
MA535440	CD59 Recombinant Rabbit Monoclonal Antibody (ARC0896)
MA535442	FAAH Recombinant Rabbit Monoclonal Antibody (ARC0903)
MA535443	MAOA Recombinant Rabbit Monoclonal Antibody (ARC0900)
MA535444	LPL Recombinant Rabbit Monoclonal Antibody (ARC0904)
MA535445	PAX7 Recombinant Rabbit Monoclonal Antibody (ARC0914)
MA535446	XRCC1 Recombinant Rabbit Monoclonal Antibody (ARC0915)
MA535447	SETD8 Recombinant Rabbit Monoclonal Antibody (ARC0916)
MA535448	GLUT3 Recombinant Rabbit Monoclonal Antibody (ARC0917)
MA535449	TrkA Recombinant Rabbit Monoclonal Antibody (ARC0906)
MA535452	gamma Catenin Recombinant Rabbit Monoclonal Antibody (ARC0909)
MA535453	DDAH2 Recombinant Rabbit Monoclonal Antibody (ARC0910)
MA535454	MSMB Recombinant Rabbit Monoclonal Antibody (ARC0912)
MA535455	PGAM1 Recombinant Rabbit Monoclonal Antibody (ARC0913)
MA535456	CD138 Recombinant Rabbit Monoclonal Antibody (ARC0918)
MA535457	ESRRA Recombinant Rabbit Monoclonal Antibody (ARC0919)
MA535459	Emerin Recombinant Rabbit Monoclonal Antibody (ARC0921)
MA535458	Apolipoprotein B Recombinant Rabbit Monoclonal Antibody (ARC0920)
MA535460	SFRP4 Recombinant Rabbit Monoclonal Antibody (ARC0923)
MA535461	Aquaporin 1 Recombinant Rabbit Monoclonal Antibody (ARC0925)
MA535462	tPA Recombinant Rabbit Monoclonal Antibody (ARC0928)
MA535463	CPS1 Recombinant Rabbit Monoclonal Antibody (ARC0929)
MA535467	RUVBL2 Recombinant Rabbit Monoclonal Antibody (ARC0933)
MA535468	CD83 Recombinant Rabbit Monoclonal Antibody (ARC0935)
MA535469	VAMP2 Recombinant Rabbit Monoclonal Antibody (ARC0936)
MA535470	CD26 Recombinant Rabbit Monoclonal Antibody (ARC0939)
MA535471	hnRNP U Recombinant Rabbit Monoclonal Antibody (ARC0942)
MA535472	RRM1 Recombinant Rabbit Monoclonal Antibody (ARC0943)
MA535474	TRPM8 Recombinant Rabbit Monoclonal Antibody (ARC0947)
MA535475	SMARCC1 Recombinant Rabbit Monoclonal Antibody (ARC0948)
MA535477	DDR2 Recombinant Rabbit Monoclonal Antibody (ARC0958)
MA535478	ALY Recombinant Rabbit Monoclonal Antibody (ARC0959)
MA535479	COX1 Recombinant Rabbit Monoclonal Antibody (ARC0960)
MA535480	PRDX2 Recombinant Rabbit Monoclonal Antibody (ARC0963)
MA535482	GABARAP Recombinant Rabbit Monoclonal Antibody (ARC1057)
MA535484	Histone H1.0 Recombinant Rabbit Monoclonal Antibody (ARC1059)
MA535485	HMGA1 Recombinant Rabbit Monoclonal Antibody (ARC1060)

MA535486	AMPK beta-1 Recombinant Rabbit Monoclonal Antibody (ARC1061)
MA535487	Calcineurin A Recombinant Rabbit Monoclonal Antibody (ARC0970)
MA535489	LYVE1 Recombinant Rabbit Monoclonal Antibody (ARC0974)
MA535490	DYNLL1 Recombinant Rabbit Monoclonal Antibody (ARC0975)
MA535491	FKBP4 Recombinant Rabbit Monoclonal Antibody (ARC0977)
MA535492	PP1 beta Recombinant Rabbit Monoclonal Antibody (ARC0981)
MA535493	UQCRC2 Recombinant Rabbit Monoclonal Antibody (ARC0982)
MA535494	DCLK1 Recombinant Rabbit Monoclonal Antibody (ARC0984)
MA535495	Laminin beta-1 Recombinant Rabbit Monoclonal Antibody (ARC0985)
MA535496	RAN Recombinant Rabbit Monoclonal Antibody (ARC0986)
MA535498	14-3-3 sigma Recombinant Rabbit Monoclonal Antibody (ARC0988)
MA535499	Stathmin 1 Recombinant Rabbit Monoclonal Antibody (ARC0989)
MA535501	NEDD4 Recombinant Rabbit Monoclonal Antibody (ARC0992)
MA535502	Tetranectin Recombinant Rabbit Monoclonal Antibody (ARC0993)
MA535503	TOP2A Recombinant Rabbit Monoclonal Antibody (ARC0994)
MA535505	UBC9 Recombinant Rabbit Monoclonal Antibody (ARC0996)
MA535507	alpha Adaptin Recombinant Rabbit Monoclonal Antibody (ARC0998)
MA535508	Carbonic Anhydrase I Recombinant Rabbit Monoclonal Antibody (ARC1062)
MA535509	DKC1 Recombinant Rabbit Monoclonal Antibody (ARC1063)
MA535510	CRMP2 Recombinant Rabbit Monoclonal Antibody (ARC1123)
MA535511	Apolipoprotein L1 Recombinant Rabbit Monoclonal Antibody (ARC1124)
MA535512	eIF5A Recombinant Rabbit Monoclonal Antibody (ARC0999)
MA535513	NuMA Recombinant Rabbit Monoclonal Antibody (ARC1000)
MA535514	Thioredoxin 2 Recombinant Rabbit Monoclonal Antibody (ARC1003)
MA535516	IP3 Receptor 1 Recombinant Rabbit Monoclonal Antibody (ARC1007)
MA535517	GRK2 Recombinant Rabbit Monoclonal Antibody (ARC1013)
MA535518	USP10 Recombinant Rabbit Monoclonal Antibody (ARC1015)
MA535520	WARS Recombinant Rabbit Monoclonal Antibody (ARC1018)
MA535521	LEF1 Recombinant Rabbit Monoclonal Antibody (ARC1019)
MA535522	GCLC Recombinant Rabbit Monoclonal Antibody (ARC1065)
MA535523	PRMT1 Recombinant Rabbit Monoclonal Antibody (ARC1066)
MA535524	Arp3 Recombinant Rabbit Monoclonal Antibody (ARC1067)
MA535525	CK2 beta Recombinant Rabbit Monoclonal Antibody (ARC1069)
MA535527	Aconitase 2 Recombinant Rabbit Monoclonal Antibody (ARC1072)
MA535528	Protein C Recombinant Rabbit Monoclonal Antibody (ARC1023)
MA535530	BMP6 Recombinant Rabbit Monoclonal Antibody (ARC1025)
MA535531	FMRP Recombinant Rabbit Monoclonal Antibody (ARC1026)
MA535532	U2AF2 Recombinant Rabbit Monoclonal Antibody (ARC1029)
MA535534	DHX9 Recombinant Rabbit Monoclonal Antibody (ARC1033)
MA535535	GSR Recombinant Rabbit Monoclonal Antibody (ARC1034)
MA535536	ACADM Recombinant Rabbit Monoclonal Antibody (ARC1035)
MA535538	CDC37 Recombinant Rabbit Monoclonal Antibody (ARC1040)
MA535539	DDAH1 Recombinant Rabbit Monoclonal Antibody (ARC1043)
MA535541	CD104 Recombinant Rabbit Monoclonal Antibody (ARC1046)
MA535542	Histone H2A.Z Recombinant Rabbit Monoclonal Antibody (ARC1048)
MA535545	IGFBP7 Recombinant Rabbit Monoclonal Antibody (ARC1052)
MA535546	Staufen Recombinant Rabbit Monoclonal Antibody (ARC1054)
MA535547	MTA2 Recombinant Rabbit Monoclonal Antibody (ARC1056)
MA535548	MAD2L2 Recombinant Rabbit Monoclonal Antibody (ARC1126)
MA535549	GST Omega 1 Recombinant Rabbit Monoclonal Antibody (ARC1128)
MA535550	TPPP Recombinant Rabbit Monoclonal Antibody (ARC1129)
MA535551	PDHB Recombinant Rabbit Monoclonal Antibody (ARC1074)
MA535552	FDFT1 Recombinant Rabbit Monoclonal Antibody (ARC1077)
MA535554	KIF11 Recombinant Rabbit Monoclonal Antibody (ARC1080)
MA535555	TACR1 Recombinant Rabbit Monoclonal Antibody (ARC1088)
MA535556	K-cadherin Recombinant Rabbit Monoclonal Antibody (ARC1093)
MA535558	SMC1 Recombinant Rabbit Monoclonal Antibody (ARC1095)
MA535559	RSK1 Recombinant Rabbit Monoclonal Antibody (ARC1096)
MA535560	TCP-1 beta Recombinant Rabbit Monoclonal Antibody (ARC1097)

MA535561	BMAL1 Recombinant Rabbit Monoclonal Antibody (ARC1101)
MA535562	NUDC Recombinant Rabbit Monoclonal Antibody (ARC1103)
MA535564	TrxR1 Recombinant Rabbit Monoclonal Antibody (ARC1106)
MA535565	Gephyrin Recombinant Rabbit Monoclonal Antibody (ARC1109)
MA535567	Pan Ras Recombinant Rabbit Monoclonal Antibody (ARC1112)
MA535568	PKD2 Recombinant Rabbit Monoclonal Antibody (ARC1113)
MA535569	TNNI2 Recombinant Rabbit Monoclonal Antibody (ARC1114)
MA535570	GTPBP4 Recombinant Rabbit Monoclonal Antibody (ARC1115)
MA535571	PIAS1 Recombinant Rabbit Monoclonal Antibody (ARC1116)
MA535572	KDM5A Recombinant Rabbit Monoclonal Antibody (ARC1120)
MA535573	BDKRB2 Recombinant Rabbit Monoclonal Antibody (ARC1121)
MA535574	Perilipin 1 Recombinant Rabbit Monoclonal Antibody (ARC1122)
MA535575	Cytokeratin 7 Recombinant Rabbit Monoclonal Antibody (ARC1267)
MA535576	MUC2 Recombinant Rabbit Monoclonal Antibody (ARC1012)
MA535577	p38 MAPK Recombinant Rabbit Monoclonal Antibody (ARC0201)
MA535581	beta Tubulin 2 Recombinant Rabbit Monoclonal Antibody (ARC0227)
MA535578	HSP70 Recombinant Rabbit Monoclonal Antibody (ARC0208)
MA535583	HAPLN1 Recombinant Rabbit Monoclonal Antibody (ARC0234)
MA535585	DNAJC15 Recombinant Rabbit Monoclonal Antibody (ARC0249)
MA535587	Histone H3.3 Recombinant Rabbit Monoclonal Antibody (ARC0255)
MA535588	Calnexin Recombinant Rabbit Monoclonal Antibody (ARC0648)
MA535591	Cyclophilin B Recombinant Rabbit Monoclonal Antibody (ARC0290)
MA535592	CD81 Recombinant Rabbit Monoclonal Antibody (ARC0615)
MA535593	JNK1/JNK2/JNK3 Recombinant Rabbit Monoclonal Antibody (ARC0291)
MA535594	RPS3 Recombinant Rabbit Monoclonal Antibody (ARC0302)
MA535595	BDNF Recombinant Rabbit Monoclonal Antibody (ARC0303)
MA535596	TSPO Recombinant Rabbit Monoclonal Antibody (ARC0308)
MA535597	Vitronectin Recombinant Rabbit Monoclonal Antibody (ARC0313)
MA535598	Calmodulin Recombinant Rabbit Monoclonal Antibody (ARC0316)
MA535601	ASH2L Recombinant Rabbit Monoclonal Antibody (ARC0326)
MA535602	RAB9 Recombinant Rabbit Monoclonal Antibody (ARC0332)
MA535604	Synuclein alpha/beta Recombinant Rabbit Monoclonal Antibody (ARC1181)
MA535630	Cullin 4A Recombinant Rabbit Monoclonal Antibody (ARC1237)
MA535605	Pan-cadherin Recombinant Rabbit Monoclonal Antibody (ARC1135)
MA535606	GRP78 Recombinant Rabbit Monoclonal Antibody (ARC1141)
MA535607	Cytochrome C Recombinant Rabbit Monoclonal Antibody (ARC1153)
MA535609	Moesin Recombinant Rabbit Monoclonal Antibody (ARC1258)
MA535610	Clathrin Heavy Chain Recombinant Rabbit Monoclonal Antibody (ARC1228)
MA535611	PRDX1 Recombinant Rabbit Monoclonal Antibody (ARC1274)
MA535612	TIMP1 Recombinant Rabbit Monoclonal Antibody (ARC1195)
MA535613	Myosin Heavy Chain Recombinant Rabbit Monoclonal Antibody (ARC1170)
MA535614	Cathepsin H/K/L/V Recombinant Rabbit Monoclonal Antibody (ARC1276)
MA535615	PI3K p85 alpha Recombinant Rabbit Monoclonal Antibody (ARC1174)
MA535616	OGG1 Recombinant Rabbit Monoclonal Antibody (ARC1233)
MA535618	Lamin B2 Recombinant Rabbit Monoclonal Antibody (ARC1252)
MA535619	MCM5 Recombinant Rabbit Monoclonal Antibody (ARC1239)
MA535620	KPNA2 Recombinant Rabbit Monoclonal Antibody (ARC1227)
MA535622	DDB1 Recombinant Rabbit Monoclonal Antibody (ARC1278)
MA535623	p14ARF Recombinant Rabbit Monoclonal Antibody (ARC1144)
MA535624	HSP90 Recombinant Rabbit Monoclonal Antibody (ARC1169)
MA535625	STAT5 alpha/beta Recombinant Rabbit Monoclonal Antibody (ARC1215)
MA535626	EGFR Recombinant Rabbit Monoclonal Antibody (ARC1139)
MA535627	Furin Recombinant Rabbit Monoclonal Antibody (ARC1221)
MA535628	JNK1/JNK3 Recombinant Rabbit Monoclonal Antibody (ARC1163)
MA535629	EEA1 Recombinant Rabbit Monoclonal Antibody (ARC1250)
MA535631	SUMO2/SUMO3 Recombinant Rabbit Monoclonal Antibody (ARC1186)
MA535632	CYP17A1 Recombinant Rabbit Monoclonal Antibody (ARC1257)
MA535633	PSME1 Recombinant Rabbit Monoclonal Antibody (ARC1254)
MA535634	UPF1 Recombinant Rabbit Monoclonal Antibody (ARC1268)

MA535636	RAC1/RAC2/RAC3 Recombinant Rabbit Monoclonal Antibody (ARC1165)
MA535666	WIF1 Recombinant Rabbit Monoclonal Antibody (ARC1161)
MA535639	PSMA3 Recombinant Rabbit Monoclonal Antibody (ARC1234)
MA535640	Cyclophilin 40 Recombinant Rabbit Monoclonal Antibody (ARC1260)
MA535641	MAD1 Recombinant Rabbit Monoclonal Antibody (ARC1188)
MA535642	GTH Recombinant Rabbit Monoclonal Antibody (ARC1200)
MA535643	PAR2 Recombinant Rabbit Monoclonal Antibody (ARC1246)
MA535644	PER2 Recombinant Rabbit Monoclonal Antibody (ARC1134)
MA535645	CBX4 Recombinant Rabbit Monoclonal Antibody (ARC1201)
MA535646	Alpha A Crystallin Recombinant Rabbit Monoclonal Antibody (ARC1245)
MA535647	ATP1B1 Recombinant Rabbit Monoclonal Antibody (ARC1232)
MA535648	EBP1 Recombinant Rabbit Monoclonal Antibody (ARC1281)
MA535651	ROR2 Recombinant Rabbit Monoclonal Antibody (ARC1197)
MA535654	ICOS Recombinant Rabbit Monoclonal Antibody (ARC1148)
MA535655	LDHB Recombinant Rabbit Monoclonal Antibody (ARC1217)
MA535656	WASP Recombinant Rabbit Monoclonal Antibody (ARC1204)
MA535657	BHMT Recombinant Rabbit Monoclonal Antibody (ARC1249)
MA535658	p47phox Recombinant Rabbit Monoclonal Antibody (ARC1206)
MA535659	Apolipoprotein M Recombinant Rabbit Monoclonal Antibody (ARC1210)
MA535661	MBNL1 Recombinant Rabbit Monoclonal Antibody (ARC1199)
MA535662	CHRM2 Recombinant Rabbit Monoclonal Antibody (ARC1225)
MA535664	Pannexin 1 Recombinant Rabbit Monoclonal Antibody (ARC1207)
MA535665	GCK Recombinant Rabbit Monoclonal Antibody (ARC1226)
MA535668	GluD1 Recombinant Rabbit Monoclonal Antibody (ARC1216)
MA535669	RUVBL1 Recombinant Rabbit Monoclonal Antibody (ARC1247)
MA535670	PLCG2 Recombinant Rabbit Monoclonal Antibody (ARC1176)
MA535671	RAD52 Recombinant Rabbit Monoclonal Antibody (ARC1183)
MA535674	UBE2C Recombinant Rabbit Monoclonal Antibody (ARC1166)
MA535675	ABO Recombinant Rabbit Monoclonal Antibody (ARC1251)
MA535676	DP1 Recombinant Rabbit Monoclonal Antibody (ARC1202)
MA535677	IRF4 Recombinant Rabbit Monoclonal Antibody (ARC1261)
MA535678	KAP3 Recombinant Rabbit Monoclonal Antibody (ARC1211)
MA535679	TEF1 Recombinant Rabbit Monoclonal Antibody (ARC1157)
MA535680	DLD Recombinant Rabbit Monoclonal Antibody (ARC1224)
MA535681	ANGPTL3 Recombinant Rabbit Monoclonal Antibody (ARC1214)
MA535682	Desmoglein 3 Recombinant Rabbit Monoclonal Antibody (ARC1190)
MA535685	Syntenin 1 Recombinant Rabbit Monoclonal Antibody (ARC1410)
MA535686	HSPE1 Recombinant Rabbit Monoclonal Antibody (ARC1411)
MA535687	NME2 Recombinant Rabbit Monoclonal Antibody (ARC1412)
MA535688	DDX17 Recombinant Rabbit Monoclonal Antibody (ARC1413)
MA535689	TSG101 Recombinant Rabbit Monoclonal Antibody (ARC0853)
MA535690	FGF1 Recombinant Rabbit Monoclonal Antibody (ARC1414)
MA535691	HSPA2 Recombinant Rabbit Monoclonal Antibody (ARC1415)
MA535692	JMJD2B Recombinant Rabbit Monoclonal Antibody (ARC1416)
MA535693	NEK2 Recombinant Rabbit Monoclonal Antibody (ARC1418)
MA535694	ADAM15 Recombinant Rabbit Monoclonal Antibody (ARC1419)
MA535695	PHOX2B Recombinant Rabbit Monoclonal Antibody (ARC1420)
MA535696	DCP1A Recombinant Rabbit Monoclonal Antibody (ARC1421)
MA535697	p57 Kip2 Recombinant Rabbit Monoclonal Antibody (ARC1424)
MA535698	PKC delta Recombinant Rabbit Monoclonal Antibody (ARC1434)
MA535699	ILF3 Recombinant Rabbit Monoclonal Antibody (ARC1438)
MA535700	SSTR5 Recombinant Rabbit Monoclonal Antibody (ARC1282)
MA535701	C2orf33 Recombinant Rabbit Monoclonal Antibody (ARC1243)
MA535702	FKBP8 Recombinant Rabbit Monoclonal Antibody (ARC1259)
MA535703	TFPI Recombinant Rabbit Monoclonal Antibody (ARC1154)
MA535704	Oncostatin M Recombinant Rabbit Monoclonal Antibody (ARC1194)
MA535705	Calpain 1 Recombinant Rabbit Monoclonal Antibody (ARC1272)
MA535706	TPM1 Recombinant Rabbit Monoclonal Antibody (ARC1283)
MA535707	MYL9 Recombinant Rabbit Monoclonal Antibody (ARC1284)

MA535708	MSH2 Recombinant Rabbit Monoclonal Antibody (ARC1285)
MA535709	MYL2 Recombinant Rabbit Monoclonal Antibody (ARC1286)
MA535710	EEF1G Recombinant Rabbit Monoclonal Antibody (ARC1287)
MA535711	CDKN2C Recombinant Rabbit Monoclonal Antibody (ARC1289)
MA535713	RANGAP1 Recombinant Rabbit Monoclonal Antibody (ARC1292)
MA535715	TLX Recombinant Rabbit Monoclonal Antibody (ARC1295)
MA535716	Haptoglobin Recombinant Rabbit Monoclonal Antibody (ARC1298)
MA535718	HPRT1 Recombinant Rabbit Monoclonal Antibody (ARC1300)
MA535719	FABP3 Recombinant Rabbit Monoclonal Antibody (ARC1302)
MA535720	ENDOG Recombinant Rabbit Monoclonal Antibody (ARC1308)
MA535721	VAMP3 Recombinant Rabbit Monoclonal Antibody (ARC1312)
MA535722	AMHR2 Recombinant Rabbit Monoclonal Antibody (ARC1313)
MA535723	DLAT Recombinant Rabbit Monoclonal Antibody (ARC1314)
MA535725	Bub3 Recombinant Rabbit Monoclonal Antibody (ARC1318)
MA535726	BST2 Recombinant Rabbit Monoclonal Antibody (ARC1321)
MA535727	Calpain S1 Recombinant Rabbit Monoclonal Antibody (ARC1326)
MA535729	HDAC8 Recombinant Rabbit Monoclonal Antibody (ARC1331)
MA535730	TrxR2 Recombinant Rabbit Monoclonal Antibody (ARC1339)
MA535732	PLCG1 Recombinant Rabbit Monoclonal Antibody (ARC1345)
MA535733	MBD3 Recombinant Rabbit Monoclonal Antibody (ARC1349)
MA535734	GYS1 Recombinant Rabbit Monoclonal Antibody (ARC1350)
MA535738	Cystatin C Recombinant Rabbit Monoclonal Antibody (ARC1357)
MA535739	14-3-3 theta Recombinant Rabbit Monoclonal Antibody (ARC1358)
MA535740	hnRNP C Recombinant Rabbit Monoclonal Antibody (ARC1366)
MA535741	Pyruvate Carboxylase Recombinant Rabbit Monoclonal Antibody (ARC1371)
MA535742	Raptor Recombinant Rabbit Monoclonal Antibody (ARC1375)
MA535743	FEN1 Recombinant Rabbit Monoclonal Antibody (ARC1376)
MA535744	UBD Recombinant Rabbit Monoclonal Antibody (ARC1379)
MA535746	SKP1 Recombinant Rabbit Monoclonal Antibody (ARC1384)
MA535747	XDH Recombinant Rabbit Monoclonal Antibody (0M8E7)
MA535749	PROX1 Recombinant Rabbit Monoclonal Antibody (ARC1391)
MA535750	PTGER2 Recombinant Rabbit Monoclonal Antibody (ARC1393)
MA535751	SRSF3 Recombinant Rabbit Monoclonal Antibody (ARC1394)
MA535753	FHIT Recombinant Rabbit Monoclonal Antibody (ARC1401)
MA535754	KHSRP Recombinant Rabbit Monoclonal Antibody (ARC1402)
MA535755	ASS1 Recombinant Rabbit Monoclonal Antibody (ARC1407)
MA535756	TMPRSS2 Recombinant Rabbit Monoclonal Antibody (ARC1439)
MA535757	LGR6 Recombinant Rabbit Monoclonal Antibody (ARC1440)
MA535758	NCK1 Recombinant Rabbit Monoclonal Antibody (ARC1441)
MA535759	PRDX4 Recombinant Rabbit Monoclonal Antibody (ARC1442)
MA535760	Hemopexin Recombinant Rabbit Monoclonal Antibody (ARC1444)
MA535761	ABCE1 Recombinant Rabbit Monoclonal Antibody (ARC1445)
MA535762	MCAK Recombinant Rabbit Monoclonal Antibody (ARC1447)
MA535763	eIF3b Recombinant Rabbit Monoclonal Antibody (ARC1448)
MA535764	REA Recombinant Rabbit Monoclonal Antibody (ARC1449)
MA535765	ABAT Recombinant Rabbit Monoclonal Antibody (ARC1450)
MA535766	CD74 Recombinant Rabbit Monoclonal Antibody (ARC1452)
MA535770	Fyn Recombinant Rabbit Monoclonal Antibody (ARC1459)
MA535771	HMGB2 Recombinant Rabbit Monoclonal Antibody (ARC1460)
MA535772	SET Recombinant Rabbit Monoclonal Antibody (ARC1462)
MA535774	PSMA2 Recombinant Rabbit Monoclonal Antibody (ARC1465)
MA535776	SOCS2 Recombinant Rabbit Monoclonal Antibody (ARC1470)
MA535778	CaMKII delta Recombinant Rabbit Monoclonal Antibody (ARC1473)
MA535780	IMPDH2 Recombinant Rabbit Monoclonal Antibody (ARC1479)
MA535782	SETD7 Recombinant Rabbit Monoclonal Antibody (ARC1483)
MA535784	Podoplanin Recombinant Rabbit Monoclonal Antibody (ARC1488)
MA535785	MCM4 Recombinant Rabbit Monoclonal Antibody (ARC1491)
MA535786	UBE1 Recombinant Rabbit Monoclonal Antibody (ARC1493)
MA535787	UBC13 Recombinant Rabbit Monoclonal Antibody (ARC1495)

MA535791	CRALBP Recombinant Rabbit Monoclonal Antibody (ARC1502)
MA535793	VPS35 Recombinant Rabbit Monoclonal Antibody (ARC1509)
MA535794	Hemoglobin alpha Recombinant Rabbit Monoclonal Antibody (ARC1510)
MA535795	DHFR Recombinant Rabbit Monoclonal Antibody (ARC1513)
MA535797	UQCRFS1 Recombinant Rabbit Monoclonal Antibody (ARC1612)
MA535799	LIN28B Recombinant Rabbit Monoclonal Antibody (ARC1616)
MA535802	AHA1 Recombinant Rabbit Monoclonal Antibody (ARC1623)
MA535803	FOXP1 Recombinant Rabbit Monoclonal Antibody (ARC1624)
MA535804	SIRP alpha Recombinant Rabbit Monoclonal Antibody (ARC1626)
MA535805	PSMB9 Recombinant Rabbit Monoclonal Antibody (ARC1629)
MA535806	RPS20 Recombinant Rabbit Monoclonal Antibody (ARC1630)
MA535809	HOXA13 Recombinant Rabbit Monoclonal Antibody (ARC1636)
MA535812	HIP Recombinant Rabbit Monoclonal Antibody (ARC1639)
MA535814	GABARAPL2 Recombinant Rabbit Monoclonal Antibody (ARC1648)
MA535815	Spectrin alpha-1 Recombinant Rabbit Monoclonal Antibody (ARC1650)
MA535817	WIPI1 Recombinant Rabbit Monoclonal Antibody (ARC1652)
MA535819	HnRNP Q Recombinant Rabbit Monoclonal Antibody (ARC1656)
MA535821	RPE65 Recombinant Rabbit Monoclonal Antibody (ARC1659)
MA535822	Nectin 2 Recombinant Rabbit Monoclonal Antibody (ARC1665)
MA535823	SF3B3 Recombinant Rabbit Monoclonal Antibody (ARC1667)
MA535825	DDX6 Recombinant Rabbit Monoclonal Antibody (ARC1673)
MA535826	EWSR1 Recombinant Rabbit Monoclonal Antibody (ARC1674)
MA535828	DC-SIGN Recombinant Rabbit Monoclonal Antibody (ARC1679)
MA535829	CARD11 Recombinant Rabbit Monoclonal Antibody (ARC1680)
MA535830	CDK8 Recombinant Rabbit Monoclonal Antibody (ARC1682)
MA535832	Phospho-mTOR (Ser2448) Recombinant Rabbit Monoclonal Antibody (ARC0094)
MA535833	Phospho-ERK1/ERK2 (Thr202,Thr185) Recombinant Rabbit Monoclonal Antibody (ARC0100)
MA535834	Phospho-JNK1/JNK2/JNK3 (Thr183, Thr221) Recombinant Rabbit Monoclonal Antibody (ARC0193)
MA535836	Phospho-STAT3 (Ser727) Recombinant Rabbit Monoclonal Antibody (ARC0150)
MA535837	Phospho-SIRT1 (Ser47) Recombinant Rabbit Monoclonal Antibody (ARC1520)
MA535838	Phospho-mTOR (Ser2481) Recombinant Rabbit Monoclonal Antibody (ARC1522)
MA535839	Phospho-beta Catenin (Ser33, Ser37) Recombinant Rabbit Monoclonal Antibody (ARC1523)
MA535840	Phospho-AKT1 (Thr450) Recombinant Rabbit Monoclonal Antibody (ARC1524)
MA535841	Phospho-AKT1 (Ser129) Recombinant Rabbit Monoclonal Antibody (ARC1525)
MA535842	Phospho-AKT1 (Ser124) Recombinant Rabbit Monoclonal Antibody (ARC1526)
MA535850	Phospho-IRF3 (Ser386) Recombinant Rabbit Monoclonal Antibody (ARC1539)
MA535845	Phospho-p53 (Thr55) Recombinant Rabbit Monoclonal Antibody (ARC1530)
MA535847	Phospho-c-Myc (Thr58, Ser62) Recombinant Rabbit Monoclonal Antibody (ARC1532)
MA535848	Phospho-c-Myc (Thr58) Recombinant Rabbit Monoclonal Antibody (ARC1534)
MA535849	Phospho-EGFR (Tyr1068) Recombinant Rabbit Monoclonal Antibody (ARC1538)
MA535853	Phospho-IkB alpha (Ser36) Recombinant Rabbit Monoclonal Antibody (ARC1543)
MA535854	Phospho-CDK1/CDK2 (Thr14) Recombinant Rabbit Monoclonal Antibody (ARC1546)
MA535855	Phospho-AMPK alpha-1 (Ser496) Recombinant Rabbit Monoclonal Antibody (ARC1547)
MA535857	Phospho-Estrogen Receptor alpha (Ser118) Recombinant Rabbit Monoclonal Antibody (ARC1549)
MA535858	Phospho-CDK2 (Tyr15) Recombinant Rabbit Monoclonal Antibody (ARC1550)
MA535859	Phospho-SMAD2 (Ser250) Recombinant Rabbit Monoclonal Antibody (ARC1552)
MA535861	Phospho-c-Raf (Ser621) Recombinant Rabbit Monoclonal Antibody (ARC1556)
MA535862	Phospho-c-Raf (Ser259) Recombinant Rabbit Monoclonal Antibody (ARC1557)
MA535864	Phospho-MEK1 (Ser298) Recombinant Rabbit Monoclonal Antibody (ARC1565)
MA535865	Phospho-MEK1 (Thr292) Recombinant Rabbit Monoclonal Antibody (ARC1566)
MA535866	Phospho-PDHA1 (Ser293) Recombinant Rabbit Monoclonal Antibody (9S3I4)
MA535867	Phospho-Smad5 (Ser463,Ser465) Recombinant Rabbit Monoclonal Antibody (ARC1568)
MA535868	Phospho-PLK1 (Thr210) Recombinant Rabbit Monoclonal Antibody (ARC1570)
MA535870	Phospho-SRC (Tyr419) Recombinant Rabbit Monoclonal Antibody (ARC1572)
MA535871	Phospho-Tau (Ser396) Recombinant Rabbit Monoclonal Antibody (ARC1573)
MA535872	Phospho-ATM (Ser1981) Recombinant Rabbit Monoclonal Antibody (ARC1575)
MA535873	Phospho-HSP27 (Ser82) Recombinant Rabbit Monoclonal Antibody (ARC1576)
MA535875	Phospho-PRKAR2A (Ser99) Recombinant Rabbit Monoclonal Antibody (ARC1579)
MA535876	Phospho-RPA2 (Thr21) Recombinant Rabbit Monoclonal Antibody (ARC1586)

MA535877	Phospho-BRAF35 (Thr401) Recombinant Rabbit Monoclonal Antibody (ARC1587)
MA535878	Phospho-PP2A alpha (Tyr307) Recombinant Rabbit Monoclonal Antibody (ARC1589)
MA535880	Phospho-Chk1 (Ser296) Recombinant Rabbit Monoclonal Antibody (ARC1593)
MA535881	Phospho-PAK4/6/7 (Ser474, Ser602, Ser560) Recombinant Rabbit Monoclonal Antibody (ARC1596)
MA535882	Phospho-Lyn (Tyr396) Recombinant Rabbit Monoclonal Antibody (ARC1597)
MA535883	Phospho-ErbB3 (Tyr1222) Recombinant Rabbit Monoclonal Antibody (ARC1599)
MA535888	Phospho-CDK1/CDK2/CDK3 (Thr14) Recombinant Rabbit Monoclonal Antibody (ARC1545)
MA535890	Phospho-Cdc37 (Ser13) Recombinant Rabbit Monoclonal Antibody (ARC1602)
MA535892	Phospho-SMC1 (Ser957) Recombinant Rabbit Monoclonal Antibody (ARC1606)
MA535893	Phospho-IRE1 alpha (Ser724) Recombinant Rabbit Monoclonal Antibody (ARC1607)
MA535895	Phospho-Crk (Tyr221) Recombinant Rabbit Monoclonal Antibody (ARC1640)
MA535896	Phospho-HSL (Ser853) Recombinant Rabbit Monoclonal Antibody (ARC1657)
MA535897	Phospho-Histone H3 (Thr3) Recombinant Rabbit Monoclonal Antibody (ARC1662)
MA538132	15-PGDH Recombinant Rabbit Monoclonal Antibody (ARC1280)
MA537914	5T4 Recombinant Rabbit Monoclonal Antibody (ARC2167)
MA538211	A2M Recombinant Rabbit Monoclonal Antibody (ARC1734)
MA537931	ABCF1 Recombinant Rabbit Monoclonal Antibody (ARC2199)
MA538176	ABHD5 Recombinant Rabbit Monoclonal Antibody (ARC1777)
MA537932	Abl2 Recombinant Rabbit Monoclonal Antibody (ARC2200)
MA538031	Acetyl-p53 (Lys370) Recombinant Rabbit Monoclonal Antibody (ARC2368)
MA537901	ACVR1 Recombinant Rabbit Monoclonal Antibody (ARC2449)
MA538039	ACVR1B Recombinant Rabbit Monoclonal Antibody (ARC1899)
MA538212	ACVRL1 Recombinant Rabbit Monoclonal Antibody (ARC1735)
MA538137	Adenosine Deaminase Recombinant Rabbit Monoclonal Antibody (ARC1152)
MA538140	ADORA1 Recombinant Rabbit Monoclonal Antibody (ARC1987)
MA538084	AGTR2 Recombinant Rabbit Monoclonal Antibody (ARC2064)
MA537960	AGO2 Recombinant Rabbit Monoclonal Antibody (ARC0189)
MA538020	AKAP14 Recombinant Rabbit Monoclonal Antibody (ARC2344)
MA538186	ALDH1A2 Recombinant Rabbit Monoclonal Antibody (ARC1797)
MA538050	ALOX5 Recombinant Rabbit Monoclonal Antibody (ARC1926)
MA538061	alpha Actinin 4 Recombinant Rabbit Monoclonal Antibody (ARC1959)
MA537942	alpha COP Recombinant Rabbit Monoclonal Antibody (ARC2215)
MA538082	alpha Internexin Recombinant Rabbit Monoclonal Antibody (ARC2054)
MA537982	Alpha Sarcoglycan Recombinant Rabbit Monoclonal Antibody (ARC2280)
MA537926	Androgen Receptor Recombinant Rabbit Monoclonal Antibody (ARC0090)
MA538041	AP2 alpha Recombinant Rabbit Monoclonal Antibody (ARC1905)
MA538043	ARHGAP5 Recombinant Rabbit Monoclonal Antibody (ARC1916)
MA538154	ARTS1 Recombinant Rabbit Monoclonal Antibody (ARC2120)
MA538002	ATAD2 Recombinant Rabbit Monoclonal Antibody (ARC2317)
MA538135	ATGL Recombinant Rabbit Monoclonal Antibody (ARC1192)
MA538164	ATG10 Recombinant Rabbit Monoclonal Antibody (ARC1425)
MA537905	ATP1A2 Recombinant Rabbit Monoclonal Antibody (ARC2458)
MA537906	ATRIP Recombinant Rabbit Monoclonal Antibody (ARC2460)
MA538094	BANF1 Recombinant Rabbit Monoclonal Antibody (ARC2085)
MA538190	BAZ1B Recombinant Rabbit Monoclonal Antibody (ARC1614)
MA537965	BCAM Recombinant Rabbit Monoclonal Antibody (ARC2252)
MA538009	BCL11B Recombinant Rabbit Monoclonal Antibody (ARC2329)
MA538012	BCL2L12 Recombinant Rabbit Monoclonal Antibody (ARC2334)
MA537912	BCL2L14 Recombinant Rabbit Monoclonal Antibody (ARC2474)
MA537831	beta Arrestin 1 Recombinant Rabbit Monoclonal Antibody (ARC2370)
MA538116	BNIP1 Recombinant Rabbit Monoclonal Antibody (ARC2137)
MA538134	Brachyury Recombinant Rabbit Monoclonal Antibody (ARC1983)
MA537889	C1 inhibitor Recombinant Rabbit Monoclonal Antibody (ARC2427)
MA538095	Calponin 1 Recombinant Rabbit Monoclonal Antibody (ARC0232)
MA537936	Calsequestrin Recombinant Rabbit Monoclonal Antibody (ARC2209)
MA538117	CAPG Recombinant Rabbit Monoclonal Antibody (ARC2138)
MA538188	Carbonic Anhydrase II Recombinant Rabbit Monoclonal Antibody (ARC1451)
MA537829	CAS Recombinant Rabbit Monoclonal Antibody (ARC1856)
MA537917	Caspase 6 Recombinant Rabbit Monoclonal Antibody (ARC0031)

MA537891	Caspr2 Recombinant Rabbit Monoclonal Antibody (ARC2430)
MA537892	CBR1 Recombinant Rabbit Monoclonal Antibody (ARC2432)
MA538201	CCR9 Recombinant Rabbit Monoclonal Antibody (ARC1706)
MA538113	CD147 Recombinant Rabbit Monoclonal Antibody (ARC0964)
MA538136	CD155 Recombinant Rabbit Monoclonal Antibody (ARC1178)
MA538206	CD166 (ALCAM) Recombinant Rabbit Monoclonal Antibody (ARC1720)
MA538151	CD266 Recombinant Rabbit Monoclonal Antibody (ARC2119)
MA537838	CD27 Recombinant Rabbit Monoclonal Antibody (ARC0625)
MA537934	CD298 Recombinant Rabbit Monoclonal Antibody (ARC2205)
MA538158	Cdc23 Recombinant Rabbit Monoclonal Antibody (ARC2133)
MA538213	Cdc25B Recombinant Rabbit Monoclonal Antibody (ARC1736)
MA538215	Cdc34 Recombinant Rabbit Monoclonal Antibody (ARC1742)
MA538004	Cdc40 Recombinant Rabbit Monoclonal Antibody (ARC2321)
MA537968	Cdc42EP1 Recombinant Rabbit Monoclonal Antibody (ARC2298)
MA538104	CDCA5 Recombinant Rabbit Monoclonal Antibody (ARC2096)
MA537939	CDH16 Recombinant Rabbit Monoclonal Antibody (ARC2212)
MA537944	CTRL Recombinant Rabbit Monoclonal Antibody (ARC2218)
MA538146	CEND Recombinant Rabbit Monoclonal Antibody (ARC2092)
MA538099	CENPC Recombinant Rabbit Monoclonal Antibody (ARC2097)
MA537879	CHMP2B Recombinant Rabbit Monoclonal Antibody (ARC2404)
MA537823	CHMP3 Recombinant Rabbit Monoclonal Antibody (ARC1849)
MA537815	CHRNA4 Recombinant Rabbit Monoclonal Antibody (ARC1822)
MA538025	CHX10 Recombinant Rabbit Monoclonal Antibody (ARC2354)
MA537911	CIRBP Recombinant Rabbit Monoclonal Antibody (ARC2472)
MA537938	CLGN Recombinant Rabbit Monoclonal Antibody (ARC2213)
MA537941	CLGN Recombinant Rabbit Monoclonal Antibody (ARC2214)
MA538168	CLOCK Recombinant Rabbit Monoclonal Antibody (ARC1819)
MA538058	c-MAF Recombinant Rabbit Monoclonal Antibody (ARC1954)
MA537967	c-Mpl Recombinant Rabbit Monoclonal Antibody (ARC2257)
MA538093	CNTNAP1 Recombinant Rabbit Monoclonal Antibody (ARC2084)
MA538207	Collagen VI Recombinant Rabbit Monoclonal Antibody (ARC1725)
MA537922	COPS3 Recombinant Rabbit Monoclonal Antibody (ARC2182)
MA538126	Creatine Kinase BB Recombinant Rabbit Monoclonal Antibody (ARC1990)
MA537908	CTNNA3 Recombinant Rabbit Monoclonal Antibody (ARC2463)
MA537837	CUGBP1 Recombinant Rabbit Monoclonal Antibody (ARC1866)
MA538155	CYP26A1 Recombinant Rabbit Monoclonal Antibody (ARC2126)
MA537946	CYP4A11/CYP4A22 Recombinant Rabbit Monoclonal Antibody (ARC2219)
MA538216	Cytokeratin 1 Recombinant Rabbit Monoclonal Antibody (ARC1745)
MA537810	Cytokeratin 4 Recombinant Rabbit Monoclonal Antibody (ARC1804)
MA538024	Cytokeratin 6 Recombinant Rabbit Monoclonal Antibody (ARC2352)
MA538120	Cytokeratin 6 Recombinant Rabbit Monoclonal Antibody (ARC2142)
MA538115	Cytokeratin 7 Recombinant Rabbit Monoclonal Antibody (ARC0978)
MA538103	DARPP-32 Recombinant Rabbit Monoclonal Antibody (ARC2095)
MA537945	DAXX Recombinant Rabbit Monoclonal Antibody (ARC2220)
MA537884	DCTN5 Recombinant Rabbit Monoclonal Antibody (ARC2412)
MA538096	Desmin Recombinant Rabbit Monoclonal Antibody (ARC0235)
MA537990	Dextrin Recombinant Rabbit Monoclonal Antibody (ARC2297)
MA537907	DMAP1 Recombinant Rabbit Monoclonal Antibody (ARC2462)
MA537883	DOK1 Recombinant Rabbit Monoclonal Antibody (ARC2410)
MA538119	DRD3 Recombinant Rabbit Monoclonal Antibody (ARC1041)
MA538045	DRP1 Recombinant Rabbit Monoclonal Antibody (ARC1918)
MA537890	Dynamin 3 Recombinant Rabbit Monoclonal Antibody (ARC2428)
MA537920	Dysferlin Recombinant Rabbit Monoclonal Antibody (ARC2178)
MA538056	E2F2 Recombinant Rabbit Monoclonal Antibody (ARC1940)
MA537840	E2F6 Recombinant Rabbit Monoclonal Antibody (ARC0628)
MA538202	EIF2A Recombinant Rabbit Monoclonal Antibody (ARC1712)
MA538067	eIF3e Recombinant Rabbit Monoclonal Antibody (ARC1997)
MA537947	EN1 Recombinant Rabbit Monoclonal Antibody (ARC2225)
MA537937	ENTPD5 Recombinant Rabbit Monoclonal Antibody (ARC2210)

MA537940	ENTPD5 Recombinant Rabbit Monoclonal Antibody (ARC2211)
MA537867	EphB3 Recombinant Rabbit Monoclonal Antibody (ARC2384)
MA538129	ERCC1 Recombinant Rabbit Monoclonal Antibody (ARC1241)
MA537989	ERP29 Recombinant Rabbit Monoclonal Antibody (ARC2295)
MA537893	Fetuin A Recombinant Rabbit Monoclonal Antibody (ARC2433)
MA538085	FHL2 Recombinant Rabbit Monoclonal Antibody (ARC2068)
MA537943	Fibrinogen alpha chain Recombinant Rabbit Monoclonal Antibody (ARC2227)
MA537933	FLAP Recombinant Rabbit Monoclonal Antibody (ARC2204)
MA537863	FMO3 Recombinant Rabbit Monoclonal Antibody (ARC2377)
MA538049	FNTB Recombinant Rabbit Monoclonal Antibody (ARC1923)
MA538083	FOXL2 Recombinant Rabbit Monoclonal Antibody (ARC2055)
MA537870	GAA Recombinant Rabbit Monoclonal Antibody (ARC2392)
MA537949	GABRB1 Recombinant Rabbit Monoclonal Antibody (ARC2229)
MA537849	GAD67 Recombinant Rabbit Monoclonal Antibody (ARC1879)
MA537896	gamma Adaptin Recombinant Rabbit Monoclonal Antibody (ARC2440)
MA538194	gamma Tubulin Recombinant Rabbit Monoclonal Antibody (ARC1684)
MA538112	GATA4 Recombinant Rabbit Monoclonal Antibody (ARC0962)
MA537836	GCLM Recombinant Rabbit Monoclonal Antibody (ARC0597)
MA537886	Gemin 2 Recombinant Rabbit Monoclonal Antibody (ARC2419)
MA538203	GLAST Recombinant Rabbit Monoclonal Antibody (ARC1714)
MA538204	GluR6 Recombinant Rabbit Monoclonal Antibody (ARC1718)
MA537843	GLUT1 Recombinant Rabbit Monoclonal Antibody (ARC0304)
MA537862	Glycerol kinase Recombinant Rabbit Monoclonal Antibody (ARC2376)
MA538001	GOLIM4 Recombinant Rabbit Monoclonal Antibody (ARC2315)
MA537839	GOLPH2 Recombinant Rabbit Monoclonal Antibody (ARC0697)
MA537880	GOT2 Recombinant Rabbit Monoclonal Antibody (ARC2405)
MA538019	GSC Recombinant Rabbit Monoclonal Antibody (ARC2343)
MA538198	GSTK1 Recombinant Rabbit Monoclonal Antibody (ARC1689)
MA537952	HB9 Recombinant Rabbit Monoclonal Antibody (ARC2233)
MA538054	H-cadherin Recombinant Rabbit Monoclonal Antibody (ARC1931)
MA537951	Hex Recombinant Rabbit Monoclonal Antibody (ARC2232)
MA537953	Hippocalcin Recombinant Rabbit Monoclonal Antibody (ARC2235)
MA538088	Histone H2A Recombinant Rabbit Monoclonal Antibody (ARC2072)
MA538030	HLA-DRB1 Recombinant Rabbit Monoclonal Antibody (ARC2360)
MA537869	HMGCS2 Recombinant Rabbit Monoclonal Antibody (ARC2390)
MA537888	HOXA9 Recombinant Rabbit Monoclonal Antibody (ARC2423)
MA538110	HP1 gamma Recombinant Rabbit Monoclonal Antibody (ARC0952)
MA538074	HRH3 Recombinant Rabbit Monoclonal Antibody (ARC2016)
MA538121	HSP105 Recombinant Rabbit Monoclonal Antibody (ARC1092)
MA538182	HSP20 Recombinant Rabbit Monoclonal Antibody (ARC1787)
MA538131	HSP90 alpha Recombinant Rabbit Monoclonal Antibody (ARC1167)
MA538023	HSPA12A Recombinant Rabbit Monoclonal Antibody (ARC2351)
MA537957	HSPA4 Recombinant Rabbit Monoclonal Antibody (ARC2237)
MA537963	Human IgM Recombinant Rabbit Monoclonal Antibody (ARC2245)
MA537992	IBA1 Recombinant Rabbit Monoclonal Antibody (ARC2301)
MA538051	IRF6 Recombinant Rabbit Monoclonal Antibody (ARC1928)
MA538170	ITCH Recombinant Rabbit Monoclonal Antibody (ARC1759)
MA537895	Jagged 2 Recombinant Rabbit Monoclonal Antibody (ARC2438)
MA538098	Kaiso Recombinant Rabbit Monoclonal Antibody (ARC2086)
MA537964	Kappa Light Chain Recombinant Rabbit Monoclonal Antibody (ARC2247)
MA538171	KARS Recombinant Rabbit Monoclonal Antibody (ARC1765)
MA538021	KDM1B Recombinant Rabbit Monoclonal Antibody (ARC2347)
MA538152	KGF Recombinant Rabbit Monoclonal Antibody (ARC2107)
MA537882	KLC1 Recombinant Rabbit Monoclonal Antibody (ARC2409)
MA538047	KLF10 Recombinant Rabbit Monoclonal Antibody (ARC1921)
MA537844	Ku80 Recombinant Rabbit Monoclonal Antibody (ARC0706)
MA538027	Lambda Light Chain Recombinant Rabbit Monoclonal Antibody (ARC2357)
MA537866	LCAT Recombinant Rabbit Monoclonal Antibody (ARC2383)
MA538044	Lefty1/Lefty2 Recombinant Rabbit Monoclonal Antibody (ARC2149)

MA537868	LGI1 Recombinant Rabbit Monoclonal Antibody (ARC2388)
MA537894	LHX2 Recombinant Rabbit Monoclonal Antibody (ARC2435)
MA537924	LHX3 Recombinant Rabbit Monoclonal Antibody (ARC2150)
MA538128	LMAN1 Recombinant Rabbit Monoclonal Antibody (ARC1978)
MA538157	LMO4 Recombinant Rabbit Monoclonal Antibody (ARC2114)
MA538208	LRP8 Recombinant Rabbit Monoclonal Antibody (ARC1730)
MA538011	MAD3 Recombinant Rabbit Monoclonal Antibody (ARC2333)
MA537821	MAGE1 Recombinant Rabbit Monoclonal Antibody (ARC1839)
MA537966	Mammaglobin A Recombinant Rabbit Monoclonal Antibody (ARC2255)
MA537988	MAPRE3 Recombinant Rabbit Monoclonal Antibody (ARC2305)
MA538161	MARK2 Recombinant Rabbit Monoclonal Antibody (ARC1862)
MA537899	MAT2A Recombinant Rabbit Monoclonal Antibody (ARC2447)
MA537827	MATR3 Recombinant Rabbit Monoclonal Antibody (ARC1854)
MA538200	MDH1 Recombinant Rabbit Monoclonal Antibody (ARC1692)
MA538123	MEK1/MEK2 Recombinant Rabbit Monoclonal Antibody (ARC0292)
MA538026	MEK3/MEK6 Recombinant Rabbit Monoclonal Antibody (ARC2356)
MA537973	MEKK2 Recombinant Rabbit Monoclonal Antibody (ARC2291)
MA538076	MELK Recombinant Rabbit Monoclonal Antibody (ARC2032)
MA538071	METAP2 Recombinant Rabbit Monoclonal Antibody (ARC2022)
MA538097	mGluR5 Recombinant Rabbit Monoclonal Antibody (ARC0250)
MA538013	MINA53 Recombinant Rabbit Monoclonal Antibody (ARC2336)
MA538166	MOV10 Recombinant Rabbit Monoclonal Antibody (ARC1817)
MA538185	MSI1 Recombinant Rabbit Monoclonal Antibody (ARC1796)
MA538016	MSI2 Recombinant Rabbit Monoclonal Antibody (ARC2341)
MA537812	MSP Recombinant Rabbit Monoclonal Antibody (ARC1809)
MA537828	MVP Recombinant Rabbit Monoclonal Antibody (ARC1855)
MA537909	NAT10 Recombinant Rabbit Monoclonal Antibody (ARC2468)
MA537860	N-cadherin Recombinant Rabbit Monoclonal Antibody (ARC0371)
MA538101	NDUFAF1 Recombinant Rabbit Monoclonal Antibody (ARC2091)
MA538178	NDUFS4 Recombinant Rabbit Monoclonal Antibody (ARC1784)
MA538077	NEK6 Recombinant Rabbit Monoclonal Antibody (ARC2034)
MA538018	NEK7 Recombinant Rabbit Monoclonal Antibody (ARC2342)
MA537970	Nkx3.1 Recombinant Rabbit Monoclonal Antibody (ARC2261)
MA537904	NR5A1 Recombinant Rabbit Monoclonal Antibody (ARC2457)
MA537986	NRG3 Recombinant Rabbit Monoclonal Antibody (ARC2289)
MA537853	NUF2 Recombinant Rabbit Monoclonal Antibody (ARC2157)
MA537898	NUP62 Recombinant Rabbit Monoclonal Antibody (ARC2444)
MA537824	OLIG1 Recombinant Rabbit Monoclonal Antibody (ARC1850)
MA538179	OPRD1 Recombinant Rabbit Monoclonal Antibody (ARC1786)
MA537971	ORM1 Recombinant Rabbit Monoclonal Antibody (ARC2263)
MA537974	Oxytocin Receptor Recombinant Rabbit Monoclonal Antibody (ARC2266)
MA537841	p16INK4a Recombinant Rabbit Monoclonal Antibody (ARC0661)
MA537928	p40 Recombinant Rabbit Monoclonal Antibody (ARC2194)
MA538068	PAF1 Recombinant Rabbit Monoclonal Antibody (ARC2019)
MA537900	PAG1 Recombinant Rabbit Monoclonal Antibody (ARC2448)
MA537994	Pan Ras Recombinant Rabbit Monoclonal Antibody (ARC2369)
MA538130	PAR4 Recombinant Rabbit Monoclonal Antibody (ARC1980)
MA538149	PARP2 Recombinant Rabbit Monoclonal Antibody (ARC2099)
MA537814	PBRM1 Recombinant Rabbit Monoclonal Antibody (ARC1820)
MA537903	PCBP1 Recombinant Rabbit Monoclonal Antibody (ARC2456)
MA538147	PCYT1A Recombinant Rabbit Monoclonal Antibody (ARC2098)
MA537864	PDLIM1 Recombinant Rabbit Monoclonal Antibody (ARC2379)
MA538177	PEN2 Recombinant Rabbit Monoclonal Antibody (ARC1779)
MA538036	PER3 Recombinant Rabbit Monoclonal Antibody (ARC1893)
MA538148	PERP Recombinant Rabbit Monoclonal Antibody (ARC2094)
MA537873	PEX19 Recombinant Rabbit Monoclonal Antibody (ARC2395)
MA538089	PF4 Recombinant Rabbit Monoclonal Antibody (ARC2076)
MA537948	PGC1 alpha/beta Recombinant Rabbit Monoclonal Antibody (ARC0153)
MA537995	PHF8 Recombinant Rabbit Monoclonal Antibody (ARC2306)

MA538237	Phospho-Cdc6 (Ser54) Recombinant Rabbit Monoclonal Antibody (ARC1693)
MA538223	Phospho-c-Jun (Ser63) Recombinant Rabbit Monoclonal Antibody (ARC0051)
MA538229	Phospho-c-Myc (Ser62) Recombinant Rabbit Monoclonal Antibody (ARC1533)
MA538222	Phospho-EGFR (Tyr1068) Recombinant Rabbit Monoclonal Antibody (ARC0013)
MA538235	Phospho-GSK3B (Ser9) Recombinant Rabbit Monoclonal Antibody (ARC0069)
MA538225	Phospho-Histone H2A.X (Ser139) Recombinant Rabbit Monoclonal Antibody (ARC0110)
MA538224	Phospho-p70 S6 Kinase (Thr421) Recombinant Rabbit Monoclonal Antibody (ARC0107)
MA538239	Phospho-PAK1/2/3 (Ser144) Recombinant Rabbit Monoclonal Antibody (ARC1698)
MA538236	Phospho-PKR (Thr446) Recombinant Rabbit Monoclonal Antibody (ARC0293)
MA538230	Phospho-POLR2A (Ser2) Recombinant Rabbit Monoclonal Antibody (ARC1540)
MA538226	Phospho-SMAD3 (Ser423) Recombinant Rabbit Monoclonal Antibody (ARC0170)
MA538231	Phospho-STAT1 (Ser727) Recombinant Rabbit Monoclonal Antibody (ARC1544)
MA538240	Phosphotyrosine Recombinant Rabbit Monoclonal Antibody (ARC1703)
MA538199	PIAS1/PIAS2 Recombinant Rabbit Monoclonal Antibody (ARC1690)
MA538144	PIK3R4 Recombinant Rabbit Monoclonal Antibody (ARC2090)
MA537910	Pirh2 Recombinant Rabbit Monoclonal Antibody (ARC2469)
MA537819	PIST Recombinant Rabbit Monoclonal Antibody (ARC1831)
MA537834	PKC alpha Recombinant Rabbit Monoclonal Antibody (ARC0197)
MA537978	PKN2 Recombinant Rabbit Monoclonal Antibody (ARC2274)
MA537885	PLA2G2A Recombinant Rabbit Monoclonal Antibody (ARC2416)
MA538111	Placental Alkaline Phosphatase Recombinant Rabbit Monoclonal Antibody (ARC0961)
MA537972	PLCB3 Recombinant Rabbit Monoclonal Antibody (ARC2269)
MA538109	PMP70 Recombinant Rabbit Monoclonal Antibody (ARC2131)
MA537998	Pokemon Recombinant Rabbit Monoclonal Antibody (ARC2320)
MA537977	PRKAR2B Recombinant Rabbit Monoclonal Antibody (ARC2272)
MA538107	Prolactin Receptor Recombinant Rabbit Monoclonal Antibody (ARC2123)
MA538108	Prostate Specific Acid Phosphatase Recombinant Rabbit Monoclonal Antibody (ARC2124)
MA538195	PRP19 Recombinant Rabbit Monoclonal Antibody (ARC1685)
MA538118	PRP8 Recombinant Rabbit Monoclonal Antibody (ARC1038)
MA538187	PTBP2 Recombinant Rabbit Monoclonal Antibody (ARC1798)
MA537997	PUM2 Recombinant Rabbit Monoclonal Antibody (ARC2308)
MA537872	RACGAP1 Recombinant Rabbit Monoclonal Antibody (ARC2394)
MA537979	RAD21 Recombinant Rabbit Monoclonal Antibody (ARC2276)
MA537981	Radixin Recombinant Rabbit Monoclonal Antibody (ARC2277)
MA538173	RAIDD Recombinant Rabbit Monoclonal Antibody (ARC1770)
MA538174	RAIDD Recombinant Rabbit Monoclonal Antibody (ARC1771)
MA537874	RANBP9 Recombinant Rabbit Monoclonal Antibody (ARC2397)
MA538205	RAP1A/RAP1B Recombinant Rabbit Monoclonal Antibody (ARC1719)
MA537925	RBM8A Recombinant Rabbit Monoclonal Antibody (ARC2189)
MA537820	RCC1 Recombinant Rabbit Monoclonal Antibody (ARC1834)
MA538006	Renalase Recombinant Rabbit Monoclonal Antibody (ARC2324)
MA538167	Rhodopsin Recombinant Rabbit Monoclonal Antibody (ARC1818)
MA537980	RING1 Recombinant Rabbit Monoclonal Antibody (ARC2278)
MA537835	ROCK1 Recombinant Rabbit Monoclonal Antibody (ARC2371)
MA537881	RPL7A Recombinant Rabbit Monoclonal Antibody (ARC2407)
MA538183	RPS6KB2 Recombinant Rabbit Monoclonal Antibody (ARC1790)
MA538153	RPSA Recombinant Rabbit Monoclonal Antibody (ARC2109)
MA538070	S100A6 Recombinant Rabbit Monoclonal Antibody (ARC2005)
MA537927	SCGN Recombinant Rabbit Monoclonal Antibody (ARC2196)
MA538005	SEC61 alpha Recombinant Rabbit Monoclonal Antibody (ARC2323)
MA538165	Septin 2 Recombinant Rabbit Monoclonal Antibody (ARC1810)
MA537935	SERCA1 ATPase Recombinant Rabbit Monoclonal Antibody (ARC2207)
MA537842	SERPINH1 Recombinant Rabbit Monoclonal Antibody (ARC0681)
MA538000	SESN1 Recombinant Rabbit Monoclonal Antibody (ARC2313)
MA538193	SFRP1 Recombinant Rabbit Monoclonal Antibody (ARC1683)
MA537983	SLC22A3 Recombinant Rabbit Monoclonal Antibody (ARC2285)
MA537865	SNAP alpha Recombinant Rabbit Monoclonal Antibody (ARC2382)
MA538150	SOX1 Recombinant Rabbit Monoclonal Antibody (ARC2100)
MA537921	SPOP Recombinant Rabbit Monoclonal Antibody (ARC2181)

MA537985	SRD5A2 Recombinant Rabbit Monoclonal Antibody (ARC2287)
MA537811	STIP1 Recombinant Rabbit Monoclonal Antibody (ARC1805)
MA538008	STK33 Recombinant Rabbit Monoclonal Antibody (ARC2330)
MA538038	STK39 Recombinant Rabbit Monoclonal Antibody (ARC1896)
MA538015	Syntaxin 1 Recombinant Rabbit Monoclonal Antibody (ARC2339)
MA537915	Synapsin II Recombinant Rabbit Monoclonal Antibody (ARC2166)
MA538156	Syntaxin 4 Recombinant Rabbit Monoclonal Antibody (ARC2113)
MA537984	Syntrophin alpha-1 Recombinant Rabbit Monoclonal Antibody (ARC2286)
MA538022	TAB3 Recombinant Rabbit Monoclonal Antibody (ARC2350)
MA537929	TACC3 Recombinant Rabbit Monoclonal Antibody (ARC2195)
MA538003	TAX1BP3 Recombinant Rabbit Monoclonal Antibody (ARC2319)
MA537916	TBR1 Recombinant Rabbit Monoclonal Antibody (ARC2198)
MA537999	TC10 Recombinant Rabbit Monoclonal Antibody (ARC2309)
MA537818	TFE3 Recombinant Rabbit Monoclonal Antibody (ARC1829)
MA537930	TGN46 Recombinant Rabbit Monoclonal Antibody (ARC2197)
MA538210	TIN2 Recombinant Rabbit Monoclonal Antibody (ARC1732)
MA538196	TIRAP Recombinant Rabbit Monoclonal Antibody (ARC1686)
MA538106	TMEFF2 Recombinant Rabbit Monoclonal Antibody (ARC2117)
MA537991	TNPO3 Recombinant Rabbit Monoclonal Antibody (ARC2310)
MA538184	TRBP Recombinant Rabbit Monoclonal Antibody (ARC1792)
MA537850	TRK fused gene Recombinant Rabbit Monoclonal Antibody (ARC1882)
MA538028	TrkB Recombinant Rabbit Monoclonal Antibody (ARC2358)
MA537871	TROY Recombinant Rabbit Monoclonal Antibody (ARC2393)
MA537902	Trypsin 2 Recombinant Rabbit Monoclonal Antibody (ARC2453)
MA538172	UAP1 Recombinant Rabbit Monoclonal Antibody (ARC1769)
MA538209	UAP56 Recombinant Rabbit Monoclonal Antibody (ARC1731)
MA537847	UCHL3 Recombinant Rabbit Monoclonal Antibody (ARC1874)
MA537913	UCP1/UCP3 Recombinant Rabbit Monoclonal Antibody (ARC2304)
MA537996	UNC84B Recombinant Rabbit Monoclonal Antibody (ARC2311)
MA537919	USP11 Recombinant Rabbit Monoclonal Antibody (ARC2176)
MA538217	USP9X Recombinant Rabbit Monoclonal Antibody (ARC1748)
MA538078	VAV2 Recombinant Rabbit Monoclonal Antibody (ARC2041)
MA537878	VEGFD Recombinant Rabbit Monoclonal Antibody (ARC2402)
MA537854	VPS18 Recombinant Rabbit Monoclonal Antibody (ARC2156)
MA537876	VPS28 Recombinant Rabbit Monoclonal Antibody (ARC2399)
MA538145	WAPL Recombinant Rabbit Monoclonal Antibody (ARC2088)
MA537987	WDR4 Recombinant Rabbit Monoclonal Antibody (ARC2292)
MA537877	XPD Recombinant Rabbit Monoclonal Antibody (ARC2401)
MA537918	ZIC1 Recombinant Rabbit Monoclonal Antibody (ARC2173)
MA542582	ELP3 Recombinant Rabbit Monoclonal Antibody (3R3U5)
MA542602	HK1 Recombinant Rabbit Monoclonal Antibody (6B7Q2)
MA542618	CSK Recombinant Rabbit Monoclonal Antibody (1A2F5)
MA542632	SAE1 Recombinant Rabbit Monoclonal Antibody (0L2N2)
MA542639	ERAB Recombinant Rabbit Monoclonal Antibody (3O4W10)
MA542640	MCL-1 Recombinant Rabbit Monoclonal Antibody (4F3W1)
MA542642	CDK9 Recombinant Rabbit Monoclonal Antibody (1C4E3)
MA542643	NGAL Recombinant Rabbit Monoclonal Antibody (1M2R3)
MA542645	Sodium Potassium ATPase Recombinant Rabbit Monoclonal Antibody (3L6F0)
MA542646	Aggrecan Recombinant Rabbit Monoclonal Antibody (6L4T2)
MA542659	BRD4 Recombinant Rabbit Monoclonal Antibody (7L7M8)
MA542660	Estrogen Receptor alpha Recombinant Rabbit Monoclonal Antibody (4L8K5)
MA542673	CHD1L Recombinant Rabbit Monoclonal Antibody (3C1O9)
MA542674	CRP Recombinant Rabbit Monoclonal Antibody (10J2R8)
MA542675	c-Myc Recombinant Rabbit Monoclonal Antibody (2O3U10)
MA542676	Cyclin D1 Recombinant Rabbit Monoclonal Antibody (7C8Z4)
MA542677	FASN Recombinant Rabbit Monoclonal Antibody (4U9S3)
MA542678	FOXA2 Recombinant Rabbit Monoclonal Antibody (9O4M2)
MA542679	Insulin Recombinant Rabbit Monoclonal Antibody (4C3Y9)
MA542681	PAI1 Recombinant Rabbit Monoclonal Antibody (6L6Y7)

MA542684	SMAD2 Recombinant Rabbit Monoclonal Antibody (2E3D8)
MA542685	SRC Recombinant Rabbit Monoclonal Antibody (7G6M9)
MA542686	SPT5 Recombinant Rabbit Monoclonal Antibody (9V9G7)
MA542687	LPP Recombinant Rabbit Monoclonal Antibody (8M2R8)
MA542688	BAG6 Recombinant Rabbit Monoclonal Antibody (8P4D9)
MA542689	Arginase 2 Recombinant Rabbit Monoclonal Antibody (1Y5U3)
MA542690	Syntaxin 1 Recombinant Rabbit Monoclonal Antibody (2B2F10)
MA542691	MYH7B Recombinant Rabbit Monoclonal Antibody (7W8J3)
MA542692	GUSB Recombinant Rabbit Monoclonal Antibody (4U0E2)
MA542693	SMG1 Recombinant Rabbit Monoclonal Antibody (5N10C7)
MA542694	MICA Recombinant Rabbit Monoclonal Antibody (5A1P3)
MA542695	CSH1 Recombinant Rabbit Monoclonal Antibody (7U9T2)
MA542696	PITX3 Recombinant Rabbit Monoclonal Antibody (3P4U8)
MA542698	PROK1 Recombinant Rabbit Monoclonal Antibody (1B4G2)
MA542699	ATP4B Recombinant Rabbit Monoclonal Antibody (3Y1Y1)
MA542700	CC2D1A Recombinant Rabbit Monoclonal Antibody (5Y6S6)
MA542701	Cytokeratin 12 Recombinant Rabbit Monoclonal Antibody (2M1K10)
MA542702	EHMT2 Recombinant Rabbit Monoclonal Antibody (2I0Q4)
MA542703	HDAC3 Recombinant Rabbit Monoclonal Antibody (1M1O7)
MA542704	Thyroid Peroxidase Recombinant Rabbit Monoclonal Antibody (3C4A4)
MA542707	STAT5 beta Recombinant Rabbit Monoclonal Antibody (6P5T0)
MA542708	YY1 Recombinant Rabbit Monoclonal Antibody (6J6L1)
MA542709	ARID1A Recombinant Rabbit Monoclonal Antibody (9D4J6)
MA542710	SynGAP Recombinant Rabbit Monoclonal Antibody (8M2S9)
MA542711	USP14 Recombinant Rabbit Monoclonal Antibody (3A3J1)
MA542713	WAVE2 Recombinant Rabbit Monoclonal Antibody (7A6C10)
MA542714	ATG12 Recombinant Rabbit Monoclonal Antibody (9R8W2)
MA542715	EBI3 Recombinant Rabbit Monoclonal Antibody (9J1H4)
MA542716	Alcohol dehydrogenase Pan Recombinant Rabbit Monoclonal Antibody (1X3C0)
MA542718	alpha Actinin 3 Recombinant Rabbit Monoclonal Antibody (4H7E6)
MA542721	Versican Recombinant Rabbit Monoclonal Antibody (6D8K9)
MA542723	E2F4 Recombinant Rabbit Monoclonal Antibody (3O6Y9)
MA542724	EDG2 Recombinant Rabbit Monoclonal Antibody (3U1Q0)
MA542725	ENT2 Recombinant Rabbit Monoclonal Antibody (7B2J3)
MA542726	SQSTM1 Recombinant Rabbit Monoclonal Antibody (5X4R1)
MA542727	Human IgA Recombinant Rabbit Monoclonal Antibody (10Q6M3)
MA542728	SOX9 Recombinant Rabbit Monoclonal Antibody (5Y4J8)
MA542730	Human IgG3 Recombinant Rabbit Monoclonal Antibody (6Q4J8)
MA542732	OTX1 Recombinant Rabbit Monoclonal Antibody (6H7J9)
MA542735	RORB Recombinant Rabbit Monoclonal Antibody (7J9T7)
MA542737	TMOD2 Recombinant Rabbit Monoclonal Antibody (5E7G5)
MA542738	Nesprin 1 Recombinant Rabbit Monoclonal Antibody (1D4D9)
MA542739	VPS11 Recombinant Rabbit Monoclonal Antibody (8L3D8)
MA542740	Talin 2 Recombinant Rabbit Monoclonal Antibody (7D1Y0)
MA542741	Argonaute 3 Recombinant Rabbit Monoclonal Antibody (8E0O6)
MA542742	GCSAM Recombinant Rabbit Monoclonal Antibody (7D1A5)
MA542743	SATB2 Recombinant Rabbit Monoclonal Antibody (6X1T3)
MA542771	SMURF2 Recombinant Rabbit Monoclonal Antibody (8O7F0)
MA542772	ACTA1 Recombinant Rabbit Monoclonal Antibody (8B9Q6)
MA542777	ACTA1 Recombinant Rabbit Monoclonal Antibody (8T8D0)
MA542787	Gelsolin Recombinant Rabbit Monoclonal Antibody (5O1O01)
MA542803	CDC27 Recombinant Rabbit Monoclonal Antibody (7U5A3)
MA542804	LSP1 Recombinant Rabbit Monoclonal Antibody (4U10N5)
MA542806	PRAS40 Recombinant Rabbit Monoclonal Antibody (2Y8Y3)
MA542812	TBK1 Recombinant Rabbit Monoclonal Antibody (4G6C7)
MA542815	CRMP5 Recombinant Rabbit Monoclonal Antibody (5L4S8)
MA542819	YBX1 Recombinant Rabbit Monoclonal Antibody (9J7M4)
MA542827	APC6 Recombinant Rabbit Monoclonal Antibody (7Y1Y10)
MA542836	MMP12 Recombinant Rabbit Monoclonal Antibody (1O3A4)

MA542842	LXR alpha Recombinant Rabbit Monoclonal Antibody (7Z4H10)
MA542846	CD3g Recombinant Rabbit Monoclonal Antibody (7J2L5)
MA542847	PFDN5 Recombinant Rabbit Monoclonal Antibody (9Y5H7)
MA542850	HAT1 Recombinant Rabbit Monoclonal Antibody (7U2O3)
MA542852	STAT4 Recombinant Rabbit Monoclonal Antibody (10H1009)
MA542859	DISC1 Recombinant Rabbit Monoclonal Antibody (4U5N1)
MA542862	IRF2 Recombinant Rabbit Monoclonal Antibody (8M2N4)
MA542863	PIK3R2 Recombinant Rabbit Monoclonal Antibody (6E2F5)
MA542865	Caveolin 3 Recombinant Rabbit Monoclonal Antibody (4G5W2)
MA542866	CDCA3 Recombinant Rabbit Monoclonal Antibody (1J5L2)
MA542867	FOXP3 Recombinant Rabbit Monoclonal Antibody (6K3D9)
MA542869	GLP-2 Recombinant Rabbit Monoclonal Antibody (3K3P5)
MA542870	Cytochrome P450 Reductase Recombinant Rabbit Monoclonal Antibody (9U1U7)
MA542871	PRMT6 Recombinant Rabbit Monoclonal Antibody (1Z10L6)
MA542874	IGF2BP2 Recombinant Rabbit Monoclonal Antibody (2O3W8)
MA542875	Activin A Recombinant Rabbit Monoclonal Antibody (1X6N9)
MA542876	JunD Recombinant Rabbit Monoclonal Antibody (5L7Y3)
MA542879	ULK3 Recombinant Rabbit Monoclonal Antibody (5Z7Z6)
MA542880	CEACAM6 Recombinant Rabbit Monoclonal Antibody (4Z3V7)
MA542881	NAP1L1 Recombinant Rabbit Monoclonal Antibody (4W8B6)
MA542890	SOX10 Recombinant Rabbit Monoclonal Antibody (2V5D6)
MA542904	MST2 Recombinant Rabbit Monoclonal Antibody (5J1I4)
MA542906	PDCD4 Recombinant Rabbit Monoclonal Antibody (3A8V10)
MA542918	JMJD2A Recombinant Rabbit Monoclonal Antibody (9A9G6)
MA542934	RALBP1 Recombinant Rabbit Monoclonal Antibody (5C10W9)
MA542937	EEF2 Recombinant Rabbit Monoclonal Antibody (4E2Z5)
MA542942	IQGAP1 Recombinant Rabbit Monoclonal Antibody (7I6X1)
MA542943	CD3d Recombinant Rabbit Monoclonal Antibody (2R1O10)
MA542950	Phospho-EGFR (Tyr1173) Recombinant Rabbit Monoclonal Antibody (4H3F8)
MA542951	Phospho-SMAD5 (Ser463, Ser465) Recombinant Rabbit Monoclonal Antibody (4D4X8)
MA542952	Phospho-FOXO3A (Ser253) Recombinant Rabbit Monoclonal Antibody (4K6R10)
MA542954	Phospho-BTK (Tyr223) Recombinant Rabbit Monoclonal Antibody (10F4W10)
MA542955	Phospho-VANGL2 (Ser5, Ser8) Recombinant Rabbit Monoclonal Antibody (2Y4L2)
MA542956	Phospho-VANGL2 (Ser11, Ser14) Recombinant Rabbit Monoclonal Antibody (4D3W2)
MA542957	Phospho-VANGL2 (Ser14, Ser17) Recombinant Rabbit Monoclonal Antibody (4C7Y1)
MA542958	Phospho-VANGL2 (Thr78, Ser79) Recombinant Rabbit Monoclonal Antibody (5K4D6)
MA535364	Flotillin 1 Recombinant Rabbit Monoclonal Antibody (ARC0765)
PA595882	PMF1 Polyclonal Antibody
PA592392	Plasma Kallikrein Polyclonal Antibody
PA590050	ZNF177 Polyclonal Antibody
PA587551	Synaptotagmin 1 Polyclonal Antibody
PA590378	Synaptotagmin 5 Polyclonal Antibody
PA587519	C/EBP alpha Polyclonal Antibody
PA589505	Dynein Polyclonal Antibody
PA589506	Dynein Polyclonal Antibody
PA590163	ZNF677 Polyclonal Antibody
PA590169	PDHA2 Polyclonal Antibody
PA587350	FGFR1 Polyclonal Antibody
PA587351	Perforin Polyclonal Antibody
PA587366	PIM1 Polyclonal Antibody
PA587363	alpha-II Spectrin Polyclonal Antibody
PA587371	BMI-1 Polyclonal Antibody
PA587381	DOK4 Polyclonal Antibody
PA587382	ERK2 Polyclonal Antibody
PA590781	NDRG3 Polyclonal Antibody
PA587387	IKK epsilon Polyclonal Antibody
PA587394	NEFL Polyclonal Antibody
PA587395	NME1 Polyclonal Antibody
PA587403	RhoA Polyclonal Antibody

PA587421	CD79a Polyclonal Antibody
PA587425	CD4 Polyclonal Antibody
PA587426	RARA Polyclonal Antibody
PA587434	IKK alpha Polyclonal Antibody
PA587435	PPIG Polyclonal Antibody
PA587436	N-cadherin Polyclonal Antibody
PA587438	MCL-1 Polyclonal Antibody
PA587439	XRCC1 Polyclonal Antibody
PA587443	PKC delta Polyclonal Antibody
PA587446	RelB Polyclonal Antibody
PA587447	IRF4 Polyclonal Antibody
PA587450	CDC25B Polyclonal Antibody
PA587458	MAX Polyclonal Antibody
PA587459	MYPT1 Polyclonal Antibody
PA587471	DcR3 Polyclonal Antibody
PA587472	DNM1L Polyclonal Antibody
PA587474	S100B Polyclonal Antibody
PA587485	SOCS3 Polyclonal Antibody
PA587486	MMP7 Polyclonal Antibody
PA587488	BMP-7 Polyclonal Antibody
PA587499	Cyclin D3 Polyclonal Antibody
PA587491	IKK beta Polyclonal Antibody
PA587500	JAK3 Polyclonal Antibody
PA587501	ErbB4 Polyclonal Antibody
PA587503	NBS1 Polyclonal Antibody
PA587505	CDX2 Polyclonal Antibody
PA587509	cIAP2 Polyclonal Antibody
PA587515	RPA70 Polyclonal Antibody
PA587520	PRKAR1A Polyclonal Antibody
PA587525	NSF Polyclonal Antibody
PA587527	MARCKS Polyclonal Antibody
PA587540	RAP1A Polyclonal Antibody
PA587542	SNX9 Polyclonal Antibody
PA587549	SNAP25 Polyclonal Antibody
PA587555	beta Arrestin 1 Polyclonal Antibody
PA587568	YAP1 Polyclonal Antibody
PA587574	RBM22 Polyclonal Antibody
PA587578	RBP1 Polyclonal Antibody
PA587612	PIAS1 Polyclonal Antibody
PA587638	Alpha-Smooth Muscle Actin Polyclonal Antibody
PA587645	KIR3DL2 Polyclonal Antibody
PA587670	ETHE1 Polyclonal Antibody
PA587696	CNPase Polyclonal Antibody
PA587702	PLD2 Polyclonal Antibody
PA587727	14-3-3 beta Polyclonal Antibody
PA587746	SGK1 Polyclonal Antibody
PA587755	14-3-3 sigma Polyclonal Antibody
PA587776	GADD45G Polyclonal Antibody
PA587786	LRIG1 Polyclonal Antibody
PA587829	B7-H3 Polyclonal Antibody
PA587831	STAR Polyclonal Antibody
PA587838	KIAA0101 Polyclonal Antibody
PA587843	RBM4 Polyclonal Antibody
PA587854	GCLC Polyclonal Antibody
PA587879	B7-H4 Polyclonal Antibody
PA587900	LMAN1 Polyclonal Antibody
PA587901	Thymidylate Synthase Polyclonal Antibody
PA587905	CHFR Polyclonal Antibody
PA587919	PHGDH Polyclonal Antibody

PA587924	CHRND Polyclonal Antibody
PA587925	Cathepsin F Polyclonal Antibody
PA587927	CYP2R1 Polyclonal Antibody
PA587930	AKR1C2 Polyclonal Antibody
PA587949	TIP47 Polyclonal Antibody
PA587964	DNMT2 Polyclonal Antibody
PA587966	BMP9 Polyclonal Antibody
PA587968	RAP1GAP Polyclonal Antibody
PA587969	HK1 Polyclonal Antibody
PA587971	TCEA1 Polyclonal Antibody
PA587974	LGR5 Polyclonal Antibody
PA587975	TRIM24 Polyclonal Antibody
PA587978	GIPC1 Polyclonal Antibody
PA587979	COPG Polyclonal Antibody
PA587981	GPSM2 Polyclonal Antibody
PA587982	BHLHE41 Polyclonal Antibody
PA587984	CDH4 Polyclonal Antibody
PA587991	RAB13 Polyclonal Antibody
PA587992	RCN1 Polyclonal Antibody
PA587993	SCN2A Polyclonal Antibody
PA587995	14-3-3 epsilon Polyclonal Antibody
PA587997	EKLF Polyclonal Antibody
PA587999	EAP30 Polyclonal Antibody
PA588001	USP24 Polyclonal Antibody
PA588002	SEC14L2 Polyclonal Antibody
PA588003	IL17C Polyclonal Antibody
PA588007	SRPRB Polyclonal Antibody
PA588016	GPHA2 Polyclonal Antibody
PA588020	NPFF2 Polyclonal Antibody
PA588024	COL4A1 Polyclonal Antibody
PA588025	RPL18 Polyclonal Antibody
PA588027	POLR3A Polyclonal Antibody
PA588029	MANEA Polyclonal Antibody
PA588034	ZP1 Polyclonal Antibody
PA588035	NUPL2 Polyclonal Antibody
PA588037	PSMD14 Polyclonal Antibody
PA588038	eIF1 Polyclonal Antibody
PA588039	RBM5 Polyclonal Antibody
PA588041	ARPC2 Polyclonal Antibody
PA588044	LHX2 Polyclonal Antibody
PA588045	FZD1 Polyclonal Antibody
PA588047	SF3A2 Polyclonal Antibody
PA588051	SDHB Polyclonal Antibody
PA588052	MYT1 Polyclonal Antibody
PA588054	MAN1A1 Polyclonal Antibody
PA588063	Cyclin D1 Polyclonal Antibody
PA588070	DNMT3B Polyclonal Antibody
PA588073	GPX1 Polyclonal Antibody
PA588087	NGF Receptor Polyclonal Antibody
PA588089	PI3K p85 alpha Polyclonal Antibody
PA588105	PD-L1 Polyclonal Antibody
PA588108	LC3B Polyclonal Antibody
PA588110	AMACR Polyclonal Antibody
PA588116	CHOP Polyclonal Antibody
PA588111	Cyclin D1 Polyclonal Antibody
PA588115	CTCF Polyclonal Antibody
PA588117	E2F8 Polyclonal Antibody
PA588119	TRAIL-R3 Polyclonal Antibody
PA588121	RAC2 Polyclonal Antibody

PA588125	TCF4 Polyclonal Antibody
PA588132	PYCARD Polyclonal Antibody
PA588135	Caspase 8 Polyclonal Antibody
PA588141	KEAP1 Polyclonal Antibody
PA588144	SMAD2 Polyclonal Antibody
PA588156	MPST Polyclonal Antibody
PA588154	A1BG Polyclonal Antibody
PA588157	HEPACAM Polyclonal Antibody
PA588159	ICT1 Polyclonal Antibody
PA588160	NAPRT1 Polyclonal Antibody
PA588161	RSU1 Polyclonal Antibody
PA588162	GREM1 Polyclonal Antibody
PA588163	CDK20 Polyclonal Antibody
PA588169	MAGEA4 Polyclonal Antibody
PA588165	HDJ2 Polyclonal Antibody
PA588167	FABP7 Polyclonal Antibody
PA588168	JMJD5 Polyclonal Antibody
PA588170	MAT2B Polyclonal Antibody
PA588171	MC3R Polyclonal Antibody
PA588172	Calmodulin 3 Polyclonal Antibody
PA588173	NAP1L1 Polyclonal Antibody
PA588174	SEC13 Polyclonal Antibody
PA588175	SEC61A1 Polyclonal Antibody
PA588177	EAP2 Polyclonal Antibody
PA588176	ST3GAL1 Polyclonal Antibody
PA588178	Aldolase C Polyclonal Antibody
PA588181	CHMP1A Polyclonal Antibody
PA588182	EXOSC3 Polyclonal Antibody
PA588184	NALP7 Polyclonal Antibody
PA588185	NT5C2 Polyclonal Antibody
PA588186	NUP85 Polyclonal Antibody
PA588188	PCDH10 Polyclonal Antibody
PA588189	PRIM2 Polyclonal Antibody
PA588190	RAB2A Polyclonal Antibody
PA588191	THOC5 Polyclonal Antibody
PA588192	GABRA3 Polyclonal Antibody
PA588194	SLC30A6 Polyclonal Antibody
PA588195	APOBEC3D Polyclonal Antibody
PA588196	MC3R Polyclonal Antibody
PA588199	SMPD2 Polyclonal Antibody
PA588200	TRIB2 Polyclonal Antibody
PA588201	CSTF3 Polyclonal Antibody
PA588202	MRPS27 Polyclonal Antibody
PA588204	GTPBP9 Polyclonal Antibody
PA588205	PRIM2 Polyclonal Antibody
PA588209	KCNK10 Polyclonal Antibody
PA588211	RPS18 Polyclonal Antibody
PA588212	SLC25A12 Polyclonal Antibody
PA588213	CHAMP1 Polyclonal Antibody
PA588215	SPR Polyclonal Antibody
PA588216	CDC123 Polyclonal Antibody
PA588217	KCNK10 Polyclonal Antibody
PA588219	DDX3Y Polyclonal Antibody
PA588220	HES1 Polyclonal Antibody
PA588221	HGK Polyclonal Antibody
PA588227	PP5 Polyclonal Antibody
PA588231	SPTLC2 Polyclonal Antibody
PA588232	THUMPD2 Polyclonal Antibody
PA588234	Calpain 1 Polyclonal Antibody

PA588235	NAC1 Polyclonal Antibody
PA588236	BRINP1 Polyclonal Antibody
PA588237	CADM2 Polyclonal Antibody
PA588238	AKAP7 Polyclonal Antibody
PA588239	MRPS9 Polyclonal Antibody
PA588240	SERPINI2 Polyclonal Antibody
PA588241	CrkL Polyclonal Antibody
PA588243	LCAT Polyclonal Antibody
PA588250	Catalase Polyclonal Antibody
PA588266	CD301 Polyclonal Antibody
PA588252	STAT5 alpha Polyclonal Antibody
PA588254	PEDF Polyclonal Antibody
PA588255	KLRB1 Polyclonal Antibody
PA588257	TGFBR2 Polyclonal Antibody
PA588260	RAB5 Polyclonal Antibody
PA588276	Cytokeratin 1 Polyclonal Antibody
PA588277	LNPEP Polyclonal Antibody
PA588281	ATG16L1 Polyclonal Antibody
PA588285	BYSL Polyclonal Antibody
PA588287	Hairless Polyclonal Antibody
PA588290	LTBR Polyclonal Antibody
PA588292	HEF1 Polyclonal Antibody
PA588293	TrkB Polyclonal Antibody
PA588295	ZSCAN4 Polyclonal Antibody
PA588296	ZSCAN4 Polyclonal Antibody
PA588300	GLG1 Polyclonal Antibody
PA588301	GLMN Polyclonal Antibody
PA588302	LRIG1 Polyclonal Antibody
PA588304	MYH10 Polyclonal Antibody
PA588307	NPTX2 Polyclonal Antibody
PA588306	NOVA1 Polyclonal Antibody
PA588309	POLR1B Polyclonal Antibody
PA588313	RPL38 Polyclonal Antibody
PA588314	SEMA7A Polyclonal Antibody
PA588315	LCN6 Polyclonal Antibody
PA588316	SPRR3 Polyclonal Antibody
PA588317	XPOT Polyclonal Antibody
PA588319	GABPA Polyclonal Antibody
PA588322	ANO6 Polyclonal Antibody
PA588323	POLK Polyclonal Antibody
PA588324	Kallikrein 6 Polyclonal Antibody
PA588327	PGCP Polyclonal Antibody
PA588329	DTNA Polyclonal Antibody
PA588330	DUSP3 Polyclonal Antibody
PA588331	GIGYF2 Polyclonal Antibody
PA588332	GSR Polyclonal Antibody
PA588334	HPS5 Polyclonal Antibody
PA588336	METTL2B Polyclonal Antibody
PA588337	PKLR Polyclonal Antibody
PA588341	RAB6A Polyclonal Antibody
PA588338	MRPS18B Polyclonal Antibody
PA588344	RAB21 Polyclonal Antibody
PA588345	RANBP9 Polyclonal Antibody
PA588346	MXI1 Polyclonal Antibody
PA588348	SEC23A Polyclonal Antibody
PA588349	PAI1 Polyclonal Antibody
PA588351	TATDN1 Polyclonal Antibody
PA588352	LILRA1 Polyclonal Antibody
PA588357	NDUFA8 Polyclonal Antibody

PA588354	PCPTP1 Polyclonal Antibody
PA588356	Muskelin Polyclonal Antibody
PA588360	RIOK2 Polyclonal Antibody
PA588361	RPS24 Polyclonal Antibody
PA588362	SART3 Polyclonal Antibody
PA588364	ACAP2 Polyclonal Antibody
PA588366	DNAI1 Polyclonal Antibody
PA588367	NGRN Polyclonal Antibody
PA588368	NALP11 Polyclonal Antibody
PA588369	PDCD7 Polyclonal Antibody
PA588372	SGK2 Polyclonal Antibody
PA588373	SLC25A24 Polyclonal Antibody
PA588377	DNAJC10 Polyclonal Antibody
PA588378	TCP1 eta Polyclonal Antibody
PA588379	CHD9 Polyclonal Antibody
PA588380	DTL Polyclonal Antibody
PA588383	LILRB2 Polyclonal Antibody
PA588384	NMB Polyclonal Antibody
PA588387	SNRPA1 Polyclonal Antibody
PA588388	SNRPF Polyclonal Antibody
PA588389	WDR1 Polyclonal Antibody
PA588392	GPN1 Polyclonal Antibody
PA588393	MED21 Polyclonal Antibody
PA588397	PPIH Polyclonal Antibody
PA588400	AGTRAP Polyclonal Antibody
PA588399	CDH11 Polyclonal Antibody
PA588401	CoA Synthase Polyclonal Antibody
PA588403	FAU Polyclonal Antibody
PA588404	KCNE4 Polyclonal Antibody
PA588405	PDE8A Polyclonal Antibody
PA588407	WSB1 Polyclonal Antibody
PA588408	ATL3 Polyclonal Antibody
PA588409	CNPY2 Polyclonal Antibody
PA588410	EFHA1 Polyclonal Antibody
PA588413	Cathepsin L Polyclonal Antibody
PA588420	DNAL4 Polyclonal Antibody
PA588421	ITPA Polyclonal Antibody
PA588422	ESAM Polyclonal Antibody
PA588425	HOXC10 Polyclonal Antibody
PA588424	H-Ras Polyclonal Antibody
PA588426	IP6K1 Polyclonal Antibody
PA588427	KCNK4 Polyclonal Antibody
PA588428	LAX1 Polyclonal Antibody
PA588429	MKRN2 Polyclonal Antibody
PA588433	NCKAP1 Polyclonal Antibody
PA588435	ZFP28 Polyclonal Antibody
PA588436	PADI2 Polyclonal Antibody
PA588437	RRBP1 Polyclonal Antibody
PA588439	SHARPIN Polyclonal Antibody
PA588440	VRK2 Polyclonal Antibody
PA588441	GABRA1 Polyclonal Antibody
PA588442	PCPTP1 Polyclonal Antibody
PA588443	ALDH18A1 Polyclonal Antibody
PA588445	PORCN Polyclonal Antibody
PA588446	PPT2 Polyclonal Antibody
PA588447	SIPA1 Polyclonal Antibody
PA588451	DNAJC10 Polyclonal Antibody
PA588454	HSF2 Polyclonal Antibody
PA588456	CIP2A Polyclonal Antibody

PA588458	AP2S1 Polyclonal Antibody
PA588461	TSNAX Polyclonal Antibody
PA588464	AMPK alpha-1 Polyclonal Antibody
PA588469	RNF10 Polyclonal Antibody
PA588472	KCC4 Polyclonal Antibody
PA588474	Syntaphilin Polyclonal Antibody
PA588475	GABPA Polyclonal Antibody
PA588476	Prostaglandin E Receptor EP1 Polyclonal Antibody
PA588481	CHRM1 Polyclonal Antibody
PA588482	TrkB Polyclonal Antibody
PA588484	Dot1L Polyclonal Antibody
PA588487	CD249 Polyclonal Antibody
PA588490	SBF1 Polyclonal Antibody
PA588491	SULT1A3 Polyclonal Antibody
PA588492	Tetranectin Polyclonal Antibody
PA588493	UBE1 Polyclonal Antibody
PA588495	RFPL3 Polyclonal Antibody
PA588497	MTBP Polyclonal Antibody
PA588498	TRPM5 Polyclonal Antibody
PA588499	A4GNT Polyclonal Antibody
PA588501	CPSF3 Polyclonal Antibody
PA588502	CPVL Polyclonal Antibody
PA588504	CPVL Polyclonal Antibody
PA588506	AID Polyclonal Antibody
PA588507	SLC25A19 Polyclonal Antibody
PA588508	EXOC4 Polyclonal Antibody
PA588509	TEKT5 Polyclonal Antibody
PA588510	TEKT5 Polyclonal Antibody
PA588512	BAI2 Polyclonal Antibody
PA588516	MCT2 Polyclonal Antibody
PA588534	Clathrin Heavy Chain Polyclonal Antibody
PA588526	ATP1B1 Polyclonal Antibody
PA588527	beta-2 Microglobulin Polyclonal Antibody
PA588538	Ataxin 3 Polyclonal Antibody
PA588558	ASK1 Polyclonal Antibody
PA588566	PP1 alpha Polyclonal Antibody
PA588578	SHP2 Polyclonal Antibody
PA588581	SHMT1 Polyclonal Antibody
PA588580	RAG2 Polyclonal Antibody
PA588582	VAMP2 Polyclonal Antibody
PA588585	PEX5 Polyclonal Antibody
PA588591	AKT2 Polyclonal Antibody
PA588606	COX2 Polyclonal Antibody
PA588614	FUBP1 Polyclonal Antibody
PA588670	EVC Polyclonal Antibody
PA588673	CXCL2 Polyclonal Antibody
PA588675	NUCB2 Polyclonal Antibody
PA588676	NUCB2 Polyclonal Antibody
PA588677	PRPS2 Polyclonal Antibody
PA588678	RAG1 Polyclonal Antibody
PA588680	DEGS1 Polyclonal Antibody
PA588692	HELQ Polyclonal Antibody
PA588686	WWOX Polyclonal Antibody
PA588688	BRWD1 Polyclonal Antibody
PA588698	PTX3 Polyclonal Antibody
PA588704	MEK1 Polyclonal Antibody
PA588705	Nkx2.5 Polyclonal Antibody
PA588706	PRDX5 Polyclonal Antibody
PA588707	TFPI Polyclonal Antibody

PA588708	HMGA1 Polyclonal Antibody
PA588710	Cullin 4B Polyclonal Antibody
PA588715	CD248 Polyclonal Antibody
PA588718	GUCY1A3 Polyclonal Antibody
PA588720	MDC1 Polyclonal Antibody
PA588721	TMEM27 Polyclonal Antibody
PA588722	LHX9 Polyclonal Antibody
PA588724	NTCP Polyclonal Antibody
PA588731	LSM1 Polyclonal Antibody
PA588726	GCS1 Polyclonal Antibody
PA588734	TRIM45 Polyclonal Antibody
PA588736	COPS8 Polyclonal Antibody
PA588742	Neuropilin 2 Polyclonal Antibody
PA588748	CLPS Polyclonal Antibody
PA588749	HNF6 Polyclonal Antibody
PA588750	FZD5 Polyclonal Antibody
PA588754	LIFR Polyclonal Antibody
PA588756	PDIA2 Polyclonal Antibody
PA588758	ITPR3 Polyclonal Antibody
PA588759	ACOX2 Polyclonal Antibody
PA588760	AKAP6 Polyclonal Antibody
PA588761	ANGEL2 Polyclonal Antibody
PA588762	ARL10 Polyclonal Antibody
PA588763	BAAT1 Polyclonal Antibody
PA588764	F1A alpha Polyclonal Antibody
PA588766	Glypican 4 Polyclonal Antibody
PA588768	ITIH4 Polyclonal Antibody
PA588769	LHFPL2 Polyclonal Antibody
PA588771	NMUR2 Polyclonal Antibody
PA588772	NRAP Polyclonal Antibody
PA588774	PARP10 Polyclonal Antibody
PA588776	PRPF18 Polyclonal Antibody
PA588777	SLA2 Polyclonal Antibody
PA588778	SLC19A1 Polyclonal Antibody
PA588779	SLC44A3 Polyclonal Antibody
PA588780	SNX13 Polyclonal Antibody
PA588781	TCL1B Polyclonal Antibody
PA588785	SEMG2 Polyclonal Antibody
PA588786	CD85j (ILT2) Polyclonal Antibody
PA588789	CDK11B Polyclonal Antibody
PA588792	GBA2 Polyclonal Antibody
PA588795	MBNL2 Polyclonal Antibody
PA588793	KIR2DL4 Polyclonal Antibody
PA588794	MAN1C1 Polyclonal Antibody
PA588796	MCF2 Polyclonal Antibody
PA588797	MLC1 Polyclonal Antibody
PA588799	PPM1J Polyclonal Antibody
PA588800	RCOR1 Polyclonal Antibody
PA588804	HLA-B Polyclonal Antibody
PA588806	SULT1A2 Polyclonal Antibody
PA588808	TARS2 Polyclonal Antibody
PA588810	TP53I11 Polyclonal Antibody
PA588811	Kaiso Polyclonal Antibody
PA588812	CRYZ Polyclonal Antibody
PA588813	NDUFS2 Polyclonal Antibody
PA588816	TLX1 Polyclonal Antibody
PA588821	FTMT Polyclonal Antibody
PA588822	GNAL Polyclonal Antibody
PA588823	GSTT1 Polyclonal Antibody

PA588826	LHFPL2 Polyclonal Antibody
PA588827	NRARP Polyclonal Antibody
PA588828	CRTH2 Polyclonal Antibody
PA588830	VGLUT1 Polyclonal Antibody
PA588831	LIF Polyclonal Antibody
PA588833	SYNPO2 Polyclonal Antibody
PA588834	TDRD1 Polyclonal Antibody
PA588835	C17orf48 Polyclonal Antibody
PA588836	XRCC6BP1 Polyclonal Antibody
PA588837	CLN5 Polyclonal Antibody
PA588839	FBXO31 Polyclonal Antibody
PA588844	Parp8 Polyclonal Antibody
PA588845	PNLIPRP2 Polyclonal Antibody
PA588846	RIPX Polyclonal Antibody
PA588847	SCGN Polyclonal Antibody
PA588849	THAP3 Polyclonal Antibody
PA588851	ZNF268 Polyclonal Antibody
PA588853	PLCH1 Polyclonal Antibody
PA588855	SMC1B Polyclonal Antibody
PA588857	AKT3 Polyclonal Antibody
PA588858	AVEN Polyclonal Antibody
PA588859	CHCHD6 Polyclonal Antibody
PA588861	DKC1 Polyclonal Antibody
PA588862	DNAJB11 Polyclonal Antibody
PA588863	MAZ Polyclonal Antibody
PA588865	RAP1B Polyclonal Antibody
PA588866	SRP14 Polyclonal Antibody
PA588867	YY1 Polyclonal Antibody
PA588872	FUT8 Polyclonal Antibody
PA588873	MTCH2 Polyclonal Antibody
PA588875	TRIM25 Polyclonal Antibody
PA588876	VAPA Polyclonal Antibody
PA588880	DLL4 Polyclonal Antibody
PA588882	HSD17B8 Polyclonal Antibody
PA588881	ECH1 Polyclonal Antibody
PA588883	Kir2.1 (KCNJ2) Polyclonal Antibody
PA588884	PGA5 Polyclonal Antibody
PA588887	SFXN1 Polyclonal Antibody
PA588888	GMFG Polyclonal Antibody
PA588890	LIN28B Polyclonal Antibody
PA588893	MARCH8 Polyclonal Antibody
PA588894	MSP Polyclonal Antibody
PA588895	PITPN Polyclonal Antibody
PA588896	SEMA3A Polyclonal Antibody
PA588898	WSB2 Polyclonal Antibody
PA588902	GBP2 Polyclonal Antibody
PA588904	RANK Polyclonal Antibody
PA588905	NKCC2 Polyclonal Antibody
PA588906	RPL32 Polyclonal Antibody
PA588908	CTU2 Polyclonal Antibody
PA588911	ARR3 Polyclonal Antibody
PA588912	FMNL1 Polyclonal Antibody
PA588914	CK1 gamma-3 Polyclonal Antibody
PA588915	DECR1 Polyclonal Antibody
PA588918	FOXF1 Polyclonal Antibody
PA588919	Filamin 2 Polyclonal Antibody
PA588920	Glypican 1 Polyclonal Antibody
PA588921	GYS2 Polyclonal Antibody
PA588922	HAL Polyclonal Antibody

PA588923	IGF2R Polyclonal Antibody
PA588924	IGFBP7 Polyclonal Antibody
PA588927	RAP Polyclonal Antibody
PA588928	CRHBP Polyclonal Antibody
PA588929	MTNR1A Polyclonal Antibody
PA588931	MYO1E Polyclonal Antibody
PA588933	OVGP1 Polyclonal Antibody
PA588934	PAFAH1B3 Polyclonal Antibody
PA588935	PPP1R7 Polyclonal Antibody
PA588936	PPP1R7 Polyclonal Antibody
PA588938	RPL27 Polyclonal Antibody
PA588939	S-arrestin Polyclonal Antibody
PA588940	p38 MAPK gamma Polyclonal Antibody
PA588941	SORL1 Polyclonal Antibody
PA588944	TBCA Polyclonal Antibody
PA588945	UPP1 Polyclonal Antibody
PA588948	IFT88 Polyclonal Antibody
PA588949	BFSP2 Polyclonal Antibody
PA588950	ITGA8 Polyclonal Antibody
PA588952	MPZL1 Polyclonal Antibody
PA588956	KNTC1 Polyclonal Antibody
PA588957	PTDSS1 Polyclonal Antibody
PA588958	GAB2 Polyclonal Antibody
PA588959	ACOT8 Polyclonal Antibody
PA588960	PRMT3 Polyclonal Antibody
PA588961	BET1 Polyclonal Antibody
PA588962	IFITM3 Polyclonal Antibody
PA588963	UBAC1 Polyclonal Antibody
PA588964	TRAFD1 Polyclonal Antibody
PA588965	GCN1L1 Polyclonal Antibody
PA588966	LAR4B Polyclonal Antibody
PA588967	NUP160 Polyclonal Antibody
PA588969	VPS39 Polyclonal Antibody
PA588977	CD266 Polyclonal Antibody
PA588974	AADAT Polyclonal Antibody
PA588978	PPME1 Polyclonal Antibody
PA588980	SF3B14 Polyclonal Antibody
PA588981	VPS29 Polyclonal Antibody
PA588982	Mig-6 Polyclonal Antibody
PA588984	TRIT1 Polyclonal Antibody
PA588985	MSL2 Polyclonal Antibody
PA588986	YY1AP1 Polyclonal Antibody
PA588987	RCC2 Polyclonal Antibody
PA588988	GKN1 Polyclonal Antibody
PA588989	Casc5 Polyclonal Antibody
PA588991	MURF3 Polyclonal Antibody
PA588992	NLN Polyclonal Antibody
PA588993	MKL2 Polyclonal Antibody
PA588994	LRRG7 Polyclonal Antibody
PA588995	USP36 Polyclonal Antibody
PA588996	CWC22 Polyclonal Antibody
PA588999	RBM26 Polyclonal Antibody
PA589001	TINAGL1 Polyclonal Antibody
PA589002	NCAPG Polyclonal Antibody
PA589000	GOLPH3 Polyclonal Antibody
PA589003	Mtmr9 Polyclonal Antibody
PA589005	DYNC2H1 Polyclonal Antibody
PA589006	SCD5 Polyclonal Antibody
PA589007	CDK5RAP3 Polyclonal Antibody

PA589008	REG4 Polyclonal Antibody
PA589009	KIRREL2 Polyclonal Antibody
PA589011	VPS25 Polyclonal Antibody
PA589012	IL1F10 Polyclonal Antibody
PA589013	TSLP Polyclonal Antibody
PA589015	RAVER1 Polyclonal Antibody
PA589016	CHCHD4 Polyclonal Antibody
PA589017	ZNF599 Polyclonal Antibody
PA589018	Munc13-4 Polyclonal Antibody
PA589019	MRGPRX1 Polyclonal Antibody
PA589020	FMN1 Polyclonal Antibody
PA589021	RAD26L Polyclonal Antibody
PA589022	NME1/NME2 Polyclonal Antibody
PA589027	DNASE1L1 Polyclonal Antibody
PA589029	MC1R Polyclonal Antibody
PA589032	UNG Polyclonal Antibody
PA589033	ZP3 Polyclonal Antibody
PA589036	RBM14 Polyclonal Antibody
PA589038	SNX5 Polyclonal Antibody
PA589039	PARVB Polyclonal Antibody
PA589041	PFTK2 Polyclonal Antibody
PA589042	SLC5A11 Polyclonal Antibody
PA589047	PATL1 Polyclonal Antibody
PA589048	POLR2F Polyclonal Antibody
PA589049	Cathepsin G Polyclonal Antibody
PA589050	DGCR14 Polyclonal Antibody
PA589051	SCYL2 Polyclonal Antibody
PA589052	SCYL2 Polyclonal Antibody
PA589056	TCEAL2 Polyclonal Antibody
PA589058	CAPS1 Polyclonal Antibody
PA589061	LSM3 Polyclonal Antibody
PA589062	ZC3H7A Polyclonal Antibody
PA589063	GNL2 Polyclonal Antibody
PA589067	CRISPLD2 Polyclonal Antibody
PA589068	C9orf89 Polyclonal Antibody
PA589069	MRPL45 Polyclonal Antibody
PA589070	OSBPL10 Polyclonal Antibody
PA589071	NTN5 Polyclonal Antibody
PA589072	OAF Polyclonal Antibody
PA589073	PIWIL3 Polyclonal Antibody
PA589074	ENT1 Polyclonal Antibody
PA589075	ENT1 Polyclonal Antibody
PA589083	WBP1 Polyclonal Antibody
PA589077	GLE1 Polyclonal Antibody
PA589079	GCNT3 Polyclonal Antibody
PA589084	WDSOF1 Polyclonal Antibody
PA589085	GPKOW Polyclonal Antibody
PA589086	VPS54 Polyclonal Antibody
PA589087	RBM11 Polyclonal Antibody
PA589088	SLC25A38 Polyclonal Antibody
PA589090	FKBP14 Polyclonal Antibody
PA589091	LARP2 Polyclonal Antibody
PA589100	AROS Polyclonal Antibody
PA589095	GABRQ Polyclonal Antibody
PA589097	CTAGE1 Polyclonal Antibody
PA589101	TTC5 Polyclonal Antibody
PA589102	GTF3C6 Polyclonal Antibody
PA589103	ARL11 Polyclonal Antibody
PA589104	SFRS12 Polyclonal Antibody

PA589105	SYNE3 Polyclonal Antibody
PA589106	GLD2 Polyclonal Antibody
PA589108	GPIHBP1 Polyclonal Antibody
PA589110	Complement Factor B Polyclonal Antibody
PA589113	DUSP12 Polyclonal Antibody
PA589112	IDH1 Polyclonal Antibody
PA589118	UMPS Polyclonal Antibody
PA589125	TSPAN8 Polyclonal Antibody
PA589126	YOD1 Polyclonal Antibody
PA589127	ACAT1 Polyclonal Antibody
PA589138	Cathepsin D Polyclonal Antibody
PA589139	DBI Polyclonal Antibody
PA589182	alpha Synuclein Polyclonal Antibody
PA589200	Apolipoprotein L1 Polyclonal Antibody
PA589232	PEO1 Polyclonal Antibody
PA589235	TBLR1 Polyclonal Antibody
PA589249	Angiogenin Polyclonal Antibody
PA589258	CYP51A1 Polyclonal Antibody
PA589275	Lysozyme Polyclonal Antibody
PA589271	Galectin 3 Polyclonal Antibody
PA589273	CD91 Polyclonal Antibody
PA589274	CD44 Polyclonal Antibody
PA589315	UGT1A1 Polyclonal Antibody
PA589318	APOBEC3G Polyclonal Antibody
PA589319	GBL Polyclonal Antibody
PA589322	Tyrosinase Polyclonal Antibody
PA589332	Albumin Polyclonal Antibody
PA589331	PKC beta-1,2 Polyclonal Antibody
PA589347	GLUT5 Polyclonal Antibody
PA589352	UBE2L6 Polyclonal Antibody
PA589359	NR5A2 Polyclonal Antibody
PA589360	RPL30 Polyclonal Antibody
PA589361	4E-BP2 Polyclonal Antibody
PA589362	Profilin 4 Polyclonal Antibody
PA589363	NDUFB7 Polyclonal Antibody
PA589364	CPLX3 Polyclonal Antibody
PA589365	GCSH Polyclonal Antibody
PA589366	CRYAB Polyclonal Antibody
PA589368	TNFAIP8L2 Polyclonal Antibody
PA589369	DENR Polyclonal Antibody
PA589370	SARNP Polyclonal Antibody
PA589372	CABP7 Polyclonal Antibody
PA589374	RND1 Polyclonal Antibody
PA589375	GEMIN8 Polyclonal Antibody
PA589376	MDFI Polyclonal Antibody
PA589377	IAH1 Polyclonal Antibody
PA589378	TMEM70 Polyclonal Antibody
PA589383	PIH1D1 Polyclonal Antibody
PA589379	RPL7A Polyclonal Antibody
PA589380	CELA3B Polyclonal Antibody
PA589384	LIME Polyclonal Antibody
PA589385	CAMLG Polyclonal Antibody
PA589387	RNF144B Polyclonal Antibody
PA589388	DPPA4 Polyclonal Antibody
PA589390	DOK5 Polyclonal Antibody
PA589394	SLC25A14 Polyclonal Antibody
PA589396	CHST11 Polyclonal Antibody
PA589398	Centaurin alpha-2 Polyclonal Antibody
PA589399	Centaurin alpha-2 Polyclonal Antibody

PA589404	SERPINB13 Polyclonal Antibody
PA589401	KLDC3 Polyclonal Antibody
PA589402	IDH3B Polyclonal Antibody
PA589403	KLF17 Polyclonal Antibody
PA589409	CMAS Polyclonal Antibody
PA589410	ZC3H10 Polyclonal Antibody
PA589411	SC65 Polyclonal Antibody
PA589412	FGL2 Polyclonal Antibody
PA589415	IFIT2 Polyclonal Antibody
PA589416	Adenylate Kinase 8 Polyclonal Antibody
PA589418	EDC3 Polyclonal Antibody
PA589420	HSF1 Polyclonal Antibody
PA589421	EHD4 Polyclonal Antibody
PA589422	P4HA3 Polyclonal Antibody
PA589425	KLHL3 Polyclonal Antibody
PA589426	KARS Polyclonal Antibody
PA589431	ZNF398 Polyclonal Antibody
PA589438	ELMO2 Polyclonal Antibody
PA589440	ASCC2 Polyclonal Antibody
PA589444	XAB2 Polyclonal Antibody
PA589445	SULF1 Polyclonal Antibody
PA589449	FcRL5 Polyclonal Antibody
PA589450	FNBP4 Polyclonal Antibody
PA589451	CARD6 Polyclonal Antibody
PA589452	RPS6KC1 Polyclonal Antibody
PA589455	ZMYM4 Polyclonal Antibody
PA589457	GCC2 Polyclonal Antibody
PA589458	PSME4 Polyclonal Antibody
PA589459	HYD Polyclonal Antibody
PA589460	URM1 Polyclonal Antibody
PA589461	NDUFB4 Polyclonal Antibody
PA589462	NXT2 Polyclonal Antibody
PA589463	LYZL1 Polyclonal Antibody
PA589464	Histone H3.3 Polyclonal Antibody
PA589465	MED22 Polyclonal Antibody
PA589467	STX10 Polyclonal Antibody
PA589468	GLYAT Polyclonal Antibody
PA589472	HS2ST1 Polyclonal Antibody
PA589473	EFCAB4B Polyclonal Antibody
PA589474	EFCAB4B Polyclonal Antibody
PA589476	Apolipoprotein L3 Polyclonal Antibody
PA589477	DRB1 Polyclonal Antibody
PA589480	Bora Polyclonal Antibody
PA589481	ZNF333 Polyclonal Antibody
PA589482	CPT1C Polyclonal Antibody
PA589483	CDHR5 Polyclonal Antibody
PA589484	GANAB Polyclonal Antibody
PA589485	GPAA1 Polyclonal Antibody
PA589486	HOXC4 Polyclonal Antibody
PA589489	ETV4 Polyclonal Antibody
PA589490	PLA2G7 Polyclonal Antibody
PA589492	PMVK Polyclonal Antibody
PA589493	KIF19 Polyclonal Antibody
PA589494	LETM2 Polyclonal Antibody
PA589497	VSIG1 Polyclonal Antibody
PA589498	RASEF Polyclonal Antibody
PA589499	eIF4E3 Polyclonal Antibody
PA589500	CPSF7 Polyclonal Antibody
PA589501	KRI1 Polyclonal Antibody

PA589502	BITE Polyclonal Antibody
PA589503	RNaseH2C Polyclonal Antibody
PA589504	SPSB1 Polyclonal Antibody
PA589510	NRIP2 Polyclonal Antibody
PA589513	TRIM9 Polyclonal Antibody
PA589514	SAP130 Polyclonal Antibody
PA589518	ATP13A1 Polyclonal Antibody
PA589519	STARD5 Polyclonal Antibody
PA589522	PKM2 Polyclonal Antibody
PA589523	HIC1 Polyclonal Antibody
PA589525	MIF Polyclonal Antibody
PA589528	VAMP8 Polyclonal Antibody
PA589536	CD34 Polyclonal Antibody
PA589537	MESDC2 Polyclonal Antibody
PA589539	SMC4 Polyclonal Antibody
PA589547	ASGR2 Polyclonal Antibody
PA589572	IGF2 Polyclonal Antibody
PA589574	FIBB Polyclonal Antibody
PA589580	SMAD5 Polyclonal Antibody
PA589582	MRE11 Polyclonal Antibody
PA589601	PTPN22 Polyclonal Antibody
PA589637	NOP16 Polyclonal Antibody
PA589657	SLC6A4 Polyclonal Antibody
PA589661	GAS7 Polyclonal Antibody
PA589675	ORP2 Polyclonal Antibody
PA589676	AGTR1 Polyclonal Antibody
PA589679	WDR33 Polyclonal Antibody
PA589681	PIM1 Polyclonal Antibody
PA589685	NGF Polyclonal Antibody
PA589694	EphA5 Polyclonal Antibody
PA589695	GABRG1 Polyclonal Antibody
PA589697	PGCF Polyclonal Antibody
PA589699	ITGA7 Polyclonal Antibody
PA589703	PCDH8 Polyclonal Antibody
PA589704	PPP2R5B Polyclonal Antibody
PA589705	PTPRG Polyclonal Antibody
PA589706	RPL21 Polyclonal Antibody
PA589707	RPL24 Polyclonal Antibody
PA589708	TSPAN31 Polyclonal Antibody
PA589713	DR3 Polyclonal Antibody
PA589716	GNMB Polyclonal Antibody
PA589717	MAN1A2 Polyclonal Antibody
PA589720	Contactin 6 Polyclonal Antibody
PA589721	MRPS2 Polyclonal Antibody
PA589723	KIRREL Polyclonal Antibody
PA589724	DCAF6 Polyclonal Antibody
PA589727	Neuroigin 2 Polyclonal Antibody
PA589728	ZNF323 Polyclonal Antibody
PA589730	Nop25 Polyclonal Antibody
PA589731	FLAD1 Polyclonal Antibody
PA589733	MYLK2 Polyclonal Antibody
PA589734	NDUFAF2 Polyclonal Antibody
PA589738	RXFP2 Polyclonal Antibody
PA589740	FBXO33 Polyclonal Antibody
PA589741	SUMF1 Polyclonal Antibody
PA589747	IL17RE Polyclonal Antibody
PA589748	LPCAT4 Polyclonal Antibody
PA589749	PARP12 Polyclonal Antibody
PA589752	RBMX2 Polyclonal Antibody

PA589753	ZNF157 Polyclonal Antibody
PA589757	RABL3 Polyclonal Antibody
PA589758	RRP36 Polyclonal Antibody
PA589768	LHPP Polyclonal Antibody
PA589769	MOCOS Polyclonal Antibody
PA589770	MRI1 Polyclonal Antibody
PA589771	NDUFAF4 Polyclonal Antibody
PA589772	PREPL Polyclonal Antibody
PA589776	ZMAT2 Polyclonal Antibody
PA589780	CAP2 Polyclonal Antibody
PA589782	GALM Polyclonal Antibody
PA589786	VAMP5 Polyclonal Antibody
PA589784	SLMAP Polyclonal Antibody
PA589785	TMEM25 Polyclonal Antibody
PA589791	FTO Polyclonal Antibody
PA589795	ZNF31 Polyclonal Antibody
PA589796	PDIR Polyclonal Antibody
PA589800	TRIM2 Polyclonal Antibody
PA589805	MRPL2 Polyclonal Antibody
PA589806	RNF181 Polyclonal Antibody
PA589808	G2E3 Polyclonal Antibody
PA589809	ZNF416 Polyclonal Antibody
PA589813	SCRN3 Polyclonal Antibody
PA589814	DHRS11 Polyclonal Antibody
PA589818	PIGZ Polyclonal Antibody
PA589820	B3GNT5 Polyclonal Antibody
PA589821	PSD2 Polyclonal Antibody
PA589827	CCDC124 Polyclonal Antibody
PA589830	ASCL4 Polyclonal Antibody
PA589832	V-ATPase G3 Polyclonal Antibody
PA589841	ZDHHC15 Polyclonal Antibody
PA589849	DDI1 Polyclonal Antibody
PA589851	TRH Polyclonal Antibody
PA589852	alpha TTP Polyclonal Antibody
PA589858	BXDC2 Polyclonal Antibody
PA589859	PEBP4 Polyclonal Antibody
PA589860	RNASE3 Polyclonal Antibody
PA589862	gamma Synuclein Polyclonal Antibody
PA589864	SPARC Polyclonal Antibody
PA589868	eIF2 beta Polyclonal Antibody
PA589869	CD335 (NKp46) Polyclonal Antibody
PA589871	ARPC1A Polyclonal Antibody
PA589872	Nicastrin Polyclonal Antibody
PA589873	PES1 Polyclonal Antibody
PA589874	PRAME Polyclonal Antibody
PA589875	ZAK Polyclonal Antibody
PA589876	KRBOX4 Polyclonal Antibody
PA589877	SMYD3 Polyclonal Antibody
PA589879	NALP12 Polyclonal Antibody
PA589881	TPH2 Polyclonal Antibody
PA589886	Calmodulin 3 Polyclonal Antibody
PA589891	GSS Polyclonal Antibody
PA589892	HCK Polyclonal Antibody
PA589894	SERPINE2 Polyclonal Antibody
PA589896	ASH2L Polyclonal Antibody
PA589897	NGDN Polyclonal Antibody
PA589898	TOMM22 Polyclonal Antibody
PA589901	TUBG1 Polyclonal Antibody
PA589905	SLC47A1 Polyclonal Antibody

PA589906	CHRM3 Polyclonal Antibody
PA589907	Ephrin B1 Polyclonal Antibody
PA589909	Gastrin Polyclonal Antibody
PA589910	SLC7A2 Polyclonal Antibody
PA589911	FAM13A Polyclonal Antibody
PA589912	CAP1 Polyclonal Antibody
PA589915	HIGD1A Polyclonal Antibody
PA589921	FBXL12 Polyclonal Antibody
PA589924	ECHDC2 Polyclonal Antibody
PA589925	GLUT9 Polyclonal Antibody
PA589926	Coronin 1B Polyclonal Antibody
PA589928	ABCG4 Polyclonal Antibody
PA589929	LZTS2 Polyclonal Antibody
PA589932	CDSN Polyclonal Antibody
PA589933	SRSF6 Polyclonal Antibody
PA589952	LACRT Polyclonal Antibody
PA589934	RFPL2 Polyclonal Antibody
PA589938	Factor XIII A Polyclonal Antibody
PA589939	PHC2 Polyclonal Antibody
PA589941	MTRR Polyclonal Antibody
PA589942	MYH13 Polyclonal Antibody
PA589943	eIF4G3 Polyclonal Antibody
PA589944	DRD3 Polyclonal Antibody
PA589947	DNAJB12 Polyclonal Antibody
PA589948	ESRP1 Polyclonal Antibody
PA589951	TM7SF4 Polyclonal Antibody
PA589957	GDNF Polyclonal Antibody
PA589959	UBE2D3 Polyclonal Antibody
PA589979	Importin 8 Polyclonal Antibody
PA589963	IDH3A Polyclonal Antibody
PA589973	COX15 Polyclonal Antibody
PA589977	PP14B Polyclonal Antibody
PA589978	LZTR2 Polyclonal Antibody
PA589980	Importin 11 Polyclonal Antibody
PA589987	B4GALT4 Polyclonal Antibody
PA589991	CD36 Polyclonal Antibody
PA589994	V-ATPase A Polyclonal Antibody
PA589995	V-ATPase A Polyclonal Antibody
PA589996	Carbonic Anhydrase VA Polyclonal Antibody
PA589997	CACNB3 Polyclonal Antibody
PA589998	MTG16 Polyclonal Antibody
PA590013	GPR6 Polyclonal Antibody
PA590007	FGF4 Polyclonal Antibody
PA590011	GNAQ Polyclonal Antibody
PA590019	BCAM Polyclonal Antibody
PA590020	BCAM Polyclonal Antibody
PA590023	NFE2L1 Polyclonal Antibody
PA590024	ORC1 Polyclonal Antibody
PA590025	PAI1 Polyclonal Antibody
PA590029	P-Glycoprotein Polyclonal Antibody
PA590030	PLXNB3 Polyclonal Antibody
PA590031	PPT1 Polyclonal Antibody
PA590032	MEK2 Polyclonal Antibody
PA590036	Radixin Polyclonal Antibody
PA590035	PTPRN2 Polyclonal Antibody
PA590039	SLC7A1 Polyclonal Antibody
PA590041	TAL2 Polyclonal Antibody
PA590043	TXK Polyclonal Antibody
PA590055	DDX3Y Polyclonal Antibody

PA590061	EDG4 Polyclonal Antibody
PA590070	SH2D3A Polyclonal Antibody
PA590073	G3BP1 Polyclonal Antibody
PA590076	ERLIN1 Polyclonal Antibody
PA590081	CA150 Polyclonal Antibody
PA590082	GALNT6 Polyclonal Antibody
PA590084	PDXDC1 Polyclonal Antibody
PA590090	IBTK Polyclonal Antibody
PA590085	SWAP70 Polyclonal Antibody
PA590088	Kallikrein 5 Polyclonal Antibody
PA590091	ATRNL1 Polyclonal Antibody
PA590094	PABP Polyclonal Antibody
PA590142	ZSWIM1 Polyclonal Antibody
PA590095	LYPD3 Polyclonal Antibody
PA590096	DKK2 Polyclonal Antibody
PA590097	RSK4 Polyclonal Antibody
PA590100	LGALSL Polyclonal Antibody
PA590109	TLR8 Polyclonal Antibody
PA590119	BATF3 Polyclonal Antibody
PA590126	KLHL8 Polyclonal Antibody
PA590128	EFHD2 Polyclonal Antibody
PA590130	TCTN1 Polyclonal Antibody
PA590139	C1orf57 Polyclonal Antibody
PA590143	LDHAL6B Polyclonal Antibody
PA590144	ORMDL3 Polyclonal Antibody
PA590145	PRPK Polyclonal Antibody
PA590147	PKD1L2 Polyclonal Antibody
PA590148	MRPL54 Polyclonal Antibody
PA590153	DPPA2 Polyclonal Antibody
PA590159	B3GNT6 Polyclonal Antibody
PA590161	ADAM32 Polyclonal Antibody
PA590166	PFTK1 Polyclonal Antibody
PA590167	Hemoglobin gamma A Polyclonal Antibody
PA590177	SGMS1 Polyclonal Antibody
PA590172	GP6 Polyclonal Antibody
PA590173	LDHC Polyclonal Antibody
PA590178	Septin-8 Polyclonal Antibody
PA590188	FOXN3 Polyclonal Antibody
PA590191	LTB4R Polyclonal Antibody
PA590192	CTRB1 Polyclonal Antibody
PA590193	ZIPK Polyclonal Antibody
PA590195	EGR2 Polyclonal Antibody
PA590199	GCK Polyclonal Antibody
PA590204	NFIC Polyclonal Antibody
PA590209	PLD1 Polyclonal Antibody
PA590212	TAS2R38 Polyclonal Antibody
PA590215	RHD Polyclonal Antibody
PA590217	RPL28 Polyclonal Antibody
PA590218	RPS26 Polyclonal Antibody
PA590220	GAT-1 Polyclonal Antibody
PA590222	TBL1X Polyclonal Antibody
PA590227	XK Polyclonal Antibody
PA590229	ZNF9 Polyclonal Antibody
PA590232	Cdc23 Polyclonal Antibody
PA590235	RASAL2 Polyclonal Antibody
PA590237	eIF5B Polyclonal Antibody
PA590239	OXSRI Polyclonal Antibody
PA590241	SPON1 Polyclonal Antibody
PA590243	RNAse H2A Polyclonal Antibody

PA590244	IFITM2 Polyclonal Antibody
PA590246	MEKK2 Polyclonal Antibody
PA590248	SLC22A7 Polyclonal Antibody
PA590249	TRAM1 Polyclonal Antibody
PA590250	HEY2 Polyclonal Antibody
PA590253	Dynamin 3 Polyclonal Antibody
PA590255	FBXO25 Polyclonal Antibody
PA590256	Claudin 17 Polyclonal Antibody
PA590257	FOXB1 Polyclonal Antibody
PA590258	FLVCR1 Polyclonal Antibody
PA590259	GPR171 Polyclonal Antibody
PA590260	PYCR2 Polyclonal Antibody
PA590261	TAS2R10 Polyclonal Antibody
PA590263	HAO2 Polyclonal Antibody
PA590267	TRAC-1 Polyclonal Antibody
PA590268	TRAC-1 Polyclonal Antibody
PA590271	STRBP Polyclonal Antibody
PA590280	CRELD1 Polyclonal Antibody
PA590291	DGCR6L Polyclonal Antibody
PA590293	HSPA12B Polyclonal Antibody
PA590294	PLD4 Polyclonal Antibody
PA590296	UPF1 Polyclonal Antibody
PA590298	FGD5 Polyclonal Antibody
PA590299	GIMAP8 Polyclonal Antibody
PA590300	QSOX2 Polyclonal Antibody
PA590301	TCTE1 Polyclonal Antibody
PA590302	NAIF1 Polyclonal Antibody
PA590304	ANGPTL5 Polyclonal Antibody
PA590306	OTUD4 Polyclonal Antibody
PA590307	TEKT1 Polyclonal Antibody
PA590309	TRHR Polyclonal Antibody
PA590310	CMKLR1 Polyclonal Antibody
PA590311	RPAP3 Polyclonal Antibody
PA590314	RXRA Polyclonal Antibody
PA590315	ULK2 Polyclonal Antibody
PA590316	EPRS Polyclonal Antibody
PA590317	FARS2 Polyclonal Antibody
PA590318	SASH1 Polyclonal Antibody
PA590319	UNC93B1 Polyclonal Antibody
PA590321	Afamin Polyclonal Antibody
PA590323	RERE Polyclonal Antibody
PA590342	FOXA1 Polyclonal Antibody
PA590325	ATP Synthase gamma Polyclonal Antibody
PA590335	EAR2 Polyclonal Antibody
PA590338	GALK1 Polyclonal Antibody
PA590339	GTF2E1 Polyclonal Antibody
PA590343	HOXD3 Polyclonal Antibody
PA590345	IVD Polyclonal Antibody
PA590348	KID Polyclonal Antibody
PA590351	MANBA Polyclonal Antibody
PA590357	SERPINB2 Polyclonal Antibody
PA590359	PKNOX1 Polyclonal Antibody
PA590361	PLS1 Polyclonal Antibody
PA590366	PCYT2 Polyclonal Antibody
PA590367	PRKAR2A Polyclonal Antibody
PA590374	CCL16 Polyclonal Antibody
PA590384	ZNF19 Polyclonal Antibody
PA590386	ZNF38 Polyclonal Antibody
PA590388	TRIM26 Polyclonal Antibody

PA590394	PSMG1 Polyclonal Antibody
PA590398	CLIC3 Polyclonal Antibody
PA590402	Cytohesin 1 Polyclonal Antibody
PA590399	MTMR6 Polyclonal Antibody
PA590403	ASC1 Polyclonal Antibody
PA590404	SLIT2 Polyclonal Antibody
PA590405	GST Omega 1 Polyclonal Antibody
PA590406	ICB-1 Polyclonal Antibody
PA590407	EEF1E1 Polyclonal Antibody
PA590411	FEZ1 Polyclonal Antibody
PA590412	FEZ1 Polyclonal Antibody
PA590414	BMS1 Polyclonal Antibody
PA590415	GIT2 Polyclonal Antibody
PA590416	KIAA0319 Polyclonal Antibody
PA590419	GFPT2 Polyclonal Antibody
PA590421	CERT Polyclonal Antibody
PA590422	SUGP2 Polyclonal Antibody
PA590425	TCIRG1 Polyclonal Antibody
PA590428	VAV3 Polyclonal Antibody
PA590429	SEMA3C Polyclonal Antibody
PA590431	LEFTY Polyclonal Antibody
PA590432	NFATC1 Polyclonal Antibody
PA590434	GRAP Polyclonal Antibody
PA590436	AFG3L2 Polyclonal Antibody
PA590437	VIP36 Polyclonal Antibody
PA590442	MGAT4A Polyclonal Antibody
PA590443	CPEB3 Polyclonal Antibody
PA590448	FAM62A Polyclonal Antibody
PA590449	SLC44A1 Polyclonal Antibody
PA590450	ORC3 Polyclonal Antibody
PA590452	EDC4 Polyclonal Antibody
PA590455	TSPAN17 Polyclonal Antibody
PA590456	SLC17A5 Polyclonal Antibody
PA590458	MYOF Polyclonal Antibody
PA590468	UTP18 Polyclonal Antibody
PA590459	MRPL46 Polyclonal Antibody
PA590463	SRPX2 Polyclonal Antibody
PA590467	FHOD1 Polyclonal Antibody
PA590470	SDF4 Polyclonal Antibody
PA590477	SLC25A39 Polyclonal Antibody
PA590483	TRSPAP1 Polyclonal Antibody
PA590486	CDCA8 Polyclonal Antibody
PA590491	ZNF331 Polyclonal Antibody
PA590495	WDR12 Polyclonal Antibody
PA590499	PCNP Polyclonal Antibody
PA590500	LYRM1 Polyclonal Antibody
PA590504	KCC2 Polyclonal Antibody
PA590506	CD4 Polyclonal Antibody
PA590507	LRRG4C Polyclonal Antibody
PA590509	TSCOT Polyclonal Antibody
PA590511	RBM25 Polyclonal Antibody
PA590515	NECAB3 Polyclonal Antibody
PA590518	MICAL1 Polyclonal Antibody
PA590519	WNK3 Polyclonal Antibody
PA590523	UBA5 Polyclonal Antibody
PA590526	WDR23 Polyclonal Antibody
PA590529	SLIRP Polyclonal Antibody
PA590530	RTBDN Polyclonal Antibody
PA590533	POMK Polyclonal Antibody

PA590548	DPP9 Polyclonal Antibody
PA590554	SDSL Polyclonal Antibody
PA590558	NTAN1 Polyclonal Antibody
PA590559	ACSM2A Polyclonal Antibody
PA590561	HEXIM2 Polyclonal Antibody
PA590565	OTUD7A Polyclonal Antibody
PA590566	LRWD1 Polyclonal Antibody
PA590570	FUT1 Polyclonal Antibody
PA590567	CDY1B Polyclonal Antibody
PA590569	DDX51 Polyclonal Antibody
PA590572	USP37 Polyclonal Antibody
PA590578	Acetyl-CoA Carboxylase Polyclonal Antibody
PA590580	NDUFB11 Polyclonal Antibody
PA590585	DDX4 Polyclonal Antibody
PA590592	SLC25A5 Polyclonal Antibody
PA590593	Arylsulfatase E Polyclonal Antibody
PA590595	TSPO Polyclonal Antibody
PA590606	ETV6 Polyclonal Antibody
PA590597	VPS51 Polyclonal Antibody
PA590599	ENTPD6 Polyclonal Antibody
PA590600	Cdc20 Polyclonal Antibody
PA590605	ETS1 Polyclonal Antibody
PA590609	FARSA Polyclonal Antibody
PA590610	FNTB Polyclonal Antibody
PA590611	GLRA2 Polyclonal Antibody
PA590612	GNGT1 Polyclonal Antibody
PA590615	hnRNP D Polyclonal Antibody
PA590616	IL10RB Polyclonal Antibody
PA590618	KV3.4 (KCNC4) Polyclonal Antibody
PA590624	NAB1 Polyclonal Antibody
PA590631	Neuroserpin Polyclonal Antibody
PA590634	PPP2R5A Polyclonal Antibody
PA590636	PKA beta Polyclonal Antibody
PA590638	PRPSAP1 Polyclonal Antibody
PA590640	PSMC1 Polyclonal Antibody
PA590651	SUMO-3 Polyclonal Antibody
PA590646	RSK1 Polyclonal Antibody
PA590649	GLAST Polyclonal Antibody
PA590658	TPI1 Polyclonal Antibody
PA590661	ZNF133 Polyclonal Antibody
PA590662	ZNF134 Polyclonal Antibody
PA590665	PPAP2B Polyclonal Antibody
PA590666	RIOK3 Polyclonal Antibody
PA590668	ZPR1 Polyclonal Antibody
PA590677	MED20 Polyclonal Antibody
PA590678	VPS4B Polyclonal Antibody
PA590681	HCC1 Polyclonal Antibody
PA590682	NFE2L3 Polyclonal Antibody
PA590683	LPIN2 Polyclonal Antibody
PA590686	CLSTN3 Polyclonal Antibody
PA590687	RASSF2 Polyclonal Antibody
PA590689	GOLGA5 Polyclonal Antibody
PA590692	FLASH Polyclonal Antibody
PA590693	Bim Polyclonal Antibody
PA590694	TENM1 Polyclonal Antibody
PA590695	TIP47 Polyclonal Antibody
PA590696	CD2BP2 Polyclonal Antibody
PA590697	ACAA2 Polyclonal Antibody
PA590699	UNC13B Polyclonal Antibody

PA590700	PAT1 Polyclonal Antibody
PA590703	AVIL Polyclonal Antibody
PA590702	POLR3C Polyclonal Antibody
PA590704	KIF1C Polyclonal Antibody
PA590705	FSTL1 Polyclonal Antibody
PA590706	HLA-DRA Polyclonal Antibody
PA590707	GPR45 Polyclonal Antibody
PA590709	LYPLA2 Polyclonal Antibody
PA590711	LSD1 Polyclonal Antibody
PA590713	NBEAL2 Polyclonal Antibody
PA590714	SPECC1L Polyclonal Antibody
PA590716	SF3B1 Polyclonal Antibody
PA590717	GCAT Polyclonal Antibody
PA590719	CCRN4L Polyclonal Antibody
PA590721	FAM98A Polyclonal Antibody
PA590722	REXO2 Polyclonal Antibody
PA590723	ZZZ3 Polyclonal Antibody
PA590725	Erc2 Polyclonal Antibody
PA590726	RAI14 Polyclonal Antibody
PA590731	TRUB2 Polyclonal Antibody
PA590732	PPA2 Polyclonal Antibody
PA590734	ZNF638 Polyclonal Antibody
PA590735	MCAT Polyclonal Antibody
PA590737	PACSIN1 Polyclonal Antibody
PA590739	MYEF2 Polyclonal Antibody
PA590742	ASCC1 Polyclonal Antibody
PA590743	RPS27L Polyclonal Antibody
PA590745	NDUFAF1 Polyclonal Antibody
PA590748	MRPL51 Polyclonal Antibody
PA590749	C3orf1 Polyclonal Antibody
PA590755	CDK5RAP1 Polyclonal Antibody
PA590758	TPCN1 Polyclonal Antibody
PA590759	BCMO1 Polyclonal Antibody
PA590760	UGT1A4 Polyclonal Antibody
PA590762	UCKL1 Polyclonal Antibody
PA590763	OXSM Polyclonal Antibody
PA590765	ELP2 Polyclonal Antibody
PA590766	DDX43 Polyclonal Antibody
PA590770	MAP7D1 Polyclonal Antibody
PA590771	ZFP64 Polyclonal Antibody
PA590773	LMBRD1 Polyclonal Antibody
PA590774	ACTR10 Polyclonal Antibody
PA590783	MRPS35 Polyclonal Antibody
PA590788	ATL2 Polyclonal Antibody
PA590790	MMP25 Polyclonal Antibody
PA590792	VPS16 Polyclonal Antibody
PA590797	DDX31 Polyclonal Antibody
PA590801	B3GNT4 Polyclonal Antibody
PA590804	CYB5B Polyclonal Antibody
PA590805	ITFG1 Polyclonal Antibody
PA590808	CD99L2 Polyclonal Antibody
PA590811	SLC4A11 Polyclonal Antibody
PA590821	PPIL4 Polyclonal Antibody
PA590824	TP53I13 Polyclonal Antibody
PA590825	ZNF766 Polyclonal Antibody
PA590829	TMEM44 Polyclonal Antibody
PA590832	ARMC6 Polyclonal Antibody
PA590833	MVB12A Polyclonal Antibody
PA590836	ADC Polyclonal Antibody

PA590839	RBP7 Polyclonal Antibody
PA590841	CHCHD1 Polyclonal Antibody
PA590842	ZFYVE27 Polyclonal Antibody
PA590843	ZFYVE27 Polyclonal Antibody
PA590849	PROTEIN S Polyclonal Antibody
PA590853	ACMSD Polyclonal Antibody
PA590855	STARD4 Polyclonal Antibody
PA590859	SIX5 Polyclonal Antibody
PA590861	PUS10 Polyclonal Antibody
PA590865	NSMCE1 Polyclonal Antibody
PA590867	MUC20 Polyclonal Antibody
PA590869	Opsin 5 Polyclonal Antibody
PA590876	CHRM1 Polyclonal Antibody
PA590881	GPX2 Polyclonal Antibody
PA590883	ZIP12 Polyclonal Antibody
PA590886	GABRB1 Polyclonal Antibody
PA590889	Laminin gamma-1 Polyclonal Antibody
PA590893	ATP1A3 Polyclonal Antibody
PA590894	CACNB2 Polyclonal Antibody
PA590896	CD80 Polyclonal Antibody
PA590899	DCT Polyclonal Antibody
PA590900	ETFB Polyclonal Antibody
PA590902	Fibulin 1 Polyclonal Antibody
PA590907	MAFG Polyclonal Antibody
PA590903	CD42b Polyclonal Antibody
PA590904	HNF4G Polyclonal Antibody
PA590908	MEKK3 Polyclonal Antibody
PA590911	TRIM27 Polyclonal Antibody
PA590912	RFX1 Polyclonal Antibody
PA590918	CDS2 Polyclonal Antibody
PA590922	SLC28A2 Polyclonal Antibody
PA590925	SEC24D Polyclonal Antibody
PA590923	CRSP9 Polyclonal Antibody
PA590924	GPRASP1 Polyclonal Antibody
PA590927	TSFM Polyclonal Antibody
PA590933	LILRB3 Polyclonal Antibody
PA590934	Neuroigin 1 Polyclonal Antibody
PA590935	Epsin 2 Polyclonal Antibody
PA590936	Cdc2L6 Polyclonal Antibody
PA591843	ALX4 Polyclonal Antibody
PA591865	Dynamin 1 Polyclonal Antibody
PA591872	PTGER2 Polyclonal Antibody
PA591874	EpoR Polyclonal Antibody
PA591876	GABRR1 Polyclonal Antibody
PA591889	KISS1 Polyclonal Antibody
PA591892	RAP Polyclonal Antibody
PA591896	CSF1R Polyclonal Antibody
PA591905	14-3-3 gamma Polyclonal Antibody
PA591916	PIBF1 Polyclonal Antibody
PA591907	Nephrin Polyclonal Antibody
PA591915	PDX1 Polyclonal Antibody
PA591925	TPM2 Polyclonal Antibody
PA591929	MYSM1 Polyclonal Antibody
PA591933	Neurogenin 1 Polyclonal Antibody
PA591942	SFTPA1 Polyclonal Antibody
PA591944	H3R17me1 Polyclonal Antibody
PA591971	CtBP1 Polyclonal Antibody
PA592014	CREBZF Polyclonal Antibody
PA592022	FAM117B Polyclonal Antibody

PA592030	DOCK2 Polyclonal Antibody
PA592037	SRSF2 Polyclonal Antibody
PA592039	ELOVL4 Polyclonal Antibody
PA592044	TMEM16B Polyclonal Antibody
PA592040	TGF beta-2 Polyclonal Antibody
PA592051	RhoGAP Polyclonal Antibody
PA592064	CLCNKA Polyclonal Antibody
PA592102	KRT81 Polyclonal Antibody
PA592094	GSTM3 Polyclonal Antibody
PA592110	NDUFB2 Polyclonal Antibody
PA592113	Nucleobindin 1 Polyclonal Antibody
PA592115	PFDN5 Polyclonal Antibody
PA592117	PITX1 Polyclonal Antibody
PA592124	PTGER3 Polyclonal Antibody
PA592164	UBL4A Polyclonal Antibody
PA592168	IFITM1 Polyclonal Antibody
PA592169	CGGBP1 Polyclonal Antibody
PA592174	NAE1 Polyclonal Antibody
PA592175	TIMELESS Polyclonal Antibody
PA592187	CHST10 Polyclonal Antibody
PA592194	NUP93 Polyclonal Antibody
PA592228	SLC38A3 Polyclonal Antibody
PA592252	CES3 Polyclonal Antibody
PA592258	EML2 Polyclonal Antibody
PA592261	MOB4 Polyclonal Antibody
PA592273	CRCP Polyclonal Antibody
PA592277	T-bet Polyclonal Antibody
PA592293	RhoF Polyclonal Antibody
PA592299	RNMTL1 Polyclonal Antibody
PA592304	TSR1 Polyclonal Antibody
PA592326	GGCT Polyclonal Antibody
PA592354	ARL13B Polyclonal Antibody
PA592356	TSPAN33 Polyclonal Antibody
PA592372	H3K14me3 Polyclonal Antibody
PA592379	ABAT Polyclonal Antibody
PA592383	PEO1 Polyclonal Antibody
PA592389	GLRX Polyclonal Antibody
PA592396	p23 Polyclonal Antibody
PA592401	SLC4A4 Polyclonal Antibody
PA592403	Apolipoprotein M Polyclonal Antibody
PA592438	BAF53A Polyclonal Antibody
PA592441	CARTPT Polyclonal Antibody
PA592450	FBP1 Polyclonal Antibody
PA592894	Transketolase Polyclonal Antibody
PA592904	PRDX5 Polyclonal Antibody
PA592962	Tnfaip6 Polyclonal Antibody
PA592982	DHRS2 Polyclonal Antibody
PA592994	NRBF2 Polyclonal Antibody
PA592995	DUSP13 Polyclonal Antibody
PA593110	IL9R Polyclonal Antibody
PA593138	NMNAT1 Polyclonal Antibody
PA593221	USP15 Polyclonal Antibody
PA593245	RECQ4 Polyclonal Antibody
PA593293	KLRB1 Polyclonal Antibody
PA593306	PMS2 Polyclonal Antibody
PA593318	RAD51B Polyclonal Antibody
PA593329	SOD3 Polyclonal Antibody
PA595747	EMX1 Polyclonal Antibody
PA597217	OAS3 Polyclonal Antibody

PA597235	HIC1 Polyclonal Antibody
PA597390	Phospho-Epha3 (Tyr779) Polyclonal Antibody
PA597411	Phospho-SLP76 (Tyr145) Polyclonal Antibody
PA596588	NOX5 Polyclonal Antibody
PA595694	CAB39 Polyclonal Antibody
PA596300	NANS Polyclonal Antibody
PA595810	SOGA1 Polyclonal Antibody
PA596309	NLRX1 Polyclonal Antibody
PA596295	MOXD1 Polyclonal Antibody
PA597400	Phospho-AS160 (Thr642) Polyclonal Antibody
PA595946	GluR1 Polyclonal Antibody
PA595776	GLIPR1 Polyclonal Antibody
PA597396	Phospho-Filamin A (Ser2152) Polyclonal Antibody
PA597406	Phospho-Vimentin (Ser56) Polyclonal Antibody
PA596711	ELMOD3 Polyclonal Antibody
PA596008	ADORA2B Polyclonal Antibody
PA596512	GM-CSF Polyclonal Antibody
PA596336	PLA2 Polyclonal Antibody
PA596112	CD86 Polyclonal Antibody
PA597211	NFIX Polyclonal Antibody
PA595734	PAG1 Polyclonal Antibody
PA596247	THAP2 Polyclonal Antibody
PA596778	SLC4A2 Polyclonal Antibody
PA597408	Phospho-LAT (Tyr171) Polyclonal Antibody
PA596378	SMAD4 Polyclonal Antibody
PA595702	TMEM132D Polyclonal Antibody
PA596176	betacellulin Polyclonal Antibody
PA595757	PSMD1 Polyclonal Antibody
PA597195	C17orf25 Polyclonal Antibody
PA595807	ZGPAT Polyclonal Antibody
PA597126	MIB1 Polyclonal Antibody
PA597397	Phospho-SHP-1 (Tyr564) Polyclonal Antibody
PA595754	NEFM Polyclonal Antibody
PA597198	FJX1 Polyclonal Antibody
PA596210	ACE Polyclonal Antibody
PA597242	HLA-DQA2 Polyclonal Antibody
PA597176	SIRP alpha Polyclonal Antibody
PA596535	TADA3L Polyclonal Antibody
PA595732	M-CSF Polyclonal Antibody
PA595790	MST4 Polyclonal Antibody
PA596838	ACVR2B Polyclonal Antibody
PA596501	TIAL1 Polyclonal Antibody
PA597221	ERC1 Polyclonal Antibody
PA595970	TRIM5 alpha/beta/gamma/delta/epsilon/iota Polyclonal Antibody
PA596385	IQGAP1 Polyclonal Antibody
PA596304	CDC37L1 Polyclonal Antibody
PA595931	HGS Polyclonal Antibody
PA596232	Biglycan Polyclonal Antibody
PA595706	SLC34A3 Polyclonal Antibody
PA596188	EB1 Polyclonal Antibody
PA597392	Phospho-JAK2 (Tyr1008) Polyclonal Antibody
PA595980	PDHA1 Polyclonal Antibody
PA597125	KIAA1199 Polyclonal Antibody
PA596206	MATN4 Polyclonal Antibody
PA596077	Caspase 3 Polyclonal Antibody
PA597363	Phospho-NFkB p65 (Ser311) Polyclonal Antibody
PA596866	SSX4 Polyclonal Antibody
PA596438	NEIL1 Polyclonal Antibody
PA596436	TIP30 Polyclonal Antibody

PA595729	CYP26A1 Polyclonal Antibody
PA595786	ITSN2 Polyclonal Antibody
PA596161	Granzyme B Polyclonal Antibody
PA595720	D2HGDH Polyclonal Antibody
PA597229	USP39 Polyclonal Antibody
PA595801	ELOVL5 Polyclonal Antibody
PA597405	Phospho-MSK1 (Ser376) Polyclonal Antibody
PA595781	PARS2 Polyclonal Antibody
PA597226	Spastin Polyclonal Antibody
PA597412	Phospho-EphA2 (Tyr594) Polyclonal Antibody
PA596246	SYNPO2L Polyclonal Antibody
PA597038	MRP2 Polyclonal Antibody
PA595697	BPIL1 Polyclonal Antibody
PA596413	IFNGR1 Polyclonal Antibody
PA595784	GMPPB Polyclonal Antibody
PA595883	P2Y12 Polyclonal Antibody
PA597398	Phospho-PEA15 (Ser104) Polyclonal Antibody
PA596377	PPAR delta Polyclonal Antibody
PA597219	PCCA Polyclonal Antibody
PA596317	PEPD Polyclonal Antibody
PA597227	ZRANB3 Polyclonal Antibody
PA596273	Alpha Sarcoglycan Polyclonal Antibody
PA595708	MEI1 Polyclonal Antibody
PA596751	Ceruloplasmin Polyclonal Antibody
PA597190	MYL12A Polyclonal Antibody
PA596605	IDH2 Polyclonal Antibody
PA597057	NF-H Polyclonal Antibody
PA595854	UGT2B7 Polyclonal Antibody
PA597391	Phospho-MYPT1 (Ser507) Polyclonal Antibody
PA595762	SCD Polyclonal Antibody
PA596491	HDLBP Polyclonal Antibody
PA595804	HKDC1 Polyclonal Antibody
PA595895	IGFBP5 Polyclonal Antibody
PA597399	Phospho-AS160 (Ser588) Polyclonal Antibody
PA595889	CYP27B1 Polyclonal Antibody
PA595792	CHDH Polyclonal Antibody
PA595699	HPS4 Polyclonal Antibody
PA595813	PPM1K Polyclonal Antibody
PA597402	Phospho-PRAS40 (Thr246) Polyclonal Antibody
PA595991	CD204 Polyclonal Antibody
PA595722	IFN beta Polyclonal Antibody
PA597171	CHST15 Polyclonal Antibody
PA595803	FAM65A Polyclonal Antibody
PA597169	SPAG8 Polyclonal Antibody
PA597208	FBXW5 Polyclonal Antibody
PA597138	eIF3j Polyclonal Antibody
PA595812	SPC24 Polyclonal Antibody
PA595971	ACADM Polyclonal Antibody
PA596441	ING4 Polyclonal Antibody
PA596163	DIABLO Polyclonal Antibody
PA597145	PFDN1 Polyclonal Antibody
PA596098	MYST1 Polyclonal Antibody
PA596087	MEK7 Polyclonal Antibody
PA597173	VBP1 Polyclonal Antibody
PA596284	ABCG1 Polyclonal Antibody
PA595780	HEBP2 Polyclonal Antibody
PA595805	KLHL15 Polyclonal Antibody
PA596224	CKMT1B Polyclonal Antibody
PA596296	MYCBP Polyclonal Antibody

PA596482	eIF3d Polyclonal Antibody
PA595748	IL1RA Polyclonal Antibody
PA596251	SLC16A2 Polyclonal Antibody
PA597192	RPL22 Polyclonal Antibody
PA595711	WBSCR27 Polyclonal Antibody
PA596360	L-Plastin Polyclonal Antibody
PA595869	Cofilin Polyclonal Antibody
PA596220	FAF1 Polyclonal Antibody
PA597237	EXOSC7 Polyclonal Antibody
PA595822	PCM1 Polyclonal Antibody
PA597197	ROGDI Polyclonal Antibody
PA597230	MAPRE3 Polyclonal Antibody
PA595808	Lsm10 Polyclonal Antibody
PA596926	Nkx3.2 Polyclonal Antibody
PA595718	C11orf79 Polyclonal Antibody
PA596156	TRAF5 Polyclonal Antibody
PA595917	CCL20 Polyclonal Antibody
PA595929	RGS4 Polyclonal Antibody
PA595758	RGS4 Polyclonal Antibody
PA596499	SRRM1 Polyclonal Antibody
PA596230	PP2A beta Polyclonal Antibody
PA595939	CYSLTR1 Polyclonal Antibody
PA596498	PTBP2 Polyclonal Antibody
PA596737	Bax Polyclonal Antibody
PA596742	14-3-3 zeta Polyclonal Antibody
PA596878	SOCS5 Polyclonal Antibody
PA595713	YIPF7 Polyclonal Antibody
PA595787	CKLF Polyclonal Antibody
PA595926	AKR1C3 Polyclonal Antibody
PA595802	NR5A1 Polyclonal Antibody
PA597236	TXNDC9 Polyclonal Antibody
PA596662	ATG10 Polyclonal Antibody
PA595821	CPSF4L Polyclonal Antibody
PA597194	TRMT6 Polyclonal Antibody
PA596162	14-3-3 theta Polyclonal Antibody
PA596915	PHF19 Polyclonal Antibody
PA595698	CK068 Polyclonal Antibody
PA596187	MSI1 Polyclonal Antibody
PA595710	SAMD3 Polyclonal Antibody
PA596354	GNAS Polyclonal Antibody
PA595721	CD62L Polyclonal Antibody
PA597183	14-3-3 eta Polyclonal Antibody
PA596733	ZNF785 Polyclonal Antibody
PA595760	RPS17 Polyclonal Antibody
PA595759	RPS17 Polyclonal Antibody
PA595774	NPRL2 Polyclonal Antibody
PA595818	LNX2 Polyclonal Antibody
PA595745	CTR2 Polyclonal Antibody
PA596500	SRSF3 Polyclonal Antibody
PA596034	DST Polyclonal Antibody
PA596412	GRIP1 Polyclonal Antibody
PA595723	ID4 Polyclonal Antibody
PA596076	ATF Polyclonal Antibody
PA597191	RPL26L1 Polyclonal Antibody
PA596808	GABARAPL1 Polyclonal Antibody
PA596173	DARPP-32 Polyclonal Antibody
PA595756	NIT1 Polyclonal Antibody
PA597212	IST1 Polyclonal Antibody
PA595746	EGR4 Polyclonal Antibody

PA596327	PSMB10 Polyclonal Antibody
PA595691	UNC50 Polyclonal Antibody
PA597134	SERTAD1 Polyclonal Antibody
PA597317	SAV1 Polyclonal Antibody
PA595847	S100A8 Polyclonal Antibody
PA597196	FAM65B Polyclonal Antibody
PA597193	PSPC1 Polyclonal Antibody
PA595968	PITX2 Polyclonal Antibody
PA595814	RBM46 Polyclonal Antibody
PA597206	KCNN1 Polyclonal Antibody
PA596011	HIST2H2BE Polyclonal Antibody
PA597160	DDT Polyclonal Antibody
PA596040	PRDX6 Polyclonal Antibody
PA597136	DDX50 Polyclonal Antibody
PA597155	BRF2 Polyclonal Antibody
PA596551	eIF2b beta Polyclonal Antibody
PA596128	BRCA2 Polyclonal Antibody
PA596617	RAP80 Polyclonal Antibody
PA596184	NRG4 Polyclonal Antibody
PA596409	CDK5 Polyclonal Antibody
PA595766	OR1A1 Polyclonal Antibody
PA597187	RFK Polyclonal Antibody
PA596090	p70 S6 Kinase Polyclonal Antibody
PA596710	ZFYVE1 Polyclonal Antibody
PA595782	COMMD5 Polyclonal Antibody
PA595705	S100A16 Polyclonal Antibody
PA596641	PAK2 Polyclonal Antibody
PA597177	RPSA Polyclonal Antibody
PA597239	USP19 Polyclonal Antibody
PA596133	HSC70 Polyclonal Antibody
PA596560	HAND2 Polyclonal Antibody
PA596517	SETDB1 Polyclonal Antibody
PA596278	SHOC2 Polyclonal Antibody
PA596170	N-WASP Polyclonal Antibody
PA596088	Presenilin 1 Polyclonal Antibody
PA596137	SIAH1 Polyclonal Antibody
PA595815	CCER1 Polyclonal Antibody
PA597232	RRP7A Polyclonal Antibody
PA595753	MGST1 Polyclonal Antibody
PA596455	SND1 Polyclonal Antibody
PA595806	YPEL3 Polyclonal Antibody
PA596197	MKP2 Polyclonal Antibody
PA596504	HuC Polyclonal Antibody
PA595778	FAM107A Polyclonal Antibody
PA595777	FAM107A Polyclonal Antibody
PA595816	CCDC11 Polyclonal Antibody
PA597410	Phospho-PDGFRB (Tyr1021) Polyclonal Antibody
PA597386	Phospho-Glucocorticoid Receptor (Ser211) Polyclonal Antibody
PA596021	CYP27A1 Polyclonal Antibody
PA595744	CMKLR1 Polyclonal Antibody
PA596293	NNT Polyclonal Antibody
PA597215	NUP50 Polyclonal Antibody
PA595800	ARHGAP22 Polyclonal Antibody
PA595791	CRTAC1 Polyclonal Antibody
PA597172	HOXC4 Polyclonal Antibody
PA595709	GPR155 Polyclonal Antibody
PA597213	HECTD1 Polyclonal Antibody
PA596579	NAA20 Polyclonal Antibody
PA596071	SPINLW1 Polyclonal Antibody

PA595834	Cdc25C Polyclonal Antibody
PA596281	ARHGEF1 Polyclonal Antibody
PA597149	SIX2 Polyclonal Antibody
PA597394	Phospho-Lamin A/C (Ser22) Polyclonal Antibody
PA597355	Phospho-STAT1 (Ser727) Polyclonal Antibody
PA595779	HRH3 Polyclonal Antibody
PA597205	SEC31A Polyclonal Antibody
PA596596	REST Polyclonal Antibody
PA596241	ROR1 Polyclonal Antibody
PA595772	ZNF263 Polyclonal Antibody
PA597357	Phospho-Stathmin 1 (Ser16) Polyclonal Antibody
PA596205	Kallikrein 15 Polyclonal Antibody
PA596111	Histone H3.1t Polyclonal Antibody
PA597220	PLA2G6 Polyclonal Antibody
PA596416	PCGF6 Polyclonal Antibody
PA597225	SARDH Polyclonal Antibody
PA595796	GPR108 Polyclonal Antibody
PA597361	Phospho-STAT2 (Tyr690) Polyclonal Antibody
PA597403	Phospho-beta Catenin (Ser675) Polyclonal Antibody
PA595763	SLIT1 Polyclonal Antibody
PA595795	FSTL5 Polyclonal Antibody
PA597418	Phospho-MYPT1 (Ser668) Polyclonal Antibody
PA596081	IL-1 alpha Polyclonal Antibody
PA597388	Phospho-VASP (Ser157) Polyclonal Antibody
PA595809	DHX57 Polyclonal Antibody
PA596919	TMEM166 Polyclonal Antibody
PA597170	SUCLG2 Polyclonal Antibody
PA596181	CCL28 Polyclonal Antibody
PA597178	MYO18A Polyclonal Antibody
PA597199	TOM1 Polyclonal Antibody
PA595703	GPR151 Polyclonal Antibody
PA595811	NOTUM Polyclonal Antibody
PA596103	CtBP2 Polyclonal Antibody
PA596534	Flotillin 1 Polyclonal Antibody
PA596505	FMRP Polyclonal Antibody
PA596398	NCOA4 Polyclonal Antibody
PA595799	CD319 (CRACC) Polyclonal Antibody
PA596425	CD319 (CRACC) Polyclonal Antibody
PA597401	Phospho-NDRG1 (Thr346) Polyclonal Antibody
PA595789	POP5 Polyclonal Antibody
PA596244	IL1RAPL1 Polyclonal Antibody
PA596722	CCL3 Polyclonal Antibody
PA596291	SF3A3 Polyclonal Antibody
PA595700	MYOCD Polyclonal Antibody
PA597416	Phospho-SRF (Ser103) Polyclonal Antibody
PA595794	SEPTIN 3 Polyclonal Antibody
PA596127	SQLE Polyclonal Antibody
PA597273	Tab2 Polyclonal Antibody
PA595693	RHCG Polyclonal Antibody
PA597407	Phospho-NDRG1 (Ser330) Polyclonal Antibody
PA596414	CD155 Polyclonal Antibody
PA597290	Aquaporin 5 Polyclonal Antibody
PA596473	MAP4 Polyclonal Antibody
PA595788	CD46 Polyclonal Antibody
PA597414	Phospho-AP2M1 (Thr156) Polyclonal Antibody
PA595764	SP4 Polyclonal Antibody
PA597393	Phospho-eIF4B (Ser422) Polyclonal Antibody
PA595761	RTN1 Polyclonal Antibody
PA597164	PTPRS Polyclonal Antibody

PA597404	Phospho-PKC theta (Thr538) Polyclonal Antibody
PA595797	ZNF471 Polyclonal Antibody
PA595798	ZNF471 Polyclonal Antibody
PA596155	Tenascin X Polyclonal Antibody
PA595751	KIF2A Polyclonal Antibody
PA595704	TAAR1 Polyclonal Antibody
PA595752	MEF2D Polyclonal Antibody
PA595923	IGF2BP2 Polyclonal Antibody
PA596025	Thromboxane synthase Polyclonal Antibody
PA597181	DDOST Polyclonal Antibody
PA595770	DCLK1 Polyclonal Antibody
PA596871	WNT9A Polyclonal Antibody
PA597165	ARL8B Polyclonal Antibody
PA596138	ChAT Polyclonal Antibody
PA596126	AIP2 Polyclonal Antibody
PA597072	PATJ Polyclonal Antibody
PA596298	IFT81 Polyclonal Antibody
PA595714	GPR141 Polyclonal Antibody
PA596756	GAA Polyclonal Antibody
PA596957	USF2 Polyclonal Antibody
PA596007	SH3KBP1 Polyclonal Antibody
PA597148	PUS1 Polyclonal Antibody
PA596085	PDGFRB Polyclonal Antibody
PA595819	NOMO2 Polyclonal Antibody
PA596916	ANAPC2 Polyclonal Antibody
PA595897	CHI3L1 Polyclonal Antibody
PA595735	FADS1 Polyclonal Antibody
PA595793	PI4K2A Polyclonal Antibody
PA597234	ST8SIA1 Polyclonal Antibody
PA597146	TESC Polyclonal Antibody
PA595696	TMC5 Polyclonal Antibody
PA595846	HMGCR Polyclonal Antibody
PA595838	ARL4D Polyclonal Antibody
PA595836	TFAP4 Polyclonal Antibody
PA596876	Artemin Polyclonal Antibody
PA596803	CKAP4 Polyclonal Antibody
PA596200	HSF2 Polyclonal Antibody
PA597385	Phospho-RNA pol II CTD (Ser2) Polyclonal Antibody
PA596807	CD2AP Polyclonal Antibody
PA596351	AKT1 Polyclonal Antibody
PA597413	Phospho-EGFR (Tyr1045) Polyclonal Antibody
PA596219	EDA Polyclonal Antibody
PA596774	PSEN2 Polyclonal Antibody
PA597409	Phospho-PTPRA (Tyr789) Polyclonal Antibody
PA597415	Phospho-NMDAR1 (Ser890) Polyclonal Antibody
PA596489	ALY Polyclonal Antibody
PA597389	Phospho-RSK1 (Ser363) Polyclonal Antibody
PA595843	CtBP2 Polyclonal Antibody
PA596120	H3K79me2 Polyclonal Antibody
PA596113	H3K4me2 Polyclonal Antibody
PA596116	H3K36me1 Polyclonal Antibody
PA596119	H3K79me1 Polyclonal Antibody
PA596121	H3K79me3 Polyclonal Antibody
PA596233	H3R2me2a Polyclonal Antibody
PA595862	ZEB1 Polyclonal Antibody
PA595871	MAML1 Polyclonal Antibody
PA597352	Phospho-Histone H3 (Thr11) Polyclonal Antibody
PA596229	PACAP Receptor Polyclonal Antibody
PA595767	Cullin 3 Polyclonal Antibody

PA595727	GM130 Polyclonal Antibody
PA595849	LAMP1 Polyclonal Antibody
PA595775	HPa1 Polyclonal Antibody
PA595826	PYCARD Polyclonal Antibody
PA595743	RPL37 Polyclonal Antibody
PA595850	Lamin B1 Polyclonal Antibody
PA595739	FUT4 Polyclonal Antibody
PA595741	UNC84B Polyclonal Antibody
PA597359	Phospho-ERK1 (Tyr204) Polyclonal Antibody
PA595828	RPL26 Polyclonal Antibody
PA595740	SON Polyclonal Antibody
PA596199	Heme oxygenase 2 Polyclonal Antibody
PA596225	SUMO-3 Polyclonal Antibody
PA595898	ATOH8 Polyclonal Antibody
PA595835	IQSEC1 Polyclonal Antibody
PA595716	HCN1 Polyclonal Antibody
PA596209	Parvalbumin Polyclonal Antibody
PA596748	CEL Polyclonal Antibody
PA595768	TWEAK Polyclonal Antibody
PA595737	FUNDC1 Polyclonal Antibody
PA597371	Phospho-eNOS (Ser1177) Polyclonal Antibody
PA595906	PIM3 Polyclonal Antibody
PA596201	TRAP1 Polyclonal Antibody
PA595881	RPP40 Polyclonal Antibody
PA595730	ASM Polyclonal Antibody
PA596862	PKC gamma Polyclonal Antibody
PA596055	CD31 Polyclonal Antibody
PA595827	NOTCH1 Polyclonal Antibody
PA595842	CENPB Polyclonal Antibody
PA595845	GSK3B Polyclonal Antibody
PA595848	Interferon alpha-21 Polyclonal Antibody
PA595852	MUC4 Polyclonal Antibody
PA595853	NAIP Polyclonal Antibody
PA595855	NPAS2 Polyclonal Antibody
PA595856	NPY2R Polyclonal Antibody
PA595858	RPL18A Polyclonal Antibody
PA595860	TBX6 Polyclonal Antibody
PA595865	ZNF711 Polyclonal Antibody
PA595870	ATP5J2 Polyclonal Antibody
PA595873	DMTF1 Polyclonal Antibody
PA595876	TMEM147 Polyclonal Antibody
PA595877	VAT1 Polyclonal Antibody
PA595878	RPP38 Polyclonal Antibody
PA595879	YKT6 Polyclonal Antibody
PA595880	YKT6 Polyclonal Antibody
PA595884	ARMER Polyclonal Antibody
PA595885	CNNM4 Polyclonal Antibody
PA595886	CIDEB Polyclonal Antibody
PA595887	FCF1 Polyclonal Antibody
PA595888	HSD17B12 Polyclonal Antibody
PA595891	PCDHA10 Polyclonal Antibody
PA595892	SLC2A4RG Polyclonal Antibody
PA595893	WRCH1 Polyclonal Antibody
PA595894	OSTC Polyclonal Antibody
PA595899	Lass5 Polyclonal Antibody
PA595900	Lass5 Polyclonal Antibody
PA595902	ZNF563 Polyclonal Antibody
PA595905	RKHD1 Polyclonal Antibody
PA5109954	MEIS1 Polyclonal Antibody

PA5109573	COMTD1 Polyclonal Antibody
PA5110599	ZNF277 Polyclonal Antibody
PA5109972	MPHOSPH6 Polyclonal Antibody
PA5110098	PER2 Polyclonal Antibody
PA5109354	ABLIM2 Polyclonal Antibody
PA5109582	CPLX4 Polyclonal Antibody
PA5109445	BAZ2A Polyclonal Antibody
PA5110088	PDE6G Polyclonal Antibody
PA5109427	ATP13A2 Polyclonal Antibody
PA5109941	MARVELD1 Polyclonal Antibody
PA5109460	BVES Polyclonal Antibody
PA5109951	MDK Polyclonal Antibody
PA5110499	TMEM55A Polyclonal Antibody
PA5109630	DGKB Polyclonal Antibody
PA5109899	LAMB2 Polyclonal Antibody
PA5109937	MAPKAPK2 Polyclonal Antibody
PA5110348	SEC24C Polyclonal Antibody
PA5109990	MTMR3 Polyclonal Antibody
PA5109759	GCNT1 Polyclonal Antibody
PA5109408	ARHGEF10L Polyclonal Antibody
PA5110392	Dopamine Transporter Polyclonal Antibody
PA5110057	OCIAD1 Polyclonal Antibody
PA5110205	PPP1CB Polyclonal Antibody
PA5110424	SPRTN Polyclonal Antibody
PA5110415	SOX30 Polyclonal Antibody
PA5109500	Cdc25A Polyclonal Antibody
PA5109821	HLA-DRA Polyclonal Antibody
PA5109420	ASMT Polyclonal Antibody
PA5110023	COBRA1 Polyclonal Antibody
PA5110398	SMARCD2 Polyclonal Antibody
PA5109629	PNAS4 Polyclonal Antibody
PA5109414	ARRDC3 Polyclonal Antibody
PA5109574	COPZ1 Polyclonal Antibody
PA5109685	p400 Polyclonal Antibody
PA5109368	AGMAT Polyclonal Antibody
PA5109598	CSMD3 Polyclonal Antibody
PA5109741	FUCA2 Polyclonal Antibody
PA5110192	POFUT2 Polyclonal Antibody
PA5109692	ERH Polyclonal Antibody
PA5110256	RAB25 Polyclonal Antibody
PA5110341	SARS2 Polyclonal Antibody
PA5109730	FLT4 Polyclonal Antibody
PA5110427	Spectrin alpha-1 Polyclonal Antibody
PA5109831	HOXC9 Polyclonal Antibody
PA5109565	Contactin 5 Polyclonal Antibody
PA5109908	LGR4 Polyclonal Antibody
PA5109848	CMG1 Polyclonal Antibody
PA5109404	Amphiregulin Polyclonal Antibody
PA5109691	ERCC3 Polyclonal Antibody
PA5109892	K-Ras Polyclonal Antibody
PA5109930	MAF Polyclonal Antibody
PA5110014	NAT2 Polyclonal Antibody
PA5109670	EBF1 Polyclonal Antibody
PA5109358	alpha Actinin 3 Polyclonal Antibody
PA5109710	FAM62B Polyclonal Antibody
PA5110195	POLE2 Polyclonal Antibody
PA5110387	SLC27A1 Polyclonal Antibody
PA5110394	SLC7A9 Polyclonal Antibody
PA5110422	SPOCK3 Polyclonal Antibody

PA5110482	THAP3 Polyclonal Antibody
PA5109716	FCER1G Polyclonal Antibody
PA5110071	PALD1 Polyclonal Antibody
PA5110588	Kaiso Polyclonal Antibody
PA5109702	F5 Polyclonal Antibody
PA5109931	MAFF Polyclonal Antibody
PA5109950	MCM8 Polyclonal Antibody
PA5109985	MRPS31 Polyclonal Antibody
PA5109614	DAAM1 Polyclonal Antibody
PA5109660	DRD2 Polyclonal Antibody
PA5109726	FKBP10 Polyclonal Antibody
PA5109845	IFI30 Polyclonal Antibody
PA5110351	SEPT5 Polyclonal Antibody
PA5109509	CEL Polyclonal Antibody
PA5110004	MYO6 Polyclonal Antibody
PA5110020	NDUFS8 Polyclonal Antibody
PA5110355	SERPINB8 Polyclonal Antibody
PA5109398	AP3D1 Polyclonal Antibody
PA5109643	DNAJB4 Polyclonal Antibody
PA5110324	RPRD2 Polyclonal Antibody
PA5109853	IL1RAPL2 Polyclonal Antibody
PA5109349	A1CF Polyclonal Antibody
PA5109903	LARP7 Polyclonal Antibody
PA5109927	LUC7L2 Polyclonal Antibody
PA5109807	HEATR1 Polyclonal Antibody
PA5110440	STK31 Polyclonal Antibody
PA5110531	TSPYL2 Polyclonal Antibody
PA5109996	MUL1 Polyclonal Antibody
PA5110216	PRDM2 Polyclonal Antibody
PA5109648	DNAJC9 Polyclonal Antibody
PA5109633	DHH Polyclonal Antibody
PA5110545	UFC1 Polyclonal Antibody
PA5109360	ACTRT3 Polyclonal Antibody
PA5110367	SIAE Polyclonal Antibody
PA5110168	PID1 Polyclonal Antibody
PA5110428	SRBD1 Polyclonal Antibody
PA5109696	ETNK1 Polyclonal Antibody
PA5109916	LIMS2 Polyclonal Antibody
PA5109359	ACTR6 Polyclonal Antibody
PA5110486	THOC3 Polyclonal Antibody
PA5110260	RAB40A Polyclonal Antibody
PA5110470	TCP11L2 Polyclonal Antibody
PA5110001	MYL4 Polyclonal Antibody
PA5110328	RRP12 Polyclonal Antibody
PA5110529	TSPAN8 Polyclonal Antibody
PA5109835	HTATSF1 Polyclonal Antibody
PA5110268	RANKL Polyclonal Antibody
PA5109889	CD94 Polyclonal Antibody
PA5109550	CLEC4C Polyclonal Antibody
PA5110060	OLIG3 Polyclonal Antibody
PA5109395	Annexin A13 Polyclonal Antibody
PA5110542	UBLCP1 Polyclonal Antibody
PA5110298	RIN1 Polyclonal Antibody
PA5110299	RIN2 Polyclonal Antibody
PA5110307	RNF17 Polyclonal Antibody
PA5109477	CADPS2 Polyclonal Antibody
PA5109868	ISY1 Polyclonal Antibody
PA5109496	CD300f Polyclonal Antibody
PA5110384	UCP4 Polyclonal Antibody

PA5109956	MGAT5B Polyclonal Antibody
PA5109742	FUT10 Polyclonal Antibody
PA5109704	FA2H Polyclonal Antibody
PA5110475	TFAP2B Polyclonal Antibody
PA5110535	TYRO3 Polyclonal Antibody
PA5110502	TNFRSF19 Polyclonal Antibody
PA5110177	PLAGL2 Polyclonal Antibody
PA5109533	CHMP7 Polyclonal Antibody
PA5109745	FXYD6 Polyclonal Antibody
PA5110251	QPRT Polyclonal Antibody
PA5110291	RHBDD1 Polyclonal Antibody
PA5109466	C2CD5 Polyclonal Antibody
PA5110320	RPL36AL Polyclonal Antibody
PA5109583	CPNE6 Polyclonal Antibody
PA5109782	GPD1L Polyclonal Antibody
PA5110453	SYF2 Polyclonal Antibody
PA5109975	MRPL15 Polyclonal Antibody
PA5109949	MCEE Polyclonal Antibody
PA5109783	GPR128 Polyclonal Antibody
PA5109422	ASZ1 Polyclonal Antibody
PA5109486	CCDC83 Polyclonal Antibody
PA5109401	APOBEC3C Polyclonal Antibody
PA5110067	OTUD6B Polyclonal Antibody
PA5109999	MYCT1 Polyclonal Antibody
PA5109836	HTATSF1 Polyclonal Antibody
PA5109778	GNAI2 Polyclonal Antibody
PA5109803	HBA1 Polyclonal Antibody
PA5110218	PRKACG Polyclonal Antibody
PA5110340	SAA Polyclonal Antibody
PA5110385	SLC26A2 Polyclonal Antibody
PA5109595	CRYGB Polyclonal Antibody
PA5109977	MRPL22 Polyclonal Antibody
PA5109622	DCUN1D1 Polyclonal Antibody
PA5109375	Adenylate Kinase 7 Polyclonal Antibody
PA5109494	CCT6B Polyclonal Antibody
PA5109838	HTRA3 Polyclonal Antibody
PA5109856	GP130 Polyclonal Antibody
PA5110254	RAB1A Polyclonal Antibody
PA5109444	BATF Polyclonal Antibody
PA5110425	SPRY3 Polyclonal Antibody
PA5110278	RBM12 Polyclonal Antibody
PA5109493	CCR7 Polyclonal Antibody
PA5109636	DHX8 Polyclonal Antibody
PA5109895	KRT12 Polyclonal Antibody
PA5110232	PSMB7 Polyclonal Antibody
PA5110235	PTCH Polyclonal Antibody
PA5109377	AKAP4 Polyclonal Antibody
PA5110223	PRPF4B Polyclonal Antibody
PA5109412	JWA Polyclonal Antibody
PA5109621	DCAF4 Polyclonal Antibody
PA5109746	FXYD7 Polyclonal Antibody
PA5110002	MYNN Polyclonal Antibody
PA5110301	RMND5A Polyclonal Antibody
PA5109735	FNDCA Polyclonal Antibody
PA5109423	ATCAY Polyclonal Antibody
PA5109768	GLB1L3 Polyclonal Antibody
PA5110586	YTHDC2 Polyclonal Antibody
PA5109357	ACSS1 Polyclonal Antibody
PA5109449	BMP10 Polyclonal Antibody

PA5109529	CHKA Polyclonal Antibody
PA5109701	EYA2 Polyclonal Antibody
PA5109744	FXD1 Polyclonal Antibody
PA5110210	PPP2R5E Polyclonal Antibody
PA5110237	PTGFR Polyclonal Antibody
PA5110293	RHD Polyclonal Antibody
PA5110322	RPL6 Polyclonal Antibody
PA5110347	SCNN1G Polyclonal Antibody
PA5110511	53BP2 Polyclonal Antibody
PA5110596	ZNF124 Polyclonal Antibody
PA5110334	RTN3 Polyclonal Antibody
PA5109566	CNTRL Polyclonal Antibody
PA5109456	BTBD1 Polyclonal Antibody
PA5110375	SLC13A3 Polyclonal Antibody
PA5109928	LYPLAL1 Polyclonal Antibody
PA5109919	LIPH Polyclonal Antibody
PA5109393	ANO5 Polyclonal Antibody
PA5109387	ANGPTL5 Polyclonal Antibody
PA5110477	TFB1M Polyclonal Antibody
PA5109367	AFM Polyclonal Antibody
PA5109429	ATP1A4 Polyclonal Antibody
PA5109474	CACNA2D1 Polyclonal Antibody
PA5109548	NNT1 Polyclonal Antibody
PA5109799	GUCY1B3 Polyclonal Antibody
PA5109820	HLA-DOA Polyclonal Antibody
PA5110084	PCOLCE Polyclonal Antibody
PA5110100	PFDN4 Polyclonal Antibody
PA5110585	XPNPEP1 Polyclonal Antibody
PA5110400	SMG7 Polyclonal Antibody
PA5109470	C6orf10 Polyclonal Antibody
PA5109556	CLSTN1 Polyclonal Antibody
PA5109571	COLGALT2 Polyclonal Antibody
PA5109356	ACOT9 Polyclonal Antibody
PA5110401	SMPD4 Polyclonal Antibody
PA5109448	BIN3 Polyclonal Antibody
PA5109575	COQ9 Polyclonal Antibody
PA5110374	KCC2 Polyclonal Antibody
PA5109654	DPEP2 Polyclonal Antibody
PA5110270	RAPH1 Polyclonal Antibody
PA5110037	NOC4L Polyclonal Antibody
PA5109734	FN3KRP Polyclonal Antibody
PA5109516	CEP63 Polyclonal Antibody
PA5110580	WDR82 Polyclonal Antibody
PA5110514	NIBP Polyclonal Antibody
PA5110559	USP44 Polyclonal Antibody
PA5109921	LLPH Polyclonal Antibody
PA5109828	HOOK3 Polyclonal Antibody
PA5109717	FCHSD1 Polyclonal Antibody
PA5110509	TOP1MT Polyclonal Antibody
PA5109646	DNAJC24 Polyclonal Antibody
PA5110025	NFAM1 Polyclonal Antibody
PA5109794	GSTT2B Polyclonal Antibody
PA5109814	H2BK12ac Polyclonal Antibody
PA5110306	RNF149 Polyclonal Antibody
PA5109462	C1r Polyclonal Antibody
PA5109578	Carboxypeptidase A3 Polyclonal Antibody
PA5109600	CK1 delta Polyclonal Antibody
PA5109662	DRP2 Polyclonal Antibody
PA5109964	AF9 Polyclonal Antibody

PA5110063	ORC2 Polyclonal Antibody
PA5109843	ICAM5 Polyclonal Antibody
PA5110275	RASSF2 Polyclonal Antibody
PA5109647	DNAJC8 Polyclonal Antibody
PA5110263	RABGAP1 Polyclonal Antibody
PA5110490	TIMM8B Polyclonal Antibody
PA5110568	VGLL1 Polyclonal Antibody
PA5110344	SCLY Polyclonal Antibody
PA5110308	RNF220 Polyclonal Antibody
PA5110578	WDR18 Polyclonal Antibody
PA5109700	EXOSC5 Polyclonal Antibody
PA5110438	STAMBPL1 Polyclonal Antibody
PA5109517	CERK Polyclonal Antibody
PA5109591	CRTC3 Polyclonal Antibody
PA5110551	UPF3A Polyclonal Antibody
PA5110183	PLVAP Polyclonal Antibody
PA5109922	LONP2 Polyclonal Antibody
PA5110498	TMEM44 Polyclonal Antibody
PA5109769	GLCCI1 Polyclonal Antibody
PA5110484	THEM4 Polyclonal Antibody
PA5110069	PACRG Polyclonal Antibody
PA5110365	SGCZ Polyclonal Antibody
PA5109833	HSCB Polyclonal Antibody
PA5110522	TRIM60 Polyclonal Antibody
PA5109471	C9orf72 Polyclonal Antibody
PA5109671	EBF3 Polyclonal Antibody
PA5110052	NUP43 Polyclonal Antibody
PA5109953	MEF2B Polyclonal Antibody
PA5109394	ANT1/ANT2/ANT3/ANT4 Polyclonal Antibody
PA5110377	SLC14A1 Polyclonal Antibody
PA5110464	TBX20 Polyclonal Antibody
PA5109747	ADAP Polyclonal Antibody
PA5109555	CLNS1A Polyclonal Antibody
PA5109703	F7 Polyclonal Antibody
PA5109933	MEKK1 Polyclonal Antibody
PA5110408	SNRPD3 Polyclonal Antibody
PA5110608	ZNF70 Polyclonal Antibody
PA5109963	MLF2 Polyclonal Antibody
PA5109816	H3K36ac Polyclonal Antibody
PA5109817	H3K4ac Polyclonal Antibody
PA5109353	ABI1 Polyclonal Antibody
PA5109603	Chymotrypsin-C Polyclonal Antibody
PA5109652	TDIF2 Polyclonal Antibody
PA5109824	HN1 Polyclonal Antibody
PA5110239	PTOV1 Polyclonal Antibody
PA5109876	KCNK15 Polyclonal Antibody
PA5110493	TMC5 Polyclonal Antibody
PA5110246	Nectin 4 Polyclonal Antibody
PA5109952	MED25 Polyclonal Antibody
PA5109943	MAS1L Polyclonal Antibody
PA5109942	MAS1L Polyclonal Antibody
PA5110389	SLC5A12 Polyclonal Antibody
PA5110607	ZNF677 Polyclonal Antibody
PA5109762	C7orf20 Polyclonal Antibody
PA5110368	SIAH2 Polyclonal Antibody
PA5110103	PGRMC2 Polyclonal Antibody
PA5109872	KCNAB2 Polyclonal Antibody
MA534959	AGK Monoclonal Antibody
PA5109749	G6PC3 Polyclonal Antibody

PA5110018	NDUFA11 Polyclonal Antibody
PA5109948	MBTPS1 Polyclonal Antibody
PA5110576	WASF3 Polyclonal Antibody
PA5109862	INSIG1 Polyclonal Antibody
PA5110558	USP22 Polyclonal Antibody
PA5110325	RSK3 Polyclonal Antibody
PA5110473	TERF2 Polyclonal Antibody
PA5109885	KLF8 Polyclonal Antibody
PA5109666	DUSP15 Polyclonal Antibody
PA5110273	RASA1 Polyclonal Antibody
PA5109673	ECM1 Polyclonal Antibody
PA5109727	FKBP1B Polyclonal Antibody
PA5109846	IFI35 Polyclonal Antibody
PA5110048	NUBP1 Polyclonal Antibody
PA5110166	PI4KB Polyclonal Antibody
PA5109978	MRPL23 Polyclonal Antibody
PA5110022	NEK4 Polyclonal Antibody
PA5110598	ZNF182 Polyclonal Antibody
PA5109454	BRAP Polyclonal Antibody
PA5109832	HS3ST3A1 Polyclonal Antibody
PA5110217	PRG3 Polyclonal Antibody
PA5109376	AKAP3 Polyclonal Antibody
PA5109997	MAD4 Polyclonal Antibody
PA5109738	FSCN2 Polyclonal Antibody
PA5110362	SFTPD Polyclonal Antibody
PA5109619	DCAF12 Polyclonal Antibody
PA5110061	Opticin Polyclonal Antibody
PA5110080	PCF11 Polyclonal Antibody
PA5110414	SOX18 Polyclonal Antibody
PA5110058	OGFOD1 Polyclonal Antibody
PA5109878	KIF17 Polyclonal Antibody
PA5109655	DPP10 Polyclonal Antibody
PA5109502	CDCA5 Polyclonal Antibody
PA5109642	LRRC50 Polyclonal Antibody
PA5109419	ASB7 Polyclonal Antibody
PA5109973	MPP7 Polyclonal Antibody
PA5109737	FREM1 Polyclonal Antibody
PA5110281	RBMY1F Polyclonal Antibody
PA5109426	ATP11C Polyclonal Antibody
PA5109689	EPN2 Polyclonal Antibody
PA5110066	OTUD6A Polyclonal Antibody
PA5109524	CHCHD2 Polyclonal Antibody
PA5110556	USP21 Polyclonal Antibody
PA5110346	SCN5A Polyclonal Antibody
PA5110313	ROS1 Polyclonal Antibody
PA5110051	gp210 Polyclonal Antibody
PA5109473	CACNA1B Polyclonal Antibody
PA5110212	PRB4 Polyclonal Antibody
PA5110361	SFRP5 Polyclonal Antibody
PA5109425	ATF6B Polyclonal Antibody
PA5110450	SURF4 Polyclonal Antibody
PA5109400	APOBEC1 Polyclonal Antibody
PA5109554	CLN8 Polyclonal Antibody
PA5110455	Brachyury Polyclonal Antibody
PA5110227	PSD3 Polyclonal Antibody
PA5109488	CCL19 Polyclonal Antibody
PA5109669	DYNLT3 Polyclonal Antibody
PA5109929	MAdCAM1 Polyclonal Antibody
PA5110169	PIGQ Polyclonal Antibody

PA5109441	BAG4 Polyclonal Antibody
PA5110289	RGL1 Polyclonal Antibody
PA5110171	PI3K p55 gamma Polyclonal Antibody
PA5109932	MAGEH1 Polyclonal Antibody
PA5109525	CHD7 Polyclonal Antibody
PA5109381	ALG9 Polyclonal Antibody
PA5110501	TMX Polyclonal Antibody
PA5110207	PPP1R9B Polyclonal Antibody
PA5109415	ARX Polyclonal Antibody
PA5110286	RCOR2 Polyclonal Antibody
PA5110059	OLFM3 Polyclonal Antibody
PA5110170	PI3K p55 gamma Polyclonal Antibody
PA5110068	OTUD7B Polyclonal Antibody
PA5110312	ROR2 Polyclonal Antibody
PA5110591	ZDHHC8 Polyclonal Antibody
PA5110047	ARK5 Polyclonal Antibody
PA5109546	CIZ1 Polyclonal Antibody
PA5109480	CAPRIN2 Polyclonal Antibody
PA5109808	HERC4 Polyclonal Antibody
PA5110388	SLC4A1 Polyclonal Antibody
PA5109936	JNK3 Polyclonal Antibody
PA5110055	OAS2 Polyclonal Antibody
PA5109911	LIG4 Polyclonal Antibody
PA5110090	PDHA1 Polyclonal Antibody
PA5110095	PDXP Polyclonal Antibody
PA5109564	CNTN1 Polyclonal Antibody
PA5109570	COL9A2 Polyclonal Antibody
PA5109661	DRG2 Polyclonal Antibody
PA5109686	EPB42 Polyclonal Antibody
PA5109767	GK2 Polyclonal Antibody
PA5109842	ICA1 Polyclonal Antibody
PA5109875	Kir7.1 Polyclonal Antibody
PA5110206	PPP1R12B Polyclonal Antibody
PA5110062	OR3A1 Polyclonal Antibody
PA5110074	PAX1 Polyclonal Antibody
PA5110075	PBX3 Polyclonal Antibody
PA5109504	PCTAIRE3 Polyclonal Antibody
PA5110099	PFAS Polyclonal Antibody
PA5109699	EXOSC10 Polyclonal Antibody
PA5110186	PNLIPRP1 Polyclonal Antibody
PA5110224	PR3 Polyclonal Antibody
PA5109352	ABCD4 Polyclonal Antibody
PA5109934	MAP4K2 Polyclonal Antibody
PA5110393	SLC6A8 Polyclonal Antibody
PA5110423	SPRR1A Polyclonal Antibody
PA5110430	SRPK2 Polyclonal Antibody
PA5110460	TBCC Polyclonal Antibody
PA5110528	TSPAN6 Polyclonal Antibody
PA5110492	TM7SF2 Polyclonal Antibody
PA5110503	Tensin 1 Polyclonal Antibody
PA5110593	ZIC1 Polyclonal Antibody
PA5110597	ZNF16 Polyclonal Antibody
PA5109536	Chordin Polyclonal Antibody
PA5109440	B4GALT2 Polyclonal Antibody
PA5110496	TMEM11 Polyclonal Antibody
PA5109417	ASAP2 Polyclonal Antibody
PA5109505	p39 Polyclonal Antibody
PA5109915	LIMD1 Polyclonal Antibody
PA5110083	PCNX Polyclonal Antibody

PA5110330	RRP9 Polyclonal Antibody
PA5110395	NHERF2 Polyclonal Antibody
PA5109869	ITGBL1 Polyclonal Antibody
PA5109787	GRHPR Polyclonal Antibody
PA5110570	VPRBP Polyclonal Antibody
PA5110305	RNF144A Polyclonal Antibody
PA5110053	NUPL1 Polyclonal Antibody
PA5110505	TOM1L1 Polyclonal Antibody
PA5110359	SF3B4 Polyclonal Antibody
PA5110515	TRIM16 Polyclonal Antibody
PA5109523	CGRRF1 Polyclonal Antibody
PA5110602	ZNF460 Polyclonal Antibody
PA5109724	FGFR1OP Polyclonal Antibody
PA5110577	WBP4 Polyclonal Antibody
PA5110417	SP140 Polyclonal Antibody
PA5109966	MLXIP Polyclonal Antibody
PA5109970	MON1B Polyclonal Antibody
PA5109682	ELL2 Polyclonal Antibody
PA5109914	LIMCH1 Polyclonal Antibody
PA5110421	SPG20 Polyclonal Antibody
PA5109389	ANKLE2 Polyclonal Antibody
PA5110333	RTF1 Polyclonal Antibody
PA5109707	FAM155B Polyclonal Antibody
PA5109713	FBXO28 Polyclonal Antibody
PA5110180	PLCL2 Polyclonal Antibody
PA5110557	USP22 Polyclonal Antibody
PA5110494	TMED3 Polyclonal Antibody
PA5109541	CHST5 Polyclonal Antibody
PA5109729	FLRT2 Polyclonal Antibody
PA5110603	ZNF500 Polyclonal Antibody
PA5109410	SGEF Polyclonal Antibody
PA5109714	FBXO8 Polyclonal Antibody
PA5110303	RNF115 Polyclonal Antibody
PA5110036	NOB1 Polyclonal Antibody
PA5109380	ALG5 Polyclonal Antibody
PA5110543	UBQLN3 Polyclonal Antibody
PA5110524	TRNT1 Polyclonal Antibody
PA5110331	RSL24D1 Polyclonal Antibody
PA5109579	CPA4 Polyclonal Antibody
PA5109774	GMIP Polyclonal Antibody
PA5110508	TOMM7 Polyclonal Antibody
PA5109886	KLHL11 Polyclonal Antibody
PA5109540	CHST12 Polyclonal Antibody
PA5109785	GPRC5D Polyclonal Antibody
PA5109991	MTMR8 Polyclonal Antibody
PA5110255	RAB20 Polyclonal Antibody
PA5109644	DNAJC12 Polyclonal Antibody
PA5109860	INPP5E Polyclonal Antibody
PA5110188	PNO1 Polyclonal Antibody
PA5110193	POGLUT1 Polyclonal Antibody
PA5109811	HHATL Polyclonal Antibody
PA5110050	NUFIP2 Polyclonal Antibody
PA5110267	RANBP10 Polyclonal Antibody
PA5110225	PRX Polyclonal Antibody
PA5110215	PRDM12 Polyclonal Antibody
PA5109712	FASTKD5 Polyclonal Antibody
PA5109676	EEFSEC Polyclonal Antibody
PA5110397	SMAP1 Polyclonal Antibody
PA5109543	CHTF18 Polyclonal Antibody

PA5110292	RHBDF1 Polyclonal Antibody
PA5109882	KLC2 Polyclonal Antibody
PA5110419	SPATS2 Polyclonal Antibody
PA5110247	PYCRL Polyclonal Antibody
PA5110601	ZNF426 Polyclonal Antibody
PA5109924	LRRK1 Polyclonal Antibody
PA5110609	ZNF768 Polyclonal Antibody
PA5109620	DCAF17 Polyclonal Antibody
PA5110534	TUBB1 Polyclonal Antibody
PA5110044	NSRP1 Polyclonal Antibody
PA5110101	PGBD1 Polyclonal Antibody
PA5110521	TRIM47 Polyclonal Antibody
PA5110283	RCCD1 Polyclonal Antibody
PA5109802	H2AFV Polyclonal Antibody
PA5110102	PGRP-L Polyclonal Antibody
PA5109562	CNTD1 Polyclonal Antibody
PA5110459	TBC1D16 Polyclonal Antibody
PA5109418	ASB17 Polyclonal Antibody
PA5109708	FAM3D Polyclonal Antibody
PA5109709	FAM3D Polyclonal Antibody
PA5110049	NUDT16 Polyclonal Antibody
PA5109560	CLVS2 Polyclonal Antibody
PA5109534	Chondrolectin Polyclonal Antibody
PA5110441	STK35 Polyclonal Antibody
PA5109917	LIN9 Polyclonal Antibody
PA5110404	SNAI3 Polyclonal Antibody
PA5110262	RAB7B Polyclonal Antibody
PA5110554	USP12 Polyclonal Antibody
PA5110560	USP46 Polyclonal Antibody
PA5109722	FGF17 Polyclonal Antibody
PA5110399	SMC2 Polyclonal Antibody
PA5110451	SV2A Polyclonal Antibody
PA5110046	NTNG1 Polyclonal Antibody
PA5110590	ZC3H12D Polyclonal Antibody
PA5110530	TSPY3 Polyclonal Antibody
PA5109988	COX3 Polyclonal Antibody
PA5110319	RPL35A Polyclonal Antibody
PA5110094	PDLIM7 Polyclonal Antibody
PA5110201	PIT1 Polyclonal Antibody
PA5109867	IRF2BP2 Polyclonal Antibody
PA5109668	DUSP8 Polyclonal Antibody
PA5109993	MT-ND6 Polyclonal Antibody
PA5109844	CD278 (ICOS) Polyclonal Antibody
PA5110336	RXRβ Polyclonal Antibody
MA534960	BDNF Monoclonal Antibody
PA5109373	AJUBA Polyclonal Antibody
PA5109818	H3K23ac Polyclonal Antibody
PA5109424	ATF5 Polyclonal Antibody
PA5110245	PURB Polyclonal Antibody
PA5109976	MRPL18 Polyclonal Antibody
PA5109983	MRPL46 Polyclonal Antibody
PA5109982	MRPL40 Polyclonal Antibody
PA5110106	PHLPP2 Polyclonal Antibody
PA5109455	BRD2 Polyclonal Antibody
PA5109981	MRPL38 Polyclonal Antibody
PA5109792	GSTA1 Polyclonal Antibody
PA5109431	ATP5G2 Polyclonal Antibody
PA5109784	GPR34 Polyclonal Antibody
PA5110519	TRIM37 Polyclonal Antibody

PA5110582	GPR177 Polyclonal Antibody
PA5110243	PTRF Polyclonal Antibody
PA5110485	THOC2 Polyclonal Antibody
PA5110202	OCT6 Polyclonal Antibody
PA5110491	TLE4 Polyclonal Antibody
PA5110081	PCIF1 Polyclonal Antibody
PA5110016	NCAPD2 Polyclonal Antibody
PA5109626	PNPO Polyclonal Antibody
PA5110337	RXRG Polyclonal Antibody
PA5110070	PADI4 Polyclonal Antibody
PA5109796	GTF2F1 Polyclonal Antibody
PA5109770	PI4KA Polyclonal Antibody
PA5110407	SNRPB Polyclonal Antibody
PA5109528	CHIT1 Polyclonal Antibody
PA5109498	CD70 Polyclonal Antibody
PA5110403	SMYD4 Polyclonal Antibody
PA5109850	P2Y6 Polyclonal Antibody
PA5109650	DNase I Polyclonal Antibody
PA5109776	GNA11 Polyclonal Antibody
PA5109587	CPT2 Polyclonal Antibody
PA5109765	GFR alpha-3 Polyclonal Antibody
PA5109453	Kallikrein 2 Polyclonal Antibody
PA5109907	Luteinizing Hormone-beta Polyclonal Antibody
PA5109906	MC2R Polyclonal Antibody
PA5110379	CTRP5 Polyclonal Antibody
PA5109967	MMP25 Polyclonal Antibody
PA5109881	KIFC1 Polyclonal Antibody
PA5110249	PZP Polyclonal Antibody
PA5110376	SLC13A4 Polyclonal Antibody
PA5109801	GULP1 Polyclonal Antibody
PA5110495	TMED9 Polyclonal Antibody
PA5110520	TRIM41 Polyclonal Antibody
PA5109651	DNASE1L1 Polyclonal Antibody
PA5109725	FGL1 Polyclonal Antibody
PA5109439	B4GALNT1 Polyclonal Antibody
PA5109827	HNRNPU Polyclonal Antibody
PA5109705	FADS3 Polyclonal Antibody
PA5110370	SIGLEC5 Polyclonal Antibody
PA5109935	MAP7 Polyclonal Antibody
PA5109492	Cyclin E2 Polyclonal Antibody
PA5109627	DDX46 Polyclonal Antibody
PA5110284	RCL1 Polyclonal Antibody
PA5109894	KRR1 Polyclonal Antibody
PA5109693	ERP44 Polyclonal Antibody
PA5110173	PISD Polyclonal Antibody
PA5110371	SIGLEC9 Polyclonal Antibody
PA5109374	AK3 Polyclonal Antibody
PA5109436	AZIN1 Polyclonal Antibody
PA5109484	CG2D1A Polyclonal Antibody
PA5109984	MRPS26 Polyclonal Antibody
PA5110561	UTS2 Polyclonal Antibody
PA5109520	Properdin Polyclonal Antibody
PA5110479	DP1 Polyclonal Antibody
PA5109871	JAK1 Polyclonal Antibody
PA5110264	RAG2 Polyclonal Antibody
PA5109672	ECE1 Polyclonal Antibody
PA5110222	PRPF31 Polyclonal Antibody
PA5109616	DAPK1 Polyclonal Antibody
PA5109858	Activin A Polyclonal Antibody

PA5109606	Cathepsin K Polyclonal Antibody
PA5110054	TAP Polyclonal Antibody
PA5110323	RPL7 Polyclonal Antibody
PA5110031	NHP2 Polyclonal Antibody
PA5109403	AQR Polyclonal Antibody
PA5110405	SNRNP200 Polyclonal Antibody
PA5109437	Azurocidin Polyclonal Antibody
PA5109467	C3a Receptor Polyclonal Antibody
PA5110316	RPL11 Polyclonal Antibody
PA5110390	SLC5A6 Polyclonal Antibody
PA5110265	RAMP1 Polyclonal Antibody
PA5109680	eIF2 gamma Polyclonal Antibody
PA5109756	GATM Polyclonal Antibody
PA5110035	NNMT Polyclonal Antibody
PA5110198	POLR2F Polyclonal Antibody
PA5109399	APLP1 Polyclonal Antibody
PA5109593	CRY2 Polyclonal Antibody
PA5109458	BTNL2 Polyclonal Antibody
PA5109481	CARD11 Polyclonal Antibody
PA5110532	TSSK3 Polyclonal Antibody
PA5109405	ARF4 Polyclonal Antibody
PA5109723	FGF6 Polyclonal Antibody
PA5109793	GSTM3 Polyclonal Antibody
PA5109766	GIMAP5 Polyclonal Antibody
PA5109884	KLF7 Polyclonal Antibody
PA5110432	SSFA2 Polyclonal Antibody
PA5110572	VPS28 Polyclonal Antibody
PA5110409	SNX15 Polyclonal Antibody
PA5110039	NPL Polyclonal Antibody
PA5109804	HBD Polyclonal Antibody
PA5109870	ITIH3 Polyclonal Antibody
PA5109628	DDX56 Polyclonal Antibody
PA5109468	Complement C4 Polyclonal Antibody
PA5109896	KRT12 Polyclonal Antibody
PA5110127	Phospho-HNF4A (Ser304) Polyclonal Antibody
PA5110153	Phospho-RAC1 (Ser71) Polyclonal Antibody
PA5110121	Phospho-FOXO1 (Ser256) Polyclonal Antibody
PA5110161	Phospho-Stathmin 1 (Ser25) Polyclonal Antibody
PA5110117	Phospho-eNOS (Ser1177) Polyclonal Antibody
PA5110157	Phospho-STAT6 (Tyr641) Polyclonal Antibody
PA5110145	Phospho-PKC zeta (Thr410) Polyclonal Antibody
PA5110130	Phospho-JAK3 (Tyr980, Tyr981) Polyclonal Antibody
PA5110111	Phospho-HP1 gamma (Ser83) Polyclonal Antibody
PA5110143	Phospho-PAK2 (Ser20) Polyclonal Antibody
PA5110108	Phospho-BCAR1 (Tyr165) Polyclonal Antibody
PA5110118	Phospho-EphA2 (Tyr772) Polyclonal Antibody
PA5110113	Phospho-CRMP2 (Thr514) Polyclonal Antibody
PA5110124	Phospho-GluR1 (Ser845) Polyclonal Antibody
PA5110123	Phospho-FOXO3A (Ser294) Polyclonal Antibody
PA5110149	Phospho-POLR2A (Ser5) Polyclonal Antibody
PA5110138	Phospho-NEDD4L (Ser448) Polyclonal Antibody
PA5110129	Phospho-IRS1 (Ser1101) Polyclonal Antibody
PA5110122	Phospho-FOXO1/FOXO3/FOXO4 (Thr24, Thr32, Thr28) Polyclonal Antibody
PA5110126	Phospho-Histone H3 (Ser10, Thr11) Polyclonal Antibody
PA5110128	Phospho-IR/IGF1R (Tyr1165, Tyr1166, Tyr1189, Tyr1190) Polyclonal Antibody
PA5110150	Phospho-PP2A alpha/beta (Tyr307) Polyclonal Antibody
PA5110160	Phospho-MST1 (STK4) (Thr387) Polyclonal Antibody
PA5110163	Phospho-YAP1 (Ser397) Polyclonal Antibody
PA5110154	Phospho-Raptor (Ser792) Polyclonal Antibody

PA5110158	Phospho-MST2 (Thr117) Polyclonal Antibody
PA5110159	Phospho-MST2 (Thr384) Polyclonal Antibody
PA5110136	Phospho-MCL1 (Thr163) Polyclonal Antibody
PA5110120	Phospho-FADD (Ser194) Polyclonal Antibody
PA5110119	Phospho-EphA3 (Tyr779) Polyclonal Antibody
PA5110109	Phospho-BCAR1 (Tyr410) Polyclonal Antibody
PA5110144	Phospho-PDGFRB (Tyr771) Polyclonal Antibody
PA5110107	Phospho-AKT1/AKT2/AKT3 (Ser473) Polyclonal Antibody
PA5116837	AAR2 Polyclonal Antibody
PA5116610	A33 Polyclonal Antibody
PA5116765	ABHD12 Polyclonal Antibody
PA5116772	ABHD4 Polyclonal Antibody
PA5116761	ABLIM3 Polyclonal Antibody
PA5116686	ACBD3 Polyclonal Antibody
PA5116749	Acetyl-FOXO1 (Lys262, Lys265, Lys274) Polyclonal Antibody
PA5116567	ACSS3 Polyclonal Antibody
PA5116639	ACTR1B Polyclonal Antibody
PA5116865	ADH1A/ADH1B/ADH1C Polyclonal Antibody
PA5116693	ADARB1 Polyclonal Antibody
PA5116801	ALG13 Polyclonal Antibody
PA5116671	ALOX12B Polyclonal Antibody
PA5116959	ALOX5 Polyclonal Antibody
PA5116800	AMOTL1 Polyclonal Antibody
PA5117012	Amphiphysin Polyclonal Antibody
PA5116856	ANKRD46 Polyclonal Antibody
PA5116737	ANP32E Polyclonal Antibody
PA5116622	APC11 Polyclonal Antibody
PA5117016	APLP2 Polyclonal Antibody
PA5116647	APOA1BP Polyclonal Antibody
PA5116874	APOBEC2 Polyclonal Antibody
PA5116779	ARHGAP30 Polyclonal Antibody
PA5116948	ASB17 Polyclonal Antibody
PA5116694	ASS1 Polyclonal Antibody
PA5116695	ATP5I Polyclonal Antibody
PA5116928	ATP6V1G2 Polyclonal Antibody
PA5117139	B4GALT1 Polyclonal Antibody
PA5116624	BAALC Polyclonal Antibody
PA5116767	BBS7 Polyclonal Antibody
PA5116531	Bcl-3 Polyclonal Antibody
PA5116537	Beclin 1 Polyclonal Antibody
PA5116905	Biotin Polyclonal Antibody
PA5116904	BRD4 Polyclonal Antibody
PA5116913	BRD4-S Polyclonal Antibody
PA5116939	BRP44 Polyclonal Antibody
PA5116938	BRP44L Polyclonal Antibody
PA5116844	BRUNOL5 Polyclonal Antibody
PA5116850	BTBD9 Polyclonal Antibody
PA5116677	C11orf9 Polyclonal Antibody
PA5116895	C1orf146 Polyclonal Antibody
PA5116589	C7orf50 Polyclonal Antibody
PA5116840	C9orf114 Polyclonal Antibody
PA5116830	CABP1 Polyclonal Antibody
PA5116964	CAF1 p150 Polyclonal Antibody
PA5116762	Calpain 7 Polyclonal Antibody
PA5116734	CaMKV Polyclonal Antibody
PA5116776	CAMSAP1 Polyclonal Antibody
PA5116781	CCDC137 Polyclonal Antibody
PA5117070	CCS Polyclonal Antibody
PA5116626	CD163L1 Polyclonal Antibody

PA5116634	CD27 Polyclonal Antibody
PA5116751	CD28 Polyclonal Antibody
PA5116889	CDH8 Polyclonal Antibody
PA5116623	CDK5RAP2 Polyclonal Antibody
PA5116732	CELA2A Polyclonal Antibody
PA5116824	GEND Polyclonal Antibody
PA5116863	cGKII Polyclonal Antibody
PA5116845	CHCHD7 Polyclonal Antibody
PA5116579	CHMP7 Polyclonal Antibody
PA5116591	CIS Polyclonal Antibody
PA5116673	CK1 alpha Polyclonal Antibody
PA5116733	Claspin Polyclonal Antibody
PA5116747	CLEC2B Polyclonal Antibody
PA5116974	CLEC4D Polyclonal Antibody
PA5116924	CLEC4E Polyclonal Antibody
PA5117052	CMG2 Polyclonal Antibody
PA5116757	COG1 Polyclonal Antibody
PA5116659	COG4 Polyclonal Antibody
PA5117048	Coilin Polyclonal Antibody
PA5116871	Collagen X Polyclonal Antibody
PA5116903	COL16A1 Polyclonal Antibody
PA5117056	COX4 Polyclonal Antibody
PA5116722	CPSF1 Polyclonal Antibody
PA5116670	CXCL6 Polyclonal Antibody
PA5116601	DEXI Polyclonal Antibody
PA5116592	Cyclin A1 Polyclonal Antibody
PA5117014	Cytochrome B5 Polyclonal Antibody
PA5116566	DAPP1 Polyclonal Antibody
PA5116575	DAZAP2 Polyclonal Antibody
PA5116570	DAZL Polyclonal Antibody
PA5116921	DCP1B Polyclonal Antibody
PA5116932	DDX55 Polyclonal Antibody
PA5116674	DDX6 Polyclonal Antibody
PA5117076	DEFA1 Polyclonal Antibody
PA5116523	DFNB31 Polyclonal Antibody
PA5116788	DHX9 Polyclonal Antibody
PA5116818	Di-Methyl Lysine Polyclonal Antibody
PA5116822	DIP2C Polyclonal Antibody
PA5116777	DJC18 Polyclonal Antibody
PA5116930	DNAJC17 Polyclonal Antibody
PA5116748	DOCK1 Polyclonal Antibody
PA5116546	DPYSL3 Polyclonal Antibody
PA5116975	EGR4 Polyclonal Antibody
PA5116554	EHD2 Polyclonal Antibody
PA5116806	eIF2b delta Polyclonal Antibody
PA5116950	EIF4A1/EIF4A2/EIF4A3 Polyclonal Antibody
PA5116997	eIF4E2 Polyclonal Antibody
PA5116573	EML2 Polyclonal Antibody
PA5116778	EML3 Polyclonal Antibody
PA5116576	EMR2 Polyclonal Antibody
PA5116786	ENAH Polyclonal Antibody
PA5116513	ENO1 Polyclonal Antibody
PA5117085	ENOX2 Polyclonal Antibody
PA5116528	ERK2 Polyclonal Antibody
PA5116754	FCGR1B Polyclonal Antibody
PA5117154	EXOSC4 Polyclonal Antibody
PA5117158	EYA1 Polyclonal Antibody
PA5116858	FAAH2 Polyclonal Antibody
PA5117027	FADD Polyclonal Antibody

PA5116811	FAM5C Polyclonal Antibody
PA5116735	FAT4 Polyclonal Antibody
PA5116649	FBXL10 Polyclonal Antibody
PA5116658	FBXO21 Polyclonal Antibody
PA5116640	FBXO3 Polyclonal Antibody
PA5116803	FBXW8 Polyclonal Antibody
PA5116797	FGF2 Polyclonal Antibody
PA5116637	FHL5 Polyclonal Antibody
PA5116593	Fibrinogen alpha chain Polyclonal Antibody
PA5117029	FIS1 Polyclonal Antibody
PA5117152	FKBP1B Polyclonal Antibody
PA5116799	FNDC5 Polyclonal Antibody
PA5117050	Folliculin Polyclonal Antibody
PA5116852	GTSF1 Polyclonal Antibody
PA5116792	FOXO1/FOXO3/FOXO4 Polyclonal Antibody
PA5116608	FUT3 Polyclonal Antibody
PA5117062	Galectin 2 Polyclonal Antibody
PA5116775	GCP5 Polyclonal Antibody
PA5116646	GCP6 Polyclonal Antibody
PA5116720	Gemin 5 Polyclonal Antibody
PA5116612	Gemin 7 Polyclonal Antibody
PA5116697	GFI1 Polyclonal Antibody
PA5116669	GlyT1 Polyclonal Antibody
PA5116532	GlyT2 Polyclonal Antibody
PA5116739	GORAB Polyclonal Antibody
PA5116678	GRINA Polyclonal Antibody
PA5117041	GRK6 Polyclonal Antibody
PA5116809	GRSF1 Polyclonal Antibody
PA5116635	GSTM5 Polyclonal Antibody
PA5116828	GTPBP6 Polyclonal Antibody
PA5117102	GYPB Polyclonal Antibody
PA5116944	H2BE105me1 Polyclonal Antibody
PA5116984	H3R26me2s Polyclonal Antibody
PA5116936	H3K27ac Polyclonal Antibody
PA5116965	H3K27me1 Polyclonal Antibody
PA5116937	H3K36ac Polyclonal Antibody
PA5117092	H3K9ac Polyclonal Antibody
PA5116966	H3R2me2s Polyclonal Antibody
PA5116574	H4K12ac Polyclonal Antibody
PA5116986	H4R3me1 Polyclonal Antibody
PA5117159	HACE1 Polyclonal Antibody
PA5116766	HEYL Polyclonal Antibody
PA5117133	HIC5 Polyclonal Antibody
PA5116542	HIP1R Polyclonal Antibody
PA5116785	Histone H3ac (pan-acetyl) Polyclonal Antibody
PA5116931	HJURP Polyclonal Antibody
PA5116755	HLA-DPB1 Polyclonal Antibody
PA5116675	HMGCR Polyclonal Antibody
PA5116698	HMGCR Polyclonal Antibody
PA5116943	hnRNP A3 Polyclonal Antibody
PA5116699	HOXB4 Polyclonal Antibody
PA5116700	Interferon alpha-10 Polyclonal Antibody
PA5117025	HPa1 Polyclonal Antibody
PA5116808	HSPB11 Polyclonal Antibody
PA5116521	HuD Polyclonal Antibody
PA5116653	IGHMBP2 Polyclonal Antibody
PA5116832	Importin 7 Polyclonal Antibody
PA5116585	INMT Polyclonal Antibody
PA5116631	ITCH Polyclonal Antibody

PA5116536	ITPRIP Polyclonal Antibody
PA5116582	JAKMIP2 Polyclonal Antibody
PA5116833	JM4 Polyclonal Antibody
PA5116679	KAL1 Polyclonal Antibody
PA5116620	KANK2 Polyclonal Antibody
PA5116848	KBTBD8 Polyclonal Antibody
PA5116625	KCNH6 Polyclonal Antibody
PA5116656	KCNK6 Polyclonal Antibody
PA5116627	KIF7 Polyclonal Antibody
PA5116581	Kir2.2 (KCNJ12) Polyclonal Antibody
PA5116877	KLF11 Polyclonal Antibody
PA5117090	KLF9 Polyclonal Antibody
PA5117037	KMT2E Polyclonal Antibody
PA5117089	LACTB Polyclonal Antibody
PA5116897	Lactyl-Histone H3 (Lys14) Polyclonal Antibody
PA5116896	Lactyl-Histone H3 (Lys18) Polyclonal Antibody
PA5116900	Lactyl-Histone H4 (Lys5) Polyclonal Antibody
PA5116496	Lamin A/C Polyclonal Antibody
PA5116611	LRRC61 Polyclonal Antibody
PA5116746	Lamin A/C Polyclonal Antibody
PA5116606	LENG8 Polyclonal Antibody
PA5116661	LHX5 Polyclonal Antibody
PA5116908	LMO2 Polyclonal Antibody
PA5116584	LRP1 Polyclonal Antibody
PA5116690	LSM14A Polyclonal Antibody
PA5117001	MA2 Polyclonal Antibody
PA5116979	MAdCAM1 Polyclonal Antibody
PA5116666	MAGOHB Polyclonal Antibody
PA5116738	MAML2 Polyclonal Antibody
PA5116680	MATN2 Polyclonal Antibody
PA5117031	MATR3 Polyclonal Antibody
PA5116920	MBD6 Polyclonal Antibody
PA5116971	MEF2C Polyclonal Antibody
PA5116873	Methyl-YAP1 (Lys342) Polyclonal Antibody
PA5116820	MHC Class II (H2-Aa) Polyclonal Antibody
PA5116586	MIOX Polyclonal Antibody
PA5116854	MITD1 Polyclonal Antibody
PA5116976	MKP-1 Polyclonal Antibody
PA5116906	MLL4 Polyclonal Antibody
PA5116821	MLLT1 Polyclonal Antibody
PA5116704	MMP13 Polyclonal Antibody
PA5116980	MMP17 Polyclonal Antibody
PA5117038	MMP2 Polyclonal Antibody
PA5116819	Mono-Methyl Arginine Polyclonal Antibody
PA5116817	Mono-Methyl Lysine Polyclonal Antibody
PA5116773	MRPL43 Polyclonal Antibody
PA5116838	MRPS18C Polyclonal Antibody
PA5117153	MRPS25 Polyclonal Antibody
PA5116789	MT-ATP6 Polyclonal Antibody
PA5116535	MTSS1 Polyclonal Antibody
PA5116524	MYH11 Polyclonal Antibody
PA5116555	MYH6 Polyclonal Antibody
PA5117077	Myosin heavy chain 1 Polyclonal Antibody
PA5116688	MYPN Polyclonal Antibody
PA5116774	NAF1 Polyclonal Antibody
PA5116780	NAP1L5 Polyclonal Antibody
PA5116807	NAPB Polyclonal Antibody
PA5116533	NARS Polyclonal Antibody
PA5116981	N-cadherin Polyclonal Antibody

PA5116636	NDUFA4 Polyclonal Antibody
PA5116705	NDUFB1 Polyclonal Antibody
PA5116753	NDUFB9 Polyclonal Antibody
PA5116609	NDUFC2 Polyclonal Antibody
PA5116577	NECAP2 Polyclonal Antibody
PA5116910	NELL2 Polyclonal Antibody
PA5116500	NeuN Polyclonal Antibody
PA5117064	NFS1 Polyclonal Antibody
PA5116934	Niban-like protein Polyclonal Antibody
PA5116583	NME6 Polyclonal Antibody
PA5116890	NOP56 Polyclonal Antibody
PA5116689	NOTCH1 Polyclonal Antibody
PA5116760	NPM3 Polyclonal Antibody
PA5117047	NR2C2 Polyclonal Antibody
PA5116925	NRCAM Polyclonal Antibody
PA5116814	OSR1 Polyclonal Antibody
PA5117022	NSD3 Polyclonal Antibody
PA5116701	NSDHL Polyclonal Antibody
PA5117155	NUP160 Polyclonal Antibody
PA5116835	NUP205 Polyclonal Antibody
PA5116879	NYREN18 Polyclonal Antibody
PA5116619	OLFM4 Polyclonal Antibody
PA5116543	OPLAH Polyclonal Antibody
PA5116847	PANK2 Polyclonal Antibody
PA5116886	Pan-cadherin Polyclonal Antibody
PA5116518	PARG Polyclonal Antibody
PA5116682	PARP2 Polyclonal Antibody
PA5116812	PCQAP Polyclonal Antibody
PA5116660	PDE11A Polyclonal Antibody
PA5116741	PDE12 Polyclonal Antibody
PA5116729	PDP2 Polyclonal Antibody
PA5116550	PEG3 Polyclonal Antibody
PA5116926	Pericentrin Polyclonal Antibody
PA5116654	PEX12 Polyclonal Antibody
PA5117231	Phospho-(Ser/Thr) ATM + ATR Substrate Polyclonal Antibody
PA5117189	Phospho-Acetyl-CoA Carboxylase (Ser79) Polyclonal Antibody
PA5117195	Phospho-alpha Synuclein (Ser129) Polyclonal Antibody
PA5117221	Phospho-AMPK alpha-1 (Ser485) Polyclonal Antibody
PA5117258	Phospho-AMPK alpha-1,2 (Thr183, Thr172) Polyclonal Antibody
PA5117269	Phospho-AMPK beta-1 (Ser182) Polyclonal Antibody
PA5117273	Phospho-ASK1 (Thr845) Polyclonal Antibody
PA5117239	Phospho-ATP1A1 (Ser16) Polyclonal Antibody
PA5117241	Phospho-Bcl-2 (Ser87) Polyclonal Antibody
PA5117233	Phospho-Aurora A/B/C (Thr288, Thr232, Thr198) Polyclonal Antibody
PA5117242	Phospho-beta Catenin (Ser29, Ser33, Ser37, Thr41) Polyclonal Antibody
PA5117262	Phospho-BRD4 (Ser1070) Polyclonal Antibody
PA5117254	Phospho-BRD4 (Thr204) Polyclonal Antibody
PA5117238	Phospho-C/EBP beta (Thr235) Polyclonal Antibody
PA5117261	Phospho-cGAS (Ser291) Polyclonal Antibody
PA5117255	Phospho-Cdc25A (Ser76) Polyclonal Antibody
PA5117191	Phospho-CDK1 (Thr161) Polyclonal Antibody
PA5117192	Phospho-CDK2 (Thr160) Polyclonal Antibody
PA5117216	Phospho-CDK9 (Thr186) Polyclonal Antibody
PA5117206	Phospho-CHK2 (Ser516) Polyclonal Antibody
PA5117265	Phospho-c-Jun (Ser63) Polyclonal Antibody
PA5117274	Phospho-c-Jun (Thr91) Polyclonal Antibody
PA5117178	Phospho-c-Myc (Ser62) Polyclonal Antibody
PA5117218	Phospho-CrkL (Tyr207) Polyclonal Antibody
PA5117249	Phospho-EGFR (Thr678) Polyclonal Antibody

PA5117213	Phospho-eIF4G (Ser1108) Polyclonal Antibody
PA5117243	Phospho-EphA2 (Ser897) Polyclonal Antibody
PA5117252	Phospho-ERK1/2 (Thr202, Tyr204) Polyclonal Antibody
PA5117253	Phospho-Ezrin/Radixin/Moesin (Thr558, Thr567, Thr564) Polyclonal Antibody
PA5117208	Phospho-FOXO3A (Ser253) Polyclonal Antibody
PA5117272	Phospho-FOXO3A (Ser425) Polyclonal Antibody
PA5117224	Phospho-Histone H3 (Thr45) Polyclonal Antibody
PA5117225	Phospho-Histone H3 (Thr6) Polyclonal Antibody
PA5117187	Phospho-HSP27 (Ser78) Polyclonal Antibody
PA5117275	Phospho-IkB alpha (Ser32) Polyclonal Antibody
PA5117248	Phospho-IKK gamma (Ser376) Polyclonal Antibody
PA5117229	Phospho-JAK2 (Tyr1007, Tyr1008) Polyclonal Antibody
PA5117257	Phospho-K-Ras (Ser89) Polyclonal Antibody
PA5117181	Phospho-LCK (Tyr394) Polyclonal Antibody
PA5117223	Phospho-LDHA (Tyr10) Polyclonal Antibody
PA5117266	Phospho-LRP6 (Ser1490) Polyclonal Antibody
PA5117175	Phospho-MAP2K4 (Thr261) Polyclonal Antibody
PA5117267	Phospho-MCM2 (Ser139) Polyclonal Antibody
PA5117184	Phospho-MEK1/MEK2 (Ser217, Ser221) Polyclonal Antibody
PA5117260	Phospho-MLKL (Ser441) Polyclonal Antibody
PA5117236	Phospho-MYL9 (Thr18, Ser19) Polyclonal Antibody
PA5117228	Phospho-MYPT1 (Thr853) Polyclonal Antibody
PA5117268	Phospho-NCF4 (Thr154) Polyclonal Antibody
PA5117180	Phospho-NFkB p50 (Ser337) Polyclonal Antibody
PA5117256	Phospho-Nrf2 (Ser40) Polyclonal Antibody
PA5117234	Phospho-p53 (Ser15) Polyclonal Antibody
PA5117235	Phospho-p53 (Ser6) Polyclonal Antibody
PA5117211	Phospho-PFKFB2 (Ser483) Polyclonal Antibody
PA5117226	Phospho-Phospholamban (Ser16) Polyclonal Antibody
PA5117237	Phospho-PKC alpha (Ser657, Tyr658) Polyclonal Antibody
PA5117210	Phospho-PKC delta (Thr505) Polyclonal Antibody
PA5117245	Phospho-PKC mu (Ser916) Polyclonal Antibody
PA5117219	Phospho-PLCG1 (Ser1248) Polyclonal Antibody
PA5117214	Phospho-PLCG2 (Tyr1217) Polyclonal Antibody
PA5117212	Phospho-PP1 alpha (Thr320) Polyclonal Antibody
PA5117190	Phospho-PRMT4 (Ser228) Polyclonal Antibody
PA5117194	Phospho-Rb (Ser780) Polyclonal Antibody
PA5117263	Phospho-SHC (Tyr427) Polyclonal Antibody
PA5117188	Phospho-SMAD2 (Ser465, Ser467) Polyclonal Antibody
PA5117232	Phospho-SMO (Ser611, Ser615) Polyclonal Antibody
PA5117185	Phospho-Stathmin 1 (Ser38) Polyclonal Antibody
PA5117244	Phospho-Tau (Ser198, Ser199, Ser202) Polyclonal Antibody
PA5117246	Phospho-Tau (Ser416) Polyclonal Antibody
PA5117230	Phospho-Tau (Thr231) Polyclonal Antibody
PA5117251	Phospho-TSC2 (Ser1387) Polyclonal Antibody
PA5117220	Phospho-TSC2 (Thr1462) Polyclonal Antibody
MA538227	Phosphotyrosine Monoclonal Antibody
PA5117247	Phospho-Tyrosine Hydroxylase (Ser40) Polyclonal Antibody
PA5117209	Phospho-ULK1 (Ser555) Polyclonal Antibody
PA5117270	Phospho-ULK1 (Ser638) Polyclonal Antibody
PA5117182	Phospho-VASP (Ser239) Polyclonal Antibody
PA5117215	Phospho-Vimentin (Ser39) Polyclonal Antibody
PA5117196	Phospho-Zap-70 (Tyr493) Polyclonal Antibody
PA5116730	PHRF1 Polyclonal Antibody
PA5116768	PI4K2B Polyclonal Antibody
PA5116882	PIGO Polyclonal Antibody
PA5116853	PIGU Polyclonal Antibody
PA5116743	PILRA Polyclonal Antibody
PA5116823	PI3K p85/p55 Polyclonal Antibody

PA5116804	PIN4 Polyclonal Antibody
PA5116559	PIWIL4 Polyclonal Antibody
PA5116664	PKC mu Polyclonal Antibody
PA5117066	PNPLA3 Polyclonal Antibody
PA5116551	POLR1E Polyclonal Antibody
PA5116630	POLRMT Polyclonal Antibody
PA5116590	POMGNT2 Polyclonal Antibody
PA5117151	POMP Polyclonal Antibody
PA5116681	PON3 Polyclonal Antibody
PA5116864	POU3F2 Polyclonal Antibody
PA5116618	PPM1F Polyclonal Antibody
PA5116861	PPP2R1A/PPP2R2A/PPP2R1B Polyclonal Antibody
PA5116862	PPP2R2B/PPP2R2C Polyclonal Antibody
PA5116556	PRDX3 Polyclonal Antibody
PA5116707	PRIM1 Polyclonal Antibody
PA5116548	PRPS2 Polyclonal Antibody
PA5116740	PRKCDBP Polyclonal Antibody
PA5116708	PROP1 Polyclonal Antibody
PA5116725	PRR16 Polyclonal Antibody
PA5116872	PTGES Polyclonal Antibody
PA5116721	PUMA alpha Polyclonal Antibody
PA5116641	PXR Polyclonal Antibody
PA5116683	RAB26 Polyclonal Antibody
PA5116763	RAB3GAP2 Polyclonal Antibody
PA5116902	RAI16 Polyclonal Antibody
PA5116790	RALGDS Polyclonal Antibody
PA5116652	RAPGEF1 Polyclonal Antibody
PA5116841	RASIP1 Polyclonal Antibody
PA5116492	Rb Polyclonal Antibody
PA5116831	RCE1 Polyclonal Antibody
PA5117068	RECK Polyclonal Antibody
PA5116724	REEP2 Polyclonal Antibody
PA5116684	RM30 Polyclonal Antibody
PA5117028	RMDN3 Polyclonal Antibody
PA5116662	RNF185 Polyclonal Antibody
PA5116866	RPS21 Polyclonal Antibody
PA5116826	RPS8 Polyclonal Antibody
PA5116933	RREB1 Polyclonal Antibody
PA5117004	RT28 Polyclonal Antibody
PA5116595	RYBP Polyclonal Antibody
PA5116597	S100A5 Polyclonal Antibody
MA538034	SARS-CoV-2 Nucleocapsid Monoclonal Antibody
PA5116742	SLC26A11 Polyclonal Antibody
PA5116941	SARS-CoV-2 NSP1 Polyclonal Antibody
PA5116945	SARS-CoV-2 NSP12 Polyclonal Antibody
PA5116894	SARS-CoV-2 Nucleocapsid Polyclonal Antibody
PA5116946	SARS-CoV-2 ORF3a Polyclonal Antibody
PA5116915	SARS-CoV-2 Spike Protein (RBD) Polyclonal Antibody
PA5117084	SCO2 Polyclonal Antibody
PA5116752	SCRIB Polyclonal Antibody
PA5116868	SEC62 Polyclonal Antibody
PA5116663	SERPINA9 Polyclonal Antibody
PA5116810	SETD1A Polyclonal Antibody
PA5116923	SETD1B Polyclonal Antibody
PA5117000	SF3A1 Polyclonal Antibody
PA5116709	SFRS8 Polyclonal Antibody
PA5117156	SH3PXD2A Polyclonal Antibody
PA5116834	SHANK2 Polyclonal Antibody
PA5116878	SIGLEC15 Polyclonal Antibody

PA5117095	SIRT3 Polyclonal Antibody
PA5116867	SLC22A3 Polyclonal Antibody
PA5116714	SLC25A16 Polyclonal Antibody
PA5116600	SLC25A17 Polyclonal Antibody
PA5116578	SLC25A23 Polyclonal Antibody
PA5116522	SLC25A43 Polyclonal Antibody
PA5116851	SLC26A9 Polyclonal Antibody
PA5117157	SLC34A2 Polyclonal Antibody
PA5116512	SLC4A7 Polyclonal Antibody
PA5117124	SLC5A7 Polyclonal Antibody
PA5116629	SLC7A11 Polyclonal Antibody
PA5116607	SLCO6A1 Polyclonal Antibody
PA5116544	SPOP Polyclonal Antibody
PA5116929	SLFN12 Polyclonal Antibody
PA5116713	SLMAP Polyclonal Antibody
PA5116846	SLTM Polyclonal Antibody
PA5116728	SMEK2 Polyclonal Antibody
PA5116655	SNAPC2 Polyclonal Antibody
PA5116692	SNIP1 Polyclonal Antibody
PA5116499	SOCS3 Polyclonal Antibody
PA5116638	SOX13 Polyclonal Antibody
PA5116553	SP3 Polyclonal Antibody
PA5116883	SPAG11A Polyclonal Antibody
PA5116829	SRSF11 Polyclonal Antibody
PA5116719	SRSF8 Polyclonal Antibody
PA5116687	STARD3NL Polyclonal Antibody
PA5116764	SUMF2 Polyclonal Antibody
PA5117123	STX12 Polyclonal Antibody
PA5116784	Succinyl-Histone H3 (Lys79) Polyclonal Antibody
PA5116617	SUCLG1 Polyclonal Antibody
PA5116891	SUPT16H Polyclonal Antibody
PA5116758	SUPT7L Polyclonal Antibody
PA5116989	SUV39H1 Polyclonal Antibody
PA5116885	TAB3 Polyclonal Antibody
PA5116580	TAF6L Polyclonal Antibody
PA5116716	Taok2 Polyclonal Antibody
PA5116736	TAS1R2 Polyclonal Antibody
PA5116558	TDO2 Polyclonal Antibody
PA5116657	TEB4 Polyclonal Antibody
PA5116718	TELO2 Polyclonal Antibody
PA5117087	TES Polyclonal Antibody
PA5116565	THOC7 Polyclonal Antibody
PA5116855	TMEM18 Polyclonal Antibody
PA5116770	TMEM48 Polyclonal Antibody
PA5116922	TMEM61 Polyclonal Antibody
PA5116605	TMX Polyclonal Antibody
PA5116723	TMX2 Polyclonal Antibody
PA5116710	TOP3A Polyclonal Antibody
PA5116849	TOPORS Polyclonal Antibody
PA5116643	TRIM62 Polyclonal Antibody
PA5116816	Tri-Methyl Lysine Polyclonal Antibody
PA5116919	Tri-Methyl Lysine Polyclonal Antibody
PA5116898	TrkA/TrkB/TrkC Polyclonal Antibody
PA5116860	TRPML2 Polyclonal Antibody
PA5116615	TSEN54 Polyclonal Antibody
PA5116843	TSHZ3 Polyclonal Antibody
PA5116869	TSN Polyclonal Antibody
PA5116628	TWIST1 Polyclonal Antibody
PA5117019	UBE2C Polyclonal Antibody

PA5116759	UBE4B Polyclonal Antibody
PA5116568	UCP4 Polyclonal Antibody
PA5116888	ULK1/ULK2 Polyclonal Antibody
PA5116665	USP33 Polyclonal Antibody
PA5116839	UTP11L Polyclonal Antibody
PA5116771	UTP6 Polyclonal Antibody
PA5116892	VAMP7 Polyclonal Antibody
PA5116859	WBP2NL Polyclonal Antibody
PA5117125	Viperin Polyclonal Antibody
PA5116881	WIBG Polyclonal Antibody
PA5116880	WWC2 Polyclonal Antibody
PA5116783	ZBTB8A Polyclonal Antibody
PA5116685	ZCCHC10 Polyclonal Antibody
PA5116614	ZDHHC18 Polyclonal Antibody
PA5116793	ZDHHC2 Polyclonal Antibody
PA5116796	ZDHHC20 Polyclonal Antibody
PA5116794	ZDHHC4 Polyclonal Antibody
PA5116795	ZDHHC7 Polyclonal Antibody
PA5116588	ZMAT2 Polyclonal Antibody
PA5116870	ZMPSTE24 Polyclonal Antibody
PA5116935	ZNF207 Polyclonal Antibody
PA5116842	ZNF444 Polyclonal Antibody
PA5116621	ZNF451 Polyclonal Antibody
PA5116668	ZNF620 Polyclonal Antibody
PA5116857	ZNF645 Polyclonal Antibody
PA5116875	ZNF652 Polyclonal Antibody
PA5116712	ZNF75A Polyclonal Antibody
MA542672	Acetyl-p53 (Lys382) Recombinant Rabbit Monoclonal Antibody
PA5120020	Rb Polyclonal Antibody
PA5120022	PCAF Polyclonal Antibody
PA5120025	GDF15 Polyclonal Antibody
PA5120027	ATG5 Polyclonal Antibody
PA5120028	BAK Polyclonal Antibody
PA5120030	BCL-XL Polyclonal Antibody
PA5120035	CRP Polyclonal Antibody
PA5120037	IGF1R (CD221) Polyclonal Antibody
PA5120038	MEK1 Polyclonal Antibody
PA5120040	PPAR gamma Polyclonal Antibody
PA5120042	CDK2 Polyclonal Antibody
PA5120045	p73 Polyclonal Antibody
PA5120049	NRG1 Polyclonal Antibody
PA5120052	C/EBP beta Polyclonal Antibody
PA5120057	IKK gamma Polyclonal Antibody
PA5120063	MSH6 Polyclonal Antibody
PA5120064	YAP1 Polyclonal Antibody
PA5120073	Cytokeratin 18 Polyclonal Antibody
PA5120077	PLCB3 Polyclonal Antibody
PA5120086	MZF1 Polyclonal Antibody
PA5120090	MMP16 Polyclonal Antibody
PA5120093	SEC11A Polyclonal Antibody
PA5120095	SLUG Polyclonal Antibody
PA5120096	MRPS12 Polyclonal Antibody
PA5120097	MST3 Polyclonal Antibody
PA5120098	GCH1 Polyclonal Antibody
PA5120103	Cullin 3 Polyclonal Antibody
PA5120104	MELK Polyclonal Antibody
PA5120111	NFkB p50/p105 Polyclonal Antibody
PA5120112	SQSTM1 Polyclonal Antibody
PA5120117	IRF3 Polyclonal Antibody

PA5120119	Vimentin Polyclonal Antibody
PA5120129	EPLIN Polyclonal Antibody
PA5120123	JAK2 Polyclonal Antibody
PA5120127	LSD1 Polyclonal Antibody
PA5120128	PAK4 Polyclonal Antibody
PA5120130	OSTM1 Polyclonal Antibody
PA5120133	NQO2 Polyclonal Antibody
PA5120134	P-Glycoprotein Polyclonal Antibody
PA5120135	Catalase Polyclonal Antibody
PA5120137	c-Rel Polyclonal Antibody
PA5120139	BET1L Polyclonal Antibody
PA5120140	JunD Polyclonal Antibody
PA5120141	Calpain 3 Polyclonal Antibody
PA5120145	SHMT2 Polyclonal Antibody
PA5120146	PFKP Polyclonal Antibody
PA5120147	PMEPA1 Polyclonal Antibody
PA5120148	BUD31 Polyclonal Antibody
PA5120150	NOX1 Polyclonal Antibody
PA5120152	BST-2 Polyclonal Antibody
PA5120153	PLSCR3 Polyclonal Antibody
PA5120154	Annexin A3 Polyclonal Antibody
PA5120155	Prothrombin Polyclonal Antibody
PA5120156	MEF2C Polyclonal Antibody
PA5120157	JAM-A (CD321) Polyclonal Antibody
PA5120184	FAM46A Polyclonal Antibody
PA5120183	VEGFB Polyclonal Antibody
PA5120185	MFN2 Polyclonal Antibody
PA5120186	MRPL42 Polyclonal Antibody
PA5120187	HOXA7 Polyclonal Antibody
PA5120188	S1P1 Polyclonal Antibody
PA5120189	HSP70 Polyclonal Antibody
PA5120190	GPX3 Polyclonal Antibody
PA5120192	AKR1C1 Polyclonal Antibody
PA5120194	RBM26 Polyclonal Antibody
PA5120196	U2AF1 Polyclonal Antibody
PA5120198	CKAP1 Polyclonal Antibody
PA5120203	Plasma Kallikrein Polyclonal Antibody
PA5120209	SLUG Polyclonal Antibody
PA5120212	TLE1 Polyclonal Antibody
PA5120213	TEF1 Polyclonal Antibody
PA5120228	EGF Polyclonal Antibody
PA5120238	GNA12 Polyclonal Antibody
PA5120239	ATP1A2 Polyclonal Antibody
PA5120240	CSAD Polyclonal Antibody
PA5120241	EZH2 Polyclonal Antibody
PA5120242	Aquaporin 4 Polyclonal Antibody
PA5120243	RhoA Polyclonal Antibody
PA5120244	AP2M1 Polyclonal Antibody
PA5120257	VASP Polyclonal Antibody
PA5120261	CSMD2 Polyclonal Antibody
PA5120262	CD299 (DC-SIGN/L) Polyclonal Antibody
PA5120264	NDUFA4L2 Polyclonal Antibody
PA5120265	ILKAP Polyclonal Antibody
PA5120267	ARFRP1 Polyclonal Antibody
PA5120268	LAIR2 Polyclonal Antibody
PA5120269	NDUFB3 Polyclonal Antibody
PA5120270	ACOT2 Polyclonal Antibody
PA5120272	TXNDC12 Polyclonal Antibody
PA5120271	p38 MAPK alpha Polyclonal Antibody

PA5120273	ALDH1L2 Polyclonal Antibody
PA5120274	Carbonic Anhydrase XIII Polyclonal Antibody
PA5120275	ETS2 Polyclonal Antibody
PA5120276	NGDN Polyclonal Antibody
PA5120277	RAP80 Polyclonal Antibody
PA5120278	Folliculin Polyclonal Antibody
PA5120280	PDP1 Polyclonal Antibody
PA5120283	MED18 Polyclonal Antibody
PA5120286	SP1 Polyclonal Antibody
PA5120288	CDH9 Polyclonal Antibody
PA5120289	EphB6 Polyclonal Antibody
PA5120290	GNAZ Polyclonal Antibody
PA5120291	ISG20 Polyclonal Antibody
PA5120292	PCDH9 Polyclonal Antibody
PA5120293	CD203c Polyclonal Antibody
PA5120295	Ribulose-phosphate 3-epimerase Polyclonal Antibody
PA5120297	SRD5A1 Polyclonal Antibody
PA5120299	FSHR Polyclonal Antibody
PA5120301	KLRC4 Polyclonal Antibody
PA5120302	UNC5C Polyclonal Antibody
PA5120303	WISP3 Polyclonal Antibody
PA5120304	FAM127A Polyclonal Antibody
PA5120305	SLC4A8 Polyclonal Antibody
PA5120306	BAG3 Polyclonal Antibody
PA5120307	TLK1 Polyclonal Antibody
PA5120308	COX4NB Polyclonal Antibody
PA5120310	GABARAPL2 Polyclonal Antibody
PA5120311	C17orf81 Polyclonal Antibody
PA5120312	BRMS1 Polyclonal Antibody
PA5120313	Cytohesin 4 Polyclonal Antibody
PA5120314	PPIL1 Polyclonal Antibody
PA5120315	GPR173 Polyclonal Antibody
PA5120316	Fggy Polyclonal Antibody
PA5120319	UNC5B Polyclonal Antibody
PA5120320	TCF4 Polyclonal Antibody
PA5120322	REC8 Polyclonal Antibody
PA5120323	Adenosine Kinase Polyclonal Antibody
PA5120325	Cortactin Polyclonal Antibody
PA5120327	RABIF Polyclonal Antibody
PA5120331	KSR1 Polyclonal Antibody
PA5120332	RRAGA Polyclonal Antibody
PA5120333	ABCF3 Polyclonal Antibody
PA5120334	PRR5 Polyclonal Antibody
PA5120335	JDP2 Polyclonal Antibody
PA5120336	RED Polyclonal Antibody
PA5120337	Kir4.2 (KCNJ15) Polyclonal Antibody
PA5120338	Menin Polyclonal Antibody
PA5120339	MYF6 Polyclonal Antibody
PA5120341	Oxytocin Polyclonal Antibody
PA5120342	PSMD11 Polyclonal Antibody
PA5120343	RNASE6 Polyclonal Antibody
PA5120344	ROBO1 Polyclonal Antibody
PA5120345	RPS25 Polyclonal Antibody
PA5120346	S100A8 Polyclonal Antibody
PA5120347	TCTA Polyclonal Antibody
PA5120348	Tenascin R Polyclonal Antibody
PA5120349	PTP4A2 Polyclonal Antibody
PA5120350	MTMR4 Polyclonal Antibody
PA5120351	MPHOSPH1 Polyclonal Antibody

PA5120352	Separase Polyclonal Antibody
PA5120353	FIG4 Polyclonal Antibody
PA5120354	FARSB Polyclonal Antibody
PA5120375	FLYWCH2 Polyclonal Antibody
PA5120355	KIF20A Polyclonal Antibody
PA5120356	IQGAP2 Polyclonal Antibody
PA5120357	WDHD1 Polyclonal Antibody
PA5120358	TNFR1 Polyclonal Antibody
PA5120359	STK38L Polyclonal Antibody
PA5120360	TBL2 Polyclonal Antibody
PA5120361	SULT1C4 Polyclonal Antibody
PA5120362	VILL Polyclonal Antibody
PA5120364	RASD1 Polyclonal Antibody
PA5120365	NUP54 Polyclonal Antibody
PA5120366	PLEKHB2 Polyclonal Antibody
PA5120368	DDX24 Polyclonal Antibody
PA5120369	ARID1B Polyclonal Antibody
PA5120370	MRPS14 Polyclonal Antibody
PA5120371	ZNF541 Polyclonal Antibody
PA5120372	TPD52L3 Polyclonal Antibody
PA5120373	DNER Polyclonal Antibody
PA5120374	GPRIN1 Polyclonal Antibody
PA5120376	SLC25A26 Polyclonal Antibody
PA5120377	GPR125 Polyclonal Antibody
PA5120380	IMPDH2 Polyclonal Antibody
PA5120383	Cathelicidin Polyclonal Antibody
PA5120384	ETV4 Polyclonal Antibody
PA5120387	RPL34 Polyclonal Antibody
PA5120388	RPS13 Polyclonal Antibody
PA5120389	Sorbitol Dehydrogenase Polyclonal Antibody
PA5120391	A2M Polyclonal Antibody
PA5120394	TBRG4 Polyclonal Antibody
PA5120396	ME3 Polyclonal Antibody
PA5120398	PLLIP Polyclonal Antibody
PA5120399	USP40 Polyclonal Antibody
PA5120400	ATP10A Polyclonal Antibody
PA5120401	ELSPBP1 Polyclonal Antibody
PA5120402	METT10D Polyclonal Antibody
PA5120403	FER1L5 Polyclonal Antibody
PA5120406	TARSL2 Polyclonal Antibody
PA5120407	GPR139 Polyclonal Antibody
PA5120408	LSM12 Polyclonal Antibody
PA5120409	TPRG1L Polyclonal Antibody
PA5120410	EDARADD Polyclonal Antibody
PA5120412	UBAC2 Polyclonal Antibody
PA5120413	FREM2 Polyclonal Antibody
PA5120414	BOLA3 Polyclonal Antibody
PA5120416	Acetyl-HMGB1 (Lys29) Polyclonal Antibody
PA5120419	GNL1 Polyclonal Antibody
PA5120420	IGFBP3 Polyclonal Antibody
PA5120421	VACHT Polyclonal Antibody
PA5120422	Protein Z Polyclonal Antibody
PA5120423	DENND4A Polyclonal Antibody
PA5120424	SYCP2 Polyclonal Antibody
PA5120425	EML4 Polyclonal Antibody
PA5120426	TRPV6 Polyclonal Antibody
PA5120428	TRAPPC6A Polyclonal Antibody
PA5120429	SGOL1 Polyclonal Antibody
PA5120430	GTPBP8 Polyclonal Antibody

PA5120433	LATS2 Polyclonal Antibody
PA5120434	ACSL1 Polyclonal Antibody
PA5120435	ALDH3A2 Polyclonal Antibody
PA5120438	Opsin 4 Polyclonal Antibody
PA5120440	ARL3 Polyclonal Antibody
PA5120441	BCAT1 Polyclonal Antibody
PA5120447	IRF5 Polyclonal Antibody
PA5120449	Prostasin Polyclonal Antibody
PA5120450	SPT6 Polyclonal Antibody
PA5120451	ELL Polyclonal Antibody
PA5120452	DGKZ Polyclonal Antibody
PA5120454	HDAC3 Polyclonal Antibody
PA5120458	CAND2 Polyclonal Antibody
PA5120459	RAB38 Polyclonal Antibody
PA5120460	P2Y10 Polyclonal Antibody
PA5120461	CTNNA3 Polyclonal Antibody
PA5120462	CARD14 Polyclonal Antibody
PA5120464	SPIRE2 Polyclonal Antibody
PA5120465	ARHGAP11B Polyclonal Antibody
PA5120466	NRN1L Polyclonal Antibody
PA5120468	Neuropeptide B Polyclonal Antibody
PA5120469	YY2 Polyclonal Antibody
PA5120470	IGFBPL1 Polyclonal Antibody
PA5120471	PDPK1 Polyclonal Antibody
PA5120472	AMPK alpha-1 Polyclonal Antibody
PA5120473	UGT2B15 Polyclonal Antibody
PA5120474	CD147 Polyclonal Antibody
PA5120476	Lamin B1 Polyclonal Antibody
PA5120477	NKIAMRE Polyclonal Antibody
PA5120478	L-Ficolin Polyclonal Antibody
PA5120480	E2F1 Polyclonal Antibody
PA5120483	PEDF Polyclonal Antibody
PA5120484	ABCA2 Polyclonal Antibody
PA5120487	CaV1.3 Polyclonal Antibody
PA5120488	Caspase 3 Polyclonal Antibody
PA5120489	CHRM1 Polyclonal Antibody
PA5120491	SAP102 Polyclonal Antibody
PA5120494	EZH2 Polyclonal Antibody
PA5120495	FGF3 Polyclonal Antibody
PA5120530	CNNM1 Polyclonal Antibody
PA5120499	GNA15 Polyclonal Antibody
PA5120500	HOXD8 Polyclonal Antibody
PA5120501	TOMM20 Polyclonal Antibody
PA5120502	INSR Polyclonal Antibody
PA5120503	Lamin B1 Polyclonal Antibody
PA5120504	MAD2 Polyclonal Antibody
PA5120505	MEKK4 Polyclonal Antibody
PA5120506	MYH9 Polyclonal Antibody
PA5120509	PAK3 Polyclonal Antibody
PA5120510	PCYT1A Polyclonal Antibody
PA5120511	PEDF Polyclonal Antibody
PA5120514	AMPK alpha-1 Polyclonal Antibody
PA5120515	TGIF Polyclonal Antibody
PA5120517	UCP3 Polyclonal Antibody
PA5120518	WIPF1 Polyclonal Antibody
PA5120519	BRPF1 Polyclonal Antibody
PA5120520	H3K4ac Polyclonal Antibody
PA5120521	Beclin 1 Polyclonal Antibody
PA5120524	eIF2 beta Polyclonal Antibody

PA5120522	B4GALT4 Polyclonal Antibody
PA5120523	SOCS2 Polyclonal Antibody
PA5120525	RAB33A Polyclonal Antibody
PA5120526	LRRRC14 Polyclonal Antibody
PA5120528	YAP1 Polyclonal Antibody
PA5120529	MMP12 Polyclonal Antibody
PA5120531	MRPS7 Polyclonal Antibody
PA5120532	TXNDC3 Polyclonal Antibody
PA5120533	HSD17B7 Polyclonal Antibody
PA5120534	CYB5R1 Polyclonal Antibody
PA5120535	EDEM2 Polyclonal Antibody
PA5120536	MIXL1 Polyclonal Antibody
PA5120537	TRIM5 beta Polyclonal Antibody
PA5120538	ZCCHC3 Polyclonal Antibody
PA5120539	SYT12 Polyclonal Antibody
PA5120540	FUK Polyclonal Antibody
PA5120541	ZNF530 Polyclonal Antibody
PA5120561	NMUR1 Polyclonal Antibody
PA5120542	SIRT1 Polyclonal Antibody
PA5120543	TET3 Polyclonal Antibody
PA5120544	PGA3 Polyclonal Antibody
PA5120545	LIM1 Polyclonal Antibody
PA5120546	HMGB2 Polyclonal Antibody
PA5120548	PIK3R2 Polyclonal Antibody
PA5120549	TIMM44 Polyclonal Antibody
PA5120552	DNMT1 Polyclonal Antibody
PA5120553	GSTA3 Polyclonal Antibody
PA5120554	ITPKA Polyclonal Antibody
PA5120555	LAD1 Polyclonal Antibody
PA5120556	PIK3R2 Polyclonal Antibody
PA5120557	ZNF143 Polyclonal Antibody
PA5120565	BAZ2A Polyclonal Antibody
PA5120566	RPL35 Polyclonal Antibody
PA5120558	H2BK5me1 Polyclonal Antibody
PA5120559	B4GALT3 Polyclonal Antibody
PA5120560	AP3S2 Polyclonal Antibody
PA5120563	TMED2 Polyclonal Antibody
PA5120564	ZPBP Polyclonal Antibody
PA5120567	Carbonic Anhydrase VB Polyclonal Antibody
PA5120568	IL-16 Polyclonal Antibody
PA5120569	CLCC1 Polyclonal Antibody
PA5120571	DECR2 Polyclonal Antibody
PA5120572	FAM127B Polyclonal Antibody
PA5120573	C13orf15 Polyclonal Antibody
PA5120574	PPHLN1 Polyclonal Antibody
PA5120576	ELOVL2 Polyclonal Antibody
PA5120577	PCDHB4 Polyclonal Antibody
PA5120578	GNPNAT1 Polyclonal Antibody
PA5120582	POF1B Polyclonal Antibody
PA5120580	MRPL34 Polyclonal Antibody
PA5120581	NOL10 Polyclonal Antibody
PA5120583	ARFGAP2 Polyclonal Antibody
PA5120584	NEK9 Polyclonal Antibody
PA5120585	PIGM Polyclonal Antibody
PA5120586	SFXN2 Polyclonal Antibody
PA5120588	ATP8B3 Polyclonal Antibody
PA5120589	APLF Polyclonal Antibody
PA5120590	VWC2L Polyclonal Antibody
PA5120591	GPR153 Polyclonal Antibody

PA5120592	CaMKII alpha Polyclonal Antibody
PA5120594	ASB2 Polyclonal Antibody
PA5120596	OSR1 Polyclonal Antibody
PA5120597	H2AK119ub Polyclonal Antibody
PA5120598	PEBP4 Polyclonal Antibody
PA5120599	MT-ND1 Polyclonal Antibody
PA5120601	MT-ND3 Polyclonal Antibody
PA5120602	GADD45A Polyclonal Antibody
PA5120603	SENP2 Polyclonal Antibody
PA5120604	MCL-1 Polyclonal Antibody
PA5120605	NNMT Polyclonal Antibody
PA5120609	CDK5 Polyclonal Antibody
PA5120610	Ub-K6 Polyclonal Antibody
PA5120611	TEX101 Polyclonal Antibody
PA5120612	Fyn Polyclonal Antibody
PA5120614	TBCCD1 Polyclonal Antibody
PA5120615	PHF6 Polyclonal Antibody
PA5120617	Ub-K63 Polyclonal Antibody
PA5120618	Paxillin Polyclonal Antibody
PA5120619	ZFPM1 Polyclonal Antibody
PA5120621	Ub-K11 Polyclonal Antibody
PA5120622	Ub-K29 Polyclonal Antibody
PA5120623	Ub-K33 Polyclonal Antibody
PA5120624	Ub-M1 Polyclonal Antibody
PA5120626	ZHX2 Polyclonal Antibody
PA5120627	MRPL55 Polyclonal Antibody
PA5120628	ATM Polyclonal Antibody
PA5120629	p130Cas Polyclonal Antibody
PA5120630	JNK1/JNK2 Polyclonal Antibody
PA5120632	GPR22 Polyclonal Antibody
PA5120633	TET3 Polyclonal Antibody
PA5120634	SLC5A3 Polyclonal Antibody
PA5120635	FZD8 Polyclonal Antibody
PA5120636	GTF3C3 Polyclonal Antibody
PA5120637	MNK1 Polyclonal Antibody
PA5120638	PLD3 Polyclonal Antibody
PA5120639	PRMT8 Polyclonal Antibody
PA5120640	HSPB7 Polyclonal Antibody
PA5120641	REPIN1 Polyclonal Antibody
PA5120642	IER5 Polyclonal Antibody
PA5120643	ZNF302 Polyclonal Antibody
PA5120644	SERC1 Polyclonal Antibody
PA5120647	Histone H1T Polyclonal Antibody
PA5120649	PRELP Polyclonal Antibody
PA5120650	DJ-1 Polyclonal Antibody
PA5120651	RAB3GAP1 Polyclonal Antibody
PA5120652	IFNAR1 Polyclonal Antibody
PA5120653	CXCL9 (MIG) Polyclonal Antibody
PA5120654	Caspase 1 Polyclonal Antibody
PA5120655	ARID1A Polyclonal Antibody
PA5120657	SMAD2/SMAD3 Polyclonal Antibody
PA5120658	GPR120 Polyclonal Antibody
PA5120659	PKM2 Polyclonal Antibody
PA5120660	PKM1 Polyclonal Antibody
PA5120661	TRRAP Polyclonal Antibody
PA5120662	EPHX2 Polyclonal Antibody
PA5120664	SNRNP40 Polyclonal Antibody
PA5120665	SHH Polyclonal Antibody
PA5120666	hnRNP C Polyclonal Antibody

PA5120667	CNTF Polyclonal Antibody
PA5120669	RhoG Polyclonal Antibody
PA5120670	Caspase 4 Polyclonal Antibody
PA5120671	FANCD2 Polyclonal Antibody
PA5120672	Connexin 50 Polyclonal Antibody
PA5120673	MAN2A2 Polyclonal Antibody
PA5120683	CDw293 Polyclonal Antibody
PA5120684	IFI16 Polyclonal Antibody
PA5120685	ASS1 Polyclonal Antibody
PA5120696	SARS-CoV-2 NSP8 Polyclonal Antibody
PA5120687	TRAF2 Polyclonal Antibody
PA5120688	RABAC1 Polyclonal Antibody
PA5120689	LPAR4 Polyclonal Antibody
PA5120690	HOXB2 Polyclonal Antibody
PA5120691	HRH1 Polyclonal Antibody
PA5120692	WRAP53 Polyclonal Antibody
PA5120693	USP11 Polyclonal Antibody
PA5120695	MAK16 Polyclonal Antibody
PA5120697	LLGL1 Polyclonal Antibody
PA5120698	SARS-CoV-2 NSP2 Polyclonal Antibody
PA5120699	SARS-CoV-2 NSP4 Polyclonal Antibody
PA5120700	SARS-CoV-2 NSP15 Polyclonal Antibody
PA5120701	SARS-CoV-2 NSP16 Polyclonal Antibody
PA5120703	RSK1/RSK2/RSK3 Polyclonal Antibody
PA5120704	RAB11A/RAB11B Polyclonal Antibody
PA5120705	TMEM119 Polyclonal Antibody
PA5120706	Asymmetric Di-Methyl Arginine Polyclonal Antibody
PA5120707	Ubiquitin Remnant Motif (K-Epsilon-GG) Polyclonal Antibody
PA5120708	SARS-CoV-2 ORF7a Polyclonal Antibody
PA5120709	SARS-CoV-2 NSP9 Polyclonal Antibody
PA5120710	ALPK1 Polyclonal Antibody
PA5120712	Aurora A/B/C Polyclonal Antibody
PA5120713	GPR81 Polyclonal Antibody
PA5120716	SARS-CoV-2 NSP10 Polyclonal Antibody
PA5120726	Proteinase K Polyclonal Antibody
PA5120727	SSR2 Polyclonal Antibody
PA5120729	KLF12 Polyclonal Antibody
PA5120731	NEURL1B Polyclonal Antibody
PA5120735	TM4SF19 Polyclonal Antibody
PA5120737	DTX1 Polyclonal Antibody
PA5120738	MED26 Polyclonal Antibody
PA5120745	MYOM2 Polyclonal Antibody
PA5120746	CD300b Polyclonal Antibody
PA5120772	IKK alpha Polyclonal Antibody
PA5120773	Enterovirus 71 3C Polyclonal Antibody
PA5120774	Enterovirus 71 3D Polyclonal Antibody
PA5120776	FANCD2 Polyclonal Antibody
PA5120782	IGF2 Polyclonal Antibody
PA5120788	POLR2A Polyclonal Antibody
PA5120791	Zyxin Polyclonal Antibody
PA5120796	Stathmin 1 Polyclonal Antibody
PA5120799	PLCG2 Polyclonal Antibody
PA5120804	MTA2 Polyclonal Antibody
PA5120808	Histone H3 Polyclonal Antibody
PA5120817	SERPINH1 Polyclonal Antibody
PA5120823	Neuropilin 2 Polyclonal Antibody
PA5120825	APLP1 Polyclonal Antibody
PA5120827	p21 Polyclonal Antibody
PA5120833	VCP Polyclonal Antibody

PA5120839	CDK1 Polyclonal Antibody
PA5120841	DAP5 Polyclonal Antibody
PA5120842	DUT Polyclonal Antibody
PA5120843	mGluR8 Polyclonal Antibody
PA5120845	MCT2 Polyclonal Antibody
PA5120847	MMP10 Polyclonal Antibody
PA5120849	MYL9 Polyclonal Antibody
PA5120850	AKAP8 Polyclonal Antibody
PA5120854	Nkx2.1 Polyclonal Antibody
PA5120855	NOTCH3 Polyclonal Antibody
PA5120857	SLC22A8 Polyclonal Antibody
PA5120858	MKP3 Polyclonal Antibody
PA5120859	NPY Polyclonal Antibody
PA5120860	MNAR Polyclonal Antibody
PA5120863	TWIST1 Polyclonal Antibody
PA5120867	SLITRK6 Polyclonal Antibody
PA5120870	CD39 Polyclonal Antibody
PA5120874	GPR17 Polyclonal Antibody
PA5120877	DRG1 Polyclonal Antibody
PA5120884	UCP2 Polyclonal Antibody
PA5120886	UBE2L6 Polyclonal Antibody
PA5120895	CDK5RAP3 Polyclonal Antibody
PA5120896	CAB39L Polyclonal Antibody
PA5120897	USP48 Polyclonal Antibody
PA5120900	H3K14me1 Polyclonal Antibody
PA5120904	LUM Polyclonal Antibody
PA5120908	Fra1 Polyclonal Antibody
PA5120918	CDC34 Polyclonal Antibody
PA5120920	DNMT1 Polyclonal Antibody
PA5120921	Cathepsin A Polyclonal Antibody
PA5120927	NFkB p52/p100 Polyclonal Antibody
PA5120928	TK1 Polyclonal Antibody
PA5120930	CREM Polyclonal Antibody
PA5120933	DDX3 Polyclonal Antibody
PA5120935	LAT Polyclonal Antibody
PA5120951	GOT1 Polyclonal Antibody
PA5120953	JMJD2B Polyclonal Antibody
PA5120964	hnRNP H1 Polyclonal Antibody
PA5120972	FGF23 Polyclonal Antibody
PA5120978	PRAS40 Polyclonal Antibody
PA5120990	FAP Polyclonal Antibody
PA5120992	ASL Polyclonal Antibody
PA5120997	Glycerol kinase Polyclonal Antibody
PA5121005	RhoD Polyclonal Antibody
PA5121008	BAP1 Polyclonal Antibody
PA5121011	LZIP Polyclonal Antibody
PA5121021	Monoacylglycerol Lipase Polyclonal Antibody
PA5121025	PIK3CG Polyclonal Antibody
PA5121046	hnRNP M Polyclonal Antibody
PA5121048	MEK5 Polyclonal Antibody
PA5121071	ELF5 Polyclonal Antibody
PA5121072	FOXN2 Polyclonal Antibody
PA5121082	H3K14ac Polyclonal Antibody
PA5121083	H3K56ac Polyclonal Antibody
PA5121086	NAT10 Polyclonal Antibody
PA5121101	RAR gamma-1,2 Polyclonal Antibody
PA5121113	ATG14 Polyclonal Antibody
PA5121121	GluD1 Polyclonal Antibody
PA5121128	MEK3 Polyclonal Antibody

PA5121130	SRF Polyclonal Antibody
PA5121134	AP2 gamma Polyclonal Antibody
PA5121180	BAIAP2L1 Polyclonal Antibody
PA5121183	RSPO1 Polyclonal Antibody
PA5121185	MAP1LC3C Polyclonal Antibody
PA5121189	EphA4 Polyclonal Antibody
PA5121199	OPRD1 Polyclonal Antibody
PA5121216	Human IgM Polyclonal Antibody
PA5121217	TAX1BP3 Polyclonal Antibody
PA5121219	FCAR Polyclonal Antibody
PA5121220	APOBEC3B Polyclonal Antibody
PA5121222	EPB41L2 Polyclonal Antibody
PA5121223	SMU1 Polyclonal Antibody
PA5121225	ATP5L Polyclonal Antibody
PA5121227	RAB11FIP1 Polyclonal Antibody
PA5121230	EGFL7 Polyclonal Antibody
PA5121231	PYGM Polyclonal Antibody
PA5121232	Desmoglein 2 (DSG2) Polyclonal Antibody
PA5121234	PHLPP1 Polyclonal Antibody
PA5121235	MEPE Polyclonal Antibody
PA5121253	FRMD6 Polyclonal Antibody
PA5121255	Phospho-AKT2 (Ser474) Polyclonal Antibody
PA5121256	Phospho-ATM (Ser1981) Polyclonal Antibody
PA5121259	Phospho-STAT3 (Tyr705) Polyclonal Antibody
PA5121267	Phospho-Cofilin (Ser3) Polyclonal Antibody
PA5121265	Phospho-Tau (Thr205) Polyclonal Antibody
PA5121266	Phospho-Tau (Ser404) Polyclonal Antibody
PA5121269	Phospho-alpha Synuclein (Tyr136) Polyclonal Antibody
PA5121270	Phospho-LCK (Tyr505) Polyclonal Antibody
PA5121273	Phospho-NFkB p100 (Ser866) Polyclonal Antibody
PA5121274	Phospho-AMPK alpha-1,2 (Thr183, Thr172) Polyclonal Antibody
PA5121275	Phospho-NFkB p65 (Ser468) Polyclonal Antibody
PA5121278	Phospho-PDPK1 (Ser241) Polyclonal Antibody
PA5121277	Phospho-STAT3 (Ser727) Polyclonal Antibody
PA5121281	Phospho-Syk (Tyr323) Polyclonal Antibody
PA5121282	Phospho-IKK alpha (Ser176) Polyclonal Antibody
PA5121285	Phospho-TrkA/TrkB (Tyr674, Tyr675, Tyr706, Tyr707) Polyclonal Antibody
PA5121286	Phospho-IRS1 (Ser636, Ser639) Polyclonal Antibody
PA5121290	Phospho-Chk1 (Ser345) Polyclonal Antibody
PA5121293	Phospho-CSF1R (Tyr708) Polyclonal Antibody
PA5121296	Phospho-ULK1 (Ser757) Polyclonal Antibody
PA5121297	Phospho-Caveolin 1 (Tyr14) Polyclonal Antibody
PA5121302	Phospho-NFkB p65 (Thr434) Polyclonal Antibody
PA5121298	Phospho-Tau (Ser202) Polyclonal Antibody
PA5121299	Phospho-NMDAR2B (Tyr1474) Polyclonal Antibody
PA5121300	Phospho-c-Cbl (Tyr700) Polyclonal Antibody
PA5121301	Phospho-NPM1 (Thr199) Polyclonal Antibody
PA5121303	Phospho-EEF2K (Ser366) Polyclonal Antibody
PA5121305	Phospho-DNM1L (Ser616) Polyclonal Antibody
PA5121304	Phospho-Histone H3 (Thr3) Polyclonal Antibody
PA5121308	Phospho-MYL9 (Ser19) Polyclonal Antibody
PA5121309	Phospho-PEA15 (Ser104) Polyclonal Antibody
PA5121310	Phospho-GYS1 (Ser641) Polyclonal Antibody
PA5121312	Phospho-STAT5A/B (Tyr694) Polyclonal Antibody
PA5121313	Phospho-IKK alpha/beta (Ser176, Ser180) Polyclonal Antibody
PA5121315	Phospho-Histone H4 (Ser1) Polyclonal Antibody
PA5121316	Phosphotyrosine Polyclonal Antibody
PA5121317	Phospho-IRE1 alpha (Ser551) Polyclonal Antibody
PA5121318	Phospho-WWTR1 (Ser89) Polyclonal Antibody

PA5121319	Phospho-SMO (Ser616) Polyclonal Antibody
PA5121320	Phospho-NFkB p65 (Ser529) Polyclonal Antibody
PA5121321	Phospho-Phospholamban (Ser16, Thr17) Polyclonal Antibody
PA5121322	Phospho-NMDAR1 (Ser896, Ser897) Polyclonal Antibody
PA5121323	Phospho-NMDAR2B (Tyr1070) Polyclonal Antibody
PA5121324	Phospho-POLR2A (Ser7) Polyclonal Antibody
PA5121325	Phospho-PKC delta (Ser645) Polyclonal Antibody
PA5121326	Phospho-NFkB p105 (Ser932) Polyclonal Antibody
PA5121327	Phospho-IkB alpha (Ser36) Polyclonal Antibody
PA5121328	Phospho-MCK10 (Tyr792) Polyclonal Antibody
PA5121329	Phospho-E-cadherin (Ser838, Ser840) Polyclonal Antibody
PA5121330	Phospho-CSF1R (Tyr723) Polyclonal Antibody
PA5121331	Phospho-PSD93 (Tyr340) Polyclonal Antibody
PA5121332	Phospho-IR/IGF1R (Tyr1165, Tyr1166, Tyr1189, Tyr1190) Polyclonal Antibody
PA5121333	Phospho-FOXO1 (Ser319) Polyclonal Antibody
PA5121334	Phospho-alpha Synuclein (Tyr133) Polyclonal Antibody
PA5121335	Phospho-MST1/MST2 (Thr183, Thr180) Polyclonal Antibody
PA5121336	Phospho-VEGF Receptor 2 (Tyr1175) Polyclonal Antibody
PA5121337	Phospho-FAK (Tyr925) Polyclonal Antibody
PA5121338	Phospho-RPA2 (Ser4, Ser8) Polyclonal Antibody
PA5121339	Phospho-p70 S6 Kinase (Thr421) Polyclonal Antibody
PA5121340	Phospho-p70 S6 Kinase (Thr421, Ser424) Polyclonal Antibody
PA5121341	Phospho-SHC (Tyr317) Polyclonal Antibody
PA5121343	Phospho-DAPP1 (Tyr139) Polyclonal Antibody
PA5121344	Phospho-PPP1R12A/PPP1R12B/PPP1R12C (Thr696) Polyclonal Antibody
PA5121348	Phospho-ASK1 (Ser967) Polyclonal Antibody
PA5121345	Phospho-AKT Pan (Thr308, Thr305, Thr309) Polyclonal Antibody
PA5121346	Phospho-LAT (Tyr191) Polyclonal Antibody
PA5121347	Phospho-MSK1 (Thr581) Polyclonal Antibody
PA5121349	Phospho-TAK1 (Thr187) Polyclonal Antibody
PA5121351	Phospho-RIP1 (Ser166) Polyclonal Antibody
PA5120266	TRMT112 Polyclonal Antibody
PA5109739	FSTL4 Polyclonal Antibody
PA5120562	DRAP1 Polyclonal Antibody
PA5120734	METTL8 Polyclonal Antibody
PA5116805	OSR2 Polyclonal Antibody
PA5116813	Symmetric Di-Methyl Arginine Polyclonal Antibody
MA542647	NMDAR1 ANTIBODY
MA542935	PTH ANTIBODY
PA5109955	MFGE8 ANTIBODY
PA5110594	ZNF10 ANTIBODY
PA5116901	LACTIC ACID LYSINE ANTIBODY
PA5116911	FILAGGRIN ANTIBODY
PA5117264	Phospho-YAP1 (Ser128) Polyclonal Antibody
PA5120119	Vimentin Polyclonal Antibody
PA587536	CASPASE 1 ANTIBODY
PA588084	Nrf2 Polyclonal Antibody
PA589059	FAM208A Polyclonal Antibody
PA589622	GLYAT ANTIBODY
PA595749	CXCR1 Polyclonal Antibody