



Certificate of Analysis

**STEMPRO® Adipocyte Differentiation
Basal Medium**
+ 4.5 g/L D-Glucose
+ 110 mg/L Sodium Pyruvate
+ GlutaMAX(TM) - I

Catalog Number: A10410
Lot Number: 2093171


Storage Temperature: 2 to 8C
Storage Instructions: Protect from light
Expiration Date: 2020-04

For research use only. Not for use in diagnostic procedures.

Sterile filtered (0.1 um)

TEST	SPECIFICATION	RESULT	UNITS
¹ Endotoxin Testing	<=0.05	<0.03	EU/mL
² Osmolality	>=290 to <=370	345	mOsm/kg
³ pH	7.0 to <=7.4	7.4	
⁴ Sterility Testing	Negative	Negative	

Read SDS



Quality Systems Department

Date: 24-Apr-2019

REFERENCES:

- 1 Current United States Pharmacopeia, <85> Bacterial Endotoxins Test.
- 2 Thermo Fisher Scientific Specifications.
- 3 Thermo Fisher Scientific Specifications.
- 4 Current Edition of USP, Thermo Fisher Scientific Modified.