

**LV-MAX (TM) Production Medium  
(+) Glutamax (TM)**

**Catalog Number:** A35834  
**Lot Number:** 2071065

**Storage Temperature:** 2 to 8C  
**Storage Instructions:** Protect from light  
**Expiration Date:** 2020-03

---


For research use only. Not for use in diagnostic procedures.

---

TEST	SPECIFICATION	RESULT	UNITS
<sup>1</sup> Endotoxin Testing	<=1.0	<0.10	EU/mL
<sup>2</sup> Osmolality	>=270 to <=295	282	mOsm/kg
<sup>3</sup> pH	>=6.9 to <=7.3	7.2	
<sup>4</sup> Sterility Testing	Negative	Negative	

Read SDS

EXAMPLE



Quality Systems Department      Date: 21-Mar-2019

**REFERENCES:**

- 1 Current United States Pharmacopeia, <85> Bacterial Endotoxins Test.
- 2 Thermo Fisher Scientific Specifications.
- 3 Thermo Fisher Scientific Specifications.
- 4 Current Edition of USP, Thermo Fisher Scientific Modified.