

Other Sera: Comparison of Tests

Test	Base Cat. Number						
	Newborn Calf Serum	Newborn Calf Serum Heat-Inactivated	Bovine Serum	Bovine Serum Heat-Inactivated	Donor Bovine Serum	Donor Bovine Serum with Iron	Horse Serum
	16010	26010	16170	26170	16030	10371	16050
Cytotoxicity Assay							
Electrophoretic Pattern	✓	✓	✓	✓	✓	✓	✓
Endotoxin (EU/ml)	✓	✓	✓	✓	✓	✓	✓
Hemoglobin (mg/dl)	✓	✓	✓	✓	✓	✓	✓
Mycoplasmas	✓	✓	✓	✓	✓	✓	✓
Osmolality	✓	✓	✓	✓	✓	✓	✓
pH	✓	✓	✓	✓	✓	✓	✓
Sp2 Performance Assay							✓
Sterility Testing (Bacteria/Fungi)	✓	✓	✓	✓	✓	✓	✓
Total Protein	✓	✓	✓	✓	✓	✓	✓
Vero Performance Assay	✓	✓	✓	✓	✓	✓	
VIRUS							
Blue Tongue Virus	✓	✓	✓	✓	✓	✓	
Bovine Adenovirus FA	✓	✓	✓	✓	✓	✓	
Bovine Parvovirus FA	✓	✓	✓	✓	✓	✓	
Bovine Respiratory Syncytial Virus FA	✓	✓	✓	✓	✓	✓	
Bovine Viral Diarrhea Virus FA	✓	✓	✓	✓	✓	✓	✓
Cytopathic Agents	✓	✓	✓	✓	✓	✓	✓
EIA Coggins							✓
Equine Herpes Virus							✓
Equine Viral Arteritis							✓
Hemadsorbing Agents	✓	✓	✓	✓	✓	✓	✓
Porcine Adenovirus							
Porcine Hemagglutinating Encephalitis							
Porcine Parvovirus							
Rabies FA	✓	✓	✓	✓	✓	✓	✓
Reovirus FA	✓	✓	✓	✓	✓	✓	✓
Transmissible Gastroenteritis Virus							

Other Sera: Comparison of Tests (cont.)

Test	Base Cat. Number					
	Horse Serum Heat-Inactivated 26050	Chicken Serum 16110	Goat Serum 16210	Lamb Serum 16070	Porcine Serum 26250	Rabbit Serum 16120
Cytotoxicity Assay		✓		✓	✓	
Electrophoretic Pattern	✓	✓	✓	✓	✓	✓
Endotoxin (EU/ml)	✓					
Hemoglobin (mg/dl)	✓	✓	✓	✓	✓	✓
Mycoplasmas	✓	✓	✓	✓	✓	✓
Osmolality	✓	✓		✓	✓	
pH	✓	✓		✓	✓	
Sp2 Performance Assay	✓					
Sterility Testing (Bacteria/Fungi)	✓	✓	✓	✓	✓	✓
Total Protein	✓	✓	✓	✓	✓	✓
Vero Performance Assay						
VIRUS						
Blue Tongue Virus						
Bovine Adenovirus FA						
Bovine Parvovirus FA						
Bovine Respiratory Syncytial Virus FA						
Bovine Viral Diarrhea Virus FA	✓				✓	
Cytopathic Agents	✓				✓	
EIA Coggins	✓					
Equine Herpes Virus	✓					
Equine Viral Arteritis	✓					
Hemadsorbing Agents	✓				✓	
Porcine Adenovirus					✓	
Porcine Hemagglutinating Encephalitis					✓	
Porcine Parvovirus					✓	
Rabies FA	✓				✓	
Reovirus FA	✓				✓	
Transmissible Gastroenteritis Virus					✓	