

Material No.: 3605-09
Revision No.: 0

Product Specification

Meets F.C.C. Requirements, Meets N.F. Requirements, GMP Manufactured Product

Test	Specification
NF – Identification A	Passes Test
NF – Identification B	Passes Test
NF – Water (H ₂ O)	≤ 0.5 %
NF – Assay (Na ₂ CO ₃) (anhydrous basis)	99.5 – 100.5 %
FCC – Identification A	Passes Test
FCC – Identification B	Passes Test
FCC – Loss on Drying	≤ 0.5 %
FCC – Assay (Na ₂ CO ₃) (anhydrous basis)	99.5 – 100.5 %
FCC – Lead (Pb)	≤ 4 mg/kg

GMP Manufactured Product
Food GMP Manufactured Product
Bulk Food Chemical
Bulk Pharmaceutical Chemical

CAUTION: For Manufacturing, processing or repackaging

No Class 1,2,3 or other solvents are used or produced in the manufacturing or purification of the product.

Packaging Site: Phillipsburg Mfg Ctr & DC



Jamie Ethier
Vice President Global Quality

For questions on this Product Specification please contact Technical Services at 855.282.6867 or +1.610.386.1700

Avantor Performance Materials, LLC

100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone 610.386.1700

Page 1 of 1