

Plastic outer support containers

Large-volume rigid plastic support containers for BioProcess Containers (BPCs)

Thermo Scientific™ outer support containers are rigid plastic vessels to support our Thermo Scientific™ BioProcess Container (BPC) products in helping meet your needs for large-volume liquid storage, handling, and transportation.

Key features

- **Reusable**—robust design and sturdy construction
- **Convenient**—fold-down walls for easy storage when not in use
- **Stackable**—both empty and full
- **Bottom access**—with ample space for a tubing assembly
- **Shippable**—qualified for shipping

Applications

- Feeding, storing, and transporting media and buffer
- Receiving and storing harvest from bioreactors
- Transporting large-volume liquids



Plastic outer support containers

Size	Outer dimensions (W x D x H)	Inner dimensions (W x D x H)	Dispensing locations, material of construction	Cat. No.
250 L	80.9 x 60.9 x 96.5 cm (31.8 x 23.9 x 38 in.)	73.9 x 54.1 x 66.0 cm (29.1 x 21.3 x 26.0 in.)	Top and bottom, polypropylene	SV50139.15
500 L	120.0 x 80.0 x 104.0 cm (47.2 x 31.5 x 41.0 in.)	113.0 x 71.9 x 70.1 cm (44.5 x 28.3 x 27.6 in.)	Top and bottom, ABS/HDPE	SV50139.12
600 L	115.5 x 115.5 x 98.0 cm (45.5 x 45.5 x 38.6 in.)	106.4 x 106.4 x 62.5 cm (41.9 x 41.9 x 24.6 in.)	Top and bottom, polypropylene	SV50139.11
1,000 L	115.5 x 115.5 x 129.4 cm (45.5 x 45.5 x 50.9 in.)	106.4 x 104.4 x 93.9 cm (41.9 x 41.1 x 37.0 in.)	Top and bottom, polypropylene	SV50139.06

Stainless steel carts

Support container Cat. No. (size)	Cat. No.
SV50139.15 (250 L)	SV50139.17
SV50139.12 (500 L)	SV50139.13
SV50139.11 (600 L)	SV50139.16
SV50139.06 (1,000 L)	

BPCs for rigid plastic support containers

Size	Support container Cat. No.	Cat. No.
200 L	SV50139.15	Aegis5-14 film SH30976.01
		CX5-14 film SH31060.01
250 L	SV50139.15	Aegis5-14 film SH30976.02
		CX5-14 film SH31060.02
500 L	SV50139.12	Aegis5-14 film SH30976.03
		CX5-14 film SH31060.03
600 L	SV50139.11	Aegis5-14 film SH30976.04
		CX5-14 film SH31060.04
1,000 L	SV50139.06	Aegis5-14 film SH30976.05
		CX5-14 film SH31060.05

Find out more at

thermofisher.com/sutstoragetransport