



# Help protect your organization and people in three simple steps with in-air pathogen surveillance

Thermo Scientific™ AerosolSense™ Sampler and Applied Biosystems™ TaqMan™ PCR Assays

01

## COLLECT SAMPLE

Position AerosolSense Sampler on a flat, unobstructed surface and insert sample cartridge.



02

## REMOVE & REPLACE SAMPLE CARTRIDGE

After sampling cycle remove sample cartridge and place in biohazard bag to send to lab.



RESULTS  
WITHIN 24 HOURS

03

## SEND SAMPLE CARTRIDGE TO LAB FOR TESTING

Verified kits available for testing:

- Applied Biosystems™ TaqPath™ COVID-19 Combo Kit
- Applied Biosystems™ TaqMan™ SARS-CoV-2, FluA, FluB RT-PCR Assay Kit (RUO)
- Applied Biosystems™ TaqMan™ SARS-CoV-2, Flu A/B, RSV RT-PCR Assay Kit (RUO)

Results:

- On-site lab: ≤ 4h result availability
- Thermo Fisher Scientific testing services: ≤ 24h result availability