

## Low and Ultra-Low Positive COVID-19 Controls

### Thermo Scientific AcroMetrix Coronavirus 2019 (COVID-19) RNA Control (RUO)


The Thermo Scientific™ AcroMetrix™ Coronavirus 2019 (COVID-19) RNA Control contains full length genomic RNA from SARS-COV-2 in proprietary buffer used to validate and monitor the performance of molecular tests for detecting COVID-19. This is a positive control kit containing two concentrations: a low positive and ultra-low positive concentration to avoid false-negative results. Both the controls are quantitated using Bio-Rad™ Droplet Digital™ PCR (ddPCR™) technology.

#### Product Features


- Product contains full length genomic RNA from SARS-COV-2
- Product stable up to 3 freeze/thaws
- Open vial stability at 2-8°C for 7 days
- Product storage -20°C or below

**AcroMetrix Coronavirus 2019 (COVID-19) RNA Control sequence design is compatible with, but not limited to, the following assays:**

- TaqPath™ COVID-19 Combo kit (A47814)
- USA CDC assay
- China CDC assay
- Japan National Institute of Infectious Diseases assay
- Thailand National Institute of Health assay
- WHO Assay from Charité Institute of Virology



For Research Use only



Catalog Number	Control Name	Copies/ $\mu$ L	Quantity
954519	AcroMetrix COVID-19 Low Positive Control	500	1 x 20 $\mu$ L
	AcroMetrix COVID-19 Ultra-Low Positive Control	100	1 x 20 $\mu$ L

#### Customer Support

**Phone:** 1-800-232-3342 - Press 1 for Customer Service, then 3 for QC

**Fax:** 1-800-829-8115

**e-mail:** [mgc-masorders@thermofisher.com](mailto:mgc-masorders@thermofisher.com)

#### Sales Support

**Phone:** 1-800-232-3342

**e-mail:** [sales.diagnostics.fmt@thermofisher.com](mailto:sales.diagnostics.fmt@thermofisher.com)

#### Technical Service and Support

**Phone:** 1-800-232-3342- option 2, then 2 for QC

#### Clinical Diagnostics

##### USA

46500 Kato Road  
Fremont, CA 94538 USA  
Tel. 800-232-3342 / + 1 510-979-5000  
[sales.diagnostics.fmt@thermofisher.com](mailto:sales.diagnostics.fmt@thermofisher.com)

##### Asia Pacific

Tel: +65 6499 9999  
[cdd.asia.info2@thermofisher.com](mailto:cdd.asia.info2@thermofisher.com)

##### Australia

Tel. + 61 1800 333 110  
[auiinfo@thermofisher.com](mailto:auiinfo@thermofisher.com)

##### Canada

Tel. 905-286-4290  
[CDD.Canada@thermofisher.com](mailto:CDD.Canada@thermofisher.com)

##### China

Tel. +86 400 650 5118  
[sales.china@thermofisher.com](mailto:sales.china@thermofisher.com)

##### France, Belgium (Brussels / Wallonia) & Luxembourg

Tel. + 33 (0) 1 40 86 65 00  
[cdd.fr.info@thermofisher.com](mailto:cdd.fr.info@thermofisher.com)

##### Germany

Tel. +49 (0) 3302 883 608  
[cdd.de.info@thermofisher.com](mailto:cdd.de.info@thermofisher.com)

##### Hong Kong

Tel. + 852 3107 7600  
[sales.hk@thermofisher.com](mailto:sales.hk@thermofisher.com)

##### Italy

[distributor.cdd@thermofisher.com](mailto:distributor.cdd@thermofisher.com)

##### Japan

Tel. +81(0)120 489 211  
[JPYOK-CDD.QC@thermofisher.com](mailto:JPYOK-CDD.QC@thermofisher.com)

##### Netherlands & Belgium (Flemish)

Tel. +49 3302 883 660  
[Info.benelux.cdd@thermofisher.com](mailto:Info.benelux.cdd@thermofisher.com)

##### New Zealand

Tel. + 61 0800933 966  
[nzinfo@thermofisher.co.nz](mailto:nzinfo@thermofisher.co.nz)

##### Nordic

Tel. + 46 8 594 691 30  
[info.nordic.cdd@thermofisher.com](mailto:info.nordic.cdd@thermofisher.com)

##### Russia

Tel. + 7 495 739 76 41  
[info.btd.moscow@thermofisher.com](mailto:info.btd.moscow@thermofisher.com)

##### South Africa

Tel. +27 11 792 6790  
[support-za.idd@thermofisher.com](mailto:support-za.idd@thermofisher.com)

##### Spain & Portugal (Iberia)

Tel. + 34 93589 8338  
[CDD.Info.ES@thermofisher.com](mailto:CDD.Info.ES@thermofisher.com)

##### Switzerland & Austria

Tel. + 41 26 663 86 70  
[cdd.ch.info@thermofisher.com](mailto:cdd.ch.info@thermofisher.com)

##### United Kingdom & Ireland

Tel. + 44 1442 868 940  
[cdd.UK.info@thermofisher.com](mailto:cdd.UK.info@thermofisher.com)

##### For Countries Not Listed

[distributor.cdd@thermofisher.com](mailto:distributor.cdd@thermofisher.com)

Learn more at: [thermofisher.com/covid-19-control](https://www.thermofisher.com/covid-19-control)

Availability of product in each country depends on local regulatory marketing authorization status.

© 2020 Thermo Fisher Scientific Inc. All rights reserved. Bio-Rad, Droplet Digital, and ddPCR are trademarks of Bio-Rad. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries.

**FR-MTL-10027916-3-AcroMetrix-COVID-19-Flyer-EN**

**ThermoFisher**  
SCIENTIFIC