



Software that Provides Flexibility

## Get more done with LabLink software for the Indiko Analyzer

LabLink QC Software streamlines your daily QC workflow by enabling instant updates of your QC testing data and lets you compare your data with other Indiko peers or in affiliated labs over the cloud.



### Automate

With LabLink xL™ software, you can focus on what matters instead of tedious data entry and manual calculations. Automatically transmit results, receive tailored reports in your inbox, get notified instantly for QC violations, and sign reports electronically to enhance your productivity.



### Simplify

LabLink xL software is simple to use yet extremely powerful. Interactive Levey-Jennings charts, electronic signature verification, and audit-friendly reports let you focus on QC and less on learning software.



### Accelerate

LabLink xL software streamlines your workflow by eliminating unnecessary steps from your daily routine. Use Westgard rules to setup QC flags and dive straight into results that require attention.

## Access your QC from anywhere

With Thermo Scientific™ LabLink xL cloud based quality assurance software, you can monitor quality control testing in real-time. The smart automation technology, incredible simplicity, and built-in productivity tools will save time and improve efficiency in your lab.

## Maximize productivity with Thermo Scientific MAS Quality Controls and LabLink xL software for Indiko analyzers

When it comes to your goal of increasing productivity in your core lab, look no further than Thermo Scientific™ MAS™ Omni•CORE™ and MAS chemTRAK™•H on the Indiko™ platforms. These third-party multianalyte controls permit truly independent validation and, paired with the flexibility offered by Indiko and Indiko Plus analyzers, they will help you maximize productivity like never before.

**#MASproductivity**

continues with Indiko analyzers



Minimize down-time, maximize throughput



Minimize waste, maximize **COST** savings



Minimize hands-on time, maximize efficiency



# Automate your workflow



Instrument



Middleware



LIS



## Hands-free QC with Auto-Connectivity Software

- No more manual data entry with fully automatic QC data transfer
- QC results are uploaded to LabLink xL software
- Save time with error-free data submission and increase confidence in results



### Automatic notifications for

- Data violations
- Connectivity issues
- Product expirations
- Product updates



### Time-saving email reports

- Automatically delivered to your inbox
- Insightful metrics boost efficiency



### MAS Omni•CORE

Thermo Scientific MAS Omni•CORE Controls consolidate general chemistry and serum protein QC processing into a single product, offering three distinct levels.

Part No.	Level	Bottles & Size	Storage & Stability	Matrix
OCR-101	1	6 x 5 mL Assayed	36 month shelf life @ -25 to -15 °C from date of manufacture 30 days open vial @ 2 to 8 °C (Bilirubin is 15 days open vial @ 2 to 8 °C)	Human serum
OCR-202	2			
OCR-303	3			



### MAS chemTRAK•H

Thermo Scientific MAS chemTRAK•H is a comprehensive multianalyte control designed to monitor general chemistry and therapeutic drug monitoring (TDM) methods. The control is liquid, ready to use, and no additional preparation is required.

Part No.	Level	Bottles & Size	Storage & Stability	Matrix
CHA-1	1	6 x 5 mL Assayed	30 month shelf life @ -25 to -15 °C from date of manufacture 7 days open vial @ 2 to 8 °C	Human serum
CHA-2	2			
CHA-3	3			



## Smartify your lab

Save time and improve efficiency in your lab with Thermo Scientific LabLink xL™ Quality Assurance Software



Learn more about LabLink xL and Auto-Connectivity Software features and requirements at [maslablink.com](http://maslablink.com)

Call your local Thermo Fisher Scientific representative or Tech Services for additional information.

#### Clinical Diagnostics

**USA**  
46500 Kato Road  
Fremont, CA 94538 USA  
Tel. 800-232-3342 / + 1 510-979-5000  
sales.diagnostics.fmt@thermofisher.com

**Asia Pacific**  
Tel. + 65 6499 9999  
cdd.asia.info2@thermofisher.com

**Australia**  
Tel. + 61 1800 333 110  
auinfo@thermofisher.com

**Canada**  
Tel. + 1 905 286 4290  
CDD.Canada@thermofisher.com

**China**  
Tel. + 86 400 650 5118  
sales.china@thermofisher.com

**France, Belgium (Brussels / Wallonia) & Luxembourg**  
Tel. + 33 (0) 1 40 86 65 00  
cdx.fr.info@thermofisher.com

**Germany**  
Tel. + 49 (0) 3302 883 608  
cdx.de.info@thermofisher.com

**Hong Kong**  
Tel. + 852 3107 7600  
SDG.HK.Sales@thermofisher.com

**Italy**  
distributor.cdd@thermofisher.com

**Japan**  
Tel. + 81 (0) 120 489 211  
JPYOK-CDD.QC@thermofisher.com

**Netherlands & Belgium (Flemish)**  
Tel. + 49 3302 883 660  
Info.benelux.cdd@thermofisher.com

**New Zealand**  
Tel. + 61 0800933 966  
nzinfo@thermofisher.co.nz

**Nordic**  
Tel. + 46 8 594 691 30  
info.nordic.cdd@thermofisher.com

**Russia**  
Tel. + 7 495 739 76 41  
info.btd.moscow@thermofisher.com

**South Africa**  
Tel. + 27 11 792 6790  
support-za.idd@thermofisher.com

**Spain & Portugal (Iberia)**  
Tel. + 34 93589 8338  
CDD.Info.ES@thermofisher.com

**Switzerland & Austria**  
Tel. + 41 26 663 86 70  
cdx.ch.info@thermofisher.com

**United Kingdom & Ireland**  
Tel. + 44 1442 868 940  
cdx.UK.info@thermofisher.com

**For Countries Not Listed**  
distributor.cdd@thermofisher.com

Find out more at [thermofisher.com/MASproductivity](http://thermofisher.com/MASproductivity)

Availability of product in each country depends on local regulatory marketing authorization status.