

Have you heard?

Chromeleon CDS extends to Capillary Electrophoresis assays with NEW CE driver

What is it? The CE Driver provides an interface for controlling the Sciex / Beckman PA 800 Plus within Chromeleon CDS

The benefits of Thermo Scientific™ Chromeleon™ Chromatography Data System extends to CE assays with the new instrument driver for Sciex™ (Beckman) PA 800 Plus! Bring in all of your CE assays that use photodiode array (PDA), ultraviolet (UV), and laser-induced fluorescence (LIF) detectors. The driver is compatible with Chromeleon CDS version 7.2.10 MUd onwards (including all 7.3 versions).

Why make the purchase? Integrate workflows, optimize practices, reduce resources

Single software solution for your IC-, LC-, GC-MS(/MS), and now CE assays! Lower operating and maintenance costs and reduce training required with centralized management of one single software to acquire, process, and report all of your laboratory data. Take advantage of the well-loved benefits Chromeleon CDS already offers!

What does it mean? Laboratories with Chromeleon CDS can easily incorporate CE assays.

Seamlessly transfer methods from 32 Karat software or utilize the guided instrument method wizard. Aggregate data across the assays to correlate data and trends and improve your lab efficiency with a common user interface.



When considering software for your CE assays, choose one that's future-proof, compliant and feature-rich. Contact your sales account manager today.

Learn more at thermofisher.com/CEDrivers4CDS