

GC parts and accessories

# GC parts and accessories


## for TRACE 1300/1600 Series Gas Chromatograph and GC-MS Systems

**Achieve optimal results and maximize the reproducibility of your chromatography**












Quality parts and accessories are designed for your Thermo Scientific™ TRACE™ 1300/1600 Series Gas Chromatograph, related GC-MS systems, and modules.



### Thermo Scientific™ Syringes

Product	Description	Quantity	Cat. no.
	<b>Syringes for Thermo Scientific™ TriPlus™ RSH SMART Autosamplers</b>		
	5 µL Fixed Needle SMART Syringe, 26c/cone	1/pack	<a href="#">365C0251-SM</a>
	10 µL Fixed Needle SMART Syringe, 26c/cone	1/pack	<a href="#">365D0291-SM</a>
	10 µL Fixed Needle Gas Tight SMART Syringe, 26c/cone	1/pack	<a href="#">365D0331-SM</a>
	10 µL Fixed Needle Gas Tight SMART Syringe, 23c/cone	1/pack	<a href="#">365D0311-SM</a>
	1 mL Fixed Needle Gas Tight SMART Syringe, 22/LC	1/pack	<a href="#">365K2811-SM</a>
	2.5 mL Fixed Needle Headspace SMART Syringe, 23G	1/pack	<a href="#">365L2321-SM</a>
	<b>Syringes for Thermo Scientific™ AI/AS 1310 or AI/AS 1610 Liquid Autosamplers</b>		
	5 µL Fixed Needle Syringe, 23/cone	1/pack	<a href="#">365C3701</a>
	10 µL Fixed Needle Syringe, 25/cone	1/pack	<a href="#">36500525</a>
	10 µL Fixed Needle Syringe, 23/cone	1/pack	<a href="#">36520060</a>
	10 µL Fixed Needle Syringe, 26/cone	6/pack	<a href="#">365D1856</a>


## Thermo Scientific™ LinerGOLD™ GC Injection Port Liners

Product	Description	Quantity	Cat. no.
	<b>Direct Straight Liner</b> 1 mm ID	5/pack	<a href="#">453A1335-UI</a>
		25/pack	<a href="#">453A2332-UI</a>
	<b>SPME Arrow Liner</b> 1.7 mm ID	3/pack	<a href="#">453A0415</a>
		5/pack	<a href="#">453A0415-UI</a>
	<b>Split Straight Liner</b> 4 mm ID, quartz wool	5/pack	<a href="#">453A2265</a>
		5/pack	<a href="#">453A2265-UI</a>
		25/pack	<a href="#">453A1262</a>
		25/pack	<a href="#">453A1262-UI</a>
	<b>Split/Splitless Precision Liner</b> 4 mm ID, quartz wool	5/pack	<a href="#">453A1255</a>
		5/pack	<a href="#">453A1255-UI</a>
		25/pack	<a href="#">453A1252</a>
	<b>Split/Splitless Single Taper Liner</b> 4 mm ID, quartz wool	25/pack	<a href="#">453A1252-UI</a>
		5/pack	<a href="#">453A1925</a>
		5/pack	<a href="#">453A1925-UI</a>
	<b>Split/Splitless Single Taper Liner</b> 4 mm ID	25/pack	<a href="#">453A2295</a>
		25/pack	<a href="#">453A2922-UI</a>
		5/pack	<a href="#">453A1345</a>
	<b>Split/Splitless Double Taper Cyclo Liner</b> 4 mm ID	5/pack	<a href="#">453A1345-UI</a>
		25/pack	<a href="#">453A1345-UI</a>
		25/pack	<a href="#">453A2342</a>
	<b>PTV Liner with three baffles</b> 1 × 2.75 × 120 mm	25/pack	<a href="#">453A2342-UI</a>
		5/pack	<a href="#">45352062</a>
	<b>PTV Siltek Baffle Liner</b> 2 mm ID	5/pack	<a href="#">453T2120</a>
	<b>PTV Concentric Baffle Liner</b> 2 mm ID	5/pack	<a href="#">453T2845-UI</a>
	<b>PTV Silcosteel Liner for OC</b> 1 mm ID with restriction	2/pack	<a href="#">45322052</a>


## Thermo Scientific™ LinerGOLD™ GC Liner Sealing Rings

Description	Quantity	Cat. no.
Liner Sealing Ring for GC SSL Injector	5/pack	<a href="#">29001320</a>
Liner Sealing Ring for GC PTV Injector	1/pack	<a href="#">29001318</a>


## Thermo Scientific™ Inlet Base Seals

Product	Description	Quantity	Cat. no.
	<b>Gold Inlet Base Seals</b> 0.8 mm ID	10/pack	<a href="#">290GA081</a>
		2/pack	<a href="#">290GA082</a>
	<b>Gold Inlet Base Seals</b> 0.8 mm ID, cross version	10/pack	<a href="#">290GA084</a>
		2/pack	<a href="#">290GA083</a>
	<b>Siltek-Treated Inlet Base Seals</b> 0.8 mm ID	10/pack	<a href="#">290GA091</a>
		2/pack	<a href="#">290GA092</a>


## Thermo Scientific™ Ferrules

Product	Description	Quantity	Cat. no.
	100% Graphite Ferrule for 0.1-0.32 mm ID columns	10/pack	<a href="#">290GA139</a>
	100% Graphite Ferrule for 0.1-0.25 mm ID columns, PTV Injector	10/pack	<a href="#">29053488</a>
	15% Graphite/85% Vespal Ferrules for 0.1-0.25 mm ID columns	10/pack	<a href="#">290VA191</a>
	15% Graphite/85% Vespal Ferrules for 0.32 mm ID columns	10/pack	<a href="#">290VA192</a>
	15% Graphite/85% Vespal Ferrules for 0.53 mm ID columns	10/pack	<a href="#">290VA193</a>
	Vespal SCP-5000 Ferrules for 0.1-0.25 mm ID columns	5/pack	<a href="#">290VT221</a>
	Column Nut for SSL injector and non-MS detectors	5/pack	<a href="#">35050458</a>
	Column Nut for PTV Injector	5/pack	<a href="#">35053221</a>
	Spring-loaded Transfer Line Nut for Thermo Scientific™ ISQ™ and Thermo Scientific™ TSQ™ GC-MS Systems	1/pack	<a href="#">1R120434-0010</a>


## Thermo Scientific™ Septas

Product	Description	Quantity	Cat. no.
	<b>BTO Low Bleed Silicone Septa</b> 11 mm diameter	50/pack, glass jar	<a href="#">31303233</a>
		50/pack, blister pack	<a href="#">31303233-BP</a>
	<b>TR-Green Low Coring Silicone Septa</b> 11 mm diameter	50/pack, glass jar	<a href="#">313G3230</a>
		50/pack, blister pack	<a href="#">313G3230-BP</a>
	<b>Marathon Silicone Septa</b> 11 mm diameter	50/pack, glass jar	<a href="#">313P3233</a>
		50/pack, blister pack	<a href="#">313P3233-BP</a>


## Thermo Scientific™ Gas Filters

Product	Description	Quantity	Cat. no.
	<b>Indicating Triple Filter for carrier gas</b> He preconditioned for GC-MS	1/pack	<a href="#">60180-825</a>
	<b>Gas Filter Kit</b> Base + He pre-conditioned Triple Filter	1/pack	<a href="#">60180-829</a>
	<b>Indicating Oxygen Cartridge Filter</b>	1/pack	<a href="#">60180-820</a>
	<b>Indicating Combi Filter</b> Moisture, hydrocarbons	1/pack	<a href="#">60180-826</a>


## Thermo Scientific™ Equipment

Product	Description	Quantity	Cat. no.
	Thermo Scientific™ GFM Pro Gas Flowmeter	1/pack	<a href="#">66002-010</a>
	Thermo Scientific™ GLD Pro Gas Leak Detector	1/pack	<a href="#">66002-004</a>
	Thermo Scientific™ GC Tool Kit for Thermo Scientific instruments	1/pack	<a href="#">60180-784</a>

## Thermo Scientific™ SureSTART™ Vials and Caps

Product	Description	Quantity	Cat. no.
	<b>2 mL clear screw glass vial</b> 9 mm short thread Thermo Scientific™ SureSTOP™ GOLD cap, Thermo Scientific™ AVCS™ screw 9 mm blue PP, white silicone/red PTFE septa vial kit	100/pack	<a href="#">6PKG592W</a>
	<b>2 mL amber screw glass vial</b> 9 mm short thread SureSTOP cap, AVCS screw 9 mm blue PP, white silicone/red PTFE vial kit	100/pack	<a href="#">6AK594W</a>
	<b>2 mL clear screw glass vial</b> 9 mm short thread cap, AVCS screw 9 mm blue PP, blue silicone/clear PTFE septa soft vial kit	100/pack	<a href="#">6AK599W</a>
	<b>2 mL clear screw glass vial</b> 9 mm short thread cap, AVCS screw 9 mm blue PP, red PTFE/white silicone/Red PTFE vial kit	100/pack	<a href="#">6AK86W</a>
	<b>2 mL clear screw glass vial</b> 9 mm short thread SureSTOP cap, AVCS screw 9 mm blue PP, white silicone/red PTFE vial kit	100/pack	<a href="#">6AK592W</a>

**Thermo Scientific™ TraceGOLD™ GC Columns**

Product	Description	Quantity	Cat. no.
	Thermo Scientific™ TraceGOLD™ <b>TG-1MS</b> GC Column 30 m × 0.25 mm × 0.25 µm	1/pack	<a href="#">26099-1420</a>
	TraceGOLD <b>TG-1MS</b> GC Column 60 m × 0.25 mm × 0.25 µm	1/pack	<a href="#">26099-1540</a>
	Thermo Scientific™ TraceGOLD™ <b>TG-5SILMS</b> GC Column 30 m × 0.25mm × 0.25 µm	1/pack	<a href="#">26096-1420</a>
	Thermo Scientific™ TraceGOLD™ <b>TG-5MS</b> GC Column 30 m × 0.32 mm × 0.25 µm	1/pack	<a href="#">26098-1430</a>
	Thermo Scientific™ TraceGOLD™ <b>TG-624SILMS</b> GC Column 30 m × 0.32 mm × 1.8 µm	1/pack	<a href="#">26059-3390</a>
	Thermo Scientific™ TraceGOLD™ <b>TG-1701MS</b> GC Column 30 m × 0.25 mm × 0.25 µm	1/pack	<a href="#">26090-1420</a>
	Thermo Scientific™ TraceGOLD™ <b>TG-WaxMS</b> GC Column 30 m × 0.25 mm × 0.25 µm	1/pack	<a href="#">26088-1420</a>
	Thermo Scientific™ TRACE™ <b>TR-Pesticide II</b> GC Column 30 m × 0.25 mm × 0.25 µm + 5 m Guard	1/pack	<a href="#">26RD142F</a>
	Thermo Scientific™ TraceGOLD™ <b>TG-PBDE</b> GC Column 15 mm × 0.25 mm × 0.10 µm	1/pack	<a href="#">26061-0350</a>
	Thermo Scientific™ TraceGOLD™ <b>TG-PAH</b> GC Column 60 m × 0.25 mm × 0.1 µm	1/pack	<a href="#">26055-0120</a>
	Thermo Scientific™ TraceGOLD™ <b>TG-Dioxin</b> GC Column 60 m × 0.25 mm × 0.25 µm	1/pack	<a href="#">26066-1540</a>
	Thermo Scientific™ TRACE™ <b>TR-FAME</b> GC Column 100 m × 0.25 mm ID × 0.2 µm	1/pack	<a href="#">260M238P</a>

 Learn more at [thermofisher.com/gcparts](https://thermofisher.com/gcparts)

**General Laboratory Equipment – Not For Diagnostic Procedures.** © 2020, 2023 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. This information is presented as an example of the capabilities of Thermo Fisher Scientific products. It is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details. **FL20889-EN 0623**

**thermo** scientific