


Preventing sample failures in bioanalysis

SOLA 30 mg Solid Phase Extraction (SPE) cartridges and plates



Thermo Scientific™ SOLA™ Solid Phase Extraction (SPE) products are designed to meet bioanalytical and clinical research demands relying on very clean extracts even at low detection limits. With the unique and innovative fritless SPE technology, which eliminates the issues with traditional loose-packed SPE formats, SOLA products are available in 5 polymeric phases, and in 10 mg and 30 mg bed weight in both 1 mL and 3 mL cartridge and 96-well plate formats. A 2 mg 96-well plate product for microelution is also available for low sample volumes and where increase in sensitivity is required.


SOLA SPE fritless design



This innovative design enables high levels of reproducibility when dealing with highly complex biological samples, not only from cartridge-to-cartridge or well-to-well but also batch-to-batch when compared to conventional SPE products.

The combination of the support material and active media components into a solid, uniform sorbent bed delivers stable and controllable flow-through characteristics and provides an added advantage when dealing with viscous biological samples—it reduces sample failures, and enables high throughput processing.

Batch-to-batch reproducibility of SOLA product compared to a conventional SPE product



Ordering information

| Product | Bed weight (mg) | Cartridge/well plates volume (mL) | Quantity | Cat. no. |
|--|-----------------|-----------------------------------|----------|-----------|
| Thermo Scientific™ SOLAμ™ Solid Phase Extraction (SPE) 96-well plates | | | | |
| SOLAμ HRP 96-well plate | 2 | 1 | 1 each | 60209-001 |
| SOLAμ SCX 96-well plate | 2 | 1 | 1 each | 60209-002 |
| SOLAμ SAX 96-well plate | 2 | 1 | 1 each | 60209-003 |
| SOLAμ WCX 96-well plate | 2 | 1 | 1 each | 60209-004 |
| SOLAμ WAX 96-well plate | 2 | 1 | 1 each | 60209-005 |
| Thermo Scientific™ SOLA™ Solid Phase Extraction (SPE) cartridges | | | | |
| SOLA HRP SPE cartridge | 10 | 1 | 100 pack | 60109-001 |
| SOLA SCX SPE cartridge | 10 | 1 | 100 pack | 60109-002 |
| SOLA SAX SPE cartridge | 10 | 1 | 100 pack | 60109-003 |
| SOLA WCX SPE cartridge | 10 | 1 | 100 pack | 60109-004 |
| SOLA WAX SPE cartridge | 10 | 1 | 100 pack | 60109-005 |
| SOLA HRP SPE cartridge | 30 | 3 | 50 pack | 60409-001 |
| SOLA SCX SPE cartridge | 30 | 3 | 50 pack | 60409-002 |
| SOLA SAX SPE cartridge | 30 | 3 | 50 pack | 60409-003 |
| SOLA WCX SPE cartridge | 30 | 3 | 50 pack | 60409-004 |
| SOLA WAX SPE cartridge | 30 | 3 | 50 pack | 60409-005 |
| Thermo Scientific™ SOLA™ Solid Phase Extraction (SPE) 96-well plates | | | | |
| SOLA HRP 96-well plate | 10 | 2 | 1 each | 60309-001 |
| SOLA SCX 96-well plate | 10 | 2 | 1 each | 60309-002 |
| SOLA SAX 96-well plate | 10 | 2 | 1 each | 60309-003 |
| SOLA WCX 96-well plate | 10 | 2 | 1 each | 60309-004 |
| SOLA WAX 96-well plate | 10 | 2 | 1 each | 60309-005 |
| SOLA HRP 96-well plate | 30 | 2 | 1 each | 60509-001 |
| SOLA SCX 96-well plate | 30 | 2 | 1 each | 60509-002 |
| SOLA SAX 96-well plate | 30 | 2 | 1 each | 60509-003 |
| SOLA WCX 96-well plate | 30 | 2 | 1 each | 60509-004 |
| SOLA WAX 96-well plate | 30 | 2 | 1 each | 60509-005 |

Expect reproducible results with sample prep, columns and vials



Don't see what you need? We would be happy to discuss your specific requirements. Please contact your local sales representative for custom orders.

Find out more at thermofisher.com/solaspe

ThermoFisher
SCIENTIFIC