

Thermo Scientific Dionex Aquion Ion Chromatography (IC) System Letter of Equivalency

The Thermo Scientific™ Dionex™ Aquion IC system is currently shipping and is a direct replacement for the Thermo Scientific™ Dionex™ ICS-1100 IC system, which is the successor to the Thermo Scientific™ Dionex™ ICS-1000 system.

The Dionex Aquion IC system has the same fit, form and function as the Dionex ICS-1100 and Dionex-1000 systems. It has been designed and validated to directly address all the applications supported by the Dionex ICS-1100 and the Dionex ICS-1000 systems (except for those using Thermo Scientific™ Dionex™ Reagent-Free™ ion chromatography with eluent regeneration (RFIC-ER)).

The Dionex Aquion IC system can use the same columns as the Dionex ICS-1100 and Dionex-1000 systems, providing equivalent chromatographic performance. Refer to the consumable manual to determine appropriate operating conditions.

Please contact your Thermo Fisher Scientific representative for more information, or visit [thermofisher.com/aquion](https://www.thermofisher.com/aquion).