



**General Arrangement Drawings
Microplates for Chromatography-Certified**

Cat.No.	Description	A	B	F	G1	H	J	K	L	M	N
		BASE LENGTH	BASE WIDTH	OVERALL HEIGHT	DIA TOP	A - I LOCATION (X)	A - I LOCATION (Y)	WELL TO WELL SPACING	WELL DEPTH	BOTTOM OF WELL DISTANCE	FLANGE SHORT SIDE
60180-P205 60180-P215	Round Well-Flat bottom	128.00	85.40	14.70	7.05	14.70	11.50	9.00	13.90	3.00	3.00
60180-P206 60180-P216	Round Well-Round bottom	128.00	85.40	14.70	7.05	14.70	11.50	9.00	13.90	3.00	3.00
60180-P207 60180-P217	Round Well-V bottom	128.00	85.40	14.70	7.05	14.70	11.50	9.00	13.90	3.00	3.00



Notes.

Drawings are for instrument set-up purposes and are not manufacturing drawings

Thermo Fisher Scientific is providing information that has been supplied by the original manufacturer and is not responsible for subsequent changes or modifications.

All dimensions mm