



38<sup>th</sup> International Symposium on Halogenated Persistent Organic Pollutants

**DioXin 2018**

& 10<sup>th</sup> International PCB Workshop

26 – 31 August 2018, Kraków, Poland

# Using Magnetic Sector DFS with DualData XL in a Commercial Dioxin Lab

- David Hope, P.Chem., CEO
- Pacific Rim Laboratories Inc.
- [dave@pacificrimlabs.com](mailto:dave@pacificrimlabs.com)

# Overview

- History
- DualData – What is it?
- DualData XL – Why?
- DualData XL – Practical Applications
- Conclusions

# Our first HRMS - VG70 arrives!



**PACIFIC RIM**  
LABORATORIES INC

# Beginnings

## Dioxin 2005 - Toronto



## Thermo POPs Symposium– Venice 2007





# More business means more instruments

First DFS (2007) – 1450 samples per year;  
2250 tests – staff of 7

Second DFS (2010) – 2530 samples per year;  
3500 tests – staff of 9

TSQ8000Evo (2014) – 3500 samples per year;  
4450 tests – staff of 11

Third DFS (2015) – 4600 samples per year;  
6750 tests – staff of 15

DualDataXL installed on DFS – projecting  
6500 samples



# Pacific Rim Laboratories Inc.


## Scientific Research and Experimental Development


- We thrive on innovation
  - PBDE method in 2005
  - 209 congener PCBs in 2005
  - Sub-ppb PAH analysis food in 2006
  - Published 2009 congener PCB by SGE HT8 column (2009)
  - Cape Tech column clean-up (2012)
  - OCPs by HRMS (2012)
  - Improved clean-up methods for dioxins/PCB (2014)
  - Rocket evaporator (2015)
  - Single run PAH and alkylated PAH on TSQ8000Evo (2015)



# What is Dual Data? PCDD/F analysis

## EPA 1613b – TCDD must elute >25 minutes










# DFS DualData XL: Staggered Injection


- Both GCs are running simultaneously.
- The injection on GC 2 is performed during the acquisition of GC 1 and vice versa.



# Flow switching 1: Column-flow is directed into the MS





Example Standard inj.




Source pressure eq. Solvent inj.



# Flow switching 2: Flow directed into service vacuum (waste)



Example Standard inj.




Source pressure eq. Solvent inj.




# Alternating Flow switching with both GCs

- **A:** Acquisition GC1
  - GC1 into MS
  - GC2 into waste




- **B:** Acquisition GC2
  - GC2 into MS
  - GC1 into waste






# DFS DualData XL: Staggered Injection Sequence Timing



# Why DualData XL?



- Cheaper than buying new instrument
- No additional floor space required
- No additional electrical considerations
- Autosampler ready to inject as soon as we get ready signal
- Can double our through put with mixed chemistries



**PACIFIC RIM**  
LABORATORIES INC

# Requirements

- DFS with dual 1310 GCs
- Older DFS can be converted, but cannot use Trace GCs.
- Will add a 4-way valve to control gases
- A gas module is installed next to your injector.
- New software





\*DFS



DualData

DFS

Config Files

GC 1

Mode : MID-Method DDXL\_WPCB\_GC2\_NEW.mid

MID: C:\Xcalibur\System\DFS\Msi\DDXL\_WPCB\_GC2\_NEW.mid

GC 2

DDXL\_GC1\_DF\_5W\_FC43\_GC2.mid

MID: C:\Xcalibur\System\DFS\Msi\DDXL\_GC1\_DF\_5W\_FC43\_GC2.mid

Acquisition Times [min]

Use Solvent Delay only if no DualData wafer is installed Solvent Delay : 0

Additional Action

Scripts:

Select:

Resolution Check

☐ Activate

Report Reference

Target Resolution 10000


Resolution Validation


Warning if below 9500.00 Error if below 8000.00 ☒ Stop Sequence on Error

Lock Peak Mode: No Check










**PACIFIC RIM**  
LABORATORIES INC




# Dual PCDD/F

Runtime DFS – 62 minutes or 23 inj/day

Runtime DualData XL – 67 minutes to complete cycle for both GCs

**42 injections in 24 h!**


82% more samples!



**PACIFIC RIM**  
LABORATORIES INC




# Column Resolution - TCDD



**PACIFIC RIM**  
LABORATORIES INC




# HxCDF resolution




**PACIFIC RIM**  
LABORATORIES INC


# TCDD - CS-Lo 0.1 pg injected



TCDD – CS-LoLo! 20 fg injected



**PACIFIC RIM**  
**LABORATORIES INC**



## Dual PCB

Runtime DFS – 68 minutes (21 inj/day)

Runtime DualData XL – for dioxin-like PCBs  
74 minutes to complete cycle for both GCs

**38 inj/day**

82% more analyses

Runtime DualData XL – for 209 congener  
PCBs 94 minutes to complete cycle for both  
GCs

**30 inj/day**


43% more analyses




**PACIFIC RIM**  
LABORATORIES INC



# HxCB @ 0.1/0.5 pg injected





# Dual BDE

Runtime DFS – 41 minutes (35 inj/day)

Runtime DualData XL – 60 minutes to complete cycle for both GCs


**47 inj/day**

34% more runs!



**PACIFIC RIM**  
LABORATORIES INC

# TeBDE @ 0.5 pg injected




**PACIFIC RIM**  
LABORATORIES INC

# DecaDBE

4/29/2017 3:25:02 PM

RT: 27.41 - 28.81 SM: 5G

CS-1  
2.5 pg



NL: 4.0000

PM

MS

HR042817

05

MS

HR042817

05

MS

HR042817

05

MS

HR042817

05

MS

HR042817

05

MS

HR042817

05

MS

HR042817

05

MS

HR042817

05

MS

HR042817

05

MS

HR042817

05

MS

HR042817

05

MS

HR042817

05

MS

HR042817

05

MS

HR042817

05

MS

HR042817

05

MS

HR042817

05

MS

HR042817

05

MS

HR042817

05

MS

HR042817

05

MS

HR042817

05

MS

HR042817

05

MS

HR042817

05

MS

HR042817

05


MS

HR042817

05



**PACIFIC RIM**  
LABORATORIES INC



# Dual OCP

Runtime DFS – 51 minutes (28 inj/day)

Runtime DualData XL – 64 minutes to complete cycle for both GCs

**44 inj/day**

57% more analyses



**PACIFIC RIM**  
LABORATORIES INC

# Our workload is not just dioxins



- PCDD/F 21%
- PCB 23%
- PAH\* 35%
- OCP 6%
- TBT\* 7%
- PBDE 2%
- Other\* 6%

• \*mix of HRMS and MS/MS




# PCDD/F & dIPCB

PCDD/F Runtime DFS – 62 minutes


PCB Runtime DFS – 48 minutes

Runtime DualData XL – 57 minutes to complete cycle for both GCs

Therefore, you can run dioxins and PCBs together in less time than it takes to run one dioxin sample!!!



**PACIFIC RIM**  
LABORATORIES INC



# PCDD/F and 209 congener PCB


PCDD/F Runtime DFS – 62 minutes

PCB Runtime DFS – 69 minutes  
(EPA1668C: PCB209 cannot elute  
before 55 minutes)

Runtime DualData XL – 77 minutes to  
complete cycle for both GCs



**PACIFIC RIM**  
LABORATORIES INC



# PCDD/F & TBT

PCDD/F Runtime DFS – 62 minutes


TBT Runtime DFS – 30 minutes

Runtime DualData XL – 56 minutes\*

\* This takes less time than running dioxins by themselves!!



**PACIFIC RIM**  
LABORATORIES INC



# PCDD/F & OCP


PCDD/F Runtime DFS – 62 minutes

OCP Runtime DFS – 51 minutes

Runtime DualData XL – 66 minutes



**PACIFIC RIM**  
LABORATORIES INC



# PCDD/F and PBDE

PCDD/F Runtime DFS – 62 minutes

PBDE Runtime DFS – 41 minutes

Runtime DualData XL – 64 minutes

Note – you must run PBDE at 10,000 resolution

In our lab we use different tuning compound for PBDE v dioxin, so would never run together. This is just an example of what could be done.




**PACIFIC RIM**  
LABORATORIES INC



# Conclusions

1. DualData XL saves time, even when running in single GC mode
2. Not limited to running the same column/program in each GC
3. Source changes less frequent as “burn off” being vented to air
4. Wafers are a consumable, but are cleanable too!
5. Source does not need to be vented when changing columns





Thank you from Pacific Rim Labs.

