

BioPharma Finder 3.0 System Requirements

The following tables list the minimum and recommended hardware and software configurations for Thermo BioPharma Finder™ 3.0 operation.

Minimum System Requirements

Table 1. Minimum system requirements

System	Requirements
Data system computer	<ul style="list-style-type: none">• Intel™ Core™ i7-4770 CPU@3.40 GHz• 8 GB registered RAM• 100 GB storage hard drive – ST1000DM -003 SCSI disk device• DVD/CD-ROM drives• Resolution display 1280 ×1024 (SXGA)
Software	<ul style="list-style-type: none">• Microsoft™ Windows™ 7 Professional SP1, 64-bit (English version)^a• Microsoft .NET 4.6.2• Microsoft Office 2010• Adobe™ Acrobat™ Reader™ DC

^a Windows 10 Professional, 64-bit (English version) is also supported.

Recommended System Requirements

For enhanced performance using the BioPharma Finder application, follow these recommendations.

Table 2. Recommended system requirements

System	Requirements
Data system computer	<ul style="list-style-type: none">• Quad core Intel Xeon™ CPU• (E5-1630v3 3.7 GHz 10MB 2133 4C)• 32 GB DDR4-2133 (4 × 8 GB) registered RAM• 2 TB storage hard drive (SATA, 7200 rpm)—512 GB solid state boot drive (SATA)• DVD-RW optical drive• Resolution display 1920 ×1080
Software	<ul style="list-style-type: none">• Microsoft™ Windows™ 7 Professional SP1, 64-bit (English version)• Microsoft .NET 4.6.2• Microsoft Office 2016• Adobe Acrobat Pro DC

Trademarks

BioPharma Finder is a trademark of Thermo Fisher Scientific Inc. in the United States.

The following are registered trademarks in the United States and other countries: Microsoft and Windows are registered trademarks of Microsoft Corporation. Intel, Intel Core, and Xeon are registered trademarks of Intel Corporation. Adobe, Acrobat, and Reader are registered trademarks of Adobe Systems Inc.