

Compound Discoverer 2.1 SP1 System Requirements

The following tables list the minimum system requirements and the recommended hardware configuration for Thermo Compound Discoverer™ 2.1 SP1 operation.

Minimum System Requirements

Table 1. Minimum system requirements

System	Requirements
Hardware	<ul style="list-style-type: none">• 3.4 GHz dual-core processor• 16GB RAM• 500 GB hard drive• DVD-ROM and USB drive• Display monitor resolution of 1920×1080
Software	<ul style="list-style-type: none">• Microsoft™ Windows™ 7 Pro SP1 (64-bit) or Windows 10 (64-bit) operating system• Microsoft .NET Framework 4.6.2• Microsoft Office 2010• Adobe™ Reader™ 11• Adobe Flash™ Player 15

Recommended Hardware Configuration

For enhanced performance using the Compound Discoverer application, follow these recommendations.

Table 2. Recommended hardware configuration

System	Requirements
Hardware	<ul style="list-style-type: none">• Intel Core™ i7 3930x @ 3.30 GHz processor (or equivalent)• 3.4 GHz dual-core processor• 32 GB RAM• 1 TB SSD (solid-state disk) hard drive for OS• 2nd 3 TB (conventional disk) hard drive for data storage• DVD-ROM and USB drive• Two 27 in. HD monitors: Display monitor resolution of 1920 ×1080

Trademarks

Compound Discoverer is a trademark of Thermo Fisher Scientific Inc. in the United States.

The following are registered trademarks in the United States and other countries: Microsoft and Windows are registered trademarks of Microsoft Corporation. Intel Core is a registered trademark of Intel Corporation. Adobe, Reader, and Flash are registered trademarks of Adobe Systems Inc.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.