

Thermo Scientific Dionex Chromeleon 7.2 Chromatography Data System

Recycle Bin




Product Spotlight

The new Recycle Bin feature in the Thermo Scientific™ Dionex™ Chromeleon™ 7.2 Chromatography Data System (CDS) software allows users to easily recover items deleted from the Console such as sequences, folders and metadata outside of sequences. Chromeleon 7.2 CDS gets you from samples to results quickly and easily – boosting your overall lab productivity.



New "Deleted Items" Recycle Bin feature.



A quick right click and your data can be restored to its original or to an alternate location.

Each Data Vault has its own Recycle Bin and user privileges can be activated to control deletion of data.



For more information, visit thermoscientific.com/chromeleon

©2013 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

SP70769_E 07/13S

Thermo
SCIENTIFIC