

SPARE PART KIT RECOMMENDATIONS

Thermo Scientific™ Model 51i-HT High Temp Total Hydrocarbon Analyzer

To keep your system running and properly collecting data, it's important to have spare parts on hand to avoid costly downtime. We developed Spare Parts Kits to ensure maximum uptime of your Thermo Scientific Instrument during its life cycle.

Spare Parts Kit # 115898-00	
Part Numbers	Description
4070	1 x Fan Filter
11290	1 x Receptacle Fuse Drawer, 110V
12190	1 x Filter Inlet Element, 0.5 micron
12188	1 x Aluminum Gasket (for P/N 12190)
11032	1 x O-ring, FID Flame Sensor
105198-00	2 x Fuse, 250 VAC, 6.25 Amp
105567-00	1 x Ignitor/Cable Assembly
108091-00	1 x PUMP REBUILD KIT
103581-00	1 x Flame Sensor Assembly
104139-00	1 x Sample Frit, FID

Extended Maintenance Kit # 115891-00	
Part Numbers	Description
103803-00	1 x Fan 4", 24V
101681-00	1 x power supply w/o cables
103762-00	1 x Solenoid Valve Assembly
105554-01	1 x Screened Orifice
103950-00	1 x Air Frit, FID
103687-00	1 x Pogo Detector Assembly
104872-00	1 x Pressure Regulator
112286-00	1 x Sample Pump

USA
27 Forge Parkway
Franklin, MA 02038
Ph: (866) 282-0430
orders.aqi@thermofisher.com

India
C/327, TTC Industrial Area
MIDC Pawane
New Mumbai 400 705, India
Ph: +91 22 4157 8800
india@thermofisher.com

China
+Units 702-715, 7th Floor
Tower West, Yonghe
Beijing, China 100007
+86 10 84193588
info.eid.china@thermofisher.com

Europe
Ion Path Road Three
Winsford, Cheshire CW7 3GA
Ph: +44 (0) 1606 548700
orders.epm.uk@thermofisher.com

Find out more at thermofisher.com/EMservice

ThermoFisher
SCIENTIFIC