
Specifications

Product Description : Ammonium Persulphate ExcelsaR-ACS 500G
Product Code : Q15055

Sr.No	Description	Unit	Standards
1	Description : White / pale yellow crystals / crystalline powder / granules (Slowly decomposes with loss of oxygen)		Should Comply
2	Solution: 10% w/v (in water) : Clear & colorless to pale yellow		Should Comply
3	Assay	%	99 Min
4	Chloride (Cl)	%	0.0005 Max
5	Heavy metals (as Pb)	%	0.002 Max
6	Insoluble matter	%	0.003 Max
7	Iron (Fe)	%	0.001 Max
8	Manganese (Mn)	%	0.00005 Max
9	Non-volatile matter	%	0.05 Max
10	Acidity (as H ₂ SO ₄) (Meq/g)	%	0.04 Max

Quality Assurance

The document has been produced electronically and is valid without signature.