
Specifications

Product Description :Aluminium Ammonium Sulphate SQ 500G
Product Code :Q16055

Sr.No	Description	Unit	Standards
1	Description : Colourless crystals / crystalline powder		Should Comply
2	Solution: 10% w/v (in water) : Almost clear & colourless		Should Comply
3	Alkalis (Na)	%	0.5 Max
4	Assay (ex Al)	%	99 Min
5	Chloride (Cl)	%	0.01 Max
6	Iron (Fe)	%	0.005 Max
7	Lead (Pb)	%	0.005 Max

Quality Assurance

The document has been produced electronically and is valid without signature.