

---

*Specifications*

---

**Product Description** : Iso Butyl Alcohol SQ 2.5L  
**Product Code** : Q25707

Sr.No	Description	Unit	Standards
1	Description : Clear colourless liquid (< 10 Hazen) (characteristic odour)		Should Comply
2	Solubility : Soluble in 20 parts of water ,miscible with alcohol and ether.		Should Comply
3	Assay (by G.C.)	%	98 Min
4	Non-volatile matter	%	0.005 Max
5	Water (by K.F.)	%	0.5 Max
6	Weight per ml (at 20°C) (In gm/ml)		0.801 - 0.804

---

**Quality Assurance**

The document has been produced electronically and is valid without signature.