
Specifications

Product Description : Phosphorous Pentoxide SQ 500G
Product Code : Q26375

Sr.No	Description	Unit	Standards
1	Description : Crystalline powder (deliquescent) (hygroscopic) (corrosive)		Should Comply
2	Assay	%	98 Min
3	Chloride (Cl)	%	0.002 Max
4	Heavy metals (as Pb)	%	0.01 Max
5	Reducing substances (as P ₂ O ₃)	%	0.02 Max

Quality Assurance

The document has been produced electronically and is valid without signature.
