
Specifications

Product Description :Sodium Hydroxide Flakes SQ 50KG
Product Code :Q2780E

Sr.No	Description	Unit	Standards
1	Description : White deliquescent flakes/chips		Should Comply
2	Solution: 10% w/v (in water) : Clear & colourless		Should Comply
3	Solution: 60% w/v (in water) : Reasonably clear & free from blackish sedimentation		Should Comply
4	Assay	%	97 Min
5	Carbonate (Na ₂ CO ₃)	%	2 Max
6	Chloride (Cl)	%	0.035 Max
7	Potassium (K)	%	0.1 Max
8	Silicate (SiO ₂)	%	0.05 Max
9	Sulphate (SO ₄)	%	0.05 Max
10	Zinc (Zn)	%	0.02 Max

Quality Assurance

The document has been produced electronically and is valid without signature.