
Specifications

Product Description : Trifluoroacetic Acid Chromatography 2ML
Product Code : Q28469

Sr.No	Description	Unit	Standards
1	Description : Clear colorless liquid (< 10 Hazen)		Should Comply
2	UV absorbance (0.1% aq. Solution) Cut-off at 282nm : Passes test		Should Comply
3	Assay (Acidimetric)	%	99.8 Min
4	Density at 20°C (In gm/ml)		1.485 - 1.495
5	Residue on evaporation	%	0.005 Max
6	Water by KF	%	0.05 Max
7	Maximum absorbance (1.0 cm cell vs H2O) at 275nm	%	0.04 Max
8	Maximum absorbance (1.0cm cell vs H2O) at 300nm	%	0.03 Max

Quality Assurance

The document has been produced electronically and is valid without signature.