
Specifications

Product Description :Pet Spirit (Pet Ether) 60°-80 200L SQ
Product Code :Q3451H

Sr.No	Description	Unit	Standards
1	Description : Clear colorless volatile liquid (< 10 Hazen)(characteristic Odor)		Should Comply
2	Non-volatile matter	%	0.002 Max
3	Water (by K.F.)	%	0.02 Max
4	Weight per ml (at 20°C) (In gm/ml)		0.66 - 0.675

Quality Assurance

The document has been produced electronically and is valid without signature.

Registered Office : Thermo Fisher Scientific India Pvt. Ltd. 403-404,B-wing, Delphi, Hiranandani Business Park, Powai, Mumbai-400076, INDIA.
Tel : +91-22-67162200/66803000, CIN-U73100MH2000PTC126872.