

Product name	TrueMark Respiratory Panel 2.0, TaqMan Array Card	TrueMark Respiratory Panel 2.0, TaqMan OpenArray Plate	Respiratory Tract Microbiota (RTM) v1, TaqMan Array Card	Respiratory Tract Microbiota (RTM) v1, OpenArray Plate
Cat. No.	A49047	A49044	A41238	A41237
Pathogen type	Verified assays for the following pathogens are pre-loaded on product			
Bacterial	<i>Bordetella (Pan)</i>	<i>Bordetella (Pan)</i>	<i>Bordetella (Pan)</i>	<i>Bordetella (Pan)</i>
	<i>Bordetella pertussis</i>	<i>Bordetella pertussis</i>	<i>Bordetella holmesii</i>	<i>Bordetella pertussis</i>
	<i>Chlamydophila pneumoniae</i>	<i>Chlamydophila pneumoniae</i>	<i>Bordetella pertussis</i>	<i>Chlamydophila pneumoniae</i>
			<i>Chlamydophila pneumoniae</i>	<i>Chlamydophila pneumoniae</i>
			<i>Coxiella burnetii</i>	
	<i>Haemophilus influenzae</i>	<i>Haemophilus influenzae</i>	<i>Haemophilus influenzae</i>	<i>Haemophilus influenzae</i>
	<i>Klebsiella pneumoniae complex</i>	<i>Klebsiella pneumoniae complex</i>	<i>Klebsiella pneumoniae complex</i>	<i>Klebsiella pneumoniae complex</i>
	<i>Legionella pneumophila</i>	<i>Legionella pneumophila</i>	<i>Legionella pneumophila</i>	<i>Legionella pneumophila</i>
<i>Moraxella catarrhalis</i>		<i>Moraxella catarrhalis</i>		
<i>Mycoplasma pneumoniae</i>	<i>Mycoplasma pneumoniae</i>	<i>Mycoplasma pneumoniae</i>	<i>Mycoplasma pneumoniae</i>	
<i>Staphylococcus aureus</i>	<i>Staphylococcus aureus</i>	<i>Staphylococcus aureus</i>	<i>Staphylococcus aureus</i>	
<i>Streptococcus pneumoniae</i>	<i>Streptococcus pneumoniae</i>	<i>Streptococcus pneumoniae</i>	<i>Streptococcus pneumoniae</i>	
Fungal	<i>Pneumocystis jirovecii</i>		<i>Pneumocystis jirovecii</i>	
Viral	Adenovirus	Adenovirus	Adenovirus	Adenovirus
	Bocavirus		Bocavirus	Bocavirus
	Coronavirus 229E	Human Coronavirus 229E	Coronavirus HKU1	Coronavirus HKU1
	Coronavirus HKU1	Human Coronavirus HKU1	Coronavirus NL63	Coronavirus NL63
	Coronavirus NL63	Human Coronavirus NL63	Coronavirus 229E	Coronavirus 229E
	Coronavirus OC43	Human Coronavirus OC43	Coronavirus OC43	Coronavirus OC43
	Enterovirus (Pan)	Human Enterovirus (pan assay)	Enterovirus	Enterovirus
	Enterovirus D68	Human Enterovirus D68	Enterovirus D68	Enterovirus D68
	Human herpesvirus 6 (HHV6)	Human herpesvirus 6 (HHV6)	Human herpesvirus 6 (HHV-6)	Human herpesvirus 6 (HHV-6)
	HHV3 (varicella zoster virus)		HHV3 (varicella zoster virus)	HHV3 (varicella zoster virus)
	HHV4 (Epstein-Barr virus)	HHV4 (Epstein-Barr virus)	HHV4 (Epstein-Barr virus)	HHV4 (Epstein-Barr virus)
	HHV5 (cytomegalovirus)	HHV5 (cytomegalovirus)	HHV5 (cytomegalovirus)	HHV5 (cytomegalovirus)
	Human Metapneumovirus (hMPV)	Human Metapneumovirus (hMPV)	Human Metapneumovirus (hMPV)	Human Metapneumovirus (hMPV)
	Influenza A (Pan)	Influenza A (Pan)	Influenza A (Pan)	Influenza A (Pan)
	Influenza A/H1-2009	Influenza A/H1-2009	Influenza A/H1-2009	Influenza A/H1-2009
	Influenza A/H3	Influenza A/H3	Influenza A/H3	Influenza A/H3
	Influenza B	Influenza B	Influenza B	Influenza B
	Measles virus		Measles virus	
	MERS-CoV		MERS-CoV	
	Mumps virus		Mumps virus	
	Parainfluenza virus 1	Human Parainfluenza virus 1	Parainfluenza virus 1	Parainfluenza virus 1
	Parainfluenza virus 2	Human Parainfluenza virus 2	Parainfluenza virus 2	Parainfluenza virus 2
	Parainfluenza virus 3	Human Parainfluenza virus 3	Parainfluenza virus 3	Parainfluenza virus 3
	Parainfluenza virus 4	Human Parainfluenza virus 4	Parainfluenza virus 4	Parainfluenza virus 4
	Parechovirus		Parechovirus	
	Respiratory Syncytial Virus A (RSV A)	Respiratory Syncytial Virus A (RSV A)	Respiratory Syncytial Virus A (RSV A)	Respiratory Syncytial Virus A (RSV A)
	Respiratory Syncytial Virus B (RSV B)	Respiratory Syncytial Virus B (RSV B)	Respiratory Syncytial Virus B (RSV B)	Respiratory Syncytial Virus B (RSV B)
	Rhinovirus	Rhinovirus	Rhinovirus	Rhinovirus
	SARS-CoV (SARS)		SARS-CoV (SARS)	
	SARS-CoV-2	SARS-CoV-2		

For Research Use Only. Not for use in diagnostic procedures