

Select your surface								
Thermo Scientific		Nunclon Delta	Nunclon Vita	Collagen I	Poly-D-Lysine	Sphera	UpCell	Untreated Surface
Hepatocytes		•	•	•		•	•	
Endothelial Cells		•		•		•	•	
Neuronal Cells	Neuronal Cells & neurosphere	•			•	•	•	
Epithelial Cells		•	•	•		•	•	
Tumor Cells		•	•	•	•	•	•	
Blood Cells	Macrophages, Dendritic Cells, Neutrophils	•				•	•	
	Lymphocytes	•						•
	Platelets	•		•			•	
Stem Cells	MSC	•	•	•		•	•	
	HSC	•			•			
	ESC and EB	•*	•**			•		

\*Requires extracellular matrix coating / \*\* Requires conditioned media from feeder cells and ROCK inhibitor