


CATALOG NO.	NOMINAL SIZE (mm)	MATERIAL	
		CLOSURE	FITTING
002158	83-B	PP	ACETAL

SIZE CODE.....	0022		0032	
THREAD FINISH (mm)...	83		83	
TREAD (T.P.I.).....	2		2	
NOM. WEIGHT (gms)....	133		154	
HOSE BARB SIZE.....	3/8"		3/8"	
NUMBER OF FITTINGS...	2		3	
	INCH	MM	INCH	MM
D (DIA).....	4.00 ±0.03	101.6 ±0.8	4.00 ±0.03	101.6 ±0.8
H (HEIGHT).....	4.25	107.9	4.25	107.9
F (DIA).....	1.25	31.8	1.25	31.8

SIZE CODE.....	0021		0031	
THREAD FINISH (mm)...	83		83	
TREAD (T.P.I.).....	2		2	
NOM. WEIGHT (gms)....	133		154	
HOSE BARB SIZE.....	1/4"		1/4"	
NUMBER OF FITTINGS...	2		3	
	INCH	MM	INCH	MM
D (DIA).....	4.00 ±0.03	101.6 ±0.8	4.00 ±0.03	101.6 ±0.8
H (HEIGHT).....	4.25	107.9	4.25	107.9
F (DIA).....	1.25	31.8	1.25	31.8

2. This is not a production drawing, but is for informational purposes only, and is subject to change without notice.

1. All dimensions are nominal and in inches and millimeters unless otherwise noted.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	JDD	APP'D PVC	PSS 00-2158-00XX	REV. 0505
CK'D	PVC	5/20/05		
DESCRIPTION:		83-B CLOSURE WITH QUICK DISCONNECTS		