


CATALOG NO.	NOMINAL CYLINDER CAPACITY (mL)							MATERIAL		
003664	10	25	50	100	250	500	1000	PP		
003665	10	25	50	100	250	500	1000	PMP		
SIZE CODE.....	-0010		-0025		-0050		-0100		-0250	
NOM. WT(gms)...	5		11		20		45		68	
	IN.	MM	IN.	MM	IN.	MM	IN.	MM	IN.	MM
A (BASE DIA)...	1.32	33.5	1.82	46.2	2.31	58.7	2.94	74.7	4.02	102.1
B (I.D.).....	0.53	13.5	0.75	19.1	0.91	23.1	1.11	28.2	1.52	38.6
C DIM.....	0.87	22.1	1.11	28.2	1.32	33.5	1.55	39.4	1.99	50.5
H (HEIGHT)....	3.85	97.8	5.53	140.5	7.00	177.8	9.00	228.6	11.6	294.6
SIZE CODE.....	-0500		-1000							
NOM. WT(gms)...	115		186							
	IN.	MM	IN.	MM						
A (BASE DIA)...	4.68	118.9	5.35	135.9						
B (I.D.).....	1.93	49.0	2.46	62.5						
C DIM.....	2.51	63.8	3.13	79.5						
H (HEIGHT)....	14.3	363.2	17.0	431.8						

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MC ²	APP'D	PVC	PSS 00-366X-XXX2	REV.
CHK'D	PVC		9/3/03		0903
DESCRIPTION:					
ECONOMY GRADUATED CYLINDERS					