


CATALOG NO.....:		7142-0015		7142-0020		7142-0030	
NOM VOL (L/GAL).....:		15/4		20/5.3		30/8	
MATERIALS	CONTAINER.....:	HDPE		HDPE		HDPE	
	COVER.....:	HDPE		HDPE		HDPE	
NOM. WEIGHT W/ COVER (gms).....:		750		1120		1420	
		INCHES	MM	INCHES	MM	INCHES	MM
A (OUTSIDE DIA OF COVER):		12.3	312	12.3	312	15.8	401
B (OUTSIDE DIA AT BOTTOM):		11.0	279	10.0	254	14.0	356
C (HT. W/ COVER).....:		9.9	252	15.5	394	12.1	307
D (HT. W/O COVER).....:		9.8	249	15.4	391	12.0	305
W.T.(NOM WALL THK).....:		0.1	2.5	0.1	2.5	0.1	2.5

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher

SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MC2	APP'D	PVC	PSS 00-7142-00XX	REV. 0402
CK'D	PVC		4/9/02		
DESCRIPTION: LARGE ROUND CONTAINERS WITH COVERS					