


CATALOG NUMBER	CAP TYPE	PACKAGING	DOUBLE BAG
141744	VENTED, EASY ON/OFF	20 BTL'S PER BAG / 20 BTL'S PER BOX	YES
142720	EASY ON/OFF	20 BTL'S PER BAG / 20 BTL'S PER BOX	NO
142744	EASY ON/OFF	20 BTL'S PER BAG / 20 BTL'S PER BOX	YES

GROWTH SURFACE AREA IN cm ²:	1450	
BOTTLE MATERIAL.....:	CRYSTAL POLYSTYRENE	
CLOSURE MATERIAL.....:	HDPE	
NOM. CAP.....:	2000mL	
CLOSURE SIZE (mm).....:	48	
	INCH	MM
A (CLOSURE O.D.).....:	2.07	52.5
B (HT. TO SHOULDER).....:	9.05	229.8
C (HT. W/ CLOSURE).....:	10.9	276.9
D (DIA.).....:	4.62	117.3
W.T. (NOM. WALL THK.):	.023	.58

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	TL	APP'D	BR	PSS 14XXXX	REV. 0711
CK'D	PB	8/1/11			

DESCRIPTION: Nunc TufRol Pleated Surface Roller Bottles