


CATALOG NUMBER	181702	
CAP TYPE	VENTED, EASY ON/OFF	
PACKAGING	2 BOTTLES PER BAG/20 BOTTLES PER BOX	
DOUBLE BAG	NO	
GROWTH SURFACE AREA	850cm <sup>2</sup>	
BOTTLE MATERIAL	CRYSTAL POLYSTYRENE	
CLOSURE MATERIAL	HDPE	
NOMINAL CAPACITY	2000mL	
CLOSURE SIZE (mm)	48	
	INCH	MM
A (CLOSURE DIAMETER)	2.07	52.6
B (HEIGHT TO SHOULDER)	9.05	229.9
C (HEIGHT W/CLOSURE)	10.9	276.9
D (DIAMETER)	4.62	117.3
E (WALL THICKNESS)	.023	.58



SPECS ON BACK

NOT TO SCALE

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

**CONFIDENTIAL**

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE P. PELYCHATY	APP'D 5/27/16	PSS 181702	REV.
CK'D JC	PWC		0316
DESCRIPTION: ROLLER BOTTLE, 850 cm <sup>2</sup> , VENTED			