


CATALOG NO.	NOMINAL CAPACITY (GAL)	MATERIALS	
		JERRICAN	CLOSURE
002240	2.5	HDPE	WHITE PP

SIZE CODE.....:	-0025	
SIZE (NOM VOL L/GAL):	10/2.5	
NOM. BRIM CAP. (L)...:	12	
NECK FINISH (mm).....:	53B	
THREAD (T.P.I.).....:	4	
NOM. WEIGHT W/ CLOSURE(gms).....:	1120	
	INCH	MM
A (NECK I.D.).....:	1.48 ±0.04	37.6 ±1.0
B (HT. W/ CLOSURE)...:	14.8 ±0.2	376 ±5
C (HT. W/O CLOSURE)..:	14.5 ±0.2	368 ±5
D (WIDTH).....:	9.70 ±0.15	246.4 ±3.8
E (DEPTH).....:	7.82 ±0.13	198.6 ±3.3
W.T.(NOM WALL THK)...:	0.20	5.1



SPECS ON BACK

NOT TO SCALE

- NOTES:
1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
 2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher

SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE TL	APP'D GRP	PSS 00-2240-0025-S	REV. 0413
CK'D KOW	5/15/13		
DESCRIPTION: JERRICANS			