


CATALOG NO.	MATERIALS	
	BOTTLE	CLOSURE
002104-0048	HDPE	PP

SIZE (NOM VOL nL/oz)	1500/ 48	
NOM. BRIM CAP. (mL):	1560	
NECK FINISH (mm):...	63-415	
THREAD (T.P.I.):....	6	
NOM. WEIGHT W/CLOSURE (gms):....	145	
	INCH	MM
A (BTLE. NECK I.D.):	2.06	52.3
B (HT. W/CLOSURE):...	11.16	283.5
C (HT. W/O CLOSURE):	10.94	277.9
D (DIA):.....	3.60	91.4
W.T.(NOM. WALL THK):	0.06	1.5



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	PB	APP'D	PVC	PSS 002104-0048	REV 0602
CHK'D	PVC		6/10/02		
DESCRIPTION:					
1.5-LITER WIDE MOUTH BOTTLE					