


CATALOG NO.	NOMINAL CAPACITY (oz)		MATERIALS		
			BOTTLE	CLOSURE	DRAW TUBE
002421	8	16	FLPE	FLPP, RED	PPCO

SIZE CODE.....:	-0250		-0500	
SIZE (NOM VOL mL/oz):	250/8		500/16	
NOM. BRIM CAP. (mL):	280		550	
NECK FINISH (mm)....:	24-415		53-415	
THREAD (T.P.I.).....:	8		6	
NOM. WEIGHT W/CLOSURE (gms).....:	30		74	
	INCH	MM	INCH	MM
A (NECK I.D.).....:	0.69 ±0.02	17.5 ±0.5	1.72 ±0.03	43.7 ±0.8
C (HT. W/O CLOSURE):	5.13 ±0.08	130.3 ±2.0	6.46 ±0.09	164.1 ±2.3
D (DIA).....:	2.40 ±0.07	61.0 ±1.8	2.88 ±0.08	73.2 ±2.0
H (NOM OVERALL HT):	6.9	175	8.5	216
W.T.(NOM. WALL THK):	0.035	0.9	0.035	0.9

NOTES:

- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
- THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher		Therno Fisher Scientific	
SCIENTIFIC		75 Panorama Creek Drive	
		Rochester, NY 14625-2385, USA	
MADE	NS	APP'D	REV.
CK'D	PVC	PVC	0999
9/20/99		PSS 00-2421-OXX0-1	
DESCRIPTION:			
FLUORINATED SOLVENT WASH BOTTLES			

REV.
0999