

CATALOG NO.	NOMINAL CAPACITY (oz)		MATERIALS	
			BOTTLE	CLOSURE
002425	16	32	LDPE	PP*


* -0506 AND -1006
ARE HDPE

SIZE CODE	-0501	-0502	-0503	-0504	-0505	-0506
CLOSURE COLOR	RED	WHITE	BLUE	YELLOW	NATURAL	WHITE
CHEMICAL	ACETONE	ETHANOL	METHANOL	ISOPROPANOL	DISTILLED WATER	SODIUM HYPOCHLORITE (BLEACH)
SIZE (NOM VOL nL/oz)	500/16					
NOM BRIM CAP (nL)	535					
NECK FINISH (mm)	28-415					
THREAD (T.P.I.)	6					
NOM.WT W/CLOSURE (gms)	60					
A (NECK I.D.)	0.84 ±.02					
H (NOM OVERALL HT)	8.3					
C (HT. W/O CLOSURE)	6.59 ±.09					
D (DIA)	2.86 ±.08					

SIZE CODE	-1002	-1003	-1004	-1005	-1006
CLOSURE COLOR	WHITE	BLUE	YELLOW	NATURAL	WHITE
CHEMICAL	ETHANOL	METHANOL	ISOPROPANOL	DISTILLED WATER	SODIUM HYPOCHLORITE (BLEACH)
SIZE (NOM VOL nL/oz)	1000/32				
NOM BRIM CAP (nL)	1060				
NECK FINISH (mm)	38-430				
THREAD (T.P.I.)	6				
NOM.WT W/CLOSURE (gms)	108				
A (NECK I.D.)	1.08 ±.03				
H (NOM OVERALL HT)	10.3				
C (HT. W/O CLOSURE)	8.35 ±.10				
D (DIA)	3.60 ±.06				

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	NS	APP'D PVC	PSS 00-2425-XXXX-1	REV.	0799
CK'D	PVC	7/15/99			
DESCRIPTION: RIGHT-TO-KNOW SAFETY WASH BOTTLES					