


CATALOG NO.	MATERIAL
3126-0175	WHITE PC

	INCH	MM
A (HEIGHT).....:	1.92	48.8
B DIA.....:	2.410 ±.020	61.21 ±.51
C (INSIDE DEPTH)....:	1.80	45.7

NOTE: FITS WITH 175mL CONICAL-BOTTOM CENTRIFUGE BOTTLES  
- (CATALOG NUMBERS: 3143-0175, 3144-0175, 3145-0175)

NOTES:

- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
- THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



NOT TO SCALE

SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	JGJ	APP'D	PVC	PSS 00-3126-0175	REV.
CHK'D	PVC		4/6/95		0495
DESCRIPTION: CONICAL-BOTTOM CENTRIFUGE BOTTLE ADAPTER					