


CATALOG NO.	NOMINAL CAPACITY, (mL)	MATERIALS	
		BOTTLE	CLOSURE
3127	250	FEP	ETFE

SIZE CODE.....:	-0250	
NOM. CAP. (mL).....:	250	
CLOSURE SIZE (mm)...:	38	
	INCH	MM
A (NECK I.D.).....:	0.96	24.4
B (HT. W/CLOSURE)...:	4.98 ±.06	126.5 ±1.5
C (DIA).....:	2.360 ±.040	59.94 ±1.02



NOT TO SCALE

SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	JGJ	APP'D	PVC	PSS 00-3127-0250	REV. 0495
CK'D	PVC		4/6/95		

DESCRIPTION:

CENTRIFUGE BOTTLE

NOTES:

- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
- THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.