


CATALOG NO. ....:	3640-0010	3640-0025	3640-0050	3640-0100				
NOMINAL RESERVOIR CAPACITY (ml)....:	500	500	1,000	1,000				
WEIGHT (gms).....:	110	125	210	230				
CALIBRATION INCREMENTS (ml)....:	.01	.01	.01	.02				
	INCH	MM	INCH	MM	INCH	MM	INCH	MM
A (BODY I.D.).....:	.37	9.4	.37	9.4	.42	10.7	.59	15.0
B (BODY O.D.).....:	.50	12.7	.50	12.7	.55	14.0	.73	18.5
H (HT W/O COVER)....:	12.6	320	20.8	528	28.3	719	28.4	721
BOTTLE CAP. (ml)....:	500	500	1,000	1,000				
TUBING MATERIAL.....:	PVC	PVC	PVC	PVC				
TUBING O.D (IN.)....:	3/16	3/16	5/16	5/16				

COMPONENT	MATERIAL
BURET COVER	LOPE
BURET BODY	ACRYLIC
STOPCOCK ASSEMBLY	—
STOPCOCK & TIP	PMP
PLUG	TFE
O-RINGS	VITON
HANDLE	PP
TUBING	PVC
CLOSURE	PP
BOTTLE	LOPE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Ponorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	MC	APP'D	PVC	PSS 00-3640-0XXX-1	REV.	1100
CK'D	PVC		11/16/00			

DESCRIPTION:

AUTOMATIC SELF-ZEROING BURETS