



CATALOG NO.....:	5115-0032	
NO. OF TUBES HELD.....:	32	
SIZE OF TUBES HELD.....:	.50mL, 1.50mL	
NOM. WEIGHT (gns).....:	1760	
	INCH	MM
A (OVERALL LENGTH).....:	9.5	241
B (OVERALL HEIGHT).....:	6.0	152
DIM W.....:	5.04	128.0
DIM H.....:	4.4	112

ACCESSORIES: (16) .50mL MICROTUBE INSERTS PROVIDED

COMPONENT	MATERIAL
BASE	BLUE PC
COVER	WHITE PC
INSERTS	GRAY PP
HANDLE	WHITE POWDER-COATED STAINLESS STEEL

2. This is not a production drawing, but is for informational purposes only, and is subject to change without notice.

1. All dimensions are nominal and in inches and millimeters unless otherwise noted.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MC	APP'D	PVC	PSS 00-5115-0032	REV. 0901
CK'D	PVC		9/13/01		
DESCRIPTION:					
LABTOP COOLER					