



CATALOG NO.....:	5116-0032	
NO. OF TUBES HELD.....:	32	
SIZE OF TUBES HELD.....:	0.50mL - 1.50mL	
NOM. WEIGHT (gms)....:	1760	
	INCH	MM
A (OVERALL LENGTH)....:	9.5	241
B (OVERALL HEIGHT).....:	6.0	152
DIM W.....:	5.04	128.0
DIM H.....:	4.4	112

COMPONENT	MATERIALS
BASE	GREEN PC
COVER	WHITE PC
INSERTS	GRAY PP
HANDLE	WHITE POWDER-COATED STAINLESS STEEL



SPECS ON BACK

NOT TO SCALE

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	PB	APP'D	GRP	PSS 00-5116-0032	REV. 0113
CK'D	KOW	1/18/13			
DESCRIPTION:				LABTOP COOLER	