


CATALOG NO.....	5245-0040	5245-0050		
SIZE CODE.....	D	E		
MAX. PIPET LENGTH (MM/IN.):	406/16	610/24		
MATERIAL	BASE.....	BLUE HDPE		
	BODY.....	HDPE		
	TUBING *.....	PVC		
	INCHES	MM	INCHES	MM
A (BASE DIA.).....	9.0	229	9.0	229
B (OVERALL WIDTH).....	10.7	272	10.5	267
D (NOM. OVERALL HEIGHT)....	30	762	38	965
H (NOM. HEIGHT).....	27	686	35	889
NOM. WT. (gms).....	1800		1700	

\* INCLUDES (2) PIECES OF TUBING AS FOLLOWS:  
(1) 1" I.D. x 1 1/4" O.D. x 18" LONG  
(1) 1/2" I.D. x 3/4" O.D. x 48" LONG

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	MCC	APP'D	PVC	PSS 00-5245-00X0	REV.	0995
CK'D	PVC		9/14/95			

DESCRIPTION:  
PIPET WASHER-RINSER