


CATALOG NUMBER.....:		0057051010		0057052020	
NOMINAL CAP., (ML)...		500		750	
MATERIAL	BOX.....:	PMP		PMP	
	COVER.....:	PPCO		PPCO	
	PLUG.....:	LOPE		LOPE	
NOM. WEIGHT (gms)....:		130		295	
		INCH	MM	INCH	MM
A (LENGTH COVER).....:		5.0	127	6.1	155
B (LENGTH BOX).....:		4.9	124	5.9	150
C (LENGTH HANDLES)...		5.4	137	6.9	175
D (OAH HEIGHT).....:		1.3	33	2.6	66
E (HEIGHT BOX).....:		1.2	30	2.4	61
F (HEIGHT TO PORT)...		0.2	5	0.4	10

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	PB	APP'D	PVC	PSS 005705XXXX	REV. 0503
CHK'D	PVC	5/9/03			
DESCRIPTION:					
STAINING BOXES					