

CATALOG NUMBER	SIZES (FITS NOMINAL TUBING I.D.)						MATERIAL
	1/8	3/16	1/4	5/16	3/8	1/2	
006152							PP


SIZE CODE.....:	-0125		-0187		-0250	
NOM. TUBING SIZE I.D. (INCHES)....:	1/8		3/16		1/4	
NOM. WT. (gms)...	0.2		0.6		1.8	
	INCHES	MM	INCHES	MM	INCHES	MM
A (I.D.).....:	0.09 ±0.01	2.3 ±0.3	0.12 ±.01	3.0 ±0.3	0.15 ±.01	3.8 ±0.3
B	0.48 ±0.02	12.2 ±0.5	0.66 ±.02	16.8 ±0.5	1.03 ±.02	26.2 ±0.5
C	0.48 ±0.02	12.2 ±0.5	0.66 ±.02	16.8 ±0.5	1.03 ±.02	26.2 ±0.5
D	0.97 ±0.03	24.6 ±0.8	1.34 ±.03	34.0 ±0.8	2.04 ±.03	51.8 ±0.8
E (O.D.).....:	0.16 ±0.01	4.1 ±0.3	0.20 ±.01	5.1 ±0.3	0.27 ±.02	6.9 ±0.5
F.....	0.62 ±0.03	15.7 ±0.8	0.86 ±.03	21.8 ±0.8	1.29 ±.03	32.8 ±0.8

SIZE CODE.....:	-0312		-0375		-0500	
NOM. TUBING SIZE I.D. (INCHES)....:	5/16		3/8		1/2	
NOM. WT. (gms)...	2.4		3.3		5.8	
	INCHES	MM	INCHES	MM	INCHES	MM
A (I.D.).....:	0.23 ±.02	5.8 ±0.5	0.30 ±.02	7.6 ±0.5	0.41 ±.02	10.4 ±0.5
B	1.11 ±.02	28.2 ±0.5	1.19 ±.02	30.2 ±0.5	1.43 ±.03	36.3 ±0.8
C	1.11 ±.02	28.2 ±0.5	1.19 ±.02	30.2 ±0.5	1.43 ±.03	36.3 ±0.8
D	2.25 ±.03	57.2 ±0.8	2.44 ±.03	62.0 ±0.8	2.93 ±.04	74.4 ±1.0
E (O.D.).....:	0.34 ±.02	8.6 ±0.5	0.41 ±.02	10.4 ±0.5	0.53 ±.02	13.5 ±0.5
F.....	1.45 ±.03	36.8 ±0.8	1.61 ±.03	40.9 ±0.8	1.96 ±.03	49.8 ±0.8

THIS DRAWING IS NOT TO BE REPRODUCED
-CONFIDENTIAL-

NOTES:

- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
- THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION
SHEET CONTAINS CONFIDENTIAL
AND/OR PROPRIETARY INFORMATION
AND IS TO BE USED ONLY FOR THE
EXPRESSED PURPOSES AGREED UPON
BY AN AUTHORIZED OFFICER OF
THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE CK'D	MCC PVC	APP'D PVC 6/21/95	PSS 00-6152-XXXX-1	REV. 0695
--------------	------------	-------------------------	--------------------	--------------

DESCRIPTION:
Y-TYPE CONNECTORS